



Technical Data

Listeria Selective Supplement (PALCAM)

FD061

An antimicrobial supplement recommended for the selective isolation and identification of *Listeria monocytogenes*.

Composition

Per vial sufficient for 500 ml medium

*Ingredients

Polymyxin B sulphate
Ceftazidime
Acriflavine hydrochloride

Concentration

5000IU
10mg
2.500mg

Directions:

Rehydrate the contents of one vial aseptically with 5 ml sterile distilled water and aseptically add to 500 ml sterile, molten, cooled (45-50°C) Listeria Identification Agar Base (PALCAM) [M1064](#) / Listeria Identification HiVeg™ Agar Base (PALCAM) [MV1064](#) / Listeria Identification Broth Base (PALCAM) [M1090](#) / Listeria Identification HiVeg™ Broth Base (PALCAM) [MV1090](#). Mix well and dispense as desired.

Storage and Shelf Life

Store at 2-8°C. Use before expiry date on the label.

Revision : 1 / 2012

* Not For Medicinal Use

Disclaimer :

User must ensure suitability of the product(s) in their application prior to use. Products conform solely to the information contained in this and other related HiMedia™ publications. The information contained in this publication is based on our research and development work and is to the best of our knowledge true and accurate. HiMedia™ Laboratories Pvt Ltd reserves the right to make changes to specifications and information related to the products at any time. Products are not intended for human or animal diagnostic or therapeutic use but for laboratory, research or further manufacturing use only, unless otherwise specified. Statements contained herein should not be considered as a warranty of any kind, expressed or implied, and no liability is accepted for infringement of any patents.