

DATE	28.03.2025		
DOC.NO	FRM-KLT-BYM-069		
PAGE NO	1		
REV.NO			
REV DATE			

	PRODUCT:	Sterile disposabl	TECHNICAL DATA SHEET surgical set for angiography		
	TRODUCT.	Sterne disposabl	Surgical Section angiography		
			To Sterile disposable surgical set for angiography Angiography field - 240x330cm (+/-5cm), integral, 100% waterproof on the entire surface. a) Critical area - triple-layered reinfor high degree of absorption, minimum dimensions: 240x80cm, with 4 holes: 2 holes of 8x10cm (+/-3cm) and 2 holes of 5x7cm (+/-1) an impregnated adhesive border around the holes (double-sided glue will not be accepted). b) Absorbent semi-critical area in 100x220cm. c) Transparent side edges, dim. 70x220cm. All areas of the field will ensure 100% tightness, will be heat-sealed, will in a dry and/or liquid environment. Super absorbent field - 35x45cm (+/-5cm), extra-absorbent, multilayered, minimum 6 layers, with adhesive tape, absorption cal	/-3cm), with inimum 240x330 cm	1
			2 minimum 500ml liquids. Field with adhesive board - 50x50cm (+/-5cm). Triple layered on the entire surface: cellulose comfort layer for the patient, wate 3 layer, absorbent layer, with impregnated adhesive board (double-sided glue will not be accepted), the adhesive board applied of	erproof 35x45 cm	1
			the edge of the field to avoid the accumulation of liquids. 4 Instrument table field - 150x190cm (+/-5cm), absorbent area - minimum 190x75cm, 2 layers, waterproof.	150x190 cm	1
			5 Hand towels, made of cellulose, minimum 30x40cm. Surgical gown XL - a) Made of non-woven material, polypropylene fibers, minimum 3 types of fibers - SMS; b) Repellent, SMS m 6 thickness minimum 35 g/m²; c) Laser-sewn sleeve to ensure impermeability to microbes and liquids, needle-sewn sleeve will be d) Sleeve length 60cm (±3cm), gown length (from shoulders to hem) - 130cm (±/-5cm), gown width in the axillary area 75cm (±3	e excluded; Length: 130 cm	2
			Sterile, powder-free gloves (sizes: N7.5 - 2 pairs; N8.0 - 1 pair) - made of deproteinized rubber, hypo-allergenic/with low potent allergies, resistant to surgical interventions with the use of metal guides; non-pyrogenic, breaking strength – not less than 18N, 7 the fingers – not less than 0.27mm; coated on the inside with biogel for double protection and breaking strength, AQL – not mo (AQL – standard for the protection of medical personnel during long-term interventions (penetration level for microorganisms ochemicals – according to the EN 374-2:2014 standard)). 8 Circular screen cover - 140cm made of polyethylene, transparent, with elastic.	ntial to cause thickness at ore than 0.65 and liquids, 140 cm	2 pai 1 pai
	Description of Pro	duct:	9 C-arm radiological tube cover, circular, 75cm made of polyethylene, transparent, with elastic. 10 120 ml polypropylene bowl, graduated, transparent, with thread and lid, minimum height 7cm, Inscribed with the label "biologi	75 cm gical risk". 7 cm	1
			11 Bowl minimum 100 ml made of polypropylene, transparent.	100 mL	1
			12 Bowl minimum 250 ml made of polypropylene, colored.	250 mL	1
			13 Kidney bowl minimum 800 ml. 14 Forceps for field fixation – 15cm (+/-3cm).	800 mL 15 cm	1
			15 Clamp – 10cm (+/-3cm). Tray 2500 ml, round, made of polypropylene, colored,with "wings", minimum height 9cm (+/-1cm) – to ensure the firm position	10 cm	1
			metal guide and prevent its slipping.	2500 ML	1
			17 Angio syringe 12 ml with rotator. 3-way OFF ramp, "wing" valves, resistant to 35 bar integrated with: a) polyurethane contrast system minimum 1800mm; b) poly	yurethane 12 mL	1
			infusion system minimum 1900mm; c) polyurethane tube, minimum 1500mm.	STD E11	1
			19 E11 scalpel. 20 Angiographic guide – minimum length 175cm, J3, 0.034-0.036".	175 cm	1
			21 Gauze pad, absorbent 10x10cm, minimum 8 layers. 22 Syringe needle 18G.	10x10 cm 18G	30 1
			23 Syringe needle 21G.	21 G	2
			24 Syringe needle 22G 25 Luer Lock syringe 20 ml.	22G	1
			26 Luer Lock syringe 20 ml. 26 Luer Lock syringe 10 ml.	20 mL 10 mL	1
			27 Luer Lock syringe 5 ml.		- 4
				5 mL	1
			28 Luer Lock syringe 2 ml. 29 Tubing bowl – 30x30cm (+/-3cm).	2 mL 30x30 cm	1 1
			28 Luer Lock syringe 2 ml.	2 mL	1
	Raw Materia	als:	28 Luer Lock syringe 2 ml. 29 Tubing bowl – 30x30cm (+/-3cm).	2 mL 30x30 cm	1
	Product Colo	ur:	28 Luer Lock syringe 2 ml. 29 Tubing bowl – 30x30cm (+/-3cm). 30 Sterile field for tubing - 100 x 150cm (+/-3cm).	2 mL 30x30 cm	1
		ur:	28 Luer Lock syringe 2 ml. 29 Tubing bowl – 30x30cm (+/-3cm). 30 Sterile field for tubing - 100 x 150cm (+/-3cm). Triflex / SMS + Reinforcement / PE + NW / Biflex / SMS / Cellulose / PE Medical Blue SP-05006-161	2 mL 30x30 cm 100x150 cm	1 1 1
	Product Colo Reference Co Weight in Gra	ur: de:	28 Luer Lock syringe 2 ml. 29 Tubing bowl – 30x30cm (+/-3cm). 30 Sterile field for tubing - 100 x 150cm (+/-3cm). Triflex / SMS + Reinforcement / PE + NW / Biflex / SMS / Cellulose / PE Medical Blue SP-05006-161 77 gsm (Triflex) / 35 gsm SMS + 110 gsm (Reinforcement) / 28 gsm(SS)+65 mic.(PE) / 56 gsm (SMS) / 60 gsm (Cellulose) / 60 mic (PE)	2 mL 30x30 cm 100x150 cm	1 1 1
	Product Colo Reference Co Weight in Gra Package:	ur: de: ms:	28 Luer Lock syringe 2 ml. 29 Tubing bowl - 30x30cm (+/-3cm). 30 Sterile field for tubing - 100 x 150cm (+/-3cm). Triflex / SMS + Reinforcement / PE + NW / Biflex / SMS / Cellulose / PE Medical Blue SP-05006-161 77 gsm (Triflex) / 35 gsm SMS + 110 gsm (Reinforcement) / 28 gsm(SS)+65 mic.(PE) / 56 gsm (SMS) / 60 gsm (Cellulose) / 60 mic (PE) Flat Pouch	2 mL 30x30 cm 100x150 cm	1 1 1
	Product Colo Reference Co Weight in Gra Package: Product: Me	ur: de: ms: edical Device Direc	28 Luer Lock syringe 2 ml. 29 Tubing bowl - 30x30cm (+/-3cm). 30 Sterile field for tubing - 100 x 150cm (+/-3cm). Triflex / SMS + Reinforcement / PE + NW / Biflex / SMS / Cellulose / PE Medical Blue SP-05006-161 77 gsm (Triflex) / 35 gsm SMS + 110 gsm (Reinforcement) / 28 gsm(SS)+65 mic.(PE) / 56 gsm (SMS) / 60 gsm (Cellulose) / 60 mic (PE) Flat Pouch ive MDR 2017/745 ECC ANNEX XI A ALL VERSIONS. / Products and materials used do not co	2 mL 30x30 cm 100x150 cm	1 1 1
	Product Colo Reference Co Weight in Gra Package: Product: Me	ur: de: ms: edical Device Directerials	28 Luer Lock syringe 2 ml. 29 Tubing bowl – 30x30cm (+/-3cm). 30 Sterile field for tubing - 100 x 150cm (+/-3cm). Triflex / SMS + Reinforcement / PE + NW / Biflex / SMS / Cellulose / PE Medical Blue SP-05006-161 77 gsm (Triflex) / 35 gsm SMS + 110 gsm (Reinforcement) / 28 gsm(SS)+65 mic.(PE) / 56 gsn (SMS) / 60 gsm (Cellulose) / 60 mic (PE) Flat Pouch ive MDR 2017/745 ECC ANNEX XI A ALL VERSIONS. / Products and materials used do not co Product Shelf Life - 5 Years PROPERTIES	2 mL 30x30 cm 100x150 cm	1 1 1
1	Product Colo Reference Co Weight in Gra Package: Product: Me	ur: de: ms: edical Device Directerials	28 Luer Lock syringe 2 ml. 29 Tubing bowl - 30x30cm (+/-3cm). 30 Sterile field for tubing - 100 x 150cm (+/-3cm). Triflex / SMS + Reinforcement / PE + NW / Biflex / SMS / Cellulose / PE Medical Blue SP-05006-161 77 gsm (Triflex) / 35 gsm SMS + 110 gsm (Reinforcement) / 28 gsm(SS)+65 mic.(PE) / 56 gsm (SMS) / 60 gsm (Cellulose) / 60 mic (PE) Flat Pouch ive MDR 2017/745 ECC ANNEX XI A ALL VERSIONS. / Products and materials used do not co	2 mL 30x30 cm 100x150 cm	1 1 1
1	Product Colo Reference Co Weight in Gra Package: Product: Mo Product Mat Unit / Size	ur: de: ms: edical Device Directerials	28 Luer Lock syringe 2 ml. 29 Tubing bowl – 30x30cm (+/-3cm). 30 Sterile field for tubing - 100 x 150cm (+/-3cm). Triflex / SMS + Reinforcement / PE + NW / Biflex / SMS / Cellulose / PE Medical Blue SP-05006-161 77 gsm (Triflex) / 35 gsm SMS + 110 gsm (Reinforcement) / 28 gsm(SS)+65 mic.(PE) / 56 gsn (SMS) / 60 gsm (Cellulose) / 60 mic (PE) Flat Pouch ive MDR 2017/745 ECC ANNEX XI A ALL VERSIONS. / Products and materials used do not co Product Shelf Life - 5 Years PROPERTIES	2 mL 30x30 cm 100x150 cm	1 1 1
	Product Colo Reference Co Weight in Gra Package: Product: Me Product Mat Unit / Size Triflex	ur: de: ms: edical Device Directerials 240x370 cm	28 Luer Lock syringe 2 ml. 29 Tubing bowl = 30x30cm (+/-3cm). 30 Sterile field for tubing - 100 x 150cm (+/-3cm). Triflex / SMS + Reinforcement / PE + NW / Biflex / SMS / Cellulose / PE Medical Blue 6SP-05006-161 77 gsm (Triflex) / 35 gsm SMS + 110 gsm (Reinforcement) / 28 gsm(SS)+65 mic.(PE) / 56 gsm (SMS) / 60 gsm (Cellulose) / 60 mic (PE) Flat Pouch ive MDR 2017/745 ECC ANNEX XI A ALL VERSIONS. / Products and materials used do not co Product Shelf Life - 5 Years PROPERTIES 1 2 3 4	2 mL 30x30 cm 100x150 cm	1 1 1
1	Product Colo Reference Co Weight in Gra Package: Product: Mo Product Mat Unit / Size Triflex Reinforcement	ur: de: ms: edical Device Directerials 240x370 cm 35x45 cm	28 Luer Lock syringe 2 ml. 29 Tubing bowl = 30x30cm (+/-3cm). 30 Sterile field for tubing - 100 x 150cm (+/-3cm). Triflex / SMS + Reinforcement / PE + NW / Biflex / SMS / Cellulose / PE Medical Blue 6SP-05006-161 77 gsm (Triflex) / 35 gsm SMS + 110 gsm (Reinforcement) / 28 gsm(SS)+65 mic.(PE) / 56 gsm (SMS) / 60 gsm (Cellulose) / 60 mic (PE) Flat Pouch ive MDR 2017/745 ECC ANNEX XI A ALL VERSIONS. / Products and materials used do not co Product Shelf Life - 5 Years PROPERTIES 1 2 3 4	2 mL 30x30 cm 100x150 cm	1 1 1
1	Product Colo Reference Co Weight in Gra Package: Product: Mo Product Mat Unit / Size Triflex Reinforcement Triflex	ur: de: ms: edical Device Directerials 240x370 cm 35x45 cm 50x50 cm	28 Luer Lock syringe 2 ml. 29 Tubing bowl = 30x30cm (+/-3cm). 30 Sterile field for tubing - 100 x 150cm (+/-3cm). Triflex / SMS + Reinforcement / PE + NW / Biflex / SMS / Cellulose / PE Medical Blue SP-05006-161 77 gsm (Triflex) / 35 gsm SMS + 110 gsm (Reinforcement) / 28 gsm(SS)+65 mic.(PE) / 56 gsm (SMS) / 60 gsm (Cellulose) / 60 mic (PE) Flat Pouch ive MDR 2017/745 ECC ANNEX XI A ALL VERSIONS. / Products and materials used do not co Product Shelf Life - 5 Years PROPERTIES 1 2 3 4	2 mL 30x30 cm 100x150 cm	1 1 1
1 1 1	Product Colo Reference Co Weight in Gra Package: Product: Mo Product Mat Unit / Size Triflex Reinforcement Triflex NW+PE	ur: de: ms: edical Device Directerials 240x370 cm 35x45 cm 50x50 cm 150x190 cm 30x40 cm Length: 60 cm Height: 130 cm Width: 75 cm	28 Luer Lock syringe 2 ml. 29 Tubing bowl = 30x30cm (+/-3cm). 30 Sterile field for tubing - 100 x 150cm (+/-3cm). Triflex / SMS + Reinforcement / PE + NW / Biflex / SMS / Cellulose / PE Medical Blue SP-05006-161 77 gsm (Triflex) / 35 gsm SMS + 110 gsm (Reinforcement) / 28 gsm(SS)+65 mic.(PE) / 56 gsm (SMS) / 60 gsm (Cellulose) / 60 mic (PE) Flat Pouch ive MDR 2017/745 ECC ANNEX XI A ALL VERSIONS. / Products and materials used do not co Product Shelf Life - 5 Years PROPERTIES 1 2 3 4	2 mL 30x30 cm 100x150 cm	1 1 1
1 1 1 4 3 2 pairs	Product Colo Reference Co Weight in Gra Package: Product: Mo Product Mat Unit / Size Triflex Reinforcement Triflex NW+PE Cellulose	ur: de: ms: edical Device Directerials 240x370 cm 35x45 cm 50x50 cm 150x190 cm 30x40 cm Length: 60 cm Height: 130 cm	28 Luer Lock syringe 2 ml.	2 mL 30x30 cm 100x150 cm	1 1 1
1 1 1 4 3	Product Colo Reference Co Weight in Gra Package: Product: Mo Product Mat Unit / Size Triflex Reinforcement Triflex NW+PE Cellulose SMS	ur: de: ms: edical Device Directerials 240x370 cm 35x45 cm 50x50 cm 150x190 cm 30x40 cm Length: 60 cm Height: 130 cm Width: 75 cm N7.5	28 Luer Lock syringe 2 ml.	2 mL 30x30 cm 100x150 cm	1 1 1
1 1 4 3 2 pairs 1 pair	Product Colo Reference Co Weight in Gra Package: Product: Mo Product Mat Unit / Size Triflex Reinforcement Triflex NW+PE Cellulose SMS ***	ur: de: ms: edical Device Directerials 240x370 cm 35x45 cm 50x50 cm 150x190 cm 30x40 cm Length: 60 cm Height: 130 cm Width: 75 cm N7.5 N8.0	28 Luer Lock syringe 2 ml.	2 mL 30x30 cm 100x150 cm	1 1 1
1 1 4 3 2 pairs 1 pair 2	Product Colo Reference Co Weight in Gra Package: Product: Mi Product Mat Unit / Size Triflex Reinforcement Triflex NW+PE Cellulose SMS ***	ur: de: ms: edical Device Directerials 240x370 cm 35x45 cm 50x50 cm 150x190 cm 40x40 cm Length: 60 cm Height: 130 cm Width: 75 cm N7.5 N8.0 140 cm	28 Luer Lock syringe 2 ml.	2 mL 30x30 cm 100x150 cm	1 1 1
1 1 4 3 2 pairs 1 pair 2	Product Colo Reference Co Weight in Gra Package: Product: Mi Product Mat Unit / Size Triflex Reinforcement Triflex NW+PE Cellulose SMS *** PE ***	ur: de: ms: edical Device Directerials 240x370 cm 35x45 cm 50x50 cm 150x190 cm 40x40 cm Length: 60 cm Height: 130 cm Width: 75 cm N7.5 N8.0 140 cm 75 cm	28 Luer Lock syringe 2 ml.	2 mL 30x30 cm 100x150 cm	1 1 1

1	***	800 mL		0			
1	***	15 cm		0			
1	***	10 cm					
1	***	2500 mL	17	18	19	20	
1	***	10 ml	D	Contract of the contract of th			
1	***	STD	and the state of t			()	
1	***	E11					
1	***	175 cm	21	22	23	24	
30	***	10x10 cm		11		I i	
1	***	18G					
2	***	21 G		11.0	A A	11.0	
1	***	22G	25	26	27	28	
1	***	20 mL	41111			0	
1	***	10 mL			10 to	8 1 1 1 1 1 1 1 1 1 1 1 1 1 1 1 1 1 1 1	
1	***	5 mL					
1	***	2 mL	29	30			
1	***	30x30 cm					
1	Biflex	100x150 cm					
	olerances: leasurement:	+/- 2% cm cm	Package Information				
IVI	easurement.	Cill	The products in the sterilization bag are double-packaged to reduce all risks during transportation. Double packaged products are put into Baymed's standard sized carton; dimensions are as follows: Height = 44 cm; Length = 40 cm ve Width = 60 cm.				
Preparation Date			QUALITY CONTROL APPROVAL				