i-gel[®] supraglottic airway

For use in anaesthesia and resuscitation

i-gel is a truly unique, single use, latex and PVC free second generation supraglottic airway management device that represents the culmination of years of extensive research and development. Made from a medical grade thermoplastic elastomer, i-gel has been designed to create a non-inflatable, anatomical seal of the pharyngeal, laryngeal and perilaryngeal structures whilst avoiding compression trauma.

The i-gel[®] supraglottic airway is supplied in a fully recyclable protective cradle or cage pack. This unique packaging protects the i-gel in transit and ensures that it maintains its anatomical shape. i-gel is available in seven sizes.

www.i-gel.com

Video available at www.intersurgical.com



15mm connector

Reliable connection to any standard catheter mount or connection


Proximal end of gastric channel

Clearly displayed product information

For quick easy reference. Includes confirmation of size and weight guidance

Position guide (adult sizes only)

Easy confirmation of optimum insertion depth



Integral bite block

Gastric channel

i-de

Reduces the possibility of airway channel occlusion

Buccal cavity stabiliser Aids insertion and eliminates the potential for rotation


Epiglottic rest

Reduces the possibility of epiglottic 'down folding' and airway obstruction

The non-inflatable cuff

Made from a unique soft gel-like material allowing ease of insertion and reduced trauma

Distal end of gastric channel





Make an enquiry

Intersurgical • Quality, innovation and choice

