

## Ra K


### Knitted vascular prosthesis with collagen

Crimped vascular prosthesis is made of polyester filament yarn using knitting technology and is impregnated with bovine collagen of type I.

#### INDICATIONS:

The knitted vascular prostheses are indicated for all types of reconstruction of arteries where application of a vascular graft of corresponding lumen is suitable.

These types of prostheses can be used in aorto-femoral, ilio-femoral regions, femoro-popliteal region above the knee joint, in the region of extra-cranial arteries and for extraanatomical bypasses.



- the prosthesis is used directly, without preparation before implantation
- the prosthesis is strong, elastic, easy to pierce, and it does not fray
- the wall is provided with a colour guiding line
- when the prostheses are bent they retain a lumen of circular shape
- the product is sterilized by irradiation
- expiration period is 3 years

## Sizes and catalogue codes:

### Linear vascular prostheses Ra 1v K

*Diameter [mm]	Catalogue code for the length 250mm	Catalogue code for the length 500mm
6	022 - 0061	022 - 0062
7	022 - 0071	022 - 0072
8	022 - 0081	022 - 0082
9	022 - 0091	022 - 0092
10	022 - 0101	022 - 0102
12	022 - 0121	022 - 0122
14	022 - 0141	022 - 0142
16	022 - 0161	022 - 0162
18	022 - 0181	022 - 0182
20	022 - 0201	022 - 0202
22	022 - 0221	022 - 0222
24	022 - 0241	022 - 0242
25	022 - 0251	022 - 0252
26	022 - 0261	022 - 0262
28	022 - 0281	022 - 0282
30	022 - 0301	022 - 0302
32	022 - 0321	022 - 0322

### Linear vascular prostheses Ra 1v K

*Diameter [mm]	*Length [mm]	Catalogue code
6 - 32	100	022 - 6...32 / 100
6 - 32	200	022 - 6...32 / 200
6 - 32	300	022 - 6...32 / 300
6 - 32	400	022 - 6...32 / 400
6 - 32	600	022 - 6...32 / 600

### Bifurcated vascular prostheses Ra bv K

*Diameter [mm]	Catalogue code for the length 400 mm
6 / 12	024 - 0612
7 / 14	024 - 0714
8 / 16	024 - 0816
9 / 18	024 - 0918
10 / 20	024 - 1020
11 / 22	024 - 1122

*\*Note: all sizes (diameter, length) can be adjusted and manufactured upon the Customer's individual requests.  
Please, consult your Distributor.*

## VUP Medical, a.s.

Šujanovo nám. 3 | 602 00 Brno

Czech Republic | [www.vup.cz](http://www.vup.cz)

Phone: +420 543 160 444

e-mail: [sales@vup.cz](mailto:sales@vup.cz)