TR Band



Radial Artery Compression Device

TR Band® is a compression device to assist patent haemostasis of the radial artery after transradial procedure. TR Band® assists in maintaining radial artery patency at the time of haemostasis, in order to prevent future radial artery occlusion (RAO).1

Product Characteristics

- Selective compression of the radial artery to allow blood return and preserve patency¹
 Transparent structure for visual control of the puncture site
- Velcro² straps make application fast and simple
- Two band lengths available for a precise fit and optimal patient comfort1
- 1 Rathore S et al. Catheter Cardiovasc Interv 2010;76:660-667
- 2 Velcro® is a registered trademark of Velcro Industries B.V.

General specifications

| Maximum air injection volume | Maximum air injection volume | 18 ml |
|------------------------------|------------------------------|-------|
|------------------------------|------------------------------|-------|

Item specifications

| Size | Code |
|-------------------|-----------|
| Large 29 cm | TRB29-LRG |
| Large 29 cm | XX*RF06L |
| Standard 24 cm | XX*RF06 |
| Standard 24 cm | TRB24-REG |

Please quote above item reference codes when placing an order