



Inara[®] Cu 375 / Cu 375 Sleek

- Effective period: 5 years
- Specially designed inverted 'U' shape frame for added comfort while inserting.
- Preloaded and ready for insertion
- Inara Cu 375 Sleek has a specially designed smaller size 'U' shape frame
- Inara Cu 375 Sleek is recommended for nulliparous women.
- Both models come with Esa - disposable uterine sound in pack.



Inara[®] Cu 250 / Cu 250 Sleek

- Effective period: 3 years
- Specially designed inverted 'U' shape frame for added comfort while inserting.
- Preloaded and ready for insertion
- Low copper content
- Inara Cu 250 Sleek has a specially designed smaller size 'U' shape frame
- Inara Cu 250 Sleek is recommended for nulliparous women.
- Both models come with Esa - disposable uterine sound in pack.

Post - Partum Range of IUDs

- Specially designed preloaded long inserter for easy fundal placement with no touch technique.
- Long nylon suture which is visible immediately after the insertion.
- To be used within 48 hours of normal delivery.



Effective period: 5 years

 Cu 375



Effective period: 10 years

 T Cu 380A



Effective Period: 5 years

Updated: December 21

✉ sales@pregna.com ☎ 91 9833001010 🌐 www.pregna.com
🏠 13 Suryoday Estate, 136 Tardeo Road, Mumbai - 400 034.


PREGNA
INTERNATIONAL LTD.

Pregna is a leading women's healthcare device manufacturer, having served over 100 million women in over 145 countries since 1991.






CryoPop[®] An avant garde cryotherapy device

- A Cryotherapy device for treatment of precancerous cervical lesions
- No cables required during operation
- Solid CO₂ based system consumes upto 75% lesser gas
- Solid state conduction makes it highly cooling efficient




Dr. Burke's ESM-UBT Kit

- Uterine Balloon Tamponade
- Effective in rapid treatment of Postpartum Bleeding
- Single use device reduces the risk of cross contamination
- Easy to assemble with minimal training required
- Majority of deaths are preventable



Accurette[®] Endometrial Biopsy Curette

- Pain free collection of endometrial tissue for pathology.
- Safe and simple procedure.
- Sufficient sample quantity, eliminate need of repetition.



esā[®] Disposable Uterine Sound

- Smooth finish with special rounded tip to avoid perforation of fundus.
- Wide handle with serrations to give better grip
- Single use device to avoid cross contamination.
- Anatomical design with high flexibility to suit the uterine cavity

Pregna's Range Of IUDs

All of our products are saliently designed with:

- High quality soft nylon suture to provide utmost comfort to the patient.
- Medical grade LDPE frames blended with Barium Sulphate for visibility under X-ray.
- Depth gauge marking on insertion tube to provide accuracy and convenience while insertion.
- Reversible contraception with 99% effectivity



eloira[®] 52 mg Levonorgestrel (LNG IUS)

- Effective period: 5 years
- Constant release of 20mg of LNG per day
- Flexible arms for perfect fundal placement
- Used in treatment of HMB (AUB & DUB) and fibroids less than (<3cm).
- Comes with Esa - disposable uterine sound



etherena[®] T Cu 380A

- Effective period: 10 years
- Specially designed Intrauterine Enabler (IUE) for aseptic loading
- Single handed operation
- Precise uterine depth adjuster for added accuracy
- Comes with Esa - disposable uterine sound.



PREGNA[®] T Cu 380A

- Effective period: 10 years
- Simple and safe insertion & removal
- Preloaded solid rod for easy insertion.