## RO Membrane w/o housing, replace CDRC050HG



RO Membrane w/o housing RAR0050HG, replace Millipore CDRC050HG

RephiLe Catalog No.: RAR0050HG Millipore Catalog No.: CDRC050HG

## **Product Description:**

This product is designed to remove most of the ions, organic carbons and particulates. This is the RO membrane inside the housing. Reuse the housing.

## Suitable for the following Millipore systems:

- Elix 35/ 70/100
- Elix Clinical 35/70
- RiOs 50/100/150/200