



Non-Invasive Ventilation Face mask

## **Bluestar PLUS**

for pulmonary ventilation's devices, modality CPAP and BILEVEL assisted and supervised.

- Innovative design that guarantees an optimal air-leak control and a high patient comfort level.
- · Lightweight.
- · Soft replaceable cushion minimizes skin breakdown.
- The kit has all the accessories for positioning.



## Kit contents: Transparent PVC mask with inflatable and replaceable cushion

• DEHP free.

with inflation valve.

- Thumbscrew adjustable forehead cushion with two grips.
- Ports for O<sub>2</sub> and CO<sub>2</sub>, pressure monitoring and a port for introducing gastric tubing and with the possibility of introducing a bronchoscope.
- Two connectors 22M and 22F.
- Pump for inflating and deflating the cushion.
- · Headgear with a rapid removal system.
- Latex free.
- The mask's connector can be connected to expiratory system (KM-903 not in the kit) provided with a hole for the CO<sub>2</sub> discharge and with a port for the enrichment of O<sub>2</sub>.
- Every kit is delivered in a box with a label indicating the contents.



Inflatable and replaceable cushions for

## **BluestarPLUS**

Code	Description	Size	For mask	Items per box
KM-3101	cushion replacement	extra large	KM-310	50
KM-3111	cushion replacement	large	KM-311	50
KM-3121	cushion replacement	medium	KM-312	50
KM-3131	cushion replacement	small	KM-313	50
KM-3141	cushion replacement	pediatric	KM-314	50





## KOO MEDICAL EQUIPMENT

Web: www.koomedical.com

Responsible for EMEA (Europa Middle East Africa) Markets - Koo Europe s.r.l. viale delle Industrie 5 - 20020 Arese (Milano) ITALY Tel: 0039 02 93589495 - Email: info@kooeurope.com Responsible for Asian Pacific Market - Koo Asia Ltd Unit 23,9 Floor Goodwill Ind. Building, No.36-44, Pak Tin Pa Street, Tsuen Wan, Hong Kong Tel: +852 2110 4780 Fax: +852 2110 4783 - Email: sales@koo-asia.com