

sale@greatmade.com

	AE001-1	Adult ear clip sensor , DB9 Male,L=1M			
	AE001-3D	Adult ear clip sensor , DB9 Male,L=3M			
	AT001-1	Animal tongue clip sensor , DB9 Male,L=1M			
	AT001-3D	Animal tongue clip sensor , DB9 Male,L=3M			
	NW001-3D- Y	Multi-site Y cable with separate wrap , L=3m ,DB9 , <b>Wrap Detachable</b>			
	NW 001-3D- F	<b>Adult/ Neonate multi-site spo2 sensor ,with FOAM wrap tape, DB9 , L=3m</b>			

**Compatibilities: Compatible with BCI** 3300, 3301, 3302,3303, 3304, 3401



**Biosys (Korea) compatible spo2 sensor**

Compatible Brand	OEM P/N	Greatmade P/N	Description(Connector, Length )	1-5 pcs	6-19 pcs	above 20 PCS
Biosys		AF002	Adult finger clip sensor , MD6M , L=3m			
		PF002	Pediatric finger clip sensor , MD6M , L=3m			
		AS002	Adult silicone soft tip sensor , MD6M , L=3m			
		PS002	Pediatric silicone soft tip sensor , MD6M , L=3m			
		NW002	Neonate silicone wrap sensor , MD6M , L=3m , <b>Non-detachable</b>			
		AE002	Adult ear clip sensor , MD6M , L=3m			
		AT002	Animal tongue clip sensor , MD6M , L=3m			
		NW 002 Y	Multi-site Y cable with separate wrap ,MD6M , <b>Wrap Detachable, L=3M</b>			
	NW 002 F	<b>Adult multi-site spo2 sensor ,with FOAM wrap tape,MD6M, L=3m</b>				

compatible with **Biosys** BPM-300



**Datex-Ohmeda Compatible reusable spo2 sensor** sale@greatmade.com

Compatible Brand	OEM P/N	Greatmade P/N	Description(Connector, Length )	1-5 pcs	6-19 pcs	above 20 PCS
Datex-Ohmeda See below compatibilities	OXY-F-UN	AF051	Adult finger clip sensor , <b>8pin Male</b> , L=1M			
		PF051	pediatric finger clip sensor , <b>8pin Male</b> , L=1M.			
		AS051	Adult silicone soft tip sensor , <b>8pin , Male</b> , L=1M			
		PS051	Pediatric silicone soft tip sensor , <b>8pin Male</b> , L=1M			
	OXY-W-UN	NW051	neonate silicone wrap sensor , <b>8pin Male</b> , L=1M , <b>Non-detachable</b>			
	OXY-E-UN	AE051	Adult ear clip sensor , <b>8pin Male</b> , L=1M			
		AT051	Animal tongue clip sensor , <b>8 pin Male</b> , L=1M			
		NW051-Y	Multi-site Y cable with separate wrap,8pin male, <b>Wrap Detachable</b>			
		NW051- F	<b>Adult/ Neonate multi-site spo2 sensor,8pin male,with FOAM wrap</b>			

Compatibilities : Use with Handheld pulse oximeters (Tuftsat, 3775 )  
 Stand Alone pulse oximeter (37\*\* series, 3800, 3900,3900P),respiratory gas monitors (4700, 5259RGM)  
 M-OSAT module , N-XOSAT option for Cardiocap.5 ; OEM patterner monitors ;  
 All S/5 monitors, all AS/3,CS/3 Monitors , Cardiocap ,Capnomac Ultima and Satellite  
**Require Adapter EX043 or EX044 or EX042**



**Datex-Ohmeda Compatible reusable spo2 sensor**

Compatible Brand	OEM P/N	Greatmade P/N	Description(Connector, Length )	1-5 pcs	6-19 pcs	above 20 PCS
Datex-Ohmeda Trusdat Pulse Oximeter	OXY-F4-MC	AF063	Adult finger clip sensor , <b>DB9 Male</b> , L=3M			
		PF063	pediatric finger clip sensor , <b>DB9 Male</b> , L=3M			
		AS063	Adult silicone soft tip sensor , <b>DB9 Male</b> , L=3M			
		PS063	Pediatric silicone soft tip sensor , <b>DB9 Male</b> , L=3M			
	OXY-W4-MC	NW063	neonate silicone wrap sensor , <b>DB9 Male</b> , L=3M , <b>Non-detachable</b>			
	OXY-E4-MC	AE063	Adult ear clip sensor <b>DB9 Male</b> L=3M			
		AT063	Animal tongue clip sensor , <b>DB9 Male</b> , L=3M			
		NW063-Y	Multi-site Y cable with separate wrap,DB9 male, <b>Wrap Detachable</b>			
		NW063- F	<b>Adult/ Neonate multi-site spo2 sensor,DB9 male,with FOAM wrap</b>			



**GE-Marquette Compatible reusable spo2 sensor**

Compatible Brand	OEM P/N	Greatmade P/N	Description(Connector, Length )	1-5 pcs	6-19 pcs	above 20 PCS
GE Dash ,Solar eagle , Tram Cardiolab/Med-lab Nellcor Non-Oximax tech		AF039	Adult finger clip sensor, <b>11PIN Male blue flat</b> ,L=3m			
		PF039	Pediatric finger clip sensor, <b>11PIN Male blue flat</b> ,L=3m			
		AS039	Adult silicone soft tip sensor, <b>11PIN Male blue flat</b> ,L=3m			
		PS039	Pediatric silicone soft tip sensor, <b>11PIN Male blue flat</b> ,L=3m			
		NW039	Neonate silicone wrap sensor, <b>11PIN Male blue flat</b> ,L=3m, <b>Non-detachable</b>			
		AE039	Adult ear clip sensor, <b>11PIN Male blue flat</b> ,L=3m			
		AT039	Animal tongue clip sensor, <b>11PIN Male blue flat</b> ,L=3m			
		NW039-Y	Multi-site Y cable with separate wrap,11pin flat, <b>Wrap Detachable</b>			
	NW039- F	<b>Adult/ Neonate multi-site spo2 sensor,11pin flat,with FOAM wrap</b>				



Compatible Brand	OEM P/N	Greatmade P/N	Description(Connector, Length )	1-5 pcs	6-19 pcs	above 20 PCS
GE Dash ,Solar eagle , Tram Cardiolab/Med-lab	LNCS-DCH-LNC-10-GE	AF011M	Adult finger clip sensor, <b>11PIN Male blue flat</b> ,L=3m			
		PF011M	Pediatric finger clip sensor, <b>11PIN Male blue flat</b> ,L=3m			
		AS11M	Adult silicone soft tip sensor, <b>11PIN Male blue flat</b> ,L=3m			

**Nihon Kohden Compatible reusable spo2 sensor (14pin)** [sale@greatmade.com](mailto:sale@greatmade.com)

Compatible Brand	OEM P/N	Greatmade P/N	Description(Connector, Length )	1-5 pcs	6-19 pcs	above 20 PCS
<b>Nihon Kohden</b>	TL-201T+	AF052	Adult finger clip sensor,14pin ,L=3M (Nihon Kohden compatible )			
BSM-2300 Lifescope i/L; BSM-4100 Lifescope P; BSM-5100 Lifescope A ; BSM-9510 Lifescope M BSM-9800 Lifescope s	K931 JL-900P	PF052	Pediatric finger clip sensor,14pin ,L=3M			
		AS052	Adult silicone soft tip sensor,14pin ,L=3M			
		PS052	Pediatric silicone soft-tip sensor,14pin ,L=3M			
		NW052	neonate wrap sensor,14pin ,L=3M, <b>Non-detachable</b>			
		AE052	Adult ear clip sensor,14pin ,L=3M			
		AT052	Animal tongue clip sensor,14pin ,L=3M			
		NW052-Y	Multi-site Y cable with separate wrap,14pin,Wrap Detachable ,L=3M			
		NW052-F	<b>Adult/ Neonate multi-site spo2 sensor,14pin,with FOAM wrap</b>			

Compatibilities : Use with **Nihon Kohden** :BSM-2300 Lifescope i/L; BSM-4100 Lifescope P;  
BSM-5100 Lifescope A ;BSM-9510 Lifescope M , BSM-9800 Lifescope s  
OPV-1500 Lifescope N, WEP-4200, WEP-8400 , DDG-2001  
Direct Connect



**Nihon Kohden Compatible reusable spo2 sensor**

Compatible Brand	OEM P/N	Greatmade P/N	Description(Connector, Length )	1-5 pcs	6-19 pcs	above 20 PCS
<b>Nihon Kohden</b>	TL-201T+	AF089	Adult finger clip sensor,3M 20pin ,L=3M			
OPV-1500 OPV-3100	JL-302T	PF089	Pediatric finger clip sensor,3M 20pin ,L=3M			
		AS089	Adult silicone soft tip sensor,3M 20pin ,L=3M			
		PS089	Pediatric silicone soft-tip sensor,3M 20pin ,L=3M			
		NW089	neonate wrap sensor,14pin ,3M 20pin, <b>Non-detachable</b>			
		AE089	Adult ear clip sensor,3M 20pin ,L=3M			
		AT089	Animal tongue clip sensor,3M 20pin ,L=3M			
		NW089-Y	Multi-site Y cable with separate wrap,3M 20pin,Wrap Detachable ,L=3M			
		NW089-F	<b>Adult/ Neonate multi-site spo2 sensor,3M 20pin,with FOAM wrap</b>			






Compatible with Fukuda ME KP-500



Compatible Brand	OEM P/N	Greatmade No.	MC001-5AG	Desceiption ( Cable and connector )	1-5pcs	6-19pcs	Above 20pcs
<b>Fukuda</b>		EC027A-R	Fukuda ME KP-500D	EKG cable , <u>IEC, DB15M&gt;DIN 3.0mm, 20KΩ resistor</u>			
KP-500D		EC027AD-R	Fukuda ME KP-500D	EKG cable , <u>AHA, DB15M&gt; DIN 3.0mm , 20KΩ resistor</u>			
		EC027B-R	Fukuda ME KP-500D	EKG cable , <u>AHA, DB15M&gt; Banana 4.0 , 20KΩ resistor</u>			
		EC027BI-R	Fukuda ME KP-500D	EKG cable , <u>IEC, DB15M&gt; Banana 4.0 , 20KΩ resistor</u>			
		EC027C-R	Fukuda ME KP-500D	EKG Cable , <u>AHA, DB15M&gt;Grabber , 20KΩ resistor</u>			
		EC027CI-R	Fukuda ME KP-500D	EKG Cable , <u>IEC, DB15M&gt;Grabber , 20KΩ resistor</u>			
		EC027D-R	Fukuda ME KP-500D	EKG Cable , <u>AHA, DB15M&gt;Snap , 20KΩ resistor</u>			
		EC027DI -R	Fukuda ME KP-500D	EKG Cable , <u>IEC, DB15M&gt;Snap, 20KΩ resistor</u>			




Compatible with Fukda ME KP-500D .



Compatible Brand	OEM P/N	Greatmade No.		Desceiption ( Cable and connector )	1-5pcs	6-19pcs	Above 20pcs
<b>Fukuda</b>		EC028A	Fukuda ME FX-031	EKG cable , <u>IEC, Round 11pin &gt;DIN 3.0mm, No resistor</u>			
FX-031		EC028AD	Fukuda ME FX-031	EKG cable , <u>AHA, Round 11pin&gt; DIN 3.0mm , No resistor</u>			
		EC028B	Fukuda ME FX-031	EKG cable , <u>AHA, Round 11pin &gt; Banana 4.0 , No resistor</u>			
		EC028BI	Fukuda ME FX-031	EKG cable , <u>IEC, Round 11pin &gt; Banana 4.0 , No resistor</u>			
		EC028C	Fukuda ME FX-031	EKG Cable , <u>AHA, Round 11pin &gt;Grabber , No resistor</u>			
		EC028CI	Fukuda ME FX-031	EKG Cable , <u>IEC, Round 11pin &gt;Grabber , No resistor</u>			
		EC028D	Fukuda ME FX-031	EKG Cable , <u>AHA, Round 11pin &gt;Snap, No resistor</u>			
		EC028DI	Fukuda ME FX-031	EKG Cable , <u>IEC, Round 11pin &gt;Snap, No resistor</u>			



***Nihon Kohden compatible One piece ECG cable with leadwires***

Compatible Brand	OEM Part No.	Greatmade No.	Description (Cable & Connector )
<b>Nihon Kohden</b>		MC051-6AS	one piece ECG cable with leadwire, <b>6-lead</b> ,AHA, Snap, 12pin
BSM2300 Lifescope I/L,		MC051-6AG	one piece ECG cable with leadwire, <b>6-lead</b> ,AHA, Grabber, 12pin
BSM-4100 Lifescope P,		MC051-3AS	one piece ECG cable with leadwire, <b>3-lead</b> ,AHA, Snap, 12pin
BSM-5100 Lifescope A		MC051-3AG	one piece ECG cable with leadwire, <b>3-lead</b> ,AHA, Grabber, 12pin
BSM-9510 Lifescope M		MC051-6IS	one piece ECG cable with leadwire, <b>6-lead</b> ,IEC, Snap , 12pin
,BSM-9800 Lifescope S ,		<b>MC051-6IG</b>	one piece ECG cable with leadwire, <b>6-lead</b> ,IEC, Grabber, 12pin
BSM-1500 Lifescope N		MC051-3IS	one piece ECG cable with leadwire, <b>3-lead</b> ,IEC, Snap , 12pin
		MC051 -3IG	one piece ECG cable with leadwire, <b>3-lead</b> ,IEC, Grabber, 12pin

**Patient Monitor Cable and Leadwire**

[www.greatmade.com](http://www.greatmade.com)



***Nihon Kohden Compatible ECG cable with leadwires***



Product Service

## EC Certificate

### Full Quality Assurance System

Directive 93/42/EEC on Medical Devices (MDD), Annex II excluding (4)  
(Devices in Class IIa, IIb or III)

No. G1 17 12 02231 002

**Manufacturer:** Shenzhen Greatmade Tech limited

3rd Floor, Building B  
Baifuli Industrial Zone, Shanghenglang  
Huahui Road, Dalang Street  
Longhua New District  
518109 Shenzhen, Guangdong Province  
PEOPLE'S REPUBLIC OF CHINA



**EC-Representative:** Prolinx GmbH

Brehmstr. 56  
40239 Duesseldorf  
GERMANY

**Product Category(ies):** Spo2 sensor

The Certification Body of TÜV SÜD Product Service GmbH declares that the aforementioned manufacturer has implemented a quality assurance system for design, manufacture and final inspection of the respective devices / device categories in accordance with MDD Annex II. This quality assurance system conforms to the requirements of this Directive and is subject to periodical surveillance. For marketing of class III devices an additional Annex II (4) certificate is mandatory. See also notes overleaf.

**Report No.:** GZ1728501

**Valid from:** 2018-04-11  
**Valid until:** 2023-04-10

**Date,** 2018-04-11

*S. Preiß*  
Stefan Preiß



TÜV SÜD Product Service GmbH is Notified Body with identification no. 0123

Page 1 of 2



Product Service

## EC Certificate

### Full Quality Assurance System

Directive 93/42/EEC on Medical Devices (MDD), Annex II excluding (4)  
(Devices in Class IIa, IIb or III)

No. G1 17 12 02231 002

**Facility(ies):**

Shenzhen Greatmade Tech limited  
3rd Floor, Building B, Baifuli Industrial Zone, Shanghenglang,  
Huahui Road, Dalang Street, Longhua New District, 518109  
Shenzhen, Guangdong Province, PEOPLE'S REPUBLIC OF  
CHINA

Page 2 of 2



Product Service

## CERTIFICATE

No. Q1N 17 12 02231 001

**Holder of Certificate:** Shenzhen Greatmade Tech limited

3rd Floor, Building B  
Baifuli Industrial Zone, Shanghenglang  
Huahui Road, Dalang Street  
Longhua New District  
518109 Shenzhen, Guangdong Province  
PEOPLE'S REPUBLIC OF CHINA

**Facility(ies):**

Shenzhen Greatmade Tech limited  
3rd Floor, Building B, Baifuli Industrial Zone,  
Shanghenglang, Huahui Road, Dalang Street,  
Longhua New District, 518109 Shenzhen,  
Guangdong Province, PEOPLE'S REPUBLIC OF  
CHINA



**Certification Mark:**



**Scope of Certificate:** Design and Development, Production  
and Distribution of Spo2 sensor,  
Patient cable and leadwire, Blood pressure cuff

**Applied Standard(s):**

EN ISO 13485:2012 + AC:2012  
Medical devices - Quality management systems -  
Requirements for regulatory purposes  
(ISO 13485:2003 + Cor. 1:2009)  
DIN EN ISO 13485:2012  
Upgrade required until 2019-03-31

The Certification Body of TÜV SÜD Product Service GmbH certifies that the company mentioned above has established and is maintaining a quality management system, which meets the requirements of the listed standard(s). See also notes overleaf.

**Report No.:** GZ1728501

**Valid from:** 2018-04-11  
**Valid until:** 2021-04-10

**Date,** 2018-04-11

*S. Preiß*  
Stefan Preiß



Page 1 of 1