



Glass staining dishes

General features

Glass staining dishes to use for special stains.

They resist to each stain solution and solvent, but they have not to be used with microwave oven.



Distributor: Bio-Optica Milano S.p.A.

Code	Description	Capacity	Size (mm)	Liquid volume
10-1410	Hellendhal dish	8 slides	50 x 30 x 90	70 ml
10-1610	Schefferdecker dish	8 slides	55 x 50 x 60	80 ml
10-1710	Coplin dish	6 slides	30 x 30 x 70	50 ml
10-2000	Complete staining set	10 slides	81 x 101 x 70	90 ml
10-1820	Spare dish and rack for 10-2000 set	10 slides	Dish 81 x 101 x 70 Rack 68 x 87 x 50	90 ml
10-1910	Spare rack for 10-2000 set	10 slides	68 x 87 x 50	-
10-2010	Spare tweezer for 10-2000 set	-	-	-

Technical features

For this kind of product, special precautions for handling and storage are not necessary.

Validity date: not applicable.

The products are realized for generic lab uses, ideal for the stain of microscope slides.

The product don't request CE mark.