


# Non-Adhesion Slides

## Product properties

- Manufactured from extra-white soda-lime glass with very low iron content
- Clean, flat and ready to use immediately
- Manufactured in accordance with ISO 8037/1 dimensions: 26 x 76 x 1.0 mm (+0/-1.0 mm)
- Tolerance for thickness: ±0.05 mm
- Planarity under 30 µm measured over 60 mm length




## Superfrost slides



Epredia Superfrost slides feature a durable, easy to write, print and read tab and may be used in general laboratory applications including histology, cytology, haematology, microbiology and urinalysis. The printed tab is impervious to virtually all common laboratory chemicals and reagents. Offered in packs of 50 slides in three case sizes. Tropical pack with sealed-foil wrap available upon request.

## Ordering Information

Product Cat. No.					
<b>Epredia Superfrost slides</b>		 <p>"Superfrost" Marking</p>			
<b>90° clipped corners</b>					
White	QTY 5000	AF00008032E01MNZ10	Other colours available		
	QTY 2000	AF00008032E01MNZ20			
	QTY 1000	AF00008032E01MNZ50			
			<b>Cut</b>	<b>90° ground</b>	<b>45° ground</b>
White	QTY 5000	AA00008032E01MNZ10	AG00008032E01MNZ10	AB00008032E01MNZ10	
	QTY 2000	AA00008032E01MNZ20	AG00008032E01MNZ20	AB00008032E01MNZ20	
	QTY 1000	AA00008032E01MNZ50	AG00008032E01MNZ50	AB00008032E01MNZ50	
Blue	QTY 5000	AA00008132E01MNZ10	AG00008132E01MNZ10	AB00008132E01MNZ10	
	QTY 2000	AA00008132E01MNZ20	AG00008132E01MNZ20	AB00008132E01MNZ20	
	QTY 1000	AA00008132E01MNZ50	AG00008132E01MNZ50	AB00008132E01MNZ50	
Pink	QTY 5000	AA00008232E01MNZ10	AG00008232E01MNZ10	AB00008232E01MNZ10	
	QTY 2000	AA00008232E01MNZ20	AG00008232E01MNZ20	AB00008232E01MNZ20	
	QTY 1000	AA00008232E01MNZ50	AG00008232E01MNZ50	AB00008232E01MNZ50	
Yellow	QTY 5000	AA00008332E01MNZ10	AG00008332E01MNZ10	AB00008332E01MNZ10	
	QTY 2000	AA00008332E01MNZ20	AG00008332E01MNZ20	AB00008332E01MNZ20	
	QTY 1000	AA00008332E01MNZ50	AG00008332E01MNZ50	AB00008332E01MNZ50	
Green	QTY 5000	AA00008432E01MNZ10	AG00008432E01MNZ10	AB00008432E01MNZ10	
	QTY 2000	AA00008432E01MNZ20	AG00008432E01MNZ20	AB00008432E01MNZ20	
	QTY 1000	AA00008432E01MNZ50	AG00008432E01MNZ50	AB00008432E01MNZ50	
Orange	QTY 5000	AA00008532E01MNZ10	AG00008532E01MNZ10	AB00008532E01MNZ10	
	QTY 2000	AA00008532E01MNZ20	AG00008532E01MNZ20	AB00008532E01MNZ20	
	QTY 1000	AA00008532E01MNZ50	AG00008532E01MNZ50	AB00008532E01MNZ50	
Violet	QTY 5000	AA00008632E01MNZ10	AG00008632E01MNZ10	AB00008632E01MNZ10	
	QTY 2000	AA00008632E01MNZ20	AG00008632E01MNZ20	AB00008632E01MNZ20	
	QTY 1000	AA00008632E01MNZ50	AG00008632E01MNZ50	AB00008632E01MNZ50	