

**EKOTEKS LABORATUVAR ve GÖZETİM  
HİZMETLERİ A.Ş.**

AB-0583-T

21003762-  
İNG

02-21

**TEST RESULT**

**WATER PERMEABILITY; ISO 811:2018**

Hydrostatic Head Tester, Textest marka Fx 3000 model  
Temperature of water 20°C. Pressure increase ratio 10 mbar/min.  
Performed in the conditioned room (20±2°C-65%±4)  
Less critical area

Sample 1  
Sample 2  
Sample 3  
Sample 4  
Sample 5

**RESULT**

25,3 cmH<sub>2</sub>O  
24,3 cmH<sub>2</sub>O  
30,8 cmH<sub>2</sub>O  
29,2 cmH<sub>2</sub>O  
25,9 cmH<sub>2</sub>O

Average

27,1 cmH<sub>2</sub>O

**REQUIREMENT**

≥ 10 cmH<sub>2</sub>O

**MICROBIAL CLEANLINESS (Bioburden)**

Test Metod: Ref: EN ISO 11737-1:2018

The sample is put in extraciton liquid after shaking well, inoculated on the agar.  
After incubation at 30 ± 1 ° C for 72 hours, growth microorganisms are counted on the agar.

Critical, Less Critical area

|   | <b>RESULTS</b>              | <b>REQUIREMENT</b>           |
|---|-----------------------------|------------------------------|
| Microbial cleanliness (cfu/ 100 cm <sup>2</sup> ) | 280 cfu/100 cm <sup>2</sup> | ≤300 cfu/100 cm <sup>2</sup> |

\*cfu= Colony forming unit.

**ZENMED TIBBİ CİHAZLAR  
VE SAĞLIK ÜRÜNLERİ SAN. VE TİC. LTD.ŞTİ.**  
Mahmutbey Mah. 2567.Sk. No:1 İş Yeri No:11,21  
Bağcılar / İSTANBUL / TÜRKİYE  
Tel: 0212 445 45 25 Tic.Sic.No: 294704  
Güneşli V.D.: 997 067 1942