# **SE-1201**

### Electrocardiograph

Version 1.4

## **Technical Specifications**

### **Physical Specifications**

#### **Dimensions**

361mm×262mm×135mm

### Weight

≤4.2kg (Excluding recorder paper and battery)

### Display

7 inch, 800×480 multicolor LCD screen (touch screen as optional)

### **Power Supply**

#### Mains Supply

Operating Voltage = 100V-240V~

Operating Frequency = 50Hz/60Hz

Input Current = 0.9-0.4A

### **Internal Li-ion Battery Pack**

Rated voltage = 14.8V

Rated capacity = 2500/5000mAh

(4/8 normal work hours)

Necessary Charge time: 3 hours

### **Recording**

### Recorder

#### Thermal dot-matrix recorder

### **Recorder Paper:**

Folded thermal paper:

### 210mm×140mm×144 pages

Folded thermal paper:

215mm×280mm×100 pages (Optional)

Folded thermal paper:

210mm×295mm×100 pages (Optional)

### **Printing Density:**

8 dots per mm / 200 dots per inch (amplitude axes)

40 dots per mm / 1000 dots per inch (time axes, @ 25 mm/s)

### **Paper Speed:**

5mm/s, 6.25mm/s, 10mm/s, 12.5mm/s, 25mm/s, 50mm/s (±3%)

#### **External Printer:**

HP2010/2132/1050/2000, HPM401/2015/2035,

HP1106/1020P, HP1112/4729/3638, HPM202D/1525/M403D



### **HR Recognition**

HR Range	Accuracy
30 BPM ~300 BPM	±1 BPM

### **ECG Unit**

Leads:	Acquisition Mode:
12 standard leads	Simultaneously 12 leads
A/D Converter:	Resolution:
24bits	0.1192uV/LSB
Time Constant:	Frequency Response:
≥5s	$0.01$ Hz $\sim 500$ Hz (-3dB)
Sampling Frequency:	Gain:
64,000 Hz	1.25, 2.5, 5, 10, 20, 10/5 mm/mV, AGC
Input Impedance:	Input Circuit Current:
≥100MΩ (10Hz)	≤0.01μA
Input Voltage Range:	Calibration Voltage:
≤±5 mVpp	$1 \text{mV} \pm 2\%$
CMRR:	QTc Formulas:
≥140dB (AC on)	Bazzet, Fridericia, Framingham, Hodges, and QRS
≥123dB (AC off)	

### **Pacemaker**

Amplitude:	Width:
500μV to ±700 mV	30μs to 2.0 m

### **Filter**

AC Filter:	DFT Filter:
50/60Hz	0.01Hz/0.05Hz/0.15Hz/0.25Hz/0.32Hz/0.5Hz/0.67Hz
EMG Filter:	LOWPASS Filter:
Off/25Hz/35Hz/45Hz	350Hz/300Hz/270Hz/150Hz/100Hz/75Hz

### **Data Transmission**

Report Format:	Data Transmission:
PDF, XML, DICOM, FDA-SCP	Wi-Fi, Ethernet, RS232
Data Management System:	HIS connection:
SE-1515 Data Management System, bi-	DICOM Worklist/DICOM Storage/HL7/GDT
directional communication	

**Data Storage:** 

800 ECGs



#### Wi-Fi

**Transmitting Frequency:** 

2400-2497MHz

**Modulation Type:** 

DSSS, CCK, OFDM

**Effective Radiated Power:** 

6 - 17dBm

### **Frequency Band:**

2400-2497MHz

**Transmitting Power:** 

6 - 17dBm

### **Safety Specifications**

### **Comply with:**

IEC 60601-1:2005/A1:2012

EN 60601-1:2006/A1:2013

IEC 60601-1-2:2014

EN 60601-1-2:2015

IEC/EN 60601-2-25

### **Patient Auxiliary Current**

### $NC < 10\mu A (AC) / < 10\mu A (DC)$

 $SFC < 50\mu A (AC) / < 50\mu A (DC)$ 

### **Anti-electric-shock type:**

Class I with internal power supply

### **Patient Leakage Current:**

 $NC < 10\mu A (AC) / < 10\mu A (DC)$ 

 $SFC < \!\! 50\mu A\left(AC\right) / \! < \!\! 50\mu A\left(DC\right)$ 

### Anti-electric-shock degree:

CF type with defibrillation-proof

### **Environment Specifications**

### **Temperature:**

Transport & Storage: -20°C (-4°F) ~

+55°C (+131°F)

Working:  $+5^{\circ}$ C ( $+41^{\circ}$ F)  $\sim +40^{\circ}$ C ( $+104^{\circ}$ F)

### **Atmospheric Pressure:**

Transport & Storage: 70kPa ~106kPa

Working: 86kPa ~106kPa

### **Relative Humidity:**

Transport & Storage: 25%~93% Non-Condensing

Working: 25%~80% Non-Condensing

