

Shanghai HBM Healthcares, Inc.

#1508-1509, No. 211, Golden Port Road, Pudong, Shanghai, China. Tel: 0086-21-58445785 www.hbmchina.com

Medispo® Medical gowns

1. Description


- --- designed with classic neck and waist ties, for easily adjustable.
- --- Knitted cuffs.
- --- Fluid Resistant.
- --- Color: Blue, Green and other colors also available.
- --- Gowns up to size 3XL, regular and extra long lengths.
- --- Comply with AAMI PB70 Level 2.

2. Composition

Polypropylene 23gsm + Polyethylene 12gsm

3. Specification

REF	Fabric Weight	Length	Packing
9730S	35gsm	110cm	
9731M	35gsm	115cm	
9732L	35gsm	120cm	10Pcs /polybag
9733XL	35gsm	125cm	100pcs /Carton
9734XXL	35gsm	130cm	
9735XXXL	35gsm	135cm	
	9730S 9731M 9732L 9733XL 9734XXL	9730S 35gsm 9731M 35gsm 9732L 35gsm 9733XL 35gsm 9734XXL 35gsm 9735XXXL 35gsm	9730S 35gsm 110cm 9731M 35gsm 115cm 9732L 35gsm 120cm 9733XL 35gsm 125cm 9734XXL 35gsm 130cm 9735XXXL 35gsm 135cm



4. Usage instruction

- ---It is used for medical staff to avoid being contaminated by blood, body fluids, and other infectious substances during contact.
- ---Isolation gown is a two-way isolation to prevent medical staff from being infected or contaminated as well as patients from being infected.
- ---When entering a strictly isolated ward, you need to wear isolation gowns.
- ---Examination and nursing require special isolation of patients. When the work clothes may be contaminated by secretions, excrement, blood or body fluids, isolation clothes shall be worn Medical staff entering the patient wards of infectious diseases that are likely to cause dissemination in the hospital and patients who require special isolation (such as extensive burns, organ transplants, and premature babies) are required to wear isolation gowns.
- ---Single use only.

5. Sterilization Non-Sterile

