

# Colenta® DX Easy Imaging CR Retrofit Kit HighCap Xr

**Extremely fast processing in a compact body streamlines diagnosis**

**This compact table-top Reader unit has outstandingly high-speed processing capability to challenge the common misconception of CR.**

## Processing Capacity

Up to 47 IPs/hr in size 14"x17" (35cmx43cm)

The Colenta DX Easy Imaging CR Retrofit Kit HighCap Xr is a table-top Reader unit which boasts a world-class high processing speed of 73 IPs/hour\*. As with all Colenta equipment the DX Easy Imaging CR Retrofit Kit HighCap Xr is easy to operate. It will help to reduce patient waiting time and greatly increase the efficiency of examination workflow.

*\*) The above value is for IPs of 8" x 10" (18cmx24cm)*



## Reading Resolution - Pixel Size

254dpi (100µm) / 508dpi (50µm)

One of the reader's key components is the Imaging Plate, or IP, which is used instead of the X-ray film to easily digitize X-ray images.

Although this is a compact machine its excellent image quality is at the same level of the one featured by the rest of the Colenta range of equipment. The Colenta's proprietary image processing technology, enhances image contrast and sharpness, without any deterioration of the details. Colenta DX Easy Imaging CR Retrofit Kit HighCap Xr supports accurate diagnosis by offering stable and optimized image quality.

## Space-Saving Design

As all the Colenta equipment the Colenta DX Easy Imaging CR Retrofit Kit HighCap Xr is easy to operate. It will help to reduce patient waiting time and to greatly increase the efficiency of examination workflow.

Colenta DX Easy Imaging CR Retrofit Kit HighCap Xr, with its light and compact table top design, can be placed on a desk, shelf or anywhere the space is limited. As this is a fully digital Reader, neither a darkroom nor an automatic processor is required. The Colenta DX Easy Imaging CR Retrofit Kit HighCap Xr can always be installed in the space formally used by a chemistry based processor.




Minimum requirement of PC: Processor min. Intel Core i5, min. 8GB Memory, 320GB SSD + HDD S-ATA 1TB, 2 x Intel GBit LAN controller, Intel VGA graphic controller with 24 bit color depth, Video Mode of graphic card supports OpenGL-Version 1.1.

# Colenta<sup>®</sup> DX Easy Imaging CR Retrofit Kit HighCap Xr

## Specifications

<b>Image Plate General X-ray</b>	14" x 17" (35cm x 43 cm), 14" x 14" (35cm x 35 cm), 10" x 12" (24cm x 30cm), 8" x 10" (18cm x 24cm), 15cm x 30cm
<b>Processing Capacity</b>	14" x 17" (35cm x 43 cm) up to 47 IPs/hr 8" x 10" (18cm x 24cm) up to 73 IPs/hr
<b>Time Required for IP Feed/Load</b>	Min. 49 sec.
<b>Time to Start on Display</b>	Min. 33 sec.
<b>Grayscale resolution</b>	12 Bit
<b>Reading Specification</b>	5 pixel/mm, 10 pixel/mm, 20 pixel/mm (Mod. XMr)
<b>Pixel Size</b>	100um, 50 um (mod. XMr)
<b>Number of Stacker</b>	1
<b>Network</b>	10 Base T/100 Base TX
<b>Weight</b>	39 kg (86 lbs.)
<b>Dimensions (WxDxH)</b>	560mm x 540mm x 392 mm (22" x 21" x 15")
<b>Power Supply Conditions</b>	Single phase 50-60Hz, AC120-240V±10%
<b>Environmental Conditions</b>	Temperature: 15-30°C Humidity: 40-80%RH (No dew condensation)



This equipment is a Class 1 laser product (IEC60825-1:2001).

*Specifications are subject to change without notice. All brand names or trademarks are the property of their respective owners. In some countries, regulatory approval may be required to import medical devices.*

## **COLENTA** Labortechnik GmbH & Co. KG

Neunkirchner Str. 117

A-2700 Wiener Neustadt, Austria

Tel.: +43(0)2622-28311-0 Fax: +43(0)2622-28311-7

**E-Mail:** [office@colenta.at](mailto:office@colenta.at)