

Product-Information-Sheet

Respiratory endoscopy

Biopsy Forceps



Article number	Article Name	Quantity value	Description	Endoscope Compatibility	Area Compatibility
028987	MAJ-209	20	Single-use, sterile suction valve	Mobile airwayscopes: MAF-TM, MAF-GM, MAF-TM2, MAF-GM2, MAF-DM2 Bronchoscopes: OES 30, 40, 60, EVIS 190, 240, 260, 290, 1100, 1200 Ultrasonic bronchoscopes: BF-UC260-OL8, -UC260FW, -UC290F Pleuroscopes: LTF-240, -260, -H290 Rhinolaryngoscopes: ENF-VT, -VT2	



- Innovative, watertight construction
- Smoother device insertion
- Secure mounting on the biopsy port
- Convenient sterile packages



Article number	Article Name	Quantity value	Description	Endoscope Compatibility	Area Compatibility
028088	MAJ-210	20	Single-use, sterile biopsy valve	Mobile airwayscopes: MAF-TM, MAF-TM2 Bronchoscopes: OES 20, 20D, 30, 40, 60, EVIS 190, 200, 240, 260, 290, 1100, 1200 Ultrasonic bronchoscopes: BF-UC260-OL8, -UC260FW, -UC290F Pleuroscopes: LTF-240, -260, -H290 Rhinolaryngoscopes: ENF-VT, -T3, -VT2 Tracheal intubation scopes: LF-2, -DP, -GP, -TP	