

Name of Product : Greenday Ultra Hygienic Pad

SIZES	NORMAL	LONG	NIGHT
Width (mm)	165 ± 5	165 ± 5	165 ± 5
Length (mm)	235 ± 5	280 ± 5	320 ± 5
Pad Length (mm)	235 ± 5	280 ± 5	320 ± 5
Top Perforated (width) (mm)	85 ± 5	85 ± 5	85 ± 5
Leak Guard (width) (mm)	130 ± 5	130 ± 5	130 ± 5
Pouch NW (width) (mm)	220 ± 5	260 ± 5	300 ± 5
Disposable Tape (width) (mm)	10 ± 2	10 ± 2	10 ± 2
Airlaid with SAP(Width) (mm)	60 ± 5	60 ± 5	60 ± 5
Airlaid with SAP( Length) (mm)	190 ± 5	235 ± 5	275 ± 5
Airlaid without SAP(Width) (mm)	70 ± 5	70 ± 5	70 ± 5
Airlaid without SAP ( Length) (mm)	210 ± 5	250 ± 5	280 ± 5
Body Siliconised paper (width) (mm)	55 ± 5	55 ± 5	55 ± 5
Body Siliconised paper (Length) (mm)	190 ± 5	230 ± 5	260 ± 5
Wing Siliconised Paper (width) (mm)	55 ± 5	55 ± 5	55 ± 5
Wing Siliconised Paper (length) (mm)	55 ± 5	80 ± 5	80 ± 5
Pad Weight (gr)	4,89 ± 1	5,95 ± 1	6,80 ± 1
Max Absorption Capacity (g) (TZS250:1985)(ISO11948)	46,17 ± 10	56,89 ± 10	66,82 ± 10
Max Retention Capacity (g) (TZS250:1985)	20,81 ± 3	25,68 ± 3	30,13 ± 3

RAW MATERIALS

RAW MATERIALS	INGREDIENTS
Nonwoven Topsheet	Polypropylene
Leg guard	Polypropylene
Glue	Hot melt
Airlaid without SAP	Airlaid Nw
Airlaid with SAP	Airlaid Nw+ Salt of Cross-linked Polyacrylic Acid
Disposable Tape	Polyethylene
Pouch NW	Polyethylene
Body Siliconised paper	Polypropylene
Wing Siliconised Paper	Polypropylene
Packaging Material	Polyethylene