


# Detro<sup>®</sup> San SFC

Concentrated Surface Disinfectant



### USAGE AREAS:

- Used in floor and surface disinfection of critical areas such as operating theater, intensive care and laboratory. It is possible to use it in fixtures by wiping method.

### PRODUCT PROPERTIES:

- Concentrated product.
- Odorless feature thanks to aldehyde, phenol and chlorine-free structure.
- Thanks to being odorless, can be used in the presence of staff or patients.

### COMPOSITION:

- Bis (3-aminopropyl) dodecylamine, dedecyldimethylammonium chloride, benzalkonium chloride, auxiliary substances.




### MICROBIOLOGICAL ACTIVITY:

- Bactericidal (including TBC), fungicidal, virucidal in 15 minutes at %0,25 concentration.

### USAGE:

- Prepared by diluting with network water according to the necessary concentration ratio. 2,5 ml of **Detro<sup>®</sup> San SFC** is added to 1000 ml (1 Liter) of water to prepare 1 liter solution at 0.25% concentration, and it is applied to floor by wet mop and to the fixtures and tables by wiping.
- The prepared solution can be used maximum for 24 hours.
- Can be changed earlier according to dirtiness level.

### PACKAGE:

| QUANTITY   | PACKING TYPE | BARCODE       |
|--|--------------|---------------|
|  5L | 4 x DRUM     | 8680097103062 |
|  1L | 12 x BOTTLE  | 8680097103208 |
|  2L | 6 x BOTTLE   | 8680097103185 |

