

2F., 12th Building, ShenZhen
SoftwareParkPhase II,Keji Middle 2nd
Road, Nanshan District, Shenzhen,
Guangdong, China
Tel:86-755-26722890
Fax:86-755-26722850
E-mail:sonoscape@sonoscape.net



Statement

To whom it may concern:

We, Sonoscape Medical Corp, as the manufacturer of P15 ultrasound system, confirm that P15 system has over 150000 processing digital channels.

About how it be calculated:

$P15 \text{ digital channels} = 128 \text{ (physical channels)} * 2 \text{ (signal receive and send)} * 5 \text{ (transducer ports)} * 3 \text{ (spatial compound lines)} * 12 \text{ (focus)} * 4 \text{ (multi beam line)}$

So, its $128 * 2 * 5 * 3 * 12 * 4 = 184320$ digital processing digital channels.

Sonoscape Medical Corp.

2020-12-29

Mr. HUANG YIBO

Vice President

SonoScape Medical Corp.

