

Endotracheal tube

Specification

Type	Size
Endotracheal tube	2.5# 3.0# 3.5# 4.0# 4.5# 5.0# 5.5# 6.0# 6.5# 7.0# 7.5# 8.0# 8.5# 9.0# 9.5# 10.0#



Endotracheal tube
1.02.01.100T

Specification

Type	Size
Preformed nasal ET tube	2.5# 3.0# 3.5# 4.0# 4.5# 5.0# 5.5# 6.0# 6.5# 7.0# 7.5# 8.0# 8.5# 9.0# 9.5# 10.0#



Preformed nasal ET tube
1.02.01.601T

Specification

Type	Size
Preformed oral ET tube	6.0# 6.5# 7.0# 7.5# 8.0# 8.5# 9.0#



Preformed oral ET tube
1.02.01.701T

Injectable ET Tube

Specification

Type	Size
Injectable ET Tube	6.0# 6.5# 7.0# 7.5# 8.0# 8.5# 9.0#

Give drugs (anesthetics, antibiotics and antispasmodics etc) to mucosa of airway

Function:

- » Relieving airway spasm
- » Reducing airway injury
- » Removing sputum



Injectable ET Tube
1.02.01.501T

Reinforced Endotracheal tube

Specification

Type	Size
Traditional	2.5# 3.0# 3.5# 4.0# 4.5# 5.0# 5.5# 6.0# 6.5# 7.0# 7.5# 8.0# 8.5# 9.0# 9.5# 10.0#



Endotracheal tube
1.02.04.100T

Specification

Type	Size
Traditional	2.5# 3.0# 3.5# 4.0# 4.5# 5.0# 5.5# 6.0# 6.5# 7.0# 7.5# 8.0# 8.5# 9.0# 9.5# 10.0#



Preformed oral ET tube

Specification

Type	Size
Traditional	2.5# 3.0# 3.5# 4.0# 4.5# 5.0# 5.5# 6.0# 6.5# 7.0# 7.5# 8.0# 8.5# 9.0# 9.5# 10.0#



Preformed nasal ET tube