



HARDIK®
INTERNATIONAL

House of Ortho Tools



SPINAL
IMPLANTS & INSTRUMENTS



• **T401** | Monoaxial Screw

MATERIAL - Titanium
DIAMETER - 4.5mm, 5.5mm, 6.5mm
LENGTH - 30mm to 50mm (5mm Increment)



• **T402** | Polyaxial Screw

MATERIAL - Titanium
DIAMETER - 4.5mm, 5.5mm, 6.5mm
LENGTH - 30mm to 50mm (5mm Increment)



• **T406** | Reduction Screw

MATERIAL - Titanium
DIAMETER - 5.5mm, 6.5mm
LENGTH - 30mm to 50mm (5mm Increment)



• **T425** | Connecting Rod

MATERIAL - Titanium
DIAMETER - 6mm
LENGTH - 50mm to 100mm (10mm Increment)
120mm to 200mm (20mm Increment)
250mm to 300mm (50mm Increment)



• **T426** | Cross Connector

MATERIAL - Titanium
LENGTH - 30mm x 34mm, 34mm x 42mm,
40mm x 53mm, 50mm x 70mm

Instruments



• 675.001 | Straight Awl



• 675.002 | Curved Awl



• 675.003 | Marker Inserter



• 675.037 | Marker Left (5 pcs)
• 675.038 | Marker Right (5 pcs)



• 675.004 | Open Awl



• 675.005 | Quick-Coupling T Handle



• 675.006 | Quick-Coupling Handle



• 675.007 | Tap 4.5mm
• 675.008 | Tap 5.5mm
• 675.009 | Tap 6.5mm



• 675.010 | Straight Probe



• 675.039 | Curved Probe



• 675.040 | Hard Probe



• 675.011 | Mono-axial Screwdriver (short)



• 675.012 | Multi-axial Screwdriver (short)



• 675.013 | Screw Holder



• 675.015 | Rod Holder



• 675.016 | Rod Rotator



• 675.017 | Rod Bender



• 675.018 | Bending Device (left and Right)



• 675.019 | Rod Pusher



• 675.021 | Spreader



• 675.022 | Compressor



• 675.023 | Anti-Rotating Sleeve



• 675.024 | Rod Compressor



• 675.025 | Reduction Breaker



• 675.027 | Wrench



• 675.028 | Mono-axial Screwdriver (long)



• 675.029 | Multi-axial Screwdriver (long)



• 675.030 | Small Wrench $\phi 4.0$



• 675.031 | Template (2 pcs)



• 675.032 | Rod Persuader



• 675.033 | Rotating Device SW5.0




• 675.034 | Rod Fork



• 675.035 | Screwdriver SW3.5



• 675.036 | Hex Screwdriver SW3.5



Implants

T401	Monoaxial Screw 4.5mm x 30mm, 35mm, 40mm, 45mm, 50mm	Each 4
T401	Monoaxial Screw 5.5mm x 30mm, 35mm, 40mm, 45mm, 50mm	Each 4
T401	Monoaxial Screw 6.5mm x 30mm, 35mm, 40mm, 45mm, 50mm	Each 4
T402	Polyaxial Screw 4.5mm x 30mm, 35mm, 40mm, 45mm, 50mm	Each 4
T402	Polyaxial Screw 5.5mm x 30mm, 35mm, 40mm, 45mm, 50mm	Each 4
T402	Polyaxial Screw 6.5mm x 30mm, 35mm, 40mm, 45mm, 50mm	Each 4
T406	Reduction Screw 5.5mm x 30mm, 35mm, 40mm, 45mm, 50mm	Each 4
T406	Reduction Screw 6.5mm x 30mm, 35mm, 40mm, 45mm, 50mm	Each 4
T425	Connecting Rod 6.0mm x 50mm, 60mm, 70mm, 80mm, 90mm, 100mm, 120mm, 140mm, 160mm, 180mm, 200mm, 250mm, 300mm	Each 2
T426	Cross Connector 30mm x 34mm, 34mm x 42mm, 40mm x 53mm, 50mm x 70mm	Each 1

AGC.37 Spinal Pedicle Instrument Set in Graphic Case



Instruments

675.001	Straight Awl	1 Each
675.002	Curved Awl	1 Each
675.003	Marker Inserter	1 Each
675.037	Marker left (5 pcs)	1 Each
675.038	Marker right (5 pcs)	1 Each
675.004	Open Awl	1 Each
675.005	Quick-Coupling Handle	1 Each
675.006	Quick-Coupling Handle	1 Each
675.007	Tap 4.5mm	1 Each
675.008	Tap 5.5mm	1 Each
675.009	Tap 6.5mm	1 Each
675.010	Straight Probe	1 Each
675.039	Curved Probe	1 Each
675.040	Hard Probe	1 Each
675.011	Mono-axial Screwdriver (short)	1 Each
675.012	Multi-axial Screwdriver (short)	1 Each
675.013	Screw Holder	1 Each
675.015	Rod Holder	1 Each
675.016	Rod Rotator	1 Each

Instruments

675.017	Rod Bender	1 Each
675.018	Bending Device (left and Right)	1 Each
675.019	Rod Pusher	1 Each
675.021	Spreader	1 Each
675.022	Compressor	1 Each
675.023	Anti-Rotating Sleeve	1 Each
675.024	Rod Compressor	1 Each
675.025	Breaker	1 Each
675.027	Wrench (2 pcs)	1 Each
675.028	Mono-axial Screwdriver (long)	1 Each
675.029	Multi-axial Screwdriver (long)	1 Each
675.030	Small Wrench ϕ 4.0	1 Each
675.031	Template (2 pcs)	1 Each
675.032	Rod Persuader	1 Each
675.033	Rotating Device SW5.0	1 Each
675.034	Rod Fork	1 Each
675.035	Screwdriver SW3.5	1 Each
675.036	Hex Screwdriver SW3.5	1 Each



• **T413** | Cervical Plate

MATERIAL - Titanium
LENGTH - 21.5mm to 23mm (1.5mm Increment)
25mm to 70mm (2.5mm Increment)



• **T414** | Cervical Bone Screw

MATERIAL - Titanium
DIAMETER - 4mm
LENGTH - 12mm to 22mm (2mm Increment)



• **T415** | Cervical Locking Screw

MATERIAL - Titanium

Instruments



◦ 677.001 | Drill Bit with handle 2.5mm for Cervical Bone Screw



◦ 677.002 | Tap with Handle 4.0mm for Cervical Bone Screw



◦ 677.003 | Cervical Drill Guide Sleeve
◦ 677.004 | Cervical Tap Guide Sleeve



◦ 677.005 | Self Holding Hexagonal Screw Driver 2.5mm Tip



◦ 677.006 | Lock Screw Driver



◦ 677.007 | Cervical Plate Holder



◦ 677.008 | Spacer Spanner



◦ 677.009 | Spacer 12mm
◦ 677.010 | Spacer 14mm
◦ 677.011 | Spacer 16mm
◦ 677.012 | Spacer 18mm
◦ 677.013 | Spacer 20mm
◦ 677.014 | Spacer 22mm



◦ 677.015 | Plate Bender Pair

PHOTO NOT AVAILABLE

Implants

T413	Cervical Plate Length 21.5mm, 23mm, 25mm, 27.5mm, 30mm, 32.5mm, 35mm, 37.5mm, 40mm, 42.5mm 45mm, 47.5mm, 50mm, 52.5mm, 55mm, 57.5mm, 60mm, 62.5mm, 65mm, 67.5mm, 70mm	1 Each
T414	Cervical Bone Screw 4.0mm x 12mm, 14mm, 16mm, 18mm, 20mm, 22mm	4 Each
T415	Cervical Locking Screw	4 Each

Instruments

677.001	Drill Bit with handle 2.5mm for Cervical Bone Screw	1 Each
677.002	Tap with Handle 4.0mm for Cervical Bone Screw	1 Each
677.003	Cervical Drill Guide Sleeve	1 Each
677.004	Cervical Tap Guide Sleeve	1 Each
677.005	Self Holding Hexagonal Screw Driver 2.5mm Tip	1 Each
677.006	Lock Screw Driver	1 Each
677.007	Cervical Plate Holder	1 Each
677.008	Spacer Spanner	1 Each
677	Spacer 12mm, 14mm, 16mm, 18mm, 20mm, 22mm	1 Each
677.015	Plate Bender Pair	1 Each



• **T418** | Tlif Cage

MATERIAL - Titanium
DEGREE - 4°, 8°, 11°
DIAMETER - 07mm to 11mm (1mm Increment)



• **T419** | Plif Cage

MATERIAL - Titanium
LENGTH - 20mm, 25mm
DIAMETER - 07mm to 11mm (1mm Increment)



• **T420** | Banana Cage

MATERIAL - Titanium
LENGTH - 25mm, 30mm
DIAMETER - 8mm, 9mm, 10mm, 11mm, 12mm



• **T416** | Spinal Eco Cage

MATERIAL - Titanium
LENGTH - 12mm, 14mm, 16mm, 18mm
DIAMETER - 20mm to 50mm (5mm Increment)

Instruments



- 680.021 | Chisel 7mm
- 680.022 | Chisel 8mm
- 680.023 | Chisel 9mm
- 680.024 | Chisel 10mm
- 680.025 | Chisel 11mm



- 680.027 | Box Chisel 7mm
- 680.028 | Box Chisel 8mm
- 680.029 | Box Chisel 9mm
- 680.030 | Box Chisel 10mm
- 680.031 | Box Chisel 11mm



- 680.032 | T Handle



- 680.033 | Box Currete (Left)
- 680.034 | Box Currete (Right)



- 680.035 | Cup Currete (Left)
- 680.036 | Cup Currete (Right)



- 680.037 | Plif Cage Introducer



- 680.038 | Banana (Kidney Cage) Introducer



- 680.039 | Cup Currete Straight



- 680.040 | Box Currete Straight

PHOTO NOT AVAILABLE

Implants

T418	Tlif Cage 4° x 7mm, 8mm, 9mm, 10mm, 11mm	1 Each
T418	Tlif Cage 8° x 7mm, 8mm, 9mm, 10mm, 11mm	1 Each
T418	Tlif Cage 11° x 7mm, 8mm, 9mm, 10mm, 11mm	1 Each
T419	Plif Cage 20mm x 7mm, 8mm, 9mm, 10mm, 11mm	1 Each
T419	Plif Cage 25mm x 7mm, 8mm, 9mm, 10mm, 11mm	1 Each
T420	Banana Cage 8mm, 9mm, 10mm 11mm, 12mm x 25mm	1 Each
T420	Banana Cage 8mm, 9mm, 10mm, 11mm, 12mm x 30mm	1 Each

Instruments

680.021	Chisel 7.0mm	1 Each
680.022	Chisel 8.0mm	1 Each
680.023	Chisel 9.0mm	1 Each
680.024	Chisel 10.0mm	1 Each
680.025	Chisel 11.0mm	1 Each
680.027	Box Chisel 7.0mm	1 Each
680.028	Box Chisel 8.0mm	1 Each
680.029	Box Chisel 9.0mm	1 Each
680.030	Box Chisel 10.0mm	1 Each
680.031	Box Chisel 11.0mm	1 Each
680.032	T Handle	1 Each
680.033	Box Currete Left	1 Each
680.034	Box Currete Right	1 Each
680.035	Cup Currete Left	1 Each
680.036	Cup Currete Right	1 Each
680.037	Plif Cage Introducer	1 Each
680.038	Banana (Kidney Cage) Introducer	1 Each
680.039	Cup Currete Straight	1 Each
680.040	Box Currete Straight	1 Each



Plot No. G-2041/42/A, Kishan Gate, Lodhika G.I.D.C.,
Kalawad Road, Metoda, RAJKOT-360 021, (INDIA).

info@hardikinternational.com | www.hardikinternational.com

☎ : +91 95740 01818 | 📞 : hardikinternational

Follow Us :     