

ENDOVASCULAR

Wide range of minimally invasive diagnostics devices dedicated to endovascular procedures enabling treatment of peripheral artery disease.



PERIPHERAL BALLOON CATHETERS

PERIPHERAL ANGIOPLASTY CATHETER LOVIX™

with OTW delivery system (Over The Wire)



Balloon catheter is an effective option for peripheral artery lesions

- High pressure
- Full range of sizes
- A low profile
- 5F Improved shaft

Technical data overview

Balloon material	nylon
Balloon diameter	1.50 - 16.00 mm
Balloon length	20-250 mm
Nominal pressure (for balloon diameters 5.00-7.00 mm)	8 atm
Nominal pressure (for balloon diameters 8.00-16.00 mm)	6 atm

Deflated time	short
Pushability and trackability	excellent
Delivery system	Over The Wire (OTW)
Delivery system lengths	80 cm, 100 cm, 120 cm, 135 cm, 170 cm



ENDOVASCULAR

20

PERIPHERAL BALLOON CATHETERS LOVIX™ - PERIPHERAL ANGIOPLASTY CATHETER-PTA

	Nominal balloon diameter [mm]	Nominal balloon length [mm]	Use with guidewire	Use with introducer [F]	Catheter length [cm]	Nominal pressure [atm]
Reference:	В	С	E			
PBA B x CE 110	1.5	20, 25, 30, 35, 40, 45, 50, 55, 60, 65, 70, 75, 80, 85, 90, 95, 100 120, 150, 200, 250	0.018"	4	110	10
			0.025"	4		
			0.025"	4		
PBA B x CE 110	2	20, 25, 30, 35, 40, 45, 50, 55, 60, 65, 70, 75, 80, 85, 90, 95, 100	0.018"	4	110	10
			0.025"	4		
		120,150, 200, 250	0.025"	4		
PBA B x CE 110	2.5	20, 25, 30, 35, 40, 45, 50, 55, 60, 65, 70, 75, 80, 85, 90, 95, 100	0.018"	4	110	10
			0.025"	4		
		120, 150, 200, 250	0.025"	4		
PBA B x CE 110	3	20, 25, 30, 35, 40, 45, 50, 55, 60, 65, 70, 75, 80, 85, 90, 95, 100	0.018"	4	110	10
			0.025"	4		
		120, 150, 200, 250	0.025"	4		
PBA B x CE 110	4	20, 25, 30, 35, 40, 45, 50, 55, 60, 65, 70, 75, 80, 85, 90, 95, 100	0.018"	5	110	10
			0.025"	5		
			0.035"	5		
		120, 150, 200, 250	0.035"	5		
PBA B x CE 110	5	20, 25, 30, 35, 40, 45, 50, 55, 60, 65, 70, 75, 80, 85, 90, 95, 100, 120	0.025"	5	110	8
			0.035"	5		
		150, 200, 250	0.035"	6		
PBA B x CE 110	6	20, 25, 30, 35, 40, 45, 50, 55, 60, 65, 70, 75, 80, 85, 90, 95, 100, 120	0.025"	5	110	8
			0.035"	5		
		150, 200, 250	0.035"	6		
PBA B x CE 110	7	20, 25, 30, 35, 40, 45, 50, 55, 60, 65, 70, 75, 80, 85, 90, 95,100,	0.035"	5	110	8
		120, 150,200, 250	0.035"	6		
PBA B x CE 110	8		0.035"	5	110	6
PBA B x CE 110	9		0.035"	6	110	6
PBA B x CE 110	10	20, 25, 30, 35, 40, 45, 50, 55, 60, 65, 70, 75, 80, 85, 90, 95, 100	0.035"	6	110	6
PBA B x CE 110	12		0.035"	7	110	6
PBABxCE110	14		0.035"	8	110	6
PBABxCE110	16		0.035"	8	110	6

When ordering use Ref. Number: PBA16x10035110.

NOTE: Catheters 80, 100, 120, 135 and 170 cm in length are also available for order.

When ordering the catheter of other length, please substitute the length of 110 in the Ref. Number with the desired length e.g. PBA10x1003580

