


RIG

DATA SHEET PAGE 1:2

TRADING NAME	HAPPY COVER® RIG
DESCRIPTION	Polyethylene vomit bag, for an easy and safe seal of its opened upper side, after use. Inside the bag there is a superabsorbent molecule-powder named SGAIEL® that transforms liquids into solids. This powder is enclosed in a water-soluble film which acts in contact with liquids. Ready for use.
TEM CODE	AM003
RDM (Italy)	1570463
CND (Italy)	A060399
PACKAGING	pcs 20
DIMENSIONS	mm. 250x250
EXPIRING TIME	5 years
BELONGING CLASS OF MD	Classe I medical device - Not sterile Pursuant to 93/42/EEC
CE MARK	Pursuant to 93/42/EEC and 2007/47 EEC of 09/02/2007
INTENDED USE	Disposable vomit system for use in public and private hospitals, first aid, Emergency and transport companies
RAW MATERIALS	Polyethylene HDPE. White cross-linked sodium polyacrylate enclosed in a water-soluble film.



HAPPYcover®

RIG

DATA SHEET PAGE 2:2

INSTRUCTIONS FOR USE	Happy Cover® RIG is not used in combination with other device. SGAIEL® powder turns vomit into a gelatinous substance in about one minute. After use close the bag by adhesive liner or make a knot. Dispose as per current legislation.
CAUTION	Do not use sharp objects to open the package Do not disperse in the environment after use. Do not use if package is open or damaged.
CONTRAINDICATIONS	No contraindications by using these medical devices.
STORAGE	Store in a cool, dry place away from heat sources and sunlight.
TRANSPORTATION	Not regulated by EEC and Italia Transport Department as dangerous product. Terrestrial transport (RID/ADR): not dangerous. Maritime transport (IMG): not dangerous. Airplane transport (ICAD/IATA): not dangerous.
APPLICABLE REGULATION	93/42/EEC of 06/14/1993 2007/47 EEC of 09/05/2007
GMDN	Code GMDN: 3923.21 Description: Polyethylene sick bag with a folding and internal double-sided adhesive pocket that allows an easy and safe sealing of its open upper side after use. Its interior is covered with a molecule-powder named "SGAIEL®" with an absorbing function that turns liquids into solids. This powder is enclosed in a water-soluble film which acts upon contact with liquid. Ready for use.
	Angimed S.r.l. Corso Umberto I, 590 ed.66 - 65015 Montesilvano (PE)

CE

CATHETER GEL

KATETER JEL / GEL POUR CATHÉTER / GEL DE CATÉTER



EN

Properties & Use

Konix Sterile Catheter Gel is a sterile lidocaine contained, clear and water-soluble sort of gel which is used before applying catheter, cystoscopy and/or else medical instrument to urethra. Covering urethra by generating a slick layer between urethral mucosa and catheter or else medical instrument is the most significant function of the gel. Catheter Gel, including Lidocaine, inside, is used before the application in order to replace urethra for the purpose of dulling the pain related to urethral manipulation. Moreover, it provides an indolent catheterization with its anaesthetic effect. Antiseptic effect of the gel insulates patients from the potential infections that may occur within the upper part and urinary bladder due to iatrogenic contamination.

Storage: Store in a cool and dry environment between 5–30 °C away from direct sunlight.

FR

Propriétés & Utilisation

Cathéter stérile est une lidocaïne stérile contenant du gel clair et soluble dans l'eau qui est utilisé avant l'application de cathéter ou d'autres instruments médicaux sur l'urètre, pour des examens cystoscopiques ou autres. Recouvrir l'urètre en générant une couche lisse entre la muqueuse urétrale et le cathéter ou d'autres instruments médicaux c'est est la plus importante fonction du gel. Konix Gel pour cathéter stérile comprenant de lidocaïne, est utilisé avant l'intervention afin de protéger l'urètre pour atténuer la douleur liée à la manipulation urétrale. En outre, la pose du cathéter devient indolore grâce à son effet anesthésique. L'effet antiseptique du gel protège les patients d'une infection éventuelle qui pourrait se déclarer dans la partie supérieure et de la vessie et due à la contamination iatrogène.

Stockage: Conserver dans un endroit frais et sec entre 5 - 30 °C à l'abri des rayons soleil.

TR

Özellikler & Kullanım Şekli

Lidokain içeren Konix Steril Kateter Jeli, üretraya kateter, sistoskopi ve/veya başka bir tıbbi alet uygulamadan önce kullanılan steril, berrak, suda çözünen kaydırıcı bir jeldir. Jelin en önemli fonksiyonu kateter veya başka bir tıbbi alet ile üretral mukoza arasında kaydırıcı bir tabaka oluşturarak üretrayı kaplamasıdır. Lidokain içeren Konix Steril Kateter Jelin uygulamasından önce üretral manipülasyonla bağlantılı olan ağrının giderilmesinde yardımcı olmak amacıyla üretranın kaydırılması için kullanılır. Ayrıca anestezi etkisi ile ağrısız bir kateterizasyon sağlar. Konix Steril Kateter Jeli'nin antiseptik etkisi iatrojenik kontaminasyondan dolayı üst bölümü ve mesanede oluşabilecek enfeksiyonlardan hastayı korur.

Saklama Koşulları: 5-30°C'nin arasında ve güneş ışığından uzakta, serin bir yerde muhafaza ediniz.

ES

Características y Modo de Empleo

Catéter contiene lidocaína y es un gel lubricante de tipo estéril, transparente y soluble en agua, que se utiliza antes de aplicar el catéter, la cistoscopia y/o otro dispositivo médico a la uretra. La función más significativa del gel es cubrir la uretra mediante la generación de una capa resbaladiza entre la mucosa de la uretra y el catéter u otro dispositivo médico. Konix Gel Estéril de Catéter que contiene lidocaína se utiliza antes, para remplazar la uretra con el fin de evitar el dolor relacionado con la su manipulación. Además, proporciona un cateterismo indolente gracias a su efecto anestésico. El efecto antiséptico de Konix Gel Estéril de Catéter protege a los pacientes de las posibles infecciones que pueden originarse dentro de la parte superior de la vejiga y que es debida a la contaminación iatrogénica.

Almacenamiento: Conserver en ambiente seco y fresco entre 5-30° C fuera del alcance directo de la luz solar.

MEDICAL PACKAGING

BOLSAPACK • PACKAGING FOR STEAM/EO/FORMOL

Pouches and rolls, manufactured with 60 g/m² medical paper and a plastic film (blue coloured). Indicators for Steam/EO/Formol are also printed in packaging.

PAPER/FILM FLAT BAG


Code	Sizes		U/Box	Box/cm	Box/Kg
MBO1025	75 x 150 mm	3 x 6"	2.700	17 x 25 x 17	3,2
MBO1000	75 x 210 mm	3 x 8"	2.500	31 x 29 x 11	4,2
MBO1001	75 x 280 mm	3 x 11"	2.400	31 x 29 x 11	5,7
MBO1003	100 x 200 mm	4 x 8"	2.400	41 x 29 x 11	5,7
MBO1004	100 x 240 mm	4 x 9"	2.400	41 x 29 x 11	6,7
MBO1005	100 x 280 mm	4 x 11"	2.400	41 x 29 x 11	7,7
MBO1026	100 x 300 mm	4 x 12"	2.500	41 x 43 x 11	8,7
MBO1006	100 x 420 mm	4 x 17"	2.400	41 x 43 x 11	11,7
MBO1022	100 x 580 mm	4 x 23"	2.000	41 x 59 x 11	13,7
MBO1007	150 x 210 mm	6 x 8"	2.000	32 x 43 x 11	7,2
MBO1008	150 x 280 mm	6 x 11"	2.000	31 x 57 x 11	9,7
MBO1009	150 x 320 mm	6 x 13"	2.000	37 x 65 x 11	11,2
MBO1010	150 x 360 mm	6 x 14"	2.000	37 x 65 x 11	12,7
MBO1021	150 x 400 mm	6 x 16"	1.200	32 x 43 x 11	8,2
MBO1011	200 x 320 mm	8 x 13"	1.800	37 x 65 x 11	13,2
MBO1012	200 x 360 mm	8 x 14"	1.800	37 x 65 x 11	14,7
MBO1013	200 x 420 mm	8 x 17"	1.200	41 x 43 x 11	11,2
MBO1023	200 x 500 mm	8 x 20"	1.200	41 x 59 x 11	14,5
MBO1014	250 x 360 mm	10 x 14"	1.200	51 x 51 x 11	12,7
MBO1015	250 x 420 mm	10 x 17"	1.200	51 x 51 x 11	15,2
MBO1016	250 x 500 mm	10 x 20"	1.200	51 x 51 x 11	17,7
MBO1017	300 x 420 mm	12 x 17"	600	32 x 43 x 11	9,2
MBO1018	300 x 500 mm	12 x 20"	600	31 x 57 x 11	10,2
MBO1024	300 x 580 mm	12 x 23"	700	41 x 59 x 11	13,9
MBO1019	400 x 580 mm	16 x 23"	600	41 x 59 x 11	15,7



Paper/Film flat bag



Paper/Filmgusset bag



Paper/Film flat roll

PAPER/FILM GUSSET BAG

Code	Sizes		U/Box	Box/cm	Box/Kg
MBF1202	200 x 55 x 400 mm	8 x 2" x 16"	600	41 x 43 x 11	6,7

MEDICAL PACKAGING

BOLSAPACK • PACKAGING FOR STEAM/EO/FORMOL

Pouches and rolls, manufactured with 60 g/m² medical paper and a plastic film (blue coloured). Indicators for Steam/EO/Formol are also printed in packaging.

PAPER/FILM SELFSEALING BAG


Code	Sizes		U/Case	U/Box	Box/cm	Box/Kg
MBA1080	75 x 200 mm	3 x 8"	200	1.200	17 x 25 x 17	2,4
MBA1081	90 x 200 mm	4 x 8 "	200	1.200	22 x 28 x 17	3,2
MBA1082	90 x 230 mm	4 x 9"	200	1.200	22 x 28 x 17	4,3
MBA1084	100 x 240 mm	4 x 9"	200	1.200	23 x 29 x 17	4,2
MBA1085	130 x 250 mm	5 x 10"	200	1.200	28 x 30 x 17	5,7
MBA1086	130 x 330 mm	5 x 13"	200	1.200	28 x 38 x 17	7,3
MBA1087	140 x 280 mm	6 x 11"	200	1.200	30 x 33 x 17	6,7
MBA1088	200 x 330 mm	8 x 13"	200	800	42 x 38 x 12	7,2
MBA1089	250 x 400 mm	10 x 16"	200	800	52 x 45 x 12	10,7
MBA1090	300 x 390 mm	12 x 15"	***	400	32 x 43 x 11	6,6
MBA1091	300 x 450 mm	12 x 18"	***	400	31 x 50 x 11	6,7



Paper/Film selfsealing bag



Selfsealing tape



Paper/Film gusset roll

PAPER/FILM FLAT ROLL

Code	Sizes		U/Box	Box/cm	Box/Kg
MRP1040	50 mm x 200 m	2" x 656 Feet	6	31 x 25 x 23	7,2
MRP1041	75 mm x 200 m	3" x 656 Feet	4	31 x 25 x 23	7,2
MRP1042	100 mm x 200 m	4" x 656 Feet	4	43 x 26 x 23	9,7
MRP1043	125 mm x 200 m	5" x 656 Feet	4	44 x 32 x 23	11,7
MRP1044	150 mm x 200 m	6" x 656 Feet	4	44 x 32 x 23	13,7
MRP1045	200 mm x 200 m	8" x 656 Feet	2	43 x 26 x 23	9,2
MRP1046	250 mm x 200 m	10" x 656 Feet	2	43 x 26 x 23	11,2
MRP1047	300 mm x 200 m	12" x 656 Feet	2	44 x 32 x 23	13,7
MRP1050	350 mm x 200 m	14" x 656 Feet	1	41 x 21 x 22	8,2
MRP1048	400 mm x 200 m	16" x 656 Feet	1	41 x 21 x 22	9,2
MRP10491	500 mm x 100 m	20" x 328 Feet	1	52 x 18 x 18	5,7

PAPER/FILM GUSSET ROLL

Code	Sizes		U/Box	Box/cm	Box/Kg
MRF1059	75 x 25 mm x 100 m	3 x 1" x 328 Feet	4	31 x 25 x 23	4,7
MRF1060	100 x 40 mm x 100 m	4 x 2" x 328 Feet	4	43 x 26 x 23	6,7
MRF1061	150 x 40 mm x 100 m	6 x 2" x 328 Feet	4	44 x 32 x 23	9,2
MRF1062	200 x 60 mm x 100 m	8 x 2" x 328 Feet	2	43 x 26 x 23	6,2
MRF1063	250 x 60 mm x 100 m	8 x 2" x 328 Feet	2	43 x 26 x 23	7,2
MRF1064	300 x 80 mm x 100 m	8 x 2" x 328 Feet	1	31 x 25 x 23	4,7
MRF1065	400 x 80 mm x 100 m	16 x 3" x 328 Feet	1	41 x 21 x 22	5,7

2-Way Balloon Catheters



Silicone Elastomer

100% Silicone

2-Way Balloon Catheter

Longlife Silicone Elastomer

Nelaton tip

Art. No.	Size (Fr/CH)	Balloon (cc/ml)	Length (cm)
9385S-0603	6	3	20
9385S-0805	8	5	20
9385S-1005	10	5	20
9385S-1210	12	5 to 15	40
9385S-1410	14	5 to 15	40
9385S-1610	16	5 to 15	40
9385S-1810	18	5 to 15	40
9385S-2010	20	5 to 15	40
9385S-2210	22	5 to 15	40
9385S-2410	24	5 to 15	40
9385S-2610	26	5 to 15	40

Tiemann tip

Art. No.	Size (Fr/CH)	Balloon (cc/ml)	Length (cm)
9310S-1210	12	5 to 15	40
9310S-1410	14	5 to 15	40
9310S-1610	16	5 to 15	40
9310S-1810	18	5 to 15	40
9310S-2010	20	5 to 15	40

Suprapubic tip

9399S-12	12	5 to 15	40
9399S-14	14	5 to 15	40
9399S-16	16	5 to 15	40
9399S-18	18	5 to 15	40
9399S-20	20	5 to 15	40
9399S-22	22	5 to 15	40
9399S-24	24	5 to 15	40

100% Silicone Nelaton

Nelaton tip

9385VS-1210	12	10	40
9385VS-1410	14	10	40
9385VS-1610	16	10	40
9385VS-1810	18	10	40
9385VS-2010	20	10	40
9385VS-2210	22	10	40
9385VS-2410	24	10	40

100% Silicone Suprapubic

Suprapubic tip

9399VS-12	12	10	40
9399VS-14	14	10	40
9399VS-16	16	10	40
9399VS-18	18	10	40
9399VS-20	20	10	40
9399VS-22	22	10	40
9399VS-24	24	10	40