



FESTERMED Twain

- For the use of the patients carried by ambulance and first aid.
- Brings the two wrist together and enables the hands at the front of the body .
- Skin touch cotton surface absorbs the sweat.



- In case of need with the velcro lock at the bottom surface can be stable by fixing to the bad or litter.
- In case of long term use does not cause allergy and irritation.

CE



FESTERMED Classic

- For the use of inpatient both bed or litter.
- Skin touch cotton surface absorbs the sweat.



- In case of long term use does not cause allergy and irritation.
- Gamma sterilized.

CE₁₇₈₃

Festermmed Ordering Information

Model	Size (cm)	Piece / Pouch	Pouch/ Box	Box Size (cm)	Box Weight	Order Code
Deluxe	34	1	25	20 x 30 x 20	1215 gr.	ALH 2006
Deluxe / Çift	34	2	12	20 x 30 x 20	1175 gr.	ALH 2007
Twain	32	1	25	20 x 30 x 20	1275 gr.	TC 2009
QRB	27	1	25	20 x 30 x 20	1790 gr.	QRB 2013
QRB / Çift	27	2	12	20 x 30 x 20	1726 gr.	QRB 2014
Klasik El	25	1	25	20 x 30 x 20	850 gr.	ALH 2001
Klasik El / Çift	25	2	12	20 x 30 x 20	815 gr.	ALH 2008
Klasik Ayak	40	1	25	20 x 30 x 20	910 gr.	ALH 2002
Klasik Ayak / Çift	40	2	12	20 x 30 x 20	885 gr.	ALH 2010
Klasik Bacak	65	1	25	20 x 30 x 20	1100 gr.	ALH 2003
Metal El	25	1	25	20 x 30 x 20	1215 gr.	ALH 2011
Metal El / Çift	25	2	12	20 x 30 x 20	1130 gr.	ALH 2012