

## Hydrophilic IOL

\_Excellent quality of raw materials manufactured by



low inflammatory response and optimal biocompatibily.

\_ldeal handling characteristics with controlled unfolding within the capsular bags

\_Shape memory techology

\_Maximal, axial radial and rotational stability

\_The unmatched material trans-parency offers the patient clear post surgical vision and optimal optical outcome, providing permanent optic quality & improved safety

\_Easily inserted through a 2.0 mm incision due to its flexibility and plasticity

\_Accompanied by cartridges manufactured in Europe

\_Wide diopter range; -10.0 to +40.0 D, in half diopter steps (out-of-range diopters available upon request)

**Products are CE-certified** 

\_Company is certified for implementing a Quality Management System according to ISO 9001:2008 & ISO 13485:2012 standards



21, Stageiriti & 24, Em.Fili str. Thessaloniki, Gr-54352, Greece T:+30 2310906660 www.sidapharm.gr