

Enteral Feeding Bag Set, PuriBag[®], Twin Bag, Pump Use



ENFit[®] Female
Connector
CE802



Medication
Y Port with Male ENFit[®]
Connector
CE805



Ucomfor[®]
Ucomfor Scientific Co., Ltd.

Enteral Feeding Products, PVC, DEHP free

Enteral Feeding Bag Set, PuriBag[®], Twin Bag, Pump Use



Feeding Bag is used mainly for nutrition transportation to patient via connecting the stomach tube using a pump. This device is sterile and for single use only, the following precautions shall be observed during use:

1. Do not use when visual damage on its packaging is detected.
2. Only for single use, and should be discarded after use.
3. Use before expiration date.

When the individual packaging is opened, it shall be used immediately.

Features:

- Double Bag design
- Drip Chamber
- Precision Roller Clamp
- ENFit[®] Connector with Vented Cap
- EO Sterile
- Silicone Tube
- Compatible with Kangaroo 324 Feeding Pump

Item No.	Volume (mL)	Tube Size (mm)	Ice Bag	Medication Port
F01033	1000	4.2	✓	✓
F01533	1500	4.2	✓	✓

Активация Windows
Чтобы активировать Windows,