

TIP.DOC. 1

DATA EMITERII:17 septembrie 2021:

PLATITI: 91000-00 LEI: Nouazeci si Unu Mii lei 00 ban
i

PLATITOR: (R) "BIOSISTEM MLD" S.R.L. CONTUL DE PLATI/CODUL IBAN MD95ML00000002251429243 CODUL FISCAL :1010600028048 /

PRESTATORUL PLATITOR	CODUL BANCII:
BC"Moldindconbank"S.A. fil."Invest" Chisinau	:MOLDMD2X329:

BENEFICIAR (R) Centrul pen	CONTUL DE PLATI/CODUL IBAN
tru achizitii publice central	MD23TRPCCC518430B01859AA
izate in sanatate	CODUL FISCAL :1016601000212 /

PRESTATORUL BENEFICIAR	CODUL BANCII:
Ministerul Finantelor - Trezoreria de Stat	:TREZMD2X :

DESTINATIA PLATII:/P102/91000,00 Pentru :	TIPUL TRANSFERULUI :
garantia pentru oferta la procedura de a:	NORMAL/URGENT :N
chizi?ie publica nr. ocds-b3wdp1-MD-1628:	:
866063582 din 18.09.2021	:

L.S.

CODUL TRANZACTIEI:101:

DATA PRIMIRII:17/09/2021

DATA EXECUTARII:

SEMNATURILE

EMITENTULUI

CONDUCATOR:Web Poiata Vitalie

MIIGYwYJKoZiHvcNAQcCoIIGVDCCB1ACAQExCzAJBqUrDgMCGqUAMAsGCSqGSib:

DQEHAaCCBGwwggRoMIIDUKADAgECAhNHAACjb1r1rqFksQ0G4AAAAAKNuMA0GCSq:

SIb3DOEBCwUAMCIxIDAeBqNVBAMTF0NFU10xLUNBLU1vbGRpbmRjb25iYW5rMB4:

DTIxMDEyODExMzgWNVoXDTI0MDEyODExNDgwNVowgZ8xCzAJBgNVBAYTAkEMRA: qYDV00IEwdNb2xkb3ZhMREwDwYDV00HEwhDaGLzaw5hdTEWMB0GA1UEChMNQmI

(semnatura electronica)

CONTABIL-SEF:Web Nasedchin Alexandr

MIIGZwYJKoZiHvcNAQcCoIIGWDCCBlQCAQExCzAJBgUrDgMCGGUAMAsGCSqGSIb3:

DQEHAaCCBHAwgqRsMIIDVKADAgECAhNHAACjcahRKqbJeq8QAAAAAKNxMA0GCSqG:

SIb3DQEBCwUAMCIxIDAeBgNVBAMTF0NFU1QxLUNBLU1vbGRpbmRjb25iYW5rMB4X:

DTIxMDEyODExMzkxOFoXDTI0MDEyODExNDkxOFowgaMxCzAJBgNVBAYTAk1EMRAw:

YDVQQIEwdNb2xkb3ZhMREwDwYDVQQHEwhDaG1zaW5hdTEWMBQGA1UEChMNQm1v

(semnatura electronica)

CONDUCTOR:

(semnatura manuala)

CONTABIL-SEF:

(semnatura manuala)

SEMNATURA PRESTATORUL

L.S.

MOTIVUL REFUZULUI

L.S.



BC "MOLDINDCONBANK" S.A.

Filiala "Invest"

Republica Moldova, MD-2068
mun. Chișinău, bd. Moscovei, 14/1
Tel. : (373-22) 43-44-81, 43-46-24
Fax : (373-22) 43-44-22
cod: MOLDMD2X329

Data 14. IAN. 2016
Nr. 03/2 - 19/23

Республика Молдова, MD-2068
мун. Кишинэу, бул. Московской, 14/1
Тел. : (373-22) 43-44-81, 43-46-24
Факс : (373-22) 43-44-22
код: MOLDMD2X329

Filiala „Invest” BC „Moldindconbank” SA confirmă existența contului curent în moneda națională al **“BIOSISTEM MLD” S.R.L. (c/f 1010600028048)**, cu **IBAN MD95ML000000002251429243**.

Codul băncii MOLDMD2X329.

Director

Nina Turcan

Director financiar



Nina Balmuş

Ex. Diana Brinza
Tel. 43-45-96

REPUBLICA



MOLDOVA

CERTIFICAT DE ÎNREGISTRARE

Societatea cu Răspundere Limitată "BIOSISTEM MLD"
— ESTE ÎNREGISTRATĂ LA CAMERA ÎNREGISTRĂRII DE STAT —

Numărul de identificare de stat - codul fiscal
1010600028048

Data înregistrării

12.08.2010

Data eliberării

12.08.2010

Svirepova Ludmila, registrator

*Funcția, numele, prenumele persoanei
care a eliberat certificatul*

L. Svirepova
semnătura

MD 0101250





I.P. "AGENȚIA SERVICII PUBLICE"

Departamentul înregistrare și licențiere a unităților de drept

EXTRAS

din Registrul de stat al persoanelor juridice

nr. 8506 din 28.04.2021

Denumirea completă: **Societatea cu Răspundere Limitată «BIOSISTEM MLD».**

Denumirea prescurtată: **«BIOSISTEM MLD» S.R.L.**

Forma juridică de organizare: **Societate cu Răspundere Limitată.**

Numărul de identificare de stat și codul fiscal: **1010600028048.**

Data înregistrării de stat: **12.08.2010.**

Sediul: **MD-2001, str. Albișoara, 16/1, ap.(of.) 7, mun. Chișinău, Republica Moldova.**

Obiectul principal de activitate:

- 1 Activitatea farmaceutică;**
- 2 Importul, fabricarea, comercializarea, asistența tehnică și (sau) reparația dispozitivelor medicale și (sau) a opticii;**
- 3 Acordarea asistenței medicale de către instituțiile medico-sanitare private;**
- 4 Comerțul cu ridicata al calculatoarelor, echipamentelor periferice și software-ului;**
- 5 Întreținerea și repararea mașinilor de birou și a tehnicii de calcul;**
- 6 Consultații în domeniul sistemelor de calcul.**

Capitalul social: **5400 lei.**

Administrator: POIATA VITALIE,

Asociați:

- 1. POIATA VITALIE 33,40 %**
- 2. NASEDCHIN ALEXANDR 33,30 %**
- 3. KOJEVNIKOV DMITRII 33,30 %.**

Prezentul extras este eliberat în temeiul art. 34 al Legii nr. 220-XVI din 19 octombrie 2007 privind înregistrarea de stat a persoanelor juridice și a întreprinzătorilor individuali și confirmă datele din Registrul de stat la data de: 28.04.2021.

Specialist coordonator
tel. 022-207-840



Lazari Aliona



EB 0358735

Lista fondatorilor Biosistem-mld SRL

Nr.	Nume, Prenume	IDNP
1.	Vitalie Poiata	0983103892591
2.	Alexandr Nasedchin	2002001070747
3.	Dmitrii Kojevnikov	0972305012362

CERTIFICAT
privind lipsa sau existența restanțelor față de bugetul public național

Nr.
№ **A2115192**

din
от **08.09.2021**

1. Destinația / Назначение

Pentru participarea la proceduri de achizitii publice

2. Date despre contribuabil / Информация о налогоплательщике

Denumirea Наименование	Codul fiscal / Numărul de identificare Фискальный код / Идентификационный номер
BIOSISTEM MLD S.R.L.	1010600028048
Adresa sediului de bază (strada, numărul) Адрес основного месторасположения (улица, номер)	Codul - Denumirea localității Код - Наименование населенного пункта
Albisoara nr.16 bl.1 of.7	0150-SEC.RISCANI

3. Atestarea lipsei sau existenței restanțelor conform datelor Sistemului Informațional Automatizat /
Подтверждение отсутствия или наличия недоимки согласно данных Информационной автоматизированной системы

La data emiterii prezentului certificat restanța față de bugetul public național constituie/ На дату выдачи данной справки недоимка перед национальным публичным бюджетом составляет:
0,00 lei/лей.

4. Valabil pînă la / Действителен до 23.09.2021

5. Autentificarea Serviciului Fiscal de Stat / Подтверждение Государственной налоговой службы



Semnătura/Подпись

Viorica CĂUȘ

Numele și prenumele/Фамилия и имя

Este extras din Sistemul Informațional al SFS SIA „Contul curent al contribuabilului”// 08.09.2021 ora 14:25:51
cu aplicarea prevederilor pct. 82-83 Ordin IFPS nr.400 din 14.03.2014 (Monitorul Oficial 72-77/399, 28.03.2014)

NOTA (3,52)

Дата предоставления 11.05.2021 10:00:47

Anexe la SNC
 "Prezentarea situatiilor financiare"
 Aprobat de Ministerul Finantelor
 al Republicii Moldova

SITUAȚIILE FINANCIARE

pentru perioada 01.01.2020 - 31.12.2020

Entitatea: BIOSISTEM MLD S.R.L.
 Cod CUIIO: 40717392
 Cod IDNO: 1010600028048

Sediuul:
 MD:
 Raionul(municipiul): 106, DOF, RISCANI
 Cod CUATM: 0150, SEC, RISCANI
 Strada: SECTORUL RISCANI STR.Albisoara nr.16 bl.1 of.7

Activitatea principală: 08646, Comerț cu ridicata al produselor farmaceutice
 Forma de proprietate: 16, Proprietate colectivă
 Forma organizatorico-juridică: 530, Societăți cu răspundere limitată

Date de contact:
 Telefon: +37322808719
 WEB:
 E-mail: zmil13@gmail.ru
 Numele și coordonatele al contabilului-șef: DI (dna) Tel.

Numărul mediu al salariaților în perioada de gestiune: 3 persoane.

Persoanele responsabile de semnarea situațiilor financiare* Nasedchin Alexandr

Unitatea de măsură: leu

BILANȚUL

la

Anexa 1

Nr. cpt.	Indicatori	Cod rd.	Sold la	
			Începutul perioadei de gestiune	Sfîrșitul perioadei de gestiune
1	2	3	4	5
	A C T I V			
	ACTIVE IMOBILIZATE			
	I. Imobilizări necorporale			
	1. Imobilizări necorporale în curs de execuție	010		
	2. Imobilizări necorporale în exploatare, total	020	487	
	din care:			
	2.1. concesiuni, licențe și mărci	021	487	
	2.2. drepturi de autor și titluri de protecție	022		
	2.3. programe informatice	023		
	2.4. alte imobilizări necorporale	024		
	3. Fond comercial	030		
	4. Avansuri acordate pentru imobilizări necorporale	040		
	Total imobilizări necorporale (rd.010 + rd.020 + rd.030 + rd.040)	050	487	
	II. Imobilizări corporale			
	1. Imobilizări corporale în curs de execuție	060		
	2. Terenuri	070		
	3. Mijloace fixe, total	080	2208593	2793637
	din care:			
	3.1. clădiri	081		
	3.2. construcții speciale	082		
	3.3. mașini, utilaje și instalații tehnice	083	2204135	2791637
	3.4. mijloace de transport	084		

A.	3.5. inventar și mobilier	085		
	3.6. alte mijloace fixe	086	4458	2000
	4. Resurse minerale	090		
	5. Active biologice imobilizate	100		
	6. Investiții imobiliare	110		
	7. Avansuri acordate pentru imobilizări corporale	120		
	Total imobilizări corporale (rd.060 + rd.070 + rd.080 + rd.090 + rd.100 + rd.110 + rd.120)	130	2208593	2793637
	III. Investiții financiare pe termen lung			
	1. Investiții financiare pe termen lung în părți neafiliate	140		
	2. Investiții financiare pe termen lung în părți afiliate, total	150		
B.	din care:	151		
	2.1. acțiuni și cote de participație deținute în părțile afiliate	152		
	2.2. Împrumuturi acordate părților afiliate	153		
	2.3 Împrumuturi acordate aferente intereselor de participație	154		
	2.4 alte investiții financiare	160		
	Total investiții financiare pe termen lung (rd.140 + rd.150)	160		
	IV. Creanțe pe termen lung și alte active imobilizate			
	1. Creanțe comerciale pe termen lung	170		
	2. Creanțe ale părților afiliate pe termen lung	180		
	Inclusiv: creanțe aferente intereselor de participație	181		
	3. Alte creanțe pe termen lung	190		
	4. Cheltuieli anticipate pe termen lung	200		
	5. Alte active imobilizate	210		
	Total creanțe pe termen lung și alte active imobilizate (rd.170 + rd.180 + rd.190 + rd.200 + rd.210)	220		
	TOTAL ACTIVE IMOBILIZATE (rd.050 + rd.130 + rd.160 + rd.220)	230	2209080	2793637
	ACTIVE CIRCULANTE			
	I. Stocuri			
	1. Materiale și obiecte de mică valoare și scurtă durată	240	54051	51978
	2. Active biologice circulante	250		
	3. Producția în curs de execuție	260		
	4. Produse și mărfuri	270	5710647	7221203
	5. Avansuri acordate pentru stocuri	280		
	Total stocuri (rd.240 + rd.250 + rd.260 + rd.270 + rd.280)	290	5764698	7273181
	II. Creanțe curente și alte active circulante			
	1. Creanțe comerciale curente	300	4337729	3912218
	2. Creanțe ale părților afiliate curente	310		
	Inclusiv: creanțe aferente intereselor de participație	311		
	3. Creanțe ale bugetului	320	166486	74631
	4. Creanțele ale personalului	330		
	5. Alte creanțe curente	340		
	6. Cheltuieli anticipate curente	350	4	2
	7. Alte active circulante	360	1647908	5756117
	Total creanțe curente și alte active circulante (rd.300 + rd.310 + rd.320 + rd.330 + rd.340 + rd.350 + rd.360)	370	6152127	9742968
	III. Investiții financiare curente			
	1. Investiții financiare curente în părți neafiliate	380		
	2. Investiții financiare curente în părți afiliate, total	390		
	din care:	391		
	2.1. acțiuni și cote de participație deținute în părțile afiliate	392		
	2.2. Împrumuturi acordate părților afiliate	393		
	2.3. Împrumuturi acordate aferente intereselor de participație			

	2.4. alte investiții financiare în părți afiliate	394		
	Total investiții financiare curente (rd.380 + rd.390)	400		
	IV. Numerar și documente bănești	410	8911899	3942779
	TOTAL ACTIVE CIRCULANTE (rd.290 + rd.370 + rd.400 + rd.410)	420	20828724	20958928
	TOTAL ACTIVE (rd.230 + rd.420)	430	23037804	23752565
	P A S I V			
	CAPITAL PROPRIU			
	I. Capital social și neînregistrat			
	1. Capital social	440	5400	5400
C.	2. Capital nevărsat	450	()	()
	3. Capital neînregistrat	460		
	4. Capital retras	470	()	()
	5. Patrimoniul primit de la stat cu drept de proprietate	480		
	Total capital social și neînregistrat (rd.440 + rd.450 + rd.460 + rd.470 + rd.480)	490	5400	5400
	II. Prime de capital	500		
	III. Rezerve			
	1. Capital de rezervă	510		
	2. Rezerve statutare	520		
	3. Alte rezerve	530		
	Total rezerve (rd.510 + rd.520 + rd.530)	540		
	IV. Profit (pierdere)			
	1. Corecții ale rezultatelor anilor precedenți	550	X	
	2. Profit nerepartizat (pierdere neacoperită) al anilor precedenți	560	21021465	12085295
	3. Profit net (pierdere netă) al perioadei de gestiune	570	X	7974831
	4. Profit utilizat al perioadei de gestiune	580	X	()
	Total profit (pierdere) (rd.550 + rd.560 + rd.570 + rd.580)	590	21021465	20060126
	V. Rezerve din reevaluare	600		
	VI. Alte elemente de capital propriu	610		
	TOTAL CAPITAL PROPRIU (rd.490 + rd.500 + rd.540 + rd.590 + rd.600 + rd.610)	620	21026865	20065526
D.	DATORII PE TERMEN LUNG			
	1. Credite bancare pe termen lung	630		
	2. Împrumuturi pe termen lung	640		
	din care:			
	2.1. Împrumuturi din emisiunea de obligațiuni	641		
	Inclusiv: împrumuturi din emisiunea de obligațiuni convertibile	642		
	2.2. alte împrumuturi pe termen lung	643		
	3. Datorii comerciale pe termen lung	650		
	4. Datorii față de părțile afiliate pe termen lung	660		
	Inclusiv: datorii aferente intereselor de participație	661		
	5. Avansuri primite pe termen lung	670		
	6. Venituri anticipate pe termen lung	680		
	7. Alte datorii pe termen lung	690		
	TOTAL DATORII PE TERMEN LUNG (rd.630 + rd.640 + rd.650 + rd.660 + rd.670 + rd.680 + rd.690)	700		
	DATORII CURENTE			
	1. Credite bancare pe termen scurt	710		
	2. Împrumuturi pe termen scurt, total	720		

E.	din care:	721		
	2.1. Împrumuturi din emisiunea de obligațiuni	722		
	Inclusiv: împrumuturi din emisiunea de obligațiuni convertibile	723		
	2.2. alte împrumuturi pe termen scurt	730	1331928	3252667
	3. Datorii comerciale curente	740		
	4. Datorii față de părțile afiliate curente	741		
	Inclusiv: datorii aferente intereselor de participație	750	159545	188105
	5. Avansuri primite curente	760	2913	50
	6. Datorii față de personal	770		
	7. Datorii privind asigurările sociale și medicale	780	434590	187676
F.	8. Datorii față de buget	790		
	9. Datorii față de proprietari	800		
	10. Venituri anticipate curente	810	81963	58541
	11. Alte datorii curente	820	2010939	3687039
	TOTAL DATORII CURENTE (rd.710 + rd.720 + rd.730 + rd.740 + rd.750 + rd.760 + rd.770 + rd.780 + rd.790 + rd.800 + rd.810)			
	PROVIZIOANE			
	1. Provizioane pentru beneficiile angajaților	830		
	2. Provizioane pentru garanții acordate cumpărătorilor/clientșilor	840		
	3. Provizioane pentru impozite	850		
	4. Alte provizioane	860		
	TOTAL PROVIZIOANE (rd.830 + rd.840 + rd.850 + rd.860)	870		
	TOTAL PASIVE (rd.620 + rd.700 + rd.820 + rd.870)	880	23037804	23752565

SITUAȚIA DE PROFIT ȘI PIERDERE

de la până la

Anexa 2

Indicatori	Cod rd.	Perioada de gestiune	
		precedenta	curenta
1	2	3	4
Venituri din vânzări, total	010	27319617	25963175
din care:			
venituri din vânzarea produselor și mărfurilor	011	26856566	25044358
venituri din prestarea serviciilor și executarea lucrărilor	012	463051	918817
venituri din contracte de construcție	013		
venituri din contracte de leasing	014		
venituri din contracte de microfinanțare	015		
alte venituri din vânzări	016		
Costul vânzărilor, total	020	15672962	15186814
din care:			
valoarea contabilă a produselor și mărfurilor vândute	021	15672962	15186814
costul serviciilor prestate și lucrărilor executate terților	022		
costuri aferente contractelor de construcție	023		
costuri aferente contractelor de leasing	024		
costuri aferente contractelor de microfinanțare	025		
alte costuri aferente vânzărilor	026		
Profit brut (pierdere brută) (rd.010 - rd.020)	030	11646655	10776361
Alte venituri din activitatea operațională	040	28586	247603
Cheltuieli de distribuie	050	16306	19740
Cheltuieli administrative	060	964136	1259776
Alte cheltuieli din activitatea operațională	070	417394	640169
Rezultatul din activitatea operațională: profit (pierdere) (rd.030 + rd.040 - rd.050 - rd.060 - rd.070)	080	10277405	9104279

Venituri financiare, total	090	490609	519239
din care:			
venituri din interese de participare	091		
inclusiv: veniturile obținute de la părțile afiliate	092		
venituri din dobânzi	093		25612
inclusiv: veniturile obținute de la părțile afiliate	094		
venituri din alte investiții financiare pe termen lung	095		
inclusiv: veniturile obținute de la părțile afiliate	096		
venituri aferente ajustărilor de valoare privind investițiile financiare pe termen lung și curente	097		
venituri din ieșirea investițiilor financiare	098		
venituri aferente diferențelor de curs valutar și de sumă	099	490609	493627
Cheltuieli financiare, total	100	686605	597528
din care:			
cheltuieli privind dobânzile	101		
inclusiv: cheltuielile aferente părților afiliate	102		
cheltuieli aferente ajustărilor de valoare privind investițiile financiare pe termen lung și curente	103		
cheltuieli aferente ieșirii investițiilor financiare	104		
cheltuieli aferente diferențelor de curs valutar și de sumă	105	686605	597528
Rezultatul: profit (pierdere) financiar(ă) (rd.090 - rd.100)	110	-195996	-78289
Venituri cu active imobilizate și excepționale	120		
Cheltuieli cu active imobilizate și excepționale	130		
Rezultatul din operațiuni cu active imobilizate și excepționale: profit (pierdere) (rd.120 - rd.130)	140		
Rezultatul din alte activități: profit (pierdere) (rd.110 + rd.140)	150	-195996	-78289
Profit (pierdere) până la impozitare (rd.080 + rd.150)	160	10081409	9025990
Cheltuieli privind impozitul pe venit	170	1178993	1051159
Profit net (pierdere netă) al perioadei de gestiune (rd.160 - rd.170)	180	8902416	7974831

SITUAȚIA MODIFICĂRILOR CAPITALULUI PROPRIU

de la până la

Anexa 3

Nr. d/o	Indicatori	Cod rd	Sold la începutul perioadei de gestiune	Majorări	Diminuări	Sold la sfârșitul perioadei de gestiune
1	2	3	4	5	6	7
I.	Capital social și neînregistrat					
	1. Capital social	010				
	2. Capital nevărsat	020 ()	()	()	()	()
	3. Capital neînregistrat	030				
	4. Capital retras	040 ()	()	()	()	()
	5. Patrimoniul primit de la stat cu drept de proprietate	050				
	Total capital social și neînregistrat (rd.010 + rd.020 + rd.030 + rd.040 + rd.050)	060				
II.	Prime de capital	070				
III.	Rezerve					
	1. Capital de rezervă	080				
	2. Rezerve statutare	090				
	3. Alte rezerve	100				
	Total rezerve (rd.080 + rd.090 + rd.100)	110				
	Profit (pierdere)					
	1. Corecții ale rezultatelor anilor precedenți	120	X			

IV.	2. Profit nerepartizat (pierdere neacoperită) al anilor precedenți	130				
	3. Profit net (pierdere netă) al perioadei de gestiune	140	X			
	4. Profit utilizat al perioadei de gestiune	150	X	()	()	()
	Total profit (pierdere) (rd.120 + rd.130 + rd.140 + rd.150)	160				
V.	Rezerve din reevaluare	170				
VI.	Alte elemente de capital propriu	180				
	Total capital propriu (rd.060 + rd.070 + rd.110 + rd.160 + rd.170 + rd.180)	190				

SITUAȚIA FLUXURILOR DE NUMERAR

de la până la

Anexa 4

Indicatori	Cod rd	Perioada de gestiune	
		precedentă	curentă
1	2	3	4
Fluxuri de numerar din activitatea operațională			
Încasări din vânzări	010		
Plăți pentru stocuri și servicii procurate	020		
Plăți către angajați și organe de asigurare socială și medicală	030		
Dobânzi plătite	040		
Plata impozitului pe venit	050		
Alte încasări	060		
Alte plăți	070		
Fluxul net de numerar din activitatea operațională (rd.010 - rd.020 - rd.030 - rd.040 - rd.050 + rd.060 - rd.070)	080		
Fluxuri de numerar din activitatea de investiții			
Încasări din vânzarea activelor imobilizate	090		
Plăți aferente intrărilor de active imobilizate	100		
Dobânzi încasate	110		
Dividende încasate	120		
Inclusiv: dividende încasate din străinătate	121		
Alte încasări (plăți)	130		
Fluxul net de numerar din activitatea de investiții (rd.090 - rd.100 + rd.110 + rd.120 ± rd.130)	140		
Fluxuri de numerar din activitatea financiară			
Încasări sub formă de credite și împrumuturi	150		
Plăți aferente rambursării creditelor și împrumuturilor	160		
Dividende plătite	170		
Inclusiv: dividende plătite nerezidenților	171		
Încasări din operațiuni de capital	180		
Alte încasări (plăți)	190		
Fluxul net de numerar din activitatea financiară (rd.150 - rd.160 - rd.170 + rd.180 ± rd.190)	200		
Fluxul net de numerar total (± rd.080 ± rd.140 ± rd.200)	210		
Diferențe de curs valutar favorabile (nefavorabile)	220		
Sold de numerar la începutul perioadei de gestiune	230		
Sold de numerar la sfârșitul perioadei de gestiune (± rd.210 ± rd.220 + rd.230)	240		

Documente atașate - Notă explicativă (fișierul pdf)

Версия для печати

Сохранить

Расписка 2

Респондент
Фискальный код: 1010600028048, наименование: BIOSISTEM MLD S.R.L.
Предоставил отчет: RSF1_21
На фискальный период: A/2020
Дата предоставления: 11.05.2021
Временная метка отчёта зарегистрированного в Информационной Системе НБС : 11.05.2021 12:26:31

National Bureau of Statistics (NBS) received the electronic version of the report, sent by you. The data provided is verified by NBS.

Версия для печати

Сохранить

Расписка

Респондент
Фискальный код: 1010600028048, наименование: BIOSISTEM MLD S.R.L.
Предоставил отчет: RSF1_21
На фискальный период: A/2020
Дата предоставления: 11.05.2021
Временная метка отчёта зарегистрированного в Системе Электронной Отчётности и отправленного в Информационную Систему БНС : 11.05.2021 10:00:47



BIOSISTEM-MLD S.R.L.

c/f 1010600028048; adresa: str. Albișoara 16/1 of.7, or. Chișinău
tel.+373-22-808517, +373-22-808719, fax +373-22-808519.
Web: www.biosistem-mld.com; e-mail: biosistem.mld@gmail.com

Către Grupul de lucru pentru evaluarea

Procedurii de achiziție nr. ocds-b3wdp1-MD-1628866063582

Din 23 aug 2021, 11:00 - 18 sept 2021, 11:00

din cadrul CAPCS

Declarație

Prin prezenta, SRL „Biosistem-mld”, declara ca :

- Va înregistra în Registrul de Stat al Dispozitivelor Medicale a Agenției Medicamentului și Dispozitivelor Medicale bunurile contractate până la momentul livrării acestora.
- Termenul de valabilitate restant (la momentul livrării) va constitui minim 80% din termenul total al produsului, dar nu mai mic de 12 luni
- Va prezenta mostre, în termen de 5 zile de la solicitarea autorității contractante

_____ Vitalie Poiata

L.Ș.



Certificate

No. Q5 105557 0001 Rev. 00

Holder of Certificate: **Meril Endo Surgery Pvt. Ltd.**

Third floor, E1- E3, Meril Park
Survey No 135/2/B & 174/2, Muktanand Marg, Chala
396191 Vapi, Gujarat
INDIA

Facility(ies):

Meril Endo Surgery Pvt. Ltd.
Third floor, E1- E3, Meril Park , Survey No 135/2/B & 174/2,
Muktanand Marg, Chala, 396191 Vapi, Gujarat, INDIA

Meril Endo Surgery Pvt. Ltd
Type A-2, Shed No. 11, Survey No. 725/P, Phase – I, GIDC,
396195 Vapi, Gujarat, INDIA

Meril Endo Surgery Pvt. Ltd
Plot No 688/10 & 11, Siddivinayak Industrial Estate, Somnath
Road, 396210 Daman, INDIA

Certification Mark:



Scope of Certificate:

Design, Development, Production, Testing, Storage, Sales and Distribution of Sterile / Non-sterile Endo Surgical Products e.g., Bulk and Finished Surgical Sutures, Polytetrafluoroethylene Pledgets, Sutures with Polytetrafluoroethylene Pledgets, Contraceptive Devices, Surgical Meshes, Mesh Fixation Devices, Bone Wax, Umbilical Cotton Tape, Surgical Kits, Surgical Staplers and Staples, Tourniquet Devices, Disposable Endoscopic Trocars, Ligating Clips, Surgical Haemostats, Skin Adhesives, Surgical Adhesives, Sealants & Surgical Needles”

Applied Standard(s):

EN ISO 13485:2016
Medical devices - Quality management systems -
Requirements for regulatory purposes
(ISO 13485:2016)
DIN EN ISO 13485:2016

The Certification Body of TÜV SÜD Product Service GmbH certifies that the company mentioned above has established and is maintaining a quality management system, which meets the requirements of the listed standard(s). See also notes overleaf.

Report No.: IND2019074
Valid from: 2020-03-30
Valid until: 2023-03-29

Date, 2020-03-30

Christoph Dicks
Head of Certification/Notified Body

EC CERTIFICATE

According to Annex II of the Directive 93/42/EEC on Medical Devices

Full Quality Assurance System

Certificate Number: 2195-MED-1929401

Manufacturer:

MERIL ENDO SURGERY PVT. LTD.

Third Floor, E1-E3, Meril Park, Survey No. 135/2/B & 174/2, Muktanand Marg,
Chala, Vapi – 396 191, Gujarat, INDIA

Product(s):

- (1) Non-Absorbable, Braided Coated Poly(Ethylene Terephthalate) Surgical Suture
- (2) Non-Absorbable, Monofilament Polyamide Surgical Suture
- (3) Non-Absorbable, Monofilament Polypropylene Surgical Suture
- (4) Non-Absorbable, Braided Coated Silk Surgical Suture
- (5) Non-Absorbable, Monofilament Stainless Steel Surgical Suture

Model(s):

Product specifications are stated on the second page.

Reference Report No:

MM0755-P001-R01, MM0755-P001-R02

Szutest, Notified Body 2195, declares that the aforementioned manufacturer has implemented a quality assurance system according to Annex II (excluding section 4), Section 3 of the directive 93/42/EEC on medical devices. This quality assurance system covers those aspects of manufacturing concerned with securing and maintaining safe conditions of the respective product(s) and conforms to the provisions of this Directive. The approved quality system is subject to surveillance pursuant to Annex II, Section 5 of Directive 93/42/EEC and unannounced audits.

Szutest must be informed of any significant changes in the design and/or construction of the product(s). For class I devices with sterile conditions the quality management system evaluation is restricted to the aspects of manufacture concerned with securing and maintaining sterile conditions. For class I devices with measuring function the quality management system evaluation is restricted to the aspects of manufacture concerned with the conformity of the devices with metrological requirements.

This EC certificate is valid till 2024-05-26.

Issue Date:

2019-10-21



Rukiye BALKAN
Deputy General Manager

SZUTEST

Certificate Number: 2195-MED-1929401

Product Specifications

Product Categories	Type (Models)	Generic Name
(1) Non-Absorbable, Braided Coated Poly(Ethylene Terephthalate) Surgical Suture	MERICRON XL™	Sterile Non-Absorbable Poly(Ethylene Terephthalate) Surgical Suture
	Aspiron™	Sterile Non-Absorbable Poly(Ethylene Terephthalate) Surgical Suture
	MERICRON XL™ P	Sterile Non-Absorbable Poly(Ethylene Terephthalate) Surgical Suture with PTFE Pledget
(2) Non-Absorbable, Monofilament Polyamide Surgical Suture	FILAMIDE™	Sterile Non-Absorbable Polyamide Surgical Suture
	Aspiron™	Sterile Non-Absorbable Polyamide Surgical Suture
(3) Non-Absorbable, Monofilament Polypropylene Surgical Suture	FILAPROP™ P	Sterile Non-Absorbable Polypropylene Surgical Suture with PTFE Pledget
(4) Non-Absorbable, Braided Silk Surgical Suture	FILASILK™ REEL	Non-Sterile Non-Absorbable Silk Surgical Suture
(5) Non-Absorbable, Monofilament Stainless Steel Surgical Suture	MERISTEEL™	Sterile Non-Absorbable Stainless Steel Surgical Suture



SZUTEST UYGUNLUK DEĞERLENDİRME A.Ş.

Tatlısu Mahallesi, Akif İnan Sk. No:1 Ümraniye 34774 İSTANBUL / TÜRKİYE

Szutest.com.tr

EC Certificate

Full Quality Assurance System

Certificate No.:
245506-2017-CE-IND-NA-PS Rev. 2.0

Project No.:
PRJC-499089-2014-MSL-IND

Valid Until:
27 May 2024

This is to certify that the quality system of:

Meril Endo Surgery Pvt Ltd

Third Floor, E1- E3, Meril Park, Survey No 135/2/B & 174/2, Muktanand Marg,
Chala, Vapi, Gujarat, India - 396191

For design, production and final product inspection/testing of:

Sterile / non - sterile surgical sutures with and without needle

Has been assessed with respect to:

The conformity assessment procedure described in Annex II of Council Directive 93/42/EEC on Medical Devices, as amended

and found to comply.

Further details of the product(s) and conditions for certification are given overleaf.

Place and Date:
Høvik, 16 July 2019



For:
DNV GL PRESAFE AS



Cathrine Wisbech

The Certificate has been digitally signed.
See www.presafe.com/digital_signatures for more info

Notice: The Certificate is subject to terms and conditions as set out in the Certification Agreement. Failure to comply may render this Certificate invalid.

EC Certificate

Full Quality Assurance System

Certificate No.:
245506-2017-CE-IND-NA-PS Rev. 2.0

Project No.:
PRJC-499089-2014-MSL-IND

Valid Until:
27 May 2024

Jurisdiction

Application of Council Directive 93/42/EEC of 14 June 1993, adopted as "Forskrift om Medisinsk Utstyr" by the Norwegian Ministry of Health and Care Services.

Certificate history:

Revision	Description	Issue Date
0.0	Supersedes DNVGL (NB 0434) Certificate no: 151561-2014-CE-IND-NA Rev 3.0 following transfer to notified body functions to DNV GL Nemko Presafe AS (NB 2460)	2017-11-15
1.0	Remove of Polypropylene Mesh	2018-01-08
2.0	Recertification and reduction in scope	2019-07-16

Products covered by this Certificate:

Product Description	Product Name	Class
Absorbable Sutures	- Megasorb™ / Aspiron™ Polyglycolic acid Braided coated Polyglycolic acid suture - Mitsu™ / Aspiron™ Polyglactin 910 and Mitsu FST™ / Aspiron™ Polyglactin 910 FST. Braided coated Poly (glycolide/l-lactide) suture - Filaxyn™ / Aspiron™ Polydioxanone suture. Monofilament Poly (p-dioxanone) suture - Filapron™ / Aspiron™ Polyglecaprone 25 suture. Monofilament poly (glycolide-co-caprolactone) suture	III*
Non-Absorbable Sutures	Filaprop™ / Aspiron™ Polypropylene Blue Monofilament Polypropylene Suture	III**

EC Certificate

Full Quality Assurance System

Certificate No.:
245506-2017-CE-IND-NA-PS Rev. 2.0

Project No.:
PRJC-499089-2014-MSL-IND

Valid Until:
27 May 2024

*Design assessment is covered by a separate EC-Design Examination Certificate No.: 245507-2017-CE-IND-NA-PS Rev. 2

**Design assessment is covered by a separate EC-Design Examination Certificate No.: 245508-2017-CE-IND-NA-PS Rev. 2

Sites covered by this certificate

Site Name	Address
Meril Endo Surgery Pvt Ltd	Third Floor, E1- E3, Meril Park, Survey No 135/2/B & 174/2, Muktanand Marg, Chala, Vapi, Gujarat, India - 396191

EU Representative

OBELIS S.A

Bd. Général Wahis, 53, 1030 Brussels, Belgium. Tel: +32.2.732.59.54. Fax: +32.2.732.60.03

E-mail: mail@obelis.net, www.obelis.net

EC Certificate

Full Quality Assurance System

Certificate No.:
245506-2017-CE-IND-NA-PS Rev. 2.0

Project No.:
PRJC-499089-2014-MSL-IND

Valid Until:
27 May 2024

Terms and conditions

The certificate is subject to the following terms and conditions:

- Any producer (see 2001/95/EC for a precise definition) is liable for damage caused by a defect in his product(s), in accordance with directive 85/374/EEC, as amended, concerning liability of defective products.
- The certificate is only valid for the products and/or manufacturing premises listed above.
- The Manufacturer shall fulfil the obligations arising out of the quality system as approved and uphold it so that it remains adequate and efficient.
- The Manufacturer shall inform Presafe of any intended updating of the quality system and Presafe will assess the changes and decide if the certificate remains valid.
- Periodical audits will be held, in order to verify that the Manufacturer maintains and applies the quality system. Presafe reserves the right, on a spot basis or based on suspicion, to pay unannounced visits.

The following may render this Certificate invalid:

- Changes in the quality system affecting production.
- Periodical audits not held within the allowed time window.


Conformity declaration and marking of product

When meeting with the terms and conditions above, the producer may draw up an EC declaration of conformity and legally affix the CE mark followed by the Notified Body identification number of Presafe.

End of Certificate

Meril Endo-Surgery

Wound Closure Devices



Other Product Ranges

Staplers | Hemostats | Energy Devices | Hernia Repair Meshes | Ligation Clips | Adhesives | Access Devices
IUDs | Umbilical Cotton Tape | Bonewax

Infrastructure

Situated on 180 acres of land and built over a 300,000 sq.ft. area, Meril's ISO 13485 & cGMP certified ultra-modern manufacturing facility is constructed to satisfy the stringent needs of our cardiovascular, orthopedics, diagnostics and endo-surgery businesses, as well as our R&D endeavors.

Fully integrated manufacturing systems ensure backward integration, flawless man-material movement and complete control over processes to eliminate production errors. All manufacturing and sterilization processes are conducted in-house, in addition to the analytical and microbiological QA/QC tests required to meet the world-class production standards.

Meril Park, Vapi





Talent

Meril's 3200+ personnel are a strong, experienced team comprising designers, engineers, chemists, microbiologists, regulatory affair experts, R&D scientists, clinical affairs experts, legal, finance, sales and marketing professionals with innovative capabilities. The team is continuously striving towards improving and saving lives every day by providing revolutionary solutions for diagnosis, prevention and treatment.

Currently Meril has direct sales & representative offices in India, Germany, Turkey, USA, China and Brazil, and a strong distribution network spread over an additional 102 countries in 6 continents.

Meril Academy

While innovations and advancements in medical technology play a vital role in contributing to the betterment of a healthy society, a larger role also lies in educating and imparting practical training to doctors, clinicians and other members of the medical fraternity. This initiative will help bridge the gaps between the available technology, its adoptions and its implementation for superior patient care.



Also, to add value to medical professionals and young medical interns, Meril Academy functions as world class institution. Spread over a massive area of 200,000 sq. ft. and surrounded with lush, green natural beauty, Meril Academy is situated at the Meril headquarters in Vapi.



Meril Endo-Surgery

In line with the objective to design and develop novel, clinically relevant, state-of-art surgical devices, Meril is expanding its endo-surgery franchise. At Meril, we strive to create and develop products through innovative thinking, cutting edge technologies and collaboration with medical professionals and healthcare organizations, with an aim to identify clinical needs and translate them into market-leading solutions.

Meril's vision is to become a 'one-stop' solution for surgeons by providing them with an extensive product portfolio comprising absorbable & non-absorbable surgical sutures, tissue sealants, absorbable hemostats, hernia repair meshes, intra uterine devices, energy devices and mechanical closure devices.



Contents

Synthetic Absorbable Sutures

MITSU	15
MITSU C+	43
MITSU FST	49
MEGASORB	55
MEGASORB PLUS	65
MEGASORB RAPID	67
FILAPRON	71
FILAXYN	79

Natural Absorbable Sutures

MERISOFT Plain Catgut	89
MERISOFT Chromic Catgut	93

Natural Non-Absorbable Sutures

FILASILK	99
----------	----

Synthetic Non-Absorbable Sutures

MERISTEEL	111
FILAMIDE	115
FILAPROP / P	127
MERICRON XL / P	149



Wound Closure

Wounds which do not heal naturally are treated with wound closure devices. The affected tissues need to be kept in opposition until the healing process provides the wound with proper strength to stand against the stress without any assistance. A wound can be treated with the help of various devices such as sutures, staplers, topical adhesives or clips. Different tissues have their own unique functions and requirements. When one decides on a particular wound closure device, one must take care of the healing needs of the particular tissue which is being treated.

The most common device used for wound closure is suture. Sutures, in general term 'stitches', are sterilized surgical threads, which are used to ligate blood vessels or to approximate tissues.

Sutures are generally categorized as:

- Absorbable and Non-Absorbable
- Natural and Synthetic
- Monofilament and Braided

Absorbable sutures have a tendency to degrade as the wound heals. Non-Absorbable sutures oppose the body's efforts to dissolve them. On complete healing of the wound, surgeons remove the non-absorbable sutures. Meril Endo-Surgery has a wide range of synthetic absorbable sutures which are made from Polyglactin 910 (MITSU family), Polydioxanone (FILAXYN), Poliglecaprone 25 (FILAPRON), Polyglycolic Acid (MEGASORB family) as well as non-absorbable sutures which are made from Polypropylene (FILAPROP), Polyamide (FILAMIDE) and Polyester (MERICRON XL). In natural absorbable sutures, Meril portfolio covers Chromic and Plain Catgut (MERISOFT) and in natural non-absorbable we have silk (FILASILK).


For placement of sutures, surgical needles are used, which are made of stainless steel. Different types and sizes of needles are used depending on the area of application. We are dedicated to assist surgeons in tissue repair by giving them a broad gamut of products with our unmatched quality and service.

Meril Endo-Surgery Needles

In pursuit of offering best quality products, Meril sutures are paired with the superior stainless steel needles for an unparalleled surgical experience.


Meril sutures come with sharper and stronger needles (300 series Stainless Steel) specifically designed to meet the ever growing needs in all the surgical disciplines.

With every needle swaged using Multi Stroke Technology (MST) and tested for a needle pull check, Meril needles are distinguished with a high penetration coating and a square body design for a firm grip.



Needle Profiles	Point type	Symbol
	Round bodied/Taper point	
	Conventional cutting	
	Reverse cutting	
	Reverse cutting prime	
	Tapercutting	
	UltraPoint reverse cutting	
	UltraPoint spatula curved	
	CSU spatula	
	SBR spatula	

Needle Shapes



Needle Specifications

Type	Taper	Reverse cutting	Diamond	Cutting	Premium cutting	Spatula	Blunt	Taper cut
Design								
Point								
Shape	Typical Applications							
Straight	Skin/subcutaneous layer closure; Purse string							
1/4 Circle	Ophthalmic/Microsurgical procedures							
3/8 Circle	Skin/subcutaneous layer/abdominal wall closure; Ophthalmic/Vascular/Plastic procedures							
1/2 Circle	Peritoneum closure; Gastrointestinal/Respiratory/Urological procedures							
5/8 Circle	Pelvic & Oral cavity; Urological procedures							
J-Shape	Laparotomy closure							
Ski	Laparoscopic procedures							

Understanding the Suture Codes

Meril code logic represents the necessary information needed to identify the required suture-needle combination. The suture codes are alphanumeric in nature to indicate the suture material, the suture size, and the reference code. The alphabetical suffixes are a key to understand the additional features about the suture and/or the needle.

Material Prefixes

ABBREVIATION	MEANING
PGN	Polyglactin 910
PGP	Polyglactin 910 (with Chlorhexidine)
PGF	Polyglactin 910 Fast
MS	Polyglycolic Acid
MSP	Polyglycolic Acid (with Chlorhexidine)
MST	Polyglycolic Acid (with Triclosan)
PCL	Poliglecaprone 25
PDX	Polydioxanone
CGL	Chromic Catgut Loop
SLK	Silk
SLS	Silk Non Needled
NYL	Polyamide
NYS	Polyamide Non Needled
PPL	Polypropylene
PPG	Polypropylene with Pledgets
STC	Steel Wire

*In case of

- Polyester (Mericon XL), codes start with ME6XX, ME6XXX, ME4XXX, MEXXXX; and polyester with pledgets (Mericon XL P), codes start with PMEXXX or MMEXXXX
- Chromic Catgut (Merisoft Chromic), codes start with ME40XX, ME42XX, ME50XX
- Chromic Catgut non needled (Merisoft Chromic), codes start with SS22XX
- Plain Catgut (Merisoft Plain), codes start with ME50XX
- Plain Catgut non needled (Merisoft Plain), codes start with SS20XX


Suture Sizes

- For sutures having a single-digit size, the size is written as: zero followed by the size digit.
E.g.: Size 01 (read as size 1 according to USP), Size 03 (read as size 3 according to USP).
- For finer sutures, the size is written as: the size digit followed by a zero.
E.g.: Size 60 (read as size 6-0 according to USP), Size 70 (read as size 7-0 according to USP), Size 100 (read as size 10-0 according to USP)

Suffixes

ABBREVIATION	MEANING
B/BL	Black Needle
BX	Black Needle; Extra Length Suture
DN	Double Needle
DNL	Double Needle; Long Length Suture
DS	Doctor Special Code
G	Green Suture
ML	Multi Strands
SL	Short Length Suture
U	Undyed Suture
V	Violet Suture
W	White Suture
X	Extra Pair of Suture
XF	Extra Fine Needle
XL	Extra Length Suture

Decoding the Codes



Abbreviation for Specification

ABBREVIATION	MEANING
HC	1/2 CIRCLE
QC	1/4 CIRCLE
CU	3/8 CIRCLE
FE	5/8 CIRCLE
RC	REVERSE CUTTING
CT	CUTTING
RB	ROUND BODIED
TC	TAPER CUT
ST	STRAIGHT
SM	SPATULA
TRO	TROCAR
BP	BLUNT POINT
TP	PLEDGET
VIO	VIOLET
UND	UNDYED

Needle Point Symbol

	TAPER CUT		SPATULA
	CUTTING		STRAIGHT
	REVERSE CUTTING		TROCAR
	ROUND BODIED		

MITSU™

Polyglactin 910 Suture

Properties

- 75% of the original is retained until initial 14 Days
- 50% of the original is retained until initial 21 Days
- There is a subsequent loss between four to five weeks post implantation
- Complete absorption usually takes place between 56 to 70 Days

Indications

- Gastrointestinal procedures
- Gynaecology/Obstetric procedures
- Ophthalmic procedures
- Orthopedic procedures
- Urology procedures
- Skin closure (intracutaneous, subcutaneous)
- Ligatures

Advantages

- Smooth knot run-down/Knot slide
- Gentle knot repositioning
- Excellent knot security
- Minimal inflammatory reaction




Product Overview


Structure	: Multifilament, Braided
Colour	: Dyed, Undyed
	
Chemical composition	: Poly (glycolide-co-l-lactide 90/10)
Coating	: Poly(glycolide-co-l-lactide 30/70) + Calcium Stearate
Sizes	: Dyed - USP 8/0 (metric 0.4) - USP 2 (metric 5) : Undyed - USP 6/0 (metric 0.7) - USP 2 (metric 5)
Type of absorption	: Hydrolysis
Sterilization	: Ethylene oxide

Description	FG Code	Size		Needle Profile Shape, Type & Length	Suture Length (cm)	Pack Size	Suture Colour
MITSU (SYNTHETIC ABSORBABLE BRAIDED COATED POLYGLACTIN 910)							
PGLA VIO 0 X 45-22MM SKKI RB VB	PGN10 2901	0		22 mm SKKI Round Body	45	12	VIO
PGLA VIO 0 X 90-30MM HC RB	PGN10 2338	0		30 mm 1/2 Circle Round Body	90	12	VIO
PGLA VIO 0 X 90-40MM HC RB	PGN10 2346	0		40 mm 1/2 Circle Round Body	90	12	VIO
PGLA VIO 0 X 140-40MM HC RB TC DN	PGN10 2346DN	0		40 mm 1/2 Circle Round Body Taper Cut	140	12	VIO
PGLA VIO 0 X 180-40MM HC RB TC DN	PGN10 2346DNL	0		40 mm 1/2 Circle Round Body Taper Cut	180	12	VIO
PGLA VIO 0 X 110-40MM HC RB	PGN10 2346DS	0		40 mm 1/2 Circle Round Body	110	12	VIO
PGLA VIO 0 X 45-40MM HC RB	PGN10 2346SL	0		40 mm 1/2 Circle Round Body	45	12	VIO
PGLA VIO 0 X 90-40MM HC TC	PGN10 2518	0		40 mm 1/2 Circle Taper Cut	90	12	VIO
PGLA VIO 0 X 90-36MM HC RC	PGN10 2534	0		36 mm 1/2 Circle Reverse Cutting	90	12	VIO
PGLA VIO 0 X 45-23MM HC RC	PGN10 2825	0		23 mm 1/2 Circle Reverse Cutting	45	12	VIO
PGLA VIO 0 X 75-31MM HC RB	PGN109138	0		31 mm 1/2 Circle Round Body	75	12	VIO
PGLA VIO 0 X 75-45MM HC RB	PGN109244	0		45 mm 1/2 Circle Round Body	75	12	VIO
PGLA VIO 0 X 75-48MM HC RB	PGN109250	0		48 mm 1/2 Circle Round Body	75	12	VIO
PGLA VIO 0 X 75-48MM HC RC	PGN109295	0		48 mm 1/2 Circle Reverse Cutting	75	12	VIO
PGLA VIO 0 X 75-40MM HC RC	PGN109320	0		40 mm 1/2 Circle Reverse Cutting	75	12	VIO
PGLA VIO 0 X 75-31MM HC TC	PGN109361	0		31 mm 1/2 Circle Taper Cut	75	12	VIO

Description	FG Code	Size		Needle Profile Shape, Type & Length	Suture Length (cm)	Pack Size	Suture Colour
MITSU (SYNTHETIC ABSORBABLE BRAIDED COATED POLYGLACTIN 910)							
PGLA VIO 0 X 75-36MM HC TC	PGN109364	0		36 mm 1/2 Circle Taper Cut	75	12	VIO
PGLA VIO 0 X 75-40MM HC TC	PGN109372	0		40 mm 1/2 Circle Taper Cut	75	12	VIO
PGLA VIO 0 X 75-45MM HC TC	PGN109376	0		45 mm 1/2 Circle Taper Cut	75	12	VIO
PGLA VIO 0 X 90-40MM HC RC	PGN109420	0		40 mm 1/2 Circle Reverse Cutting	90	12	VIO
PGLA VIO 0 X 90-48MM HC RB	PGN109450	0		48 mm 1/2 Circle Round Body	90	12	VIO
PGLA VIO 0 X 90-36MM HC TC	PGN109464	0		36 mm 1/2 Circle Taper Cut	90	12	VIO
PGLA UND 0 X 70-36MM HC RB	PGN100260	0		36 mm 1/2 Circle Round Body	70	12	UND
PGLA UND 0 X 70-36MM HC CT	PGN100267	0		36 mm 1/2 Circle Conventional Cutting	70	12	UND
PGLA UND 0 X 70-26MM HC RB	PGN100270	0		26 mm 1/2 Circle Round Body	70	12	UND
PGLA VIO 0 X 90-36MM HC RB	PGN100346	0		36 mm 1/2 Circle Round Body	90	12	VIO
PGLA VIO 0 X 150-40MM HC RB LOOP	PGN1001407	0		40 mm 1/2 Circle Round Body	150	12	VIO
PGLA VIO 0 X 70-110MM FE RB	PGN100376	0		110 mm 5/8 Circle Round Body	70	12	VIO
PGLA VIO 0 X 70-40MM HC RB H	PGN100583	0		40 mm 1/2 Circle Round Body	70	12	VIO
PGLA VIO 0 X 45-21MM HC RB	PGN1080113	0		21 mm 1/2 Circle Round Body	45	12	VIO
PGLA VIO 0 X 45-26MM HC RB	PGN102713	0		26 mm 1/2 Circle Round Body	45	12	VIO
PGLA VIO 0 X 70-26MM HC RB	PGN100334	0		26 mm 1/2 Circle Round Body	70	12	VIO
PGLA VIO 0 X 70-26MM FE RB	PGN100604	0		26 mm 5/8 Circle Round Body	70	12	VIO
PGLA VIO 0 X 70-36MM FE RB	PGN100377	0		36 mm 5/8 Circle Round Body	70	12	VIO

Description	FG Code	Size		Needle Profile Shape, Type & Length	Suture Length (cm)	Pack Size	Suture Colour
MITSU (SYNTHETIC ABSORBABLE BRAIDED COATED POLYGLACTIN 910)							
PGLA VIO 0 X 70-40MM HC RC	PGN100485	0		40 mm 1/2 Circle Reverse Cutting	70	12	VIO
PGLA UND 0 X 70-40MM HC RC	PGN100268	0		40 mm 1/2 Circle Reverse Cutting	70	12	UND
PGLA VIO 0 X 70-40MM FE RB	PGN100381	0		40 mm 5/8 Circle Round Body	70	12	VIO
PGLA VIO 0 X 75-40MM HC RB BLUNT	PGN100280	0		40 mm 1/2 Circle Blunt Point	75	12	VIO
PGLA VIO 0 X 75-30MM HC TC	PGN100518	0		30 mm 1/2 Circle Taper Cut	75	12	VIO
PGLA VIO 0 X 75-30MM CU RB	PGN109030	0		30 mm 3/8 Circle Round Body	75	12	VIO
PGLA VIO 0 X 70-36MM CU RC	PGN109036	0		36 mm 3/8 Circle Reverse Cutting	70	12	VIO
PGLA UND 0 X 90-36MM HC RB	PGN10 946U	0		36 mm 1/2 Circle Round Body	90	12	UND
PGLA VIO 0 X 75-31MM J Type RB	PGN10 9138J	0		31 mm J Shape Round Body	75	12	VIO
PGLA VIO 0 X 70-30MM HC RB	PGN10 2342	0		30 mm 1/2 Circle Round Body	70	12	VIO
PGLA VIO 0 X 70-40MM HC RB	PGN10 2546	0		40 mm 1/2 Circle Round Body	70	12	VIO
PGLA VIO 0 X 90-40MM HC BP	PGN10 2996	0		40 mm 1/2 Circle Blunt Point	90	12	VIO
PGLA UND 0 X 75-40MM HC TC	PGN10 9901U	0		40 mm 1/2 Circle Taper Cut	75	12	UND
PGLA VIO 0 X 75-36MM J TYPE RB	PGN109221J	0		36 mm J Shape Round Body	75	12	VIO
PGLA VIO 0 X 75-40MM HC RB	PGN109230	0		40 mm 1/2 Circle Round Body	75	12	VIO
PGLA UND 0 X 90-36MM HC RB EX	PGN100946	0		36 mm 1/2 Circle Round Body	90	12	UND

Description	FG Code	Size		Needle Profile Shape, Type & Length	Suture Length (cm)	Pack Size	Suture Colour
MITSU (SYNTHETIC ABSORBABLE BRAIDED COATED POLYGLACTIN 910)							
PGLA VIO 0 X 70-40MM HC RB EX	PGN100352	0		40 mm 1/2 Circle Round Body	70	12	VIO
PGLA UND 0 X 75-40MM HC TC EX	PGN109901	0		40 mm 1/2 Circle Taper Cut	75	12	UND
PGLA VIO 0 X 90-40MM HC RB C1	PGN10 2358	0		40 mm 1/2 Circle Round Body	90	12	VIO
PGLA VIO 0 X 90-36MM HC TC C1	PGN10 2517	0		36 mm 1/2 Circle Taper Cut	90	12	VIO
PGLA UND 0 X 90-40MM HC RB	PGN10 2346U	0		40 mm 1/2 Circle Round Body	90	12	UND
PGLA VIO 0 X 90-36MM HC RC-3D	PGN102534Z	0		36 mm 1/2 Circle Reverse Cutting	90	36	VIO
PGLA VIO 0 X 70-36MM HC RB	PGN100324	0		36 mm 1/2 Circle Round Body	70	12	VIO
PGLA VIO 0 X 90-70MM HC RB	PGN100582	0		70 mm 1/2 Circle Round Body	90	12	VIO
PGLA VIO 0 X 70-30MM CU RC	PGN100610	0		30 mm 3/8 Circle Reverse Cutting	70	12	VIO
PGLA VIO 1 X 90-63MM CU RB BLUNT	PGN019391	1		63 mm 3/8 Circle Blunt Point	90	12	VIO
PGLA VIO 1 X 90-40MM HC RB TN	PGN01 2347	1		40 mm 1/2 Circle Round Body	90	12	VIO
PGLA VIO 1 X 140-40MM HC RB TC DN	PGN01 2347DN	1		40 mm 1/2 Circle Round Body Taper Cut	140	12	VIO
PGLA VIO 1 X 180-40MM HC RB TC DN	PGN01 2347DNL	1		40 mm 1/2 Circle Round Body Taper Cut	180	12	VIO
PGLA VIO 1 X 110-40MM HC RB	PGN01 2347DS	1		40 mm 1/2 Circle Round Body	110	12	VIO
PGLA VIO 1 X 45-40MM HC RB	PGN01 2347SL	1		40 mm 1/2 Circle Round Body	45	12	VIO
PGLA VIO 1 X 120-36MM HC RB	PGN01 2350	1		36 mm 1/2 Circle Round Body	120	12	VIO

Description	FG Code	Size		Needle Profile Shape, Type & Length	Suture Length (cm)	Pack Size	Suture Colour
MITSU (SYNTHETIC ABSORBABLE BRAIDED COATED POLYGLACTIN 910)							
PGLA VIO 1 X 70-40MM HC RB	PGN01 2359	1		40 mm 1/2 Circle Round Body	70	12	VIO
PGLA VIO 1 X 90-40MM HC RB TH	PGN01 2360	1		40 mm 1/2 Circle Round Body	90	12	VIO
PGLA VIO 1 X 90-40MM HC RC	PGN01 2421	1		40 mm 1/2 Circle Reverse Cutting	90	12	VIO
PGLA VIO 1 X 100-50MM HC RB	PGN01 2448	1		50 mm 1/2 Circle Round Body	100	12	VIO
PGLA VIO 1 X 35-23MM HC RC	PGN01 2826	1		23 mm 1/2 Circle Reverse Cutting	35	12	VIO
PGLA VIO 1 X 70-23MM HC RC	PGN01 2826XL	1		23 mm 1/2 Circle Reverse Cutting	70	12	VIO
PGLA VIO 1 X 90-40MM HC TC	PGN01 2519	1		40 mm 1/2 Circle Taper Cut	90	12	VIO
PGLA VIO 1 X 90-50MM HC RC	PGN01 2438	1		50 mm 1/2 Circle Reverse Cutting	90	12	VIO
PGLA VIO 1 X 75-31MM HC RB	PGN019213	1		31 mm 1/2 Circle Round Body	75	12	VIO
PGLA VIO 1 X 75-45MM HC RB	PGN019245	1		45 mm 1/2 Circle Round Body	75	12	VIO
PGLA VIO 1 X 100-80MM HC RB	PGN019289	1		80 mm 1/2 Circle Round Body	100	12	VIO
PGLA VIO 1 X 75-48MM HC RC	PGN019296	1		48 mm 1/2 Circle Reverse Cutting	75	12	VIO
PGLA VIO 1 X 75-40MM HC RC	PGN019321	1		40 mm 1/2 Circle Reverse Cutting	75	12	VIO
PGLA VIO 1 X 75-40MM HC RB H	PGN019340	1		40 mm 1/2 Circle Round Body	75	12	VIO
PGLA VIO 1 X 75-31MM HC TC	PGN019362	1		31 mm 1/2 Circle Taper Cut	75	12	VIO
PGLA VIO 1 X 75-36MM HC TC	PGN019365	1		36 mm 1/2 Circle Taper Cut	75	12	VIO
PGLA VIO 1 X 75-40MM HC TC	PGN019368	1		40 mm 1/2 Circle Taper Cut	75	12	VIO
PGLA VIO 1 X 75-45MM HC TC	PGN019377	1		45 mm 1/2 Circle Taper Cut	75	12	VIO


Description	FG Code	Size		Needle Profile Shape, Type & Length	Suture Length (cm)	Pack Size	Suture Colour
MITSU (SYNTHETIC ABSORBABLE BRAIDED COATED POLYGLACTIN 910)							
PGLA VIO 1 X 90-48MM HC RB H	PGN019451	1		48 mm 1/2 Circle Round Body	90	12	VIO
PGLA UND 1 X 75-40MM HC TC	PGN019906	1		40 mm 1/2 Circle Taper Cut	75	12	UND
PGLA UND 1 X 70-48MM HC RC	PGN010569	1		48 mm 1/2 Circle Reverse Cutting	70	12	UND
PGLA VIO 1 X 150-40MM HC RB LOOP	PGN011408	1		40 mm 1/2 Circle Round Body	150	12	VIO
PGLA VIO 1 X 150-48MM HC RB LOOP	PGN011428	1		48 mm 1/2 Circle Round Body	150	12	VIO
PGLA VIO 1 X 90-31MM HC RB	PGN010253	1		31 mm 1/2 Circle Round Body	90	12	VIO
PGLA VIO 1 X 70-48MM HC RB	PGN010321	1		48 mm 1/2 Circle Round Body	70	12	VIO
PGLA VIO 1 X 70-22MM HC RB H	PGN010583	1		22 mm 1/2 Circle Round Body	70	12	VIO
PGLA VIO 1 X 45-26MM HC RB	PGN0176513	1		26 mm 1/2 Circle Round Body	45	12	VIO
PGLA VIO 1 X 70-26MM HC RB	PGN010335	1		26 mm 1/2 Circle Round Body	70	12	VIO
PGLA VIO 1 X 90-26MM HC RB	PGN010372	1		26 mm 1/2 Circle Round Body	90	12	VIO
PGLA UND 1 X 70-36MM HC RC	PGN010535	1		36 mm 1/2 Circle Reverse Cutting	70	12	UND
PGLA VIO 1 X 70-36MM HC RB	PGN010341	1		36 mm 1/2 Circle Round Body	70	12	VIO
PGLA UND 1 X 70-36MM HC RB	PGN010347	1		36 mm 1/2 Circle Round Body	70	12	UND
PGLA UND 1 X 90-36MM HC RB BP	PGN019473	1		36 mm 1/2 Circle Blunt Point	90	12	UND
PGLA UND 1 X 70-40MM HC RB	PGN010281	1		40 mm 1/2 Circle Round Body	70	12	UND
PGLA UND 1 X 90-40MM HC RB	PGN010959	1		40 mm 1/2 Circle Round Body	90	12	UND
PGLA VIO 1 X 70-65MM HC RB	PGN010650	1		65 mm 1/2 Circle Round Body	70	12	VIO

Description	FG Code	Size		Needle Profile Shape, Type & Length	Suture Length (cm)	Pack Size	Suture Colour
MITSU (SYNTHETIC ABSORBABLE BRAIDED COATED POLYGLACTIN 910)							
PGLA VIO 1 X 75-40MM HC RB BLUNT	PGN010353	1		40 mm 1/2 Circle Blunt Point	75	12	VIO
PGLA VIO 1 X 90-35MM HC RB	PGN010261	1		35 mm 1/2 Circle Round Body	90	12	VIO
PGLA VIO 1 X 45-36MM HC RB (8 in 1)	PGN010741	1		36 mm 1/2 Circle Round Body	45	12	VIO
PGLA VIO 1 X 90-45MM HC BP	PGN01 9998	1		45 mm 1/2 Circle Blunt Point	90	12	VIO
PGLA VIO 1 X 100-65MM CU RB BLUNT	PGN019165	1		65 mm 3/8 Circle Blunt Point	100	12	VIO
PGLA VIO 1 X 75-30MM HC RB	PGN019214	1		30 mm 1/2 Circle Round Body	75	12	VIO
PGLA VIO 1 X 75-30MM HC RB	PGN019214	1		30 mm 1/2 Circle Round Body	75	12	VIO
PGLA VIO 1 X 150-48MM HC RC LOOP	PGN019150	1		48 mm 1/2 Circle Reverse Cutting	150	12	VIO
PGLA VIO 1 X 100-55MM HC RB	PGN01 2449	1		55 mm 1/2 Circle Round Body	100	12	VIO
PGLA UND 1 X 90-36MM HC RB	PGN01 947U	1		36 mm 1/2 Circle Round Body	90	12	UND
PGLA VIO 1 X 90-36MM HC TC	PGN01 2319	1		36 mm 1/2 Circle Taper Cut	90	12	VIO
PGLA VIO 1 X 75-31MM J TC	PGN01 9335	1		31 mm J Shape Taper Cut	75	12	VIO
PGLA VIO 1 X 70-45MM CU RB	PGN01 2326	1		45 mm 3/8 Circle Round Body	70	12	VIO
PGLA VIO 1 X 100-45MM HC RB	PGN01 2347XL	1		45 mm 1/2 Circle Round Body	100	12	VIO
PGLA VIO 1 X 70-40MM HC TC	PGN01 2363	1		40 mm 1/2 Circle Taper Cut	70	12	VIO
PGLA VIO 1 X 90-40MM CU RC	PGN01 2422	1		40 mm 3/8 Circle Reverse Cutting	90	12	VIO
PGLA VIO 1 X 75-22MM HC RC TH	PGN019287	1		22 mm 1/2 Circle Reverse Cutting	75	12	VIO
PGLA VIO 1 X 90-40MM HC RB BP	PGN019997	1		40 mm 1/2 Circle Blunt Point	90	12	VIO


Description	FG Code	Size		Needle Profile Shape, Type & Length	Suture Length (cm)	Pack Size	Suture Colour
MITSU (SYNTHETIC ABSORBABLE BRAIDED COATED POLYGLACTIN 910)							
PGLA VIO 1 X 75-110MM FE RB	PGN019111	1		110 mm 5/8 Circle Round Body	75	12	VIO
PGLA VIO 2 X 90-48MM HC TC	PGN021937	2		48 mm 1/2 Circle Taper Cut	90	12	VIO
PGLA VIO 2 X 90-40MM HC RC	PGN02 2478	2		40 mm 1/2 Circle Reverse Cutting	90	12	VIO
PGLA VIO 2 X 75-45MM HC RB	PGN02 9345N	2		45 mm 1/2 Circle Round Body	75	12	VIO
PGLA VIO 2 X 75-22MM HC RC	PGN029287	2		22 mm 1/2 Circle Reverse Cutting	75	12	VIO
PGLA VIO 2 X 75-48MM HC RC	PGN029297	2		48 mm 1/2 Circle Reverse Cutting	75	12	VIO
PGLA VIO 2 X 75-36MM HC TC	PGN029357	2		36 mm 1/2 Circle Taper Cut	75	12	VIO
PGLA VIO 2 X 75-45MM HC TC	PGN029378	2		45 mm 1/2 Circle Taper Cut	75	12	VIO
PGLA VIO 2 X 90-48MM HC RB	PGN029452	2		48 mm 1/2 Circle Round Body	90	12	VIO
PGLA VIO 2 X 90-48MM HC RC	PGN029497	2		48 mm 1/2 Circle Reverse Cutting	90	12	VIO
PGLA VIO 2 X 45-40MM HC RC	PGN020719	2		40 mm 1/2 Circle Reverse Cutting	45	12	VIO
PGLA VIO 2 X 75-40MM HC RB	PGN029240	2		40 mm 1/2 Circle Round Body	75	12	VIO
PGLA VIO 2 X 75-50MM HC RB	PGN029250	2		50 mm 1/2 Circle Round Body	75	12	VIO
PGLA VIO 2 X 75-35MM HC TC	PGN029257	2		35 mm 1/2 Circle Taper Cut	75	12	VIO
PGLA VIO 2 X 90-37MM HC RB	PGN029237	2		37 mm 1/2 Circle Round Body	90	12	VIO
PGLA VIO 2 X 75-36MM CU RC	PGN029236	2		36 mm 3/8 Circle Reverse Cutting	75	12	VIO
PGLA VIO 2 X 90-40MM HC RB	PGN02 2452	2		40 mm 1/2 Circle Round Body	90	12	VIO
PGLA VIO 2 X 75-48MM HC RB	PGN029252	2		48 mm 1/2 Circle Round Body	75	12	VIO

Description	FG Code	Size		Needle Profile Shape, Type & Length	Suture Length (cm)	Pack Size	Suture Colour
MITSU (SYNTHETIC ABSORBABLE BRAIDED COATED POLYGLACTIN 910)							
PGLA VIO 2 X 75-45MM HC RB H	PGN029246	2		45 mm 1/2 Circle Round Body	75	12	VIO
PGLA VIO 2 X 75-30MM HC RB	PGN029230	2		30 mm 1/2 Circle Round Body	75	12	VIO
PGLA VIO 10-0 X 10-6.5MM CU SM	PGN100 960	10-0		6.5 mm 3/8 Circle Spatula Micropoint	10	12	VIO
PGLA VIO 10-0 X 10-6MM CU SM DN	PGN100960DN	10-0		6 mm 3/8 Circle Spatula Micropoint Spatula Micropoint	10	12	VIO
PGLA VIO 2-0 X 90-55MM ST CT	PGN209713	2-0		55 mm Straight Conventional Cutting	90	12	VIO
PGLA VIO 2-0 X 90-45MM CU RC	PGN209282	2-0		45 mm 3/8 Circle Reverse Cutting	90	12	VIO
PGLA VIO 2-0 X 90-30MM HC RB	PGN20 2317	2-0		30 mm 1/2 Circle Round Body	90	12	VIO
PGLA VIO 2-0 X 140-30MM HC RB TC DN	PGN20 2317DN	2-0		30 mm 1/2 Circle Round Body Taper Cut	140	12	VIO
PGLA VIO 2-0 X 70-30MM HC RB	PGN20 2341	2-0		30 mm 1/2 Circle Round Body	70	12	VIO
PGLA VIO 2-0 X 90-40MM HC RB	PGN20 2345	2-0		40 mm 1/2 Circle Round Body	90	12	VIO
PGLA VIO 2-0 X 90-40MM HC RC	PGN20 2382	2-0		40 mm 1/2 Circle Reverse Cutting	90	12	VIO
PGLA UND 2-0 X 90-26MM CU CT	PGN20 2404U	2-0		26 mm 3/8 Circle Conventional Cutting	90	12	UND
PGLA VIO 2-0 X 45-26MM HC RB	PGN200775	2-0		26 mm 1/2 Circle Round Body	45	12	VIO
PGLA VIO 2-0 X 75-26MM HC RB	PGN209121	2-0		26 mm 1/2 Circle Round Body	75	12	VIO
PGLA VIO 2-0 X 70-26MM HC RB	PGN209123	2-0		26 mm 1/2 Circle Round Body	70	12	VIO
PGLA VIO 2-0 X 75-31MM HC RB	PGN209136	2-0		31 mm 1/2 Circle Round Body	75	12	VIO
PGLA VIO 2-0 X 75-36MM HC RB	PGN209140	2-0		36 mm 1/2 Circle Round Body	75	12	VIO
PGLA VIO 2-0 X 75-40MM HC RB	PGN209150	2-0		40 mm 1/2 Circle Round Body	75	12	VIO
PGLA VIO 2-0 X 75-65MM HC RB	PGN209158	2-0		65 mm 1/2 Circle Round Body	75	12	VIO

PGLA VIO 2-0 X 75-26MM HC TC	PGN20 9350	2-0		26 mm 1/2 Circle Taper Cut	75	12	VIO
PGLA VIO 2-0 X 75-40MM HC TC	PGN209359	2-0		40 mm 1/2 Circle Taper Cut	75	12	VIO
PGLA VIO 2-0 X 75-31MM HC TC	PGN209360	2-0		31 mm 1/2 Circle Taper Cut	75	12	VIO
PGLA VIO 2-0 X 75-36MM HC TC	PGN209363	2-0		36 mm 1/2 Circle Taper Cut	75	12	VIO
PGLA VIO 2-0 X 75-45MM HC TC	PGN209375	2-0		45 mm 1/2 Circle Taper Cut	75	12	VIO
PGLA VIO 2-0 X 90-36MM HC RB	PGN209440	2-0		36 mm 1/2 Circle Round Body	90	12	VIO
PGLA VIO 2-0 X 90-36MM HC TC	PGN209463	2-0		36 mm 1/2 Circle Taper Cut	90	12	VIO
PGLA VIO 2-0 X 90-40MM HC TC	PGN209466	2-0		40 mm 1/2 Circle Taper Cut	90	12	VIO
PGLA VIO 2-0 X 100-31MM HC RB	PGN209637	2-0		31 mm 1/2 Circle Round Body	100	12	VIO
PGLA UND 2-0 X 75-40MM HC TC	PGN209904	2-0		40 mm 1/2 Circle Taper Cut	75	12	UND
PGLA UND 2-0 X 75-36MM HC RB	PGN200258	2-0		36 mm 1/2 Circle Round Body	75	12	UND
PGLA UND 2-0 X 70-26MM HC RB	PGN200417	2-0		26 mm 1/2 Circle Round Body	70	12	UND
PGLA VIO 2-0 X 70-36MM HC RC	PGN200466	2-0		36 mm 1/2 Circle Reverse Cutting	70	12	VIO
PGLA UND 2-0 X 90-30MM CU RC	PGN200589	2-0		30 mm 3/8 Circle Reverse Cutting	90	12	UND
PGLA UND 2-0 X 70-26MM HC RC	PGN200869	2-0		26 mm 1/2 Circle Reverse Cutting	70	12	UND
PGLA VIO 2-0 X 45-17MM HC RB	PGN200306	2-0		17 mm 1/2 Circle Round Body	45	12	VIO
PGLA VIO 2-0 X 70-21MM HC RB	PGN200824	2-0		21 mm 1/2 Circle Round Body	70	12	VIO
PGLA VIO 2-0 X 70-22MM HC RC	PGN20461	2-0		22 mm 1/2 Circle Reverse Cutting	70	12	VIO
PGLA UND 2-0 X 70-24MM CU RC	PGN200443	2-0		24 mm 3/8 Circle Reverse Cutting	70	12	UND


Description	FG Code	Size		Needle Profile Shape, Type & Length	Suture Length (cm)	Pack Size	Suture Colour
MITSU (SYNTHETIC ABSORBABLE BRAIDED COATED POLYGLACTIN 910)							
PGLA VIO 2-0 X 45-36MM HC RC	PGN200710	2-0		36 mm 1/2 Circle Reverse Cutting	45	12	VIO
PGLA UND 2-0 X 70-36MM HC RB BP	PGN200259	2-0		36 mm 1/2 Circle Blunt Point	70	12	UND
PGLA VIO 2-0 X 70-36MM FE RB	PGN200375	2-0		36 mm 5/8 Circle Round Body	70	12	VIO
PGLA VIO 2-0 X 70-40MM HC RC	PGN200478	2-0		40 mm 1/2 Circle Reverse Cutting	70	12	VIO
PGLA VIO 2-0 X 70-40MM HC RB	PGN200351	2-0		40 mm 1/2 Circle Round Body	70	12	VIO
PGLA VIO 2-0 X 75-30MM HC RB	PGN200368	2-0		30 mm 1/2 Circle Round Body	75	12	VIO
PGLA VIO 2-0 X 75-30MM HC TC	PGN200517	2-0		30 mm 1/2 Circle Taper Cut	75	12	VIO
PGLA VIO 2-0 X 75-35MM HC RB BLUNT	PGN200738	2-0		35 mm 1/2 Circle Blunt Point	75	12	VIO
PGLA VIO 2-0 X 90-35MM CU RC	PGN209390	2-0		35 mm 3/8 Circle Reverse Cutting	90	12	VIO
PGLA VIO 2-0 X 75-25MM HC CT	PGN209225	2-0		25 mm 1/2 Circle Conventional Cutting	75	12	VIO
PGLA VIO 2-0 X 75-45MM HC CT	PGN209245	2-0		45 mm 1/2 Circle Conventional Cutting	75	12	VIO
PGLA VIO 2-0 X 75-35MM CU CT	PGN209235	2-0		35 mm 3/8 Circle Conventional Cutting	75	12	VIO
PGLA VIO 2-0 X 75-26MM CU CT	PGN209226	2-0		26 mm 3/8 Circle Conventional Cutting	75	12	VIO
PGLA VIO 2-0 X 75-24MM CU RB	PGN209224	2-0		24 mm 3/8 Circle Round Body	75	12	VIO
PGLA VIO 2-0 X 70-24MM CU RC	PGN200453	2-0		24 mm 3/8 Circle Reverse Cutting	70	12	VIO
PGLA VIO 2-0 X 90-48MM HC RB	PGN209290	2-0		48 mm 1/2 Circle Round Body	90	12	VIO
PGLA UND 2-0 X 90-36MM HC TC	PGN209464	2-0		36 mm 1/2 Circle Taper Cut	90	12	UND
PGLA VIO 2-0 X 70-31MM HC RC	PGN209231	2-0		31 mm 1/2 Circle Reverse Cutting	70	12	VIO

Description	FG Code	Size		Needle Profile Shape, Type & Length	Suture Length (cm)	Pack Size	Suture Colour
MITSU (SYNTHETIC ABSORBABLE BRAIDED COATED POLYGLACTIN 910)							
PGLA UND 2-0 X 90-36MM HC RB	PGN20 945U	2-0		36 mm 1/2 Circle Round Body	90	12	UND
PGLA VIO 2-0 X 90-30MM CU RC	PGN20 2316	2-0		30 mm 3/8 Circle Reverse Cutting	90	12	VIO
PGLA VIO 2-0 X 75-22MM HC TC	PGN20 9828	2-0		22 mm 1/2 Circle Taper Cut	75	12	VIO
PGLA VIO 2-0 X 45-30MM HC RB	PGN20 2317SL	2-0		30 mm 1/2 Circle Round Body	45	12	VIO
PGLA VIO 2-0 X 90-26MM HC RB	PGN20 2356	2-0		26 mm 1/2 Circle Round Body	90	12	VIO
PGLA VIO 2-0 X 35-26MM HC RB	PGN20 2356SL	2-0		26 mm 1/2 Circle Round Body	35	12	VIO
PGLA VIO 2-0 X 70-26MM HC RB VB	PGN20 2122	2-0		26 mm 1/2 Circle Round Body	70	12	VIO
PGLA VIO 2-0 X 90-30MM HC RB DN	PGN20 2318DN	2-0		30 mm 1/2 Circle Round Body Round Body	90	12	VIO
PGLA VIO 2-0 X 45-19MM SKKI RB	PGN20 2903	2-0		19 mm SKKI Round Body	45	12	VIO
PGLA UND 2-0 X 76-60MM ST CT	PGN20 2390U	2-0		60 mm Straight Conventional Cutting	76	12	UND
PGLA VIO 2-0 X 70-36MM HC RB	PGN20 323	2-0		36 mm 1/2 Circle Round Body	70	12	VIO
PGLA VIO 2-0 X 70-26MM FE RB	PGN20 602	2-0		26 mm 5/8 Circle Round Body	70	12	VIO
PGLA VIO 2-0 X 90-60MM ST CT	PGN209718	2-0		60 mm Straight Conventional Cutting	90	12	VIO
PGLA UND 2-0 X 75-26MM CU CT	PGN202404USL	2-0		26 mm 3/8 Circle Conventional Cutting	75	12	UND
PGLA VIO 2-0 X 75-36MM J TYPE RB	PGN209220J	2-0		36 mm J Shape Round Body	75	12	VIO
PGLA UND 2-0 X 75-36MM CU CT	PGN209532	2-0		36 mm 3/8 Circle Conventional Cutting	75	12	UND
PGLA UND 2-0 X 75-26MM CU RC	PGN209572	2-0		26 mm 3/8 Circle Reverse Cutting	75	12	UND
PGLA UND 2-0 X 45-26MM CU CT	PGN209522	2-0		26 mm 3/8 Circle Conventional Cutting	45	12	UND
PGLA VIO 2-0 X 75-60MM ST CT	PGN2013360	2-0		60 mm Straight Conventional Cutting	75	12	VIO

Description	FG Code	Size		Needle Profile Shape, Type & Length	Suture Length (cm)	Pack Size	Suture Colour
MITSU (SYNTHETIC ABSORBABLE BRAIDED COATED POLYGLACTIN 910)							
PGLA VIO 2-0 X 75-35MM CU RC	PGN209390SL	2-0		35 mm 3/8 Circle Reverse Cutting	75	12	VIO
PGLA VIO 2-0 X 90-26MM HC RB EX	PGN200369	2-0		26 mm 1/2 Circle Round Body	90	12	VIO
PGLA VIO 2-0 X 70-26MM FE RB EX	PGN200602	2-0		26 mm 5/8 Circle Round Body	70	12	VIO
PGLA UND 2-0 X 90-36MM HC RB EX	PGN200945	2-0		36 mm 1/2 Circle Round Body	90	12	UND
PGLA VIO 2-0 X 70-36MM HC RB EX	PGN200339	2-0		36 mm 1/2 Circle Round Body	70	12	VIO
PGLA VIO 2-0 X 75-31MM J Type RB	PGN20 9136J	2-0		31 mm J Shape Round Body	75	12	VIO
PGLA UND 2-0 X 70-51MM ST CT	PGN209273	2-0		51 mm Straight Conventional Cutting	70	12	UND
PGLA VIO 2-0 X 75-26MM HC RB-3D	PGN209121Z	2-0		26 mm 1/2 Circle Round Body	75	36	VIO
PGLA UND 2-0 X 70-24MM CU RC-3D	PGN200443Z	2-0		24 mm 3/8 Circle Reverse Cutting	70	36	UND
PGLA VIO 2-0 X 70-36MM HC RC-3D	PGN200466Z	2-0		36 mm 1/2 Circle Reverse Cutting	70	36	VIO
PGLA VIO 2-0 X 75-36MM HC RB-3D	PGN209140Z	2-0		36 mm 1/2 Circle Round Body	75	36	VIO
PGLA VIO 2-0 X 90-22MM HC RB-3D	PGN200312XLZ	2-0		22 mm 1/2 Circle Round Body	90	36	VIO
PGLA UND 2-0 X 70-40MM HC RC	PGN200267U	2-0		40 mm 1/2 Circle Reverse Cutting	70	12	UND
PGLA UND 2-0 X 75-55MM ST CT	PGN209713U	2-0		55 mm Straight Conventional Cutting	75	12	UND
PGLA VIO 2-0 X 75-17MM HC RB	PGN200306XL	2-0		17 mm 1/2 Circle Round Body	75	12	VIO
PGLA UND 2-0 X 70-22MM HC RC	PGN206934U	2-0		22 mm 1/2 Circle Reverse Cutting	70	12	UND
PGLA UND 2-0 X 90-40MM HC RB	PGN200957U	2-0		40 mm 1/2 Circle Round Body	90	12	UND
PGLA UND 3-0 X 90-25MM FE RB	PGN309730	3-0		25 mm 5/8 Circle Round Body	90	12	UND
PGLA VIO 3-0 X 90-25MM CU CT	PGN302643	3-0		25 mm 3/8 Circle Conventional Cutting	90	12	VIO

Meril Endo-Surgery Product Catalogue | **Wound Closure** 29

Description	FG Code	Size		Needle Profile Shape, Type & Length	Suture Length (cm)	Pack Size	Suture Colour
MITSU (SYNTHETIC ABSORBABLE BRAIDED COATED POLYGLACTIN 910)							
PGLA UND 3-0 X 90-25MM CU CT	PGN30 2328U	3-0		25 mm 3/8 Circle Conventional Cutting	90	12	UND
PGLA VIO 3-0 X 45-22MM CU CT	PGN30 2401	3-0		22 mm 3/8 Circle Conventional Cutting	45	12	VIO
PGLA UND 3-0 X 90-19MM CU CT	PGN30 2402U	3-0		19 mm 3/8 Circle Conventional Cutting	90	12	UND
PGLA VIO 3-0 X 70-20MM HC RB	PGN30 2437	3-0		20 mm 1/2 Circle Round Body	70	12	VIO
PGLA VIO 3-0 X 45-16MM CU CT	PGN30 2471	3-0		16 mm 3/8 Circle Conventional Cutting	45	12	VIO
PGLA VIO 3-0 X 90-22MM HC CT	PGN30 2472	3-0		22 mm 1/2 Circle Conventional Cutting	90	12	VIO
PGLA VIO 3-0 X 90-36MM HC TC	PGN30 2516	3-0		36 mm 1/2 Circle Taper Cut	90	12	VIO
PGLA VIO 3-0 X 75-26MM HC RB	PGN309120	3-0		26 mm 1/2 Circle Round Body	75	12	VIO
PGLA VIO 3-0 X 75-31MM HC RB	PGN30 9331N	3-0		31 mm 1/2 Circle Round Body	75	12	VIO
PGLA UND 3-0 X 45-19MM CU RC	PGN304970	3-0		19 mm 3/8 Circle Reverse Cutting	45	12	UND
PGLA UND 3-0 X 90-17MM HC TC	PGN305180	3-0		17 mm 1/2 Circle Taper Cut	90	12	UND
PGLA UND 3-0 X 45-16MM HC RC	PGN300657	3-0		16 mm 1/2 Circle Reverse Cutting	45	12	UND
PGLA VIO 3-0 X 70-17MM HC RB	PGN300305	3-0		17 mm 1/2 Circle Round Body	70	12	VIO
PGLA UND 3-0 X 45-19MM CU CT	PGN300824	3-0		19 mm 3/8 Circle Conventional Cutting	45	12	UND
PGLA UND 3-0 X 70-19MM CU RC	PGN300427	3-0		19 mm 3/8 Circle Reverse Cutting	70	12	UND
PGLA VIO 3-0 X 45-22MM HC RC	PGN3079013	3-0		22 mm 1/2 Circle Reverse Cutting	45	12	VIO
PGLA VIO 3-0 X 70-22MM HC RC	PGN300460	3-0		22 mm 1/2 Circle Reverse Cutting	70	12	VIO
PGLA VIO 3-0 X 70-22MM HC RB	PGN300311	3-0		22 mm 1/2 Circle Round Body	70	12	VIO
PGLA UND 3-0 X 90-24MM CU RC	PGN300442	3-0		24 mm 3/8 Circle Reverse Cutting	90	12	UND
PGLA VIO 3-0 X 45-26MM CU RC	PGN3078713	3-0		26 mm 3/8 Circle Reverse Cutting	45	12	VIO
PGLA UND 3-0 X 70-26MM HC RC	PGN300868	3-0		26 mm 1/2 Circle Reverse Cutting	70	12	UND
PGLA VIO 3-0 X 90-26MM HC RB	PGN300527	3-0		26 mm 1/2 Circle Round Body	90	12	VIO

Description	FG Code	Size		Needle Profile Shape, Type & Length	Suture Length (cm)	Pack Size	Suture Colour
MITSU (SYNTHETIC ABSORBABLE BRAIDED COATED POLYGLACTIN 910)							
PGLA UND 3-0 X 90-36MM HC RB	PGN300344	3-0		36 mm 1/2 Circle Round Body	90	12	UND
PGLA VIO 3-0 X 90-40MM HC RC	PGN30 9389	3-0		40 mm 1/2 Circle Reverse Cutting	90	12	VIO
PGLA VIO 3-0 X 90-40MM HC RB	PGN30 9388	3-0		40 mm 1/2 Circle Round Body	90	12	VIO
PGLA VIO 3-0 X 75-30MM HC RB	PGN30864	3-0		30 mm 1/2 Circle Round Body	75	12	VIO
PGLA VIO 3-0 X 75-25MM HC RB	PGN300528	3-0		25 mm 1/2 Circle Round Body	75	12	VIO
PGLA VIO 3-0 X 70-25MM HC RB	PGN300526	3-0		25 mm 1/2 Circle Round Body	70	12	VIO
PGLA VIO 3-0 X 75-17MM HC TC	PGN300216	3-0		17 mm 1/2 Circle Taper Cut	75	12	VIO
PGLA VIO 3-0 X 90-35MM HC RB	PGN300738	3-0		35 mm 1/2 Circle Round Body	90	12	VIO
PGLA VIO 3-0 X 90-25MM CU RB	PGN309325	3-0		25 mm 3/8 Circle Round Body	90	12	VIO
PGLA VIO 3-0 X 90-26MM HC CT	PGN309326	3-0		26 mm 1/2 Circle Conventional Cutting	90	12	VIO
PGLA VIO 3-0 X 75-26MM CU RC	PGN3078713XL	3-0		26 mm 3/8 Circle Reverse Cutting	75	12	VIO
PGLA VIO 3-0 X 70-60MM ST CT	PGN300664	3-0		60 mm Straight Conventional Cutting	70	12	VIO
PGLA VIO 3-0 X 70-19MM CU RC	PGN300398	3-0		19 mm 3/8 Circle Reverse Cutting	70	12	VIO
PGLA UND 3-0 X 75-30MM CU RC	PGN309330U	3-0		30 mm 3/8 Circle Reverse Cutting	75	12	UND
PGLA VIO 3-0 X 75-30MM CU RC	PGN309329	3-0		30 mm 3/8 Circle Reverse Cutting	75	12	VIO
PGLA VIO 3-0 X 70-30MM HC RC	PGN309330	3-0		30 mm 1/2 Circle Reverse Cutting	70	12	VIO
PGLA VIO 3-0 X 70-26MM HC RB VB	PGN30 2123	3-0		26 mm 1/2 Circle Round Body	70	12	VIO
PGLA VIO 3-0 X 45-17MM HC TC	PGN30 2515	3-0		17 mm 1/2 Circle Taper Cut	45	12	VIO
PGLA VIO 3-0 X 70-12MM HC RB	PGN30 2427	3-0		12 mm 1/2 Circle Round Body	70	12	VIO

34 **Wound Closure** | Meril Endo-Surgery Product Catalogue

Description	FG Code	Size		Needle Profile Shape, Type & Length	Suture Length (cm)	Pack Size	Suture Colour
MITSU (SYNTHETIC ABSORBABLE BRAIDED COATED POLYGLACTIN 910)							
PGLA UND 4-0 X 45-13MM CU CT	PGN4083565U	4-0		13 mm 3/8 Circle Conventional Cutting	45	12	UND
PGLA VIO 4-0 X 75-30MM HC RB	PGN409430	4-0		30 mm 1/2 Circle Round Body	75	12	VIO
PGLA UND 4-0 X 70-26MM HC RB	PGN400415U	4-0		26 mm 1/2 Circle Round Body	70	12	UND
PGLA UND 5-0 X 90 -08MM HC SM	PGN509567	5-0		8 mm 1/2 Circle Spatula Micropoint	90	12	UND
PGLA VIO 5-0 X 90-16MM CU RC	PGN509436	5-0		16 mm 3/8 Circle Reverse Cutting	90	12	VIO
PGLA VIO 5-0 X 45-16MM HC RB	PGN50 2303	5-0		16 mm 1/2 Circle Round Body	45	12	VIO
PGLA VIO 5-0 X 45-12MM CU RC	PGN50 2493	5-0		12 mm 3/8 Circle Reverse Cutting	45	12	VIO
PGLA UND 5-0 X 70-13MM HC RB	PGN500433	5-0		13 mm 1/2 Circle Round Body	70	12	UND
PGLA UND 5-0 X 45-08MM QC SM	PGN500671	5-0		8 mm 1/4 Circle Spatula Micropoint	45	12	UND
PGLA UND 5-0 X 45-11MM HC RC	PGN500511	5-0		11 mm 1/2 Circle Reverse Cutting	45	12	UND
PGLA UND 5-0 X 45-11MM CU RC	PGN500490	5-0		11 mm 3/8 Circle Reverse Cutting	45	12	UND
PGLA UND 5-0 X 45-13MM CU CT	PGN5083465	5-0		13 mm 3/8 Circle Conventional Cutting	45	12	UND
PGLA UND 5-0 X 45-16MM CU CT	PGN500844	5-0		16 mm 3/8 Circle Conventional Cutting	45	12	UND
PGLA UND 5-0 X 45-16MM HC RC	PGN500506	5-0		16 mm 1/2 Circle Reverse Cutting	45	12	UND
PGLA UND 5-0 X 70-17MM HC RB	PGN500213	5-0		17 mm 1/2 Circle Round Body	70	12	UND
PGLA UND 5-0 X 45-19MM CU RC	PGN500495	5-0		19 mm 3/8 Circle Reverse Cutting	45	12	UND
PGLA UND 5-0 X 45-16MM CU RC	PGN500500	5-0		16 mm 3/8 Circle Reverse Cutting	45	12	UND
PGLA VIO 5-0 X 45-19MM CU RC	PGN500391	5-0		19 mm 3/8 Circle Reverse Cutting	45	12	VIO
PGLA VIO 5-0 X 75-16MM CU CT	PGN500845	5-0		16 mm 3/8 Circle Conventional Cutting	75	12	VIO
PGLA VIO 5-0 X 75-26MM HC RB	PGN509526	5-0		26 mm 1/2 Circle Round Body	75	12	VIO
PGLA UND 5-0 X 45-13MM CU RC	PGN50 2872U	5-0		13 mm 3/8 Circle Reverse Cutting	45	12	UND
PGLA VIO 5-0 X 70-10MM CU RB	PGN50 2395	5-0		10 mm 3/8 Circle Round Body	70	12	VIO
PGLA VIO 5-0 X 45-16MM CU CT	PGN50 2442	5-0		16 mm 3/8 Circle Conventional Cutting	45	12	VIO

Description	FG Code	Size		Needle Profile Shape, Type & Length	Suture Length (cm)	Pack Size	Suture Colour
MITSU (SYNTHETIC ABSORBABLE BRAIDED COATED POLYGLACTIN 910)							
PGLA VIO 6-0 X 70-13MM HC RB	PGN600134	6-0		13 mm 1/2 Circle Round Body	70	12	VIO
PGLA VIO 6-0 X 45-8MM HC SM	PGN600562	6-0		8 mm 1/2 Circle Spatula Micropoint	45	12	VIO
PGLA UND 6-0 X 45-11MM HC RC	PGN600510	6-0		11 mm 1/2 Circle Reverse Cutting	45	12	UND
PGLA UND 6-0 X 45-13MM CU RC	PGN600492	6-0		13 mm 3/8 Circle Reverse Cutting	45	12	UND
PGLA UND 6-0 X 70-17MM HC RB	PGN600212	6-0		17 mm 1/2 Circle Round Body	70	12	UND
PGLA UND 6-0 X 45-13MM CU CT	PGN60 2831U	6-0		13 mm 3/8 Circle Conventional Cutting	45	12	UND
PGLA VIO 6-0 X 45-08MM QC SM DN	PGN60 9552DN	6-0		8 mm 1/4 Circle Spatula Micropoint Spatula Micropoint	45	12	VIO
PGLA VIO 6-0 X 30-10MM FE RB	PGN609610	6-0		10 mm 5/8 Circle Round Body	30	12	VIO
PGLA UND 6-0 X 45-8MM QC SM DN	PGN609752	6-0		8 mm 1/4 Circle Spatula Micropoint Spatula Micropoint	45	12	UND
PGLA UND 6-0 X 45-13MM HC RB	PGN609981	6-0		13 mm 1/2 Circle Round Body	45	12	UND
PGLA VIO 6-0 X 45-08MM QC SM DN C1	PGN600563	6-0		8 mm 1/4 Circle Spatula Micropoint Spatula Micropoint	45	12	VIO
PGLA UND 6-0 X 45-13MM CU CT EX	PGN6083365	6-0		13 mm 3/8 Circle Conventional Cutting	45	12	UND
PGLA VIO 6-0 X 45-12MM CU RC	PGN60 2671	6-0		12 mm 3/8 Circle Reverse Cutting	45	12	VIO
PGLA UND 6-0 X 45-8MM HC SM DN	PGN609756DN	6-0		8 mm 1/2 Circle Spatula Micropoint Spatula Micropoint	45	12	UND
PGLA UND 6-0 X 45-11MM CU SM DN	PGN609760DN	6-0		11 mm 3/8 Circle Spatula Micropoint Spatula Micropoint	45	12	UND
PGLA VIO 6-0 X 45-8MM HC SM DN	PGN609566DN	6-0		8 mm 1/2 Circle Spatula Micropoint Spatula Micropoint	45	12	VIO
PGLA VIO 6-0 X 45-11MM CU RC DN	PGN609537	6-0		11 mm 3/8 Circle Reverse Cutting Reverse Cutting	45	12	VIO
PGLA VIO 6-0 X 45-8MM QC SM DN-3D	PGN602670DNZ	6-0		8 mm 1/4 Circle Spatula Micropoint Spatula Micropoint	45	36	VIO
PGLA UND 6-0 X 45-16MM CU RC	PGN600499U	6-0		16 mm 3/8 Circle Reverse Cutting	45	12	UND
PGLA UND 6-0 X 45-11MM CU RC	PGN600489U	6-0		11 mm 3/8 Circle Reverse Cutting	45	12	UND

Description	FG Code	Size		Needle Profile Shape, Type & Length	Suture Length (cm)	Pack Size	Suture Colour
MITSU (SYNTHETIC ABSORBABLE BRAIDED COATED POLYGLACTIN 910)							
PGLA VIO 6-0 X 70-17MM HC RB	PGN600212V	6-0		17 mm 1/2 Circle Round Body	70	12	VIO
PGLA VIO 7-0 X 30-06MM CU SM	PGN70 2561	7-0		6 mm 3/8 Circle Spatula Micropoint	30	12	VIO
PGLA VIO 7-0 X 45-7MM HC SM	PGN709565	7-0		7 mm 1/2 Circle Spatula Micropoint	45	12	VIO
PGLA VIO 7-0 X 30-6.5MM CU SM DN	PGN70 9561DN	7-0		6.5 mm 3/8 Circle Spatula Micropoint Spatula Micropoint	30	12	VIO
PGLA VIO 7-0 X 90-6.5MM CU SM DN	PGN709561	7-0		6.5 mm 3/8 Circle Spatula Micropoint Spatula Micropoint	90	12	VIO
PGLA VIO 7-0 X 30-5MM HC SM	PGN700575	7-0		5 mm 1/2 Circle Spatula Micropoint	30	12	VIO
PGLA VIO 7-0 X 45-7MM HC SM DN	PGN709565DN	7-0		7 mm 1/2 Circle Spatula Micropoint Spatula Micropoint	45	12	VIO
PGLA VIO 8-0 X 30-9MM CU RB	PGN80 9577	8-0		9 mm 3/8 Circle Round Body	30	12	VIO
PGLA VIO 8-0 X 30-6MM CU SM DN	PGN80 2348DN	8-0		6 mm 3/8 Circle Spatula Micropoint Spatula Micropoint	30	12	VIO
PGLA VIO 8-0 X 30-08MM CU RB	PGN80 2301	8-0		8 mm 3/8 Circle Round Body	30	12	VIO
PGLA VIO 8-0 X 30-6.5MM CU SM DN	PGN80 9560DN	8-0		6.5 mm 3/8 Circle Spatula Micropoint Spatula Micropoint	30	12	VIO
PGLA VIO 8-0 X 90-7.5MM HC SM DN	PGN809564	8-0		7.5 mm 1/2 Circle Spatula Micropoint Spatula Micropoint	90	12	VIO
PGLA VIO 8-0 X 45-6MM CU SM DN	PGN809559	8-0		6 mm 3/8 Circle Spatula Micropoint Spatula Micropoint	45	12	VIO
PGLA VIO 8-0 X 30-6MM CU SM	PGN809558	8-0		6 mm 3/8 Circle Spatula Micropoint	30	12	VIO
PGLA VIO 8-0 X 45-6MM CU SM	PGN800547	8-0		6 mm 3/8 Circle Spatula Micropoint	45	12	VIO
PGLA VIO 9-0 X 10-6.5MM CU SM	PGN900950	9-0		6.5 mm 3/8 Circle Spatula Micropoint	10	12	VIO

Description	FG Code	Size	Needle Profile Shape, Type & Length	Suture Length (cm)	Pack Size	Suture Colour
MITSU (SYNTHETIC ABSORBABLE BRAIDED COATED POLYGLACTIN 910)						
PGLA VIO 0 X 350 (70X5)	PGN100208	0	Non needed	70	12	VIO
PGLA VIO 0 X 540 (12X45)	PGN100209	0	Non needed	45	12	VIO
PGLA VIO 0 X 75	PGN100675	0	Non needed	75	12	VIO
PGLA VIO 0 X 135 (45X3)	PGNS10 2601	0	Non needed	45	12	VIO
PGLA VIO 0 X 150	PGN109026	0	Non needed	150	12	VIO
PGLA VIO 0 X 60	PGN100114	0	Non needed	60	12	VIO
PGLA VIO 0 X 50 (50X3)	PGN100617	0	Non needed	50	12	VIO
PGLA VIO 0 X 250	PGN100250	0	Non needed	250	12	VIO
PGLA VIO 0 X 45 (45X3)	PGN100616	0	Non needed	45	12	VIO
PGLA VIO 0 X 60 (60X12)	PGNS100034	0	Non needed	60	12	VIO
PGLA UND 0 X 180 (60X6)	PGN100207U	0	Non needed	60	12	UND
PGLA UND 0 X 45 (12X45)	PGNS100912U	0	Non needed	45	12	UND
PGLA VIO 0 X 45 (45X6)	PGN101227	0	Non needed	45	12	VIO
PGLA VIO 1 X 45 (45X3)	PGN016173	1	Non needed	45	12	VIO
PGLA VIO 1 X 350 (70X5)	PGN016175	1	Non needed	70	12	VIO
PGLA VIO 1 X 75	PGN010908	1	Non needed	75	12	VIO
PGLA VIO 1 X 150	PGN019027	1	Non needed	150	12	VIO
PGLA VIO 1 X 60	PGN010115	1	Non needed	60	12	VIO
PGLA VIO 1 X 250-2D	PGN01250Y	1	Non needed	250	24	VIO
PGLA VIO 1 X 250 (50X5)	PGN010617	1	Non needed	50	12	VIO
PGLA VIO 1 X 540 (45X12)	PGN014512	1	Non needed	45	12	VIO
PGLA VIO 1 X 45 (45X6)	PGNS011228	1	Non needed	45	12	VIO
PGLA VIO 1 X 60 (60X13)	PGNS016013	1	Non needed	60	12	VIO
PGLA VIO 1 X 45 (45X13)	PGNS014513	1	Non needed	45	12	VIO
PGLA UND 1 X 45 (45x6)	PGNS014506U	1	Non needed	45	12	VIO
PGLA VIO 2 X 45 (45X3)	PGN026183	2	Non needed	45	12	VIO
PGLA VIO 2 X 45(12X45)	PGNS020908	2	Non needed	45	12	VIO
PGLA UND 2 X 45 (12X45)	PGNS020914U	2	Non needed	45	12	UND
PGLA VIO 2 X 150	PGN020116	2	Non needed	150	12	VIO
PGLA VIO 2 X 80	PGNS02080	2	Non needed	80	12	VIO
PGLA VIO 1 X 250	PGN01250	1	Non needed	250	12	VIO
PGLA VIO 2-0 X 45	PGN200105	2-0	Non needed	45	12	VIO
PGLA VIO 2-0 X 150 (50X3)	PGN200616	2-0	Non needed	50	12	VIO
PGLA VIO 2-0 X 540 (45X12)	PGN200617	2-0	Non needed	45	12	VIO
PGLA VIO 2-0 X 250 (50X5)	PGN200613	2-0	Non needed	50	12	VIO
PGLA VIO 2-0 X 350 (70X5)	PGN201216	2-0	Non needed	35	12	VIO
PGLA VIO 2-0 X 270 (45X6)	PGN201226	2-0	Non needed	45	12	VIO
PGLA VIO 2-0 X 135 (45X3)	PGNS20 2602	2-0	Non needed	45	12	VIO
PGLA VIO 2-0 X 150	PGN209025	2-0	Non needed	150	12	VIO
PGLA VIO 2-0 X 45 (45X3)	PGN200615	2-0	Non needed	45	12	VIO
PGLA VIO 2-0 X 60 (60X12)	PGNS20033	2-0	Non needed	60	12	VIO
PGLA VIO 2-0 X 75 (75X6)	PGNS200039	2-0	Non needed	75	12	VIO
PGLA VIO 2-0 X 250-2D	PGN200223Y	2-0	Non needed	250	24	VIO
PGLA VIO 2-0 X 250	PGN200223	2-0	Non needed	250	12	VIO

Description	FG Code	Size	Needle Profile Shape, Type & Length	Suture Length (cm)	Pack Size	Suture Colour
MITSU (SYNTHETIC ABSORBABLE BRAIDED COATED POLYGLACTIN 910)						
PGLA UND 2-0 X 135 (45X3)	PGN200286U	2-0	Non needled	45	12	UND
PGLA VIO 2-0 X 80	PGNS20080	2-0	Non needled	80	12	VIO
PGLA VIO 2-0 X 60 (60X13)	PGNS206013	2-0	Non needled	60	12	VIO
PGLA UND 2-0 X 45 (45x12)	PGNS204512U	2-0	Non needled	45	12	VIO
PGLA VIO 3-0 X 45 (45X3)	PGN300414	3-0	Non needled	45	12	VIO
PGLA VIO 3-0 X 45 (45X6)	PGN300180	3-0	Non needled	45	12	VIO
PGLA VIO 3-0 X 135 (45X3)	PGNS30 2603	3-0	Non needled	45	12	VIO
PGLA VIO 3-0 X 250 (50X5)	PGN300614	3-0	Non needled	50	12	VIO
PGLA VIO 3-0 X 540 (45X12)	PGN300904	3-0	Non needled	45	12	VIO
PGLA VIO 3-0 X 150	PGN309024	3-0	Non needled	150	12	VIO
PGLA VIO 3-0 X 70 (70X5)	PGNS308445	3-0	Non needled	70	12	VIO
PGLA VIO 3-0 X 60 (60X12)	PGNS300032	3-0	Non needled	60	12	VIO
PGLA VIO 3-0 X 80	PGNS30080	3-0	Non needled	80	12	VIO
PGLA VIO 4-0 X 45 (45X3)	PGN400515	4-0	Non needled	45	12	VIO
PGLA VIO 4-0 X 135 (45X3)	PGNS40 2604	4-0	Non needled	45	12	VIO
PGLA VIO 4-0 X 45 (45X12)	PGN400903	4-0	Non needled	45	12	VIO
PGLA VIO 4-0 X 45(45X6)	PGNS401224	4-0	Non needled	45	12	VIO
PGLA UND 4-0 X 45 (45X6)	PGNS400977U	4-0	Non needled	45	12	UND
PGLA UND 4-0 X 135 (45X3)	PGN400284U	4-0	Non needled	45	12	UND
PGLA VIO 4-0 X 80	PGNS40080	4-0	Non needled	80	12	VIO
PGLA VIO 5-0 X 150	PGNS509022	5-0	Non needled	150	12	VIO

MITSU C+™

Polyglactin 910 Suture with Chlorhexidine

Properties

- Retention of 75% knot tensile strength after 14 days
- Retention of 50% knot tensile strength after 21 days
- Complete mass absorption within 56 to 70 days

Indications


- Gastrointestinal procedures
- Gynaecology/Obstetric procedures
- Orthopedic procedures
- Urology procedures
- Skin closure (intracutaneous, subcutaneous)
- Ligatures

Advantages




- Smooth knot run-down/Knot slide
- Gentle knot repositioning
- Excellent knot security
- Minimal inflammatory reaction



Product Overview

Structure	: Multifilament, Braided
Colour	: Dyed, Undyed
	
Chemical composition	: Poly(glycolide-co-L-lactide 90/10)
Coating	: Poly(glycolide-co-L-lactide 30/70) + Calcium Stearate + antibacterial agent Chlorhexidine
Sizes	: USP 5/0 (metric 1) - USP 2 (metric 5)
Type of absorption	: Hydrolysis
Sterilization	: Ethylene oxide

Description	FG Code	Size		Needle Profile Shape, Type & Length	Suture Length (cm)	Pack Size	Suture Colour
MITSU C+ (SYNTHETIC ABSORBABLE BRAIDED COATED POLYGLACTIN 910 WITH CHLORHEXIDINE)							
PGLA C PLUS VIO 0 X 90-30MM HC RB	PGP10 2338	0		30 mm 1/2 Circle Round Body	90	12	VIO
PGLA C PLUS VIO 0 X 90-40MM HC RB TN	PGP10 2346	0		40 mm 1/2 Circle Round Body	90	12	VIO
PGLA C PLUS VIO 0 X 90-40MM HC RB TH	PGP10 2358	0		40 mm 1/2 Circle Round Body	90	12	VIO
PGLA C PLUS VIO 0 X 90-40MM HC TC	PGP10 2518	0		40 mm 1/2 Circle Taper Cut	90	12	VIO
PGLA C PLUS VIO 0 X 90-36MM HC RC	PGP10 2534	0		36 mm 1/2 Circle Reverse Cutting	90	12	VIO
PGLA C PLUS VIO 0 X 70-26MM HC RB	PGP100334	0		26 mm 1/2 Circle Round Body	70	12	VIO
PGLA C PLUS VIO 0 X 70-36MM HC RB	PGP100341	0		36 mm 1/2 Circle Round Body	70	12	VIO
PGLA C PLUS VIO 0 X 45-48MM HC RB	PGP100764	0		48 mm 1/2 Circle Round Body	45	12	VIO
PGLA C PLUS UND 0 X 90-40MM HC TC	PGP109901	0		40 mm 1/2 Circle Taper Cut	90	12	UND
PGLA C PLUS VIO 0 X 90-36MM HC TC	PGP10 2517	0		36 mm 1/2 Circle Taper Cut	90	12	VIO
PGLA C PLUS UND 0 X 90-36MM HC RB	PGP10 946U	0		36 mm 1/2 Circle Round Body	90	12	UND
PGLA C PLUS VIO 0 X 70-36MM FE RB	PGP100376	0		36 mm 5/8 Circle Round Body	70	12	VIO
PGLA C PLUS VIO 0 X 70-26MM FE RB	PGP100603	0		26 mm 5/8 Circle Round Body	70	12	VIO
PGLA C PLUS VIO 1 X 90-40MM HC RB TN	PGP01 2347	1		40 mm 1/2 Circle Round Body	90	12	VIO
PGLA C PLUS VIO 1 X 120-36MM HC RB	PGP01 2350	1		36 mm 1/2 Circle Round Body	120	12	VIO
PGLA C PLUS VIO 1 X 90-40MM HC RC	PGP01 2421	1		40 mm 1/2 Circle Reverse Cutting	90	12	VIO

Description	FG Code	Size		Needle Profile Shape, Type & Length	Suture Length (cm)	Pack Size	Suture Colour
MITSU C+ (SYNTHETIC ABSORBABLE BRAIDED COATED POLYGLACTIN 910 WITH CHLORHEXIDINE)							
PGLA C PLUS VIO 1 X 35-23MM HC RC	PGP01 2826	1		23 mm 1/2 Circle Reverse Cutting	35	12	VIO
PGLA C PLUS VIO 1 X 90-40MM HC RB TH	PGP01 2360	1		40 mm 1/2 Circle Round Body	90	12	VIO
PGLA C PLUS VIO 1 X 70-31MM HC RB	PGP010080	1		31 mm 1/2 Circle Round Body	70	12	VIO
PGLA C PLUS VIO 1 X 70-36MM HC RB	PGP010342	1		36 mm 1/2 Circle Round Body	70	12	VIO
PGLA C PLUS VIO 1 X 45-48MM HC RB	PGP010765	1		48 mm 1/2 Circle Round Body	45	12	VIO
PGLA C PLUS VIO 1 X 70-48MM HC RB	PGP010365	1		48 mm 1/2 Circle Round Body	70	12	VIO
PGLA C PLUS UND 1 X 90-48MM HC RB	PGP010977	1		48 mm 1/2 Circle Round Body	90	12	UND
PGLA C PLUS VIO 1 X 70-23MM HC RC	PGP01 2826XL	1		23 mm 1/2 Circle Reverse Cutting	70	12	VIO
PGLA C PLUS VIO 1 X 70-40MM HC RB	PGP01 2359	1		40 mm 1/2 Circle Round Body	70	12	VIO
PGLA C PLUS VIO 1 X 90-40MM HC TC	PGP01 2519	1		40 mm 1/2 Circle Taper Cut	90	12	VIO
PGLA C PLUS UND 1 X 90-40MM HC TC	PGP019374	1		40 mm 1/2 Circle Taper Cut	90	12	UND
PGLA C PLUS VIO 1 X 90-40MM HC TC EX	PGP019373	1		40 mm 1/2 Circle Taper Cut	90	12	VIO
PGLA C PLUS VIO 2 X 90-48MM HC RB	PGP020372	2		48 mm 1/2 Circle Round Body	90	12	VIO
PGLA C PLUS VIO 2 X 90-40MM HC RC	PGP02 2478	2		40 mm 1/2 Circle Reverse Cutting	90	12	VIO
PGLA C PLUS VIO 2-0 X 90-30MM HC RB	PGP20 2317	2-0		30 mm 1/2 Circle Round Body	90	12	VIO
PGLA C PLUS VIO 2-0 X 90-40MM HC RB	PGP20 2357	2-0		40 mm 1/2 Circle Round Body	90	12	VIO

Description	FG Code	Size		Needle Profile Shape, Type & Length	Suture Length (cm)	Pack Size	Suture Colour
MITSU C+ (SYNTHETIC ABSORBABLE BRAIDED COATED POLYGLACTIN 910 WITH CHLORHEXIDINE)							
PGLA C PLUS VIO 2-0 X 90-40MM HC RC	PGP20 2382	2-0		40 mm 1/2 Circle Reverse Cutting	90	12	VIO
PGLA C PLUS VIO 2-0 X 70-22MM HC RB	PGP200328	2-0		22 mm 1/2 Circle Round Body	70	12	VIO
PGLA C PLUS VIO 2-0 X 70-26MM HC RB	PGP200333	2-0		26 mm 1/2 Circle Round Body	70	12	VIO
PGLA C PLUS VIO 2-0 X 70-36MM HC RB	PGP200339	2-0		36 mm 1/2 Circle Round Body	70	12	VIO
PGLA C PLUS VIO 2-0 X 90-36MM HC TC	PGP200517	2-0		36 mm 1/2 Circle Taper Cut	90	12	VIO
PGLA C PLUS UND 2-0 X 90-26MM CU CT	PGP20 2404U	2-0		26 mm 3/8 Circle Conventional Cutting	90	12	UND
PGLA C PLUS UND 2-0 X 90-36MM HC RB	PGP20 945U	2-0		36 mm 1/2 Circle Round Body	90	12	UND
PGLA C PLUS VIO 2-0 X 70-36MM FE RB	PGP20375	2-0		36 mm 5/8 Circle Round Body	70	12	VIO
PGLA C PLUS VIO 2-0 X 70-26MM FE RB	PGP206023	2-0		26 mm 5/8 Circle Round Body	70	12	VIO
PGLA C PLUS UND 3-0 X 90-25MM CU CT	PGP30 2328U	3-0		25 mm 3/8 Circle Conventional Cutting	90	12	UND
PGLA C PLUS VIO 3-0 X 70-20MM HC RB	PGP30 2437	3-0		20 mm 1/2 Circle Round Body	70	12	VIO
PGLA C PLUS VIO 3-0 X 90-22MM HC CT	PGP30 2472	3-0		22 mm 1/2 Circle Conventional Cutting	90	12	VIO
PGLA C PLUS UND 3-0 X 70-17MM HC RB	PGP300215	3-0		17 mm 1/2 Circle Round Body	70	12	UND
PGLA C PLUS VIO 3-0 X 70-22MM HC RB	PGP300327	3-0		22 mm 1/2 Circle Round Body	70	12	VIO
PGLA C PLUS UND 3-0 X 70-24MM CU RC	PGP304423	3-0		24 mm 3/8 Circle Reverse Cutting	70	12	UND
PGLA C PLUS VIO 3-0 X 45-26MM HC RB	PGP300773	3-0		26 mm 1/2 Circle Round Body	45	12	VIO
PGLA C PLUS VIO 3-0 X 70-26MM HC RB	PGP300332	3-0		26 mm 1/2 Circle Round Body	70	12	VIO
PGLA C PLUS VIO 3-0 X 75-31MM HC RB	PGP300245	3-0		31 mm 1/2 Circle Round Body	75	12	VIO
PGLA C PLUS VIO 3-0 X 90-36MM HC TC	PGP30 2516	3-0		36 mm 1/2 Circle Taper Cut	90	12	VIO

Description	FG Code	Size		Needle Profile Shape, Type & Length	Suture Length (cm)	Pack Size	Suture Colour
MITSU C+ (SYNTHETIC ABSORBABLE BRAIDED COATED POLYGLACTIN 910 WITH CHLORHEXIDINE)							
PGLA C PLUS UND 3-0 X 70-24MM CU CT	PGP30 936U	3-0		24 mm 3/8 Circle Conventional Cutting	70	12	UND
PGLA C PLUS VIO 3-0 X 45-17MM HC RB	PGP300306	3-0		17 mm 1/2 Circle Round Body	45	12	VIO
PGLA C PLUS VIO 3-0 X 70-36MM HC RB	PGP300338	3-0		36 mm 1/2 Circle Round Body	70	12	VIO
PGLA C PLUS UND 3-0 X 70-60MM ST CT	PGP300663U	3-0		60 mm Straight Conventional Cutting	70	12	UND
PGLA C PLUS VIO 4-0 X 70-20MM HC RB	PGP40 2304	4-0		20 mm 1/2 Circle Round Body	70	12	VIO
PGLA C PLUS VIO 4-0 X 70-17MM HC RB	PGP400304	4-0		17 mm 1/2 Circle Round Body	70	12	VIO
PGLA C PLUS UND 4-0 X 70-60MM ST CT	PGP400662U	4-0		60 mm Straight Conventional Cutting	70	12	UND
PGLA C PLUS VIO 4-0 X 70-26MM HC RB	PGP400774	4-0		26 mm 1/2 Circle Round Body	70	12	VIO
PGLA C PLUS VIO 5-0 X 70-17MM HC RB	PGP500303	5-0		17 mm 1/2 Circle Round Body	70	12	VIO
PGLA C PLUS VIO 5-0 X 45-16MM HC RB	PGP50 2303	5-0		16 mm 1/2 Circle Round Body	45	12	VIO
PGLA C PLUS VIO 0 X 350 (70X5)	PGP101217	0		Non needled	70	12	VIO

MITSU FST™

Polyglactin 910 Fast Suture

Properties

- 50% of the original is retained until initial 5 Days
- approximately 10-14 Days post implantation, the original tensile strength is completely lost
- Complete absorption usually takes place by 42 Days

Indications

- Gynaecology/Obstetric procedures (e.g.: episiotomy repair)
- Urology procedures
- Ophthalmic procedures (e.g.: conjunctiva suturing)
- Odontology procedures (e.g.: oral mucosa)
- Skin closure (intracutaneous, subcutaneous)
- Ligatures

Advantages

- High initial tensile strength
- Predictable and constant degradation rate
- Excellent knot security
- Excellent handling properties
- Quick mass absorption



Product Overview

Structure	: Multifilament, Braided
Colour	: Undyed
Chemical composition	: Poly(glycolide-co-l-lactide 90/10)
Sizes	: USP 8/0 (metric 0.7) - USP 2 (metric 5)
Type of absorption	: Hydrolysis
Sterilization	: Gamma radiation

Description	FG Code	Size		Needle Profile Shape, Type & Length	Suture Length (cm)	Pack Size	Suture Colour
MITSU FST (SYNTHETIC ABSORBABLE BRAIDED COATED POLYGLACTIN 910 FAST)							
PGLA FAST UND 0 X 110-36MM HC TC	PGF10 2763	0		36 mm 1/2 Circle Taper Cut	110	12	UND
PGLA FAST UND 0 X 90-40MM HC RB	PGF109977	0		40 mm 1/2 Circle Round Body	90	12	UND
PGLA FAST UND 0 X 90-48MM HC RB	PGF1016951	0		48 mm 1/2 Circle Round Body	90	12	UND
PGLA FAST UND 0 X 90-36MM HC RB	PGF109036	0		36 mm 1/2 Circle Round Body	90	12	UND
PGLA FAST UND 0 X 60-36MM CU RC	PGF10 2721	0		36 mm 3/8 Circle Reverse Cutting	60	12	UND
PGLA FAST UND 0 X 110-40MM HC RB	PGF10 2763RB	0		40 mm 1/2 Circle Round Body	110	12	UND
PGLA FAST UND 0 X 90-36MM HC TC	PGF10 9963	0		36 mm 1/2 Circle Taper Cut	90	12	UND
PGLA FAST UND 0 X 90-30MM CU RC-3D	PGF102287Z	0		30 mm 3/8 Circle Reverse Cutting	90	36	UND
PGLA FAST UND 1 X 90-40MM HC TC	PGF016190	1		40 mm 1/2 Circle Taper Cut	90	12	UND
PGLA FAST UND 1 X 90-40MM HC RB	PGF016191	1		40 mm 1/2 Circle Round Body	90	12	UND
PGLA FAST UND 1 X 90-36MM HC RB	PGF016192	1		36 mm 1/2 Circle Round Body	90	12	UND
PGLA FAST UND 1 X 90-48MM HC RB	PGF016193	1		48 mm 1/2 Circle Round Body	90	12	UND
PGLA FAST UND 1 X 90-36MM HC TC	PGF019190	1		36 mm 1/2 Circle Taper Cut	90	12	UND
PGLA FAST UND 0 X 90-45MM HC TC	PGF109947	0		45 mm 1/2 Circle Taper Cut	90	12	UND
PGLA FAST UND 2-0 X 90-40MM HC TC	PGF209946	2-0		40 mm 1/2 Circle Taper Cut	90	12	UND
PGLA FAST UND 2-0 X 90-35MM CU RC	PGF209938	2-0		35 mm 3/8 Circle Reverse Cutting	90	12	UND
PGLA FAST UND 2-0 X 100-36MM HC TC	PGF20 2761	2-0		36 mm 1/2 Circle Taper Cut	100	12	UND

Description	FG Code	Size		Needle Profile Shape, Type & Length	Suture Length (cm)	Pack Size	Suture Colour
MITSU FST (SYNTHETIC ABSORBABLE BRAIDED COATED POLYGLACTIN 910 FAST)							
PGLA FAST UND 2-0 X 90-36MM HC TC	PGF20 2762	2-0		36 mm 1/2 Circle Taper Cut	90	12	UND
PGLA FAST UND 2-0 X 140-36MM HC RC RB DN	PGF20 2777DN	2-0		36 mm 1/2 Circle Reverse Cutting Round Body	140	12	UND
PGLA FAST UND 2-0 X 75-31MM HC RB	PGF209975	2-0		31 mm 1/2 Circle Round Body	75	12	UND
PGLA FAST UND 2-0 X 90-48MM HC RB	PGF209978	2-0		48 mm 1/2 Circle Round Body	90	12	UND
PGLA FAST UND 2-0 X 70-26MM HC RB	PGF204171	2-0		26 mm 1/2 Circle Round Body	70	12	UND
PGLA FAST UND 2-0 X 90-36MM HC RB	PGF209451	2-0		36 mm 1/2 Circle Round Body	90	12	UND
PGLA FAST UND 2-0 X 135-36MM HC RB RC DN	PGF209321	2-0		36 mm 1/2 Circle Round Body Reverse Cutting	135	12	UND
PGLA FAST UND 2-0 X 90-36MM HC RC	PGF209236	2-0		36 mm 1/2 Circle Reverse Cutting	90	12	UND
PGLA FAST UND 2-0 X 70-19MM CU RC	PGF200428	2-0		19 mm 3/8 Circle Reverse Cutting	70	12	UND
PGLA FAST UND 2-0 X 60-36MM CU RC	PGF20 2720	2-0		36 mm 3/8 Circle Reverse Cutting	60	12	UND
PGLA FAST UND 2-0 X 110-36MM HC TC	PGF20 2762XL	2-0		36 mm 1/2 Circle Taper Cut	110	12	UND
PGLA FAST UND 2-0 X 120-36MM HC TC	PGF209948	2-0		36 mm 1/2 Circle Taper Cut	120	12	UND
PGLA FAST UND 2-0 X 75-26MM CU RC	PGF20 9933	2-0		26 mm 3/8 Circle Reverse Cutting	75	12	UND
PGLA FAST UND 2-0 X 100-36MM HC TC SP	PGF20 2761SP	2-0		36 mm 1/2 Circle Taper Cut	100	12	UND
PGLA FAST UND 2-0 X 90-36MM HC TC SP	PGF20 2762SP	2-0		36 mm 1/2 Circle Taper Cut	90	12	UND

Description	FG Code	Size		Needle Profile Shape, Type & Length	Suture Length (cm)	Pack Size	Suture Colour
MITSU FST (SYNTHETIC ABSORBABLE BRAIDED COATED POLYGLACTIN 910 FAST)							
PGLA FAST UND 2-0 X 90-26MM CU RC	PGF209941	2-0		26 mm 3/8 Circle Reverse Cutting	90	12	UND
PGLA FAST UND 2-0 X 90-26MM CU RC-3D	PGF209941Z	2-0		26 mm 3/8 Circle Reverse Cutting	90	36	UND
PGLA FAST UND 2-0 X 75-24MM CU RC-3D	PGF209275Z	2-0		24 mm 3/8 Circle Reverse Cutting	75	36	UND
PGLA FAST UND 2-0 X 75-60MM ST CT	PGF209944	2-0		60 mm Straight Conventional Cutting	75	12	UND
PGLA FAST UND 3-0 X 75-16MM CU CT	PGF30 9919	3-0		16 mm 3/8 Circle Conventional Cutting	75	12	UND
PGLA FAST UND 3-0 X 90-60MM ST CT	PGF309943	3-0		60 mm Straight Conventional Cutting	90	12	UND
PGLA FAST UND 3-0 X 90-25MM FE RB	PGF309929	3-0		25 mm 5/8 Circle Round Body	90	12	UND
PGLA FAST UND 3-0 X 90-35MM CU RC	PGF309937	3-0		35 mm 3/8 Circle Reverse Cutting	90	12	UND
PGLA FAST UND 3-0 X 70-26MM CU RC	PGF30 2732	3-0		26 mm 3/8 Circle Reverse Cutting	70	12	UND
PGLA FAST UND 3-0 X 45-26MM CU RC	PGF30 2732SL	3-0		26 mm 3/8 Circle Reverse Cutting	45	12	UND
PGLA FAST UND 3-0 X 75-19MM CU RC	PGF309923	3-0		19 mm 3/8 Circle Reverse Cutting	75	12	UND
PGLA FAST UND 3-0 X 75-22MM HC RB	PGF309974	3-0		22 mm 1/2 Circle Round Body	75	12	UND
PGLA FAST UND 3-0 X 90-36MM HC RB	PGF309441	3-0		36 mm 1/2 Circle Round Body	90	12	UND
PGLA FAST UND 3-0 X 45-22MM HC CT	PGF30 9935	3-0		22 mm 1/2 Circle Conventional Cutting	45	12	UND
PGLA FAST UND 3-0 X 70-18MM HC RB	PGF300318	3-0		18 mm 1/2 Circle Round Body	70	12	UND
PGLA FAST UND 3-0 X 45-16MM CU CT	PGF30 2719	3-0		16 mm 3/8 Circle Conventional Cutting	45	12	UND
PGLA FAST UND 3-0 X 70-22MM HC CT	PGF30 2735	3-0		22 mm 1/2 Circle Conventional Cutting	70	12	UND
PGLA FAST UND 3-0 X 70-26MM HC RB	PGF30 2764	3-0		26 mm 1/2 Circle Round Body	70	12	UND
PGLA FAST UND 3-0 X 45-60MM ST CT	PGF309943SL	3-0		60 mm Straight Conventional Cutting	45	12	UND
PGLA FAST UND 3-0 X 75-26MM CU RC	PGF309932	3-0		26 mm 3/8 Circle Reverse Cutting	75	12	UND

Description	FG Code	Size		Needle Profile Shape, Type & Length	Suture Length (cm)	Pack Size	Suture Colour
MITSU FST (SYNTHETIC ABSORBABLE BRAIDED COATED POLYGLACTIN 910 FAST)							
PGLA FAST UND 5-0 X 45-11MM CU RC	PGF50 9915	5-0		11 mm 3/8 Circle Reverse Cutting	45	12	UND
PGLA FAST UND 5-0 X 75-16MM CU RC	PGF502296	5-0		16 mm 3/8 Circle Reverse Cutting	75	12	UND
PGLA FAST UND 5-0 X 70-12MM CU RC	PGF502493	5-0		12 mm 3/8 Circle Reverse Cutting	70	12	UND
PGLA FAST UND 5-0 X 45-16MM CU CT	PGF508441	5-0		16 mm 3/8 Circle Conventional Cutting	45	12	UND
PGLA FAST UND 6-0 X 45-11MM CU RC	PGF60 9913	6-0		11 mm 3/8 Circle Reverse Cutting	45	12	UND
PGLA FAST UND 2-0 X 150 (50X3)	PGF200516	2-0		Non needled	50	12	UND

MEGASORB™

Polyglycolic Acid Suture

Properties

- 70%-80% of the original is retained until initial 14 Days
- 50% of the original is retained until initial 21 Days
- There is a subsequent loss between four to five weeks post implantation
- Complete absorption usually takes place between 60 and 90 Days

Indications


- Gastrointestinal procedures
- Gynaecology/Obstetric procedures
- Ophthalmic procedures
- Orthopedic procedures
- Urology procedures
- Skin closure (intracutaneous, subcutaneous)
- Ligatures

Advantages


- High tensile strength
- Excellent knotting ability
- Smooth passage through tissue
- Easy handling



Product Overview

Structure	: Multifilament, Braided
Colour	: Dyed, Undyed
	
Chemical composition	: Polyglycolic acid
Coating	: Polycaprolactone and Calcium Stearate
Sizes	: Violet - USP 8/0 (metric 0.4) - USP 2 (metric 5) : Undyed - USP 6/0 (metric 0.7) - USP 2 (metric 5)
Type of absorption	: Hydrolysis
Sterilisation	: Ethylene oxide

Description	FG Code	Size		Needle Profile Shape, Type & Length	Suture Length (cm)	Pack Size	Suture Colour
MEGASORB (SYNTHETIC ABSORBABLE BRAIDED COATED POLYGLYCOLIC ACID)							
PGA VIO 0 X 90-30MM HC RB	MS2338	0		30 mm 1/2 Circle Round Body	90	12	VIO
PGA VIO 0 X 90-40MM HC RB	MS2346	0		40 mm 1/2 Circle Round Body	90	12	VIO
PGA VIO 0 X 140-40MM HC RB TC DN	MS2346DN	0		40 mm 1/2 Circle Round Body Taper Cut	140	12	VIO
PGA VIO 0 X 180-40MM HC RB TC DN	MS2346DNL	0		40 mm 1/2 Circle Round Body Taper Cut	180	12	VIO
PGA VIO 0 X 45-40MM HC RB	MS2346S	0		40 mm 1/2 Circle Round Body	45	12	VIO
PGA VIO 0 X 90-40MM HC TC	MS2518	0		40 mm 1/2 Circle Taper Cut	90	12	VIO
PGA VIO 0 X 90-36MM HC RC	MS2534	0		36 mm 1/2 Circle Reverse Cutting	90	12	VIO
PGA VIO 0 X 75-40MM HC RB	MS100352	0		40 mm 1/2 Circle Round Body	75	12	VIO
PGA VIO 0 X 70-30MM HC RB EX	MS100808	0		30 mm 1/2 Circle Round Body	70	12	VIO
PGA VIO 0 X 70-22MM HC RB	MS100329	0		22 mm 1/2 Circle Round Body	70	12	VIO
PGA VIO 0 X 70-30MM CU RC	MS100587	0		30 mm 3/8 Circle Reverse Cutting	70	12	VIO
PGA VIO 0 X 100-40MM HC RB	MS2346XL	0		40 mm 1/2 Circle Round Body	100	12	VIO
PGA VIO 0 X 70-36MM HC RB	MS100260	0		36 mm 1/2 Circle Round Body	70	12	VIO
PGA VIO 0 X 75-26MM HC RB	MS100334	0		26 mm 1/2 Circle Round Body	75	12	VIO
PGA VIO 0 X 100-48MM HC RB	MS109250XL	0		48 mm 1/2 Circle Round Body	100	12	VIO

Description	FG Code	Size		Needle Profile Shape, Type & Length	Suture Length (cm)	Pack Size	Suture Colour
MEGASORB (SYNTHETIC ABSORBABLE BRAIDED COATED POLYGLYCOLIC ACID)							
PGA VIO 0 X 100-48MM HC RB	MS109250XL	0		48 mm 1/2 Circle Round Body	100	12	VIO
PGA VIO 0 X 75-50MM HC RB	MS109251	0		50 mm 1/2 Circle Round Body	75	12	VIO
PGA VIO 0 X 45-40MM HC TC	MS2518SL	0		40 mm 1/2 Circle Taper Cut	45	12	VIO
PGA UND 0 X 90-40MM HC RB	MS2346U	0		40 mm 1/2 Circle Round Body	90	12	UND
PGA VIO 0 X 70-30MM HC RB	MS2342	0		30 mm 1/2 Circle Round Body	70	12	VIO
PGA VIO 0 X 75-45MM HC RB	MS109244	0		45 mm 1/2 Circle Round Body	75	12	VIO
PGA VIO 0 X 75-48MM HC RB	MS109250	0		48 mm 1/2 Circle Round Body	75	12	VIO
PGA VIO 1 X 90-40MM HC RB TN	MS2347	1		40 mm 1/2 Circle Round Body	90	12	VIO
PGA VIO 1 X 140-40MM HC RB TC DN	MS2347DN	1		40 mm 1/2 Circle Round Body Taper Cut	140	12	VIO
PGA VIO 1 X 180-40MM HC RB TC DN	MS2347DNL	1		40 mm 1/2 Circle Round Body Taper Cut	180	12	VIO
PGA VIO 1 X 45-40MM HC RB	MS2347S	1		40 mm 1/2 Circle Round Body	45	12	VIO
PGA VIO 1 X 70-40MM HC RB	MS2359	1		40 mm 1/2 Circle Round Body	70	12	VIO
PGA VIO 1 X 90-40MM HC RB TH	MS2360	1		40 mm 1/2 Circle Round Body	90	12	VIO
PGA VIO 1 X 90-40MM HC RC	MS2421	1		40 mm 1/2 Circle Reverse Cutting	90	12	VIO

Description	FG Code	Size		Needle Profile Shape, Type & Length	Suture Length (cm)	Pack Size	Suture Colour
MEGASORB (SYNTHETIC ABSORBABLE BRAIDED COATED POLYGLYCOLIC ACID)							
PGA VIO 1 X 70-37MM HC RB	MS010947	1		37 mm 1/2 Circle Round Body	70	12	VIO
PGA VIO 1 X 90-48MM HC RC	MS011058	1		48 mm 1/2 Circle Reverse Cutting	90	12	VIO
PGA VIO 1 X 90-48MM HC RB	MS011428	1		48 mm 1/2 Circle Round Body	90	12	VIO
PGA VIO 1 X 90-31MM HC RB	MS010253	1		31 mm 1/2 Circle Round Body	90	12	VIO
PGA VIO 1 X 150-45MM HC RC	MS011059	1		45 mm 1/2 Circle Reverse Cutting	150	12	VIO
PGA VIO 1 X 90-30MM CU RC	MS019130	1		30 mm 3/8 Circle Reverse Cutting	90	12	VIO
PGA VIO 1 X 180-50/40MM CU RB HC TC DN	MS2160DN	1		50/40 mm 3/8 Circle 1/2 Circle Round Body Taper Cut	180	12	VIO
PGA VIO 1 X 110-40MM HC RB	MS2347DS	1		40 mm 1/2 Circle Round Body	110	12	VIO
PGA VIO 1 X 75-30MM HC RB	MS019214	1		30 mm 1/2 Circle Round Body	75	12	VIO
PGA VIO 1 X 75-35MM HC RB	MS010261	1		35 mm 1/2 Circle Round Body	75	12	VIO
PGA VIO 1 X 75-45MM HC RB	MS019245	1		45 mm 1/2 Circle Round Body	75	12	VIO
PGA VIO 1 X 100-48MM HC RB	MS011428XL	1		48 mm 1/2 Circle Round Body	100	12	VIO
PGA VIO 1 X 90-50MM HC RB	MS012448	1		50 mm 1/2 Circle Round Body	90	12	VIO
PGA VIO 1 X 90-36MM HC RC	MS010468	1		36 mm 1/2 Circle Reverse Cutting	90	12	VIO
PGA VIO 1 X 75-37MM HC TC	MS012137	1		37 mm 1/2 Circle Taper Cut	75	12	VIO
PGA VIO 1 X 70-26MM HC RB	MS012335	1		26 mm 1/2 Circle Round Body	70	12	VIO

Description	FG Code	Size		Needle Profile Shape, Type & Length	Suture Length (cm)	Pack Size	Suture Colour
MEGASORB (SYNTHETIC ABSORBABLE BRAIDED COATED POLYGLYCOLIC ACID)							
PGA VIO 4-0 X 75-22MM HC RB	MS409422	4-0		22 mm 1/2 Circle Round Body	75	12	VIO
PGA VIO 4-0 X 70-16MM HC RB	MS2305	4-0		16 mm 1/2 Circle Round Body	70	12	VIO
PGA VIO 4-0 X 45-12MM CU RC	MS2494	4-0		12 mm 3/8 Circle Reverse Cutting	45	12	VIO
PGA UND 4-0 X 75-26MM CU CT	MS9918U	4-0		26 mm 3/8 Circle Conventional Cutting	75	12	UND
PGA VIO 5-0 X 45-16MM HC RB	MS2303	5-0		16 mm 1/2 Circle Round Body	45	12	VIO
PGA VIO 5-0 X 70-17MM HC RB	MS500303	5-0		17 mm 1/2 Circle Round Body	70	12	VIO
PGA VIO 5-0 X 70-13MM HC RB	MS502433	5-0		13 mm 1/2 Circle Round Body	70	12	VIO
PGA VIO 5-0 X 70-36MM HC RB	MS502536	5-0		36 mm 1/2 Circle Round Body	70	12	VIO
PGA VIO 5-0 X 45-13MM CU RB	MS502513	5-0		13 mm 3/8 Circle Round Body	45	12	VIO
PGA VIO 5-0 X 75-16MM HC RB	MS2303XL	5-0		16 mm 1/2 Circle Round Body	75	12	VIO
PGA VIO 5-0 X 75-16MM CU RC	MS500500XL	5-0		16 mm 3/8 Circle Reverse Cutting	75	12	VIO
PGA VIO 5-0 X 45-12MM CU RC	MS2493	5-0		12 mm 3/8 Circle Reverse Cutting	45	12	VIO
PGA VIO 5-0 X 45-15MM HC CT	MS500557	5-0		15 mm 1/2 Circle Conventional Cutting	45	12	VIO
PGA VIO 6-0 X 70-13MM HC RB	MS602134	6-0		13 mm 1/2 Circle Round Body	70	12	VIO
PGA VIO 6-0 X 45-13MM CU RB	MS602613	6-0		13 mm 3/8 Circle Round Body	45	12	VIO
PGA VIO 6-0 X 75-16MM CU RC	MS602612	6-0		16 mm 3/8 Circle Reverse Cutting	75	12	VIO
PGA VIO 6-0 X 45-8MM QC SM DN	MS2670DN	6-0		8 mm 1/4 Circle Spatula Micropoint Spatula Micropoint	45	12	VIO
PGA VIO 6-0 X 45-12MM CU RC	MS602645	6-0		12 mm 3/8 Circle Reverse Cutting	45	12	VIO
PGA VIO 7-0 X 30-6MM QC SM DN	MS2349DN	7-0		6 mm 1/4 Circle Spatula Micropoint Spatula Micropoint	30	12	VIO
PGA VIO 8-0 X 30-6MM CU SM DN	MS2348DN	8-0		6 mm 3/8 Circle Spatula Micropoint Spatula Micropoint	30	12	VIO

Description	FG Code	Size	Needle Profile Shape, Type & Length	Suture Length (cm)	Pack Size	Suture Colour
MEGASORB (SYNTHETIC ABSORBABLE BRAIDED COATED POLYGLYCOLIC ACID)						
PGA VIO 0 X 45 (45 X3)	MS100616	0	Non needled	45	12	VIO
PGA VIO 0 X 180	MS2614NS	0	Non needled	180	12	VIO
PGA VIO 0 X 270 (45X6)	MS100062	0	Non needled	45	12	VIO
PGA VIO 0 X 150	MS100065	0	Non needled	150	12	VIO
PGA VIO 1 X 45 (45X3)	MS016173	1	Non needled	45	12	VIO
PGA VIO 1 X 45 (45X6)	MS014506	1	Non needled	45	12	VIO
PGA VIO 2 X 270 (45X6)	MS020082	2	Non needled	45	12	VIO
PGA VIO 2-0 X 45 (45X3)	MS200615	2-0	Non needled	45	12	VIO
PGA VIO 2-0 X 180	MS2615NS	2-0	Non needled	180	12	VIO
PGA VIO 2-0 X 150	MS300055	2-0	Non needled	150	12	VIO
PGA VIO 3-0 X 45 (45X3)	MS300414	3-0	Non needled	45	12	VIO
PGA VIO 3-0 X 50 (50X10)	MS30500	3-0	Non needled	50	12	VIO
PGA VIO 3-0 X 270 (45X6)	MS300042	3-0	Non needled	45	12	VIO
PGA VIO 3-0 X 150	MS300045	3-0	Non needled	150	12	VIO
PGA VIO 4-0 X 135 (45X3)	Ms400284	4-0	Non needled	45	12	VIO

MEGASORB PLUS™

Polyglycolic Acid Suture with Chlorhexidine

Properties

- Retention of 70% knot tensile strength after 14 days
- Complete mass absorption within 60 to 90 days

Indications

- Gastrointestinal procedures
- Gynaecology/Obstetric procedures
- Orthopedic procedures
- Urology procedures
- Skin closure (intracutaneous, subcutaneous)
- Ligatures


Advantages

- Smooth knot run-down
- Gentle knot repositioning
- High knot holding capacity
- Minimized inflammatory reaction



Also available **MEGASORB T+** Polyglycolic Acid with Triclosan

Product Overview

Structure	: Multifilament, Braided
Colour	: Dyed, Undyed
	
Chemical composition	: Polyglycolic acid
Coating	: Polycaprolactone + calcium stearate + Antibacterial agent Chlorhexidine
Sizes	: USP 8/0 (metric 0.4) - USP 2 (metric 5)
Type of absorption	: Hydrolysis
Sterilization	: Ethylene oxide

Description	FG Code	Size		Needle Profile Shape, Type & Length	Suture Length (cm)	Pack Size	Suture Colour
MEGASORB PLUS (POLYGLYCOLIC ACID SUTURE WITH CHLORHEXIDINE)							
PGA C PLUS VIO 0 X 90-40MM HC RB	MSP2346	0		40 mm 1/2 Circle Round Body	90	12	VIO
PGA C PLUS VIO 1 X 90-40MM HC RB TH	MSP2360	1		40 mm 1/2 Circle Round Body	90	12	VIO
PGA C PLUS VIO 1 X 90-40MM HC RB TN	MSP2347	1		40 mm 1/2 Circle Round Body	90	12	VIO
PGA C PLUS VIO 1 X 90-40MM HC RC	MSP2421	1		40 mm 1/2 Circle Reverse Cutting	90	12	VIO
PGA C PLUS VIO 2-0 X 90-30MM HC RB	MSP2317	2-0		30 mm 1/2 Circle Round Body	90	12	VIO
PGA C PLUS VIO 3-0 X 70-20MM HC RB	MSP2437	3-0		20 mm 1/2 Circle Round Body	70	12	VIO

MEGASORB RAPID™

Polyglycolic Acid Fast Suture

Properties

- 50% of the original is retained until initial 6 to 8 Days
- There is a subsequent loss by 17 days post implantation
- Complete absorption usually takes place by 42 Days

Indications

- Gynaecology/Obstetric procedures (e.g.: episiotomy repair)
- Urology procedures
- Ophthalmic procedures (e.g.: conjunctiva suturing)
- Odontology procedures (e.g.: oral mucosa)
- Paediatric procedures
- Skin closure (intracutaneous, subcutaneous)
- Ligatures


Advantages

- High initial tensile strength
- Predictable and constant degradation rate
- Good knot security
- Excellent handling properties
- Quick mass absorption



Product Overview

Structure	: Multifilament, Braided
Colour	: Undyed
Chemical composition	: Polyglycolic acid
Sizes	: USP 6/0 (metric 0.7) - USP 2 (metric 5)
Type of absorption	: Hydrolysis
Sterilization	: Gamma radiation

Description	FG Code	Size		Needle Profile Shape, Type & Length	Suture Length (cm)	Pack Size	Suture Colour
MEGASORB RAPID (POLYGLYCOLIC ACID FAST SUTURE)							
PGA RAPID UND 0 X 90-36MM HC RB	MS109037	0		36 mm 1/2 Circle Round Body	90	12	UND
PGA RAPID UND 1 X 90-40MM HC RB	MS019140	1		40 mm 1/2 Circle Round Body	90	12	UND
PGA RPID UND 2-0 X 140-36MM HC RC RB DN	MS2777	2-0		36 mm 1/2 Circle Reverse Cutting Round Body	140	12	UND
PGA RPID UND 2-0 X 90-36MM HC TC	MS2762	2-0		36 mm 1/2 Circle Taper Cut	90	12	UND
PGA RAPID UND 2-0 X 90-26MM CU RC	MS209933	2-0		26 mm 3/8 Circle Reverse Cutting	90	12	UND
PGA RAPID UND 2-0 X 90-36MM HC RB	MS209451	2-0		36 mm 1/2 Circle Round Body	90	12	UND
PGA RAPID UND 2-0 X 75-30MM HC RB	MS209975	2-0		30 mm 1/2 Circle Round Body	75	12	UND
PGA RAPID UND 2-0 X 90-36MM HC RC	MS209236	2-0		36 mm 1/2 Circle Reverse Cutting	36	12	UND
PGA RAPID UND 3-0 X 70-26MM CU RC	MS302732	3-0		26 mm 3/8 Circle Reverse Cutting	70	12	UND
PGA RAPID UND 3-0 X 70-22MM HC RB	MS309974	3-0		22 mm 1/2 Circle Round Body	70	12	UND
PGA RAPID UND 3-0 X 75-19MM CU RC	MS309923	3-0		19 mm 3/8 Circle Reverse Cutting	75	12	UND
PGA RAP UND 3-0 X 75-26MM CU RC	MS9932	3-0		26 mm 3/8 Circle Reverse Cutting	75	12	UND
PGA RAPID UND 3-0 X 45-22MM HC RC	MS306230	3-0		22 mm 1/2 Circle Reverse Cutting	45	12	UND
PGA RAPID UND 3-0 X 70-26MM HC RB	MS304161	3-0		26 mm 1/2 Circle Round Body	70	12	UND
PGA RAPID UND 3-0 X 45-22MM HC CT	MS309927	3-0		22 mm 1/2 Circle Conventional Cutting	45	12	UND
PGA RAPID UND 4-0 X 45-16MM CU CT	MS408451	4-0		16 mm 3/8 Circle Conventional Cutting	45	12	UND
PGA RAPID UND 4-0 X 75-17MM HC RB	MS409970	4-0		17 mm 1/2 Circle Round Body	75	12	UND

Description	FG Code	Size		Needle Profile Shape, Type & Length	Suture Length (cm)	Pack Size	Suture Colour
MEGASORB RAPID (POLYGLYCOLIC ACID FAST SUTURE)							
PGA RAPID UND 4-0 X 75-19MM CU RC	MS409922	4-0		19 mm 3/8 Circle Reverse Cutting	75	12	UND
PGA RAPID UND 4-0 X 75-24MM CU RC	MS409924	4-0		24 mm 3/8 Circle Reverse Cutting	75	12	UND
PGA RAP UND 4-0 X 75-26MM CU CT	MS9918	4-0		26 mm 3/8 Circle Conventional Cutting	75	12	UND
PGA RAPID UND 4-0 X 70-19MM CU CT	MS406529	4-0		19 mm 3/8 Circle Conventional Cutting	70	12	UND
PGA RAPID UND 4-0 X 75-17MM HC RC	MS409418	4-0		17 mm 1/2 Circle Reverse Cutting	17	12	UND
PGA RAPID UND 5-0 X 45-12MM CU RC	MS9915	5-0		12 mm 3/8 Circle Reverse Cutting	45	12	UND
PGA RAPID UND 2-0 X 45 (45X3)	MS200616	2-0		Non needed	45	12	UND
PGA RAPID UND 3-0 X 45 (45X3)	MS300285	3-0		Non needed	45	12	UND

FILAPRON™

Poliglecaprone 25 Suture

Properties

- 80%-60% of the original is retained until initial 7 Days (DYED)
- 60%-50% of the original is retained until initial 7 Days (UNDYED)
- 40%-30% of the original is retained until initial 14 Days (DYED)
- 30%-20% of the original is retained until initial 7 Days (UNDYED)
- There is a subsequent loss 28 days post implantation
- Complete absorption usually takes place between 90 to 120 Days

Indications

When short wound support is desirable

- Skinclosure (intracutaneous, subcutaneous)
- Ligatures


Advantages

- Short-term absorption profile
- High knot pull tensile strength
- Superior handling features



FILAPRON

Product Overview

Structure	: Monofilament
Colour	: Dyed, Undyed
	
Chemical composition	: Poly(glycolide-co-caprolactone)PGCL
Sizes	: USP 6/0 (metric 0.7) - USP 1 (metric 4)
Type of absorption	: Hydrolysis
Sterilization	: Ethylene oxide

Description	FG Code	Size		Needle Profile Shape, Type & Length	Suture Length (cm)	Pack Size	Suture Colour
FILAPRON (POLIGLECAPRONE 25 SUTURE)							
PGCL VIO 0 X 70-31MM HC RB H	PCL103442	0		31 mm 1/2 Circle Round Body	70	12	VIO
PGCL VIO 0 X 70-36MM HC RB	PCL103443	0		36 mm 1/2 Circle Round Body	70	12	VIO
PGCL VIO 0 X 70-36MM HC TC	PCL103489	0		36 mm 1/2 Circle Taper Cut	70	12	VIO
PGCL VIO 0 X 70-36MM HC RC	PCL103601	0		36 mm 1/2 Circle Reverse Cutting	70	12	VIO
PGCL VIO 0 X 70-40MM HC RB	PCL103758	0		40 mm 1/2 Circle Round Body	70	12	VIO
PGCL VIO 0 X 70-19MM CU RC	PCL101019	0		19 mm 3/8 Circle Reverse Cutting	70	12	VIO
PGCL VIO 0 X 70-30MM HC RB	PCL10 1742V	0		30 mm 1/2 Circle Round Body	70	12	VIO
PGCL UND 0 X 70-30MM HC RB	PCL10 1642U	0		30 mm 1/2 Circle Round Body	70	12	UND
PGCL UND 0 X 45-23MM HC RB	PCL10 1245U	0		23 mm 1/2 Circle Round Body	45	12	UND
PGCL VIO 0 X 70-31MM HC RB	PCL10 2442V	0		31 mm 1/2 Circle Round Body	70	12	VIO
PGCL VIO 0 X 90-26MM HC RB	PCL104271	0		26 mm 1/2 Circle Round Body	90	12	VIO
PGCL VIO 0 X 70-26MM HC RB	PCL104270	0		26 mm 1/2 Circle Round Body	70	12	VIO
PGCL VIO 1 X 70-36MM HC RB	PCL013457	1		36 mm 1/2 Circle Round Body	70	12	VIO
PGCL VIO 1 X 90-45MM HC RB	PCL013727	1		45 mm 1/2 Circle Round Body	90	12	VIO
PGCL VIO 1 X 70-31MM HC RB	PCL013755	1		31 mm 1/2 Circle Round Body	70	12	VIO
PGCL VIO 1 X 70-40MM HC RB	PCL013759	1		40 mm 1/2 Circle Round Body	70	12	VIO
PGCL VIO 1 X 90-40MM HC RC	PCL013770	1		40 mm 1/2 Circle Reverse Cutting	90	12	VIO
PGCL VIO 1 X 70-48MM HC TC	PCL011148	1		48 mm 1/2 Circle Taper Cut	70	12	VIO

Description	FG Code	Size		Needle Profile Shape, Type & Length	Suture Length (cm)	Pack Size	Suture Colour
FILAPRON (POLIGLECAPRONE 25 SUTURE)							
PGCL VIO 1 X 90-48MM HC RB	PCL013148	1		48 mm 1/2 Circle Round Body	90	12	VIO
PGCL VIO 2 X 90-45MM HC RB	PCL023945	2		45 mm 1/2 Circle Round Body	90	12	VIO
PGCL VIO 2-0 X 90-23MM HC RC	PCL200459	2-0		23 mm 1/2 Circle Reverse Cutting	90	12	VIO
PGCL VIO 2-0 X 70-31MM CU RB	PCL203416	2-0		31 mm 3/8 Circle Round Body	70	12	VIO
PGCL UND 2-0 X 70-31MM HC RB	PCL203431	2-0		31 mm 1/2 Circle Round Body	70	12	UND
PGCL UND 2-0 X 70-26MM HC TC	PCL203432	2-0		26 mm 1/2 Circle Taper Cut	70	12	UND
PGCL VIO 2-0 X 70-26MM HC TC	PCL203440	2-0		26 mm 1/2 Circle Taper Cut	70	12	VIO
PGCL VIO 2-0 X 70-26MM HC RB	PCL203448	2-0		26 mm 1/2 Circle Round Body	70	12	VIO
PGCL VIO 2-0 X 70-36MM HC RB	PCL203463	2-0		36 mm 1/2 Circle Round Body	70	12	VIO
PGCL VIO 2-0 X 70-36MM HC RC	PCL203600	2-0		36 mm 1/2 Circle Reverse Cutting	70	12	VIO
PGCL UND 2-0 X 70-60MM ST RB	PCL203651	2-0		60 mm Straight Round Body	70	12	UND
PGCL VIO 2-0 X 70-22MM HC RB	PCL203662	2-0		22 mm 1/2 Circle Round Body	70	12	VIO
PGCL UND 2-0 X 70-24MM CU RC	PCL203327	2-0		24 mm 3/8 Circle Reverse Cutting	70	12	UND
PGCL VIO 2-0 X 70-26MM FE RB	PCL20610	2-0		26 mm 5/8 Circle Round Body	70	12	VIO
PGCL UND 2-0 X 70-26MM FE RB	PCL20 1602U	2-0		26 mm 5/8 Circle Round Body	70	12	UND
PGCL UND 2-0 X 70-26MM HC RB VB	PCL20 1665U	2-0		26 mm 1/2 Circle Round Body	70	12	UND
PGCL UND 2-0 X 70-30MM HC RB	PCL20 1641U	2-0		30 mm 1/2 Circle Round Body	70	12	UND
PGCL UND 2-0 X 70-36MM HC RB	PCL20 1666U	2-0		36 mm 1/2 Circle Round Body	70	12	UND
PGCL VIO 2-0 X 70-26MM HC RB VB	PCL20 1765V	2-0		26 mm 1/2 Circle Round Body	70	12	VIO

Description	FG Code	Size		Needle Profile Shape, Type & Length	Suture Length (cm)	Pack Size	Suture Colour
FILAPRON (POLIGLECAPRONE 25 SUTURE)							
PGCL VIO 3-0 X 70-22MM HC RB	PCL303661	3-0		22 mm 1/2 Circle Round Body	70	12	VIO
PGCL UND 3-0 X 45-19MM CU RC DN	PCL300497	3-0		19 mm 3/8 Circle Reverse Cutting Reverse Cutting	45	12	UND
PGCL VIO 3-0 X 70-26MM HC RB	PCL301326	3-0		26 mm 1/2 Circle Round Body	70	12	VIO
PGCL UND 3-0 X 45-25MM CU CT	PCL30 1326SLU	3-0		25 mm 3/8 Circle Conventional Cutting	45	12	UND
PGCL UND 3-0 X 70-22MM HC RB VB	PCL30 1663U	3-0		22 mm 1/2 Circle Round Body	70	12	UND
PGCL VIO 3-0 X 70-26MM HC RB VB	PCL30 1764V	3-0		26 mm 1/2 Circle Round Body	70	12	VIO
PGCL UND 3-0 X 70-26MM CU RC	PCL30 3326U	3-0		26 mm 3/8 Circle Reverse Cutting	70	12	UND
PGCL VIO 3-0 X 70-25MM CU CT	PCL30 1326V	3-0		25 mm 3/8 Circle Conventional Cutting	70	12	VIO
PGCL UND 3-0 X 70-60MM ST CT	PCL30 3650U	3-0		60 mm Straight Conventional Cutting	70	12	UND
PGCL UND 3-0 X 70-19MM CU RC	PCL303202	3-0		19 mm 3/8 Circle Reverse Cutting	70	12	UND
PGCL UND 3-0 X 60-26MM CU RC	PCL303123	3-0		26 mm 3/8 Circle Reverse Cutting	60	12	UND
PGCL VIO 3-0 X 70-26MM HC RB DN	PCL303416	3-0		26 mm 1/2 Circle Round Body Round Body	70	12	VIO
PGCL UND 3-0 X 70-26MM CU RC EX	PCL303213	3-0		26 mm 3/8 Circle Reverse Cutting	70	12	UND
PGCL UND 3-0 X 70-26MM HC RB	PCL303446	3-0		26 mm 1/2 Circle Round Body	70	12	UND
PGCL UND 3-0 X 70-26MM CU RC EX-3D	PCL303213Z	3-0		26 mm 3/8 Circle Reverse Cutting	70	36	UND
PGCL VIO 3-0 X 90-26MM HC RB	PCL301326XL	3-0		26 mm 1/2 Circle Round Body	90	12	VIO
PGCL VIO 3-0 X 60-22MM CU RB	PCL306022	3-0		22 mm 3/8 Circle Round Body	60	12	VIO
PGCL VIO 3-0 X 70-17MM HC RB DN	PCL303437DN	3-0		17 mm 1/2 Circle Round Body Round Body	70	12	VIO
PGCL VIO 3-0 X 70-22MM HC RB DN	PCL303661DN	3-0		22 mm 1/2 Circle Round Body Round Body	70	12	VIO
PGCL UND 4-0 X 45-19MM CU RC	PCL403206	4-0		19 mm 3/8 Circle Reverse Cutting	45	12	UND

Description	FG Code	Size		Needle Profile Shape, Type & Length	Suture Length (cm)	Pack Size	Suture Colour
FILAPRON (POLIGLECAPRONE 25 SUTURE)							
PGCL VIO 5-0 X 90-16MM CU RC	PCL503204	5-0		16 mm 3/8 Circle Reverse Cutting	90	12	VIO
PGCL VIO 5-0 X 70-13MM HC RB	PCL501513	5-0		13 mm 1/2 Circle Round Body	70	12	VIO
PGCL VIO 5-0 X 70-17MM HC RB	PCL501517	5-0		17 mm 1/2 Circle Round Body	70	12	VIO
PGCL VIO 5-0 X 90-12MM CU RC	PCL501512	5-0		12 mm 3/8 Circle Reverse Cutting	90	12	VIO
PGCL VIO 5-0 X 70-19MM HC RB	PCL501518	5-0		19 mm 1/2 Circle Round Body	70	12	VIO
PGCL UND 5-0 X 45-13MM CU RC	PCL50 3203U	5-0		13 mm 3/8 Circle Reverse Cutting	45	12	UND
PGCL UND 5-0 X 45-16MM CU RC	PCL50 3204U	5-0		16 mm 3/8 Circle Reverse Cutting	45	12	UND
PGCL UND 5-0 X 45-16MM CU CT	PCL50 844U	5-0		16 mm 3/8 Circle Conventional Cutting	45	12	UND
PGCL VIO 5-0 X 45-13MM CU CT	PCL503203	5-0		13 mm 3/8 Circle Conventional Cutting	45	12	VIO
PGCL UND 5-0 X 70-19MM CU RC	PCL503209U	5-0		19 mm 3/8 Circle Reverse Cutting	70	12	UND
PGCL VIO 6-0 X 45-13MM HC RB	PCL603224	6-0		13 mm 1/2 Circle Round Body	45	12	VIO
PGCL VIO 6-0 X 70-13MM HC RB	PCL601613	6-0		13 mm 1/2 Circle Round Body	70	12	VIO
PGCL UND 6-0 X 45-13MM SKKI CT	PCL60 833U	6-0		13 mm SKKI Conventional Cutting	45	12	UND
PGCL UND 6-0 X 45-13MM CU RC	PCL603215	6-0		13 mm 3/8 Circle Reverse Cutting	45	12	UND
PGCL UND 6-0 X 45-11MM CU RC	PCL603214	6-0		11 mm 3/8 Circle Reverse Cutting	45	12	UND

FILAXYN™

Polydioxanone Suture

Properties

- 70% of the original is retained until initial 14 Days
- 50%-60% of the original is retained until initial 28 Days
- There is a subsequent minimal absorption until about 90 th post implantation day
- Complete absorption usually takes place between 180 - 210 days

Indications

- When prolonged wound support is desirable
- Abdominal wall closure
- Orthopedic procedures
- Pediatric sternum closure
- Oncocological procedures

Advantages



- High initial tensile strength
- Predictable absorption
- Prolonged wound support
- Smooth passage through tissues
- Minimal tissue reaction
- Good handling characteristics
- Standard monofilament knot for security






Product Overview

Structure	: Monofilament
Colour	: Violet
Chemical Composition	: Poly (p-dioxanone)
Sizes	: USP 6/0 (metric 1.5) - USP 2 (metric 5)
Type of absorption	: Hydrolysis
Sterilization	: Ethylene oxide

Description	FG Code	Size		Needle Profile Shape, Type & Length	Suture Length (cm)	Pack Size	Suture Colour
FILAXYN (POLYDIOXANONE SUTURE)							
PDO VIO 0 X 90-30MM HC RB	PDX109333	0		30 mm 1/2 Circle Round Body	90	12	VIO
PDO VIO 0 X 90-36MM HC RC	PDX100467	0		36 mm 1/2 Circle Reverse Cutting	90	12	VIO
PDO VIO 0 X 90-40MM HC RB	PDX10 9371	0		40 mm 1/2 Circle Round Body	90	12	VIO
PDO VIO 0 X 70-30MM HC RB	PDX10 9254	0		30 mm 1/2 Circle Round Body	70	12	VIO
PDO VIO 0 X 150-40MM HC RB LOOP	PDX10 9261	0		40 mm 1/2 Circle Round Body	150	12	VIO
PDO VIO 0 X 90-50MM HC RB	PDX10 9210	0		50 mm 1/2 Circle Round Body	90	12	VIO
PDO VIO 0 X 90-48MM HC RB	PDX109255	0		48 mm 1/2 Circle Round Body	90	12	VIO
PDO VIO 0 X 150-48MM HC RB LOOP	PDX109262	0		48 mm 1/2 Circle Round Body	150	12	VIO
PDO VIO 0 X 90-36MM HC TC	PDX109355	0		36 mm 1/2 Circle Taper Cut	90	12	VIO
PDO VIO 0 X 70-26MM HC RB	PDX100334	0		26 mm 1/2 Circle Round Body	70	12	VIO
PDO VIO 0 X 70-36MM HC RB BLUNT	PDX100341	0		36 mm 1/2 Circle Blunt Point	70	12	VIO
PDO VIO 0 X 90-40MM HC TC	PDX10 9371TC	0		40 mm 1/2 Circle Taper Cut	90	12	VIO
PDO VIO 0 X 70-40MM HC RB	PDX10 9233	0		40 mm 1/2 Circle Round Body	70	12	VIO
PDO VIO 0 X 70-30MM J TYPE TC TH	PDX109333J	0		30 mm J Shape Taper Cut	70	12	VIO
PDO VIO 0 X 90-40MM HC TC EX	PDX109366	0		40 mm 1/2 Circle Taper Cut	90	12	VIO
PDO VIO 0 X 150-48MM HC BP LOOP	PDX109966	0		48 mm 1/2 Circle Blunt Point	150	12	VIO

Description	FG Code	Size		Needle Profile Shape, Type & Length	Suture Length (cm)	Pack Size	Suture Colour
FILAXYN (POLYDIOXANONE SUTURE)							
PDO VIO 0 X 90-45MM HC TC	PDX109356	0		45 mm 1/2 Circle Taper Cut	90	12	VIO
PDO VIO 0 X 90-48MM HC RC	PDX100469	0		48 mm 1/2 Circle Reverse Cutting	90	12	VIO
PDO VIO 1 X 90-48MM HC TC	PDX011936	1		48 mm 1/2 Circle Taper Cut	90	12	VIO
PDO VIO 1 X 90-40MM HC TC	PDX019370	1		40 mm 1/2 Circle Taper Cut	90	12	VIO
PDO VIO 1 X 90-23MM HC RC	PDX019286	1		23 mm 1/2 Circle Reverse Cutting	90	12	VIO
PDO VIO 1 X 150-50MM HC RB LOOP	PDX01 9262	1		50 mm 1/2 Circle Round Body	150	12	VIO
PDO VIO 1 X 90-40MM HC RC	PDX01 9221	1		40 mm 1/2 Circle Reverse Cutting	90	12	VIO
PDO VIO 1 X 90-45MM HC RB	PDX01 9248	1		45 mm 1/2 Circle Round Body	90	12	VIO
PDO VIO 1 X 90-50MM HC RC	PDX01 9374	1		50 mm 1/2 Circle Reverse Cutting	90	12	VIO
PDO VIO 1 X 90-48MM HC RB	PDX01 9255	1		48 mm 1/2 Circle Round Body	90	12	VIO
PDO VIO 1 X 90-40MM HC RB	PDX01 9352	1		40 mm 1/2 Circle Round Body	90	12	VIO
PDO VIO 1 X 150-48MM HC RB LOOP	PDX019263	1		48 mm 1/2 Circle Round Body	150	12	VIO
PDO VIO 1 X 70-35MM HC RB	PDX010341	1		35 mm 1/2 Circle Round Body	70	12	VIO
PDO VIO 1 X 90-35MM HC TC	PDX019382	1		35 mm 1/2 Circle Taper Cut	90	12	VIO
PDO VIO 1 X 70-30MM HC RB	PDX019211	1		30 mm 1/2 Circle Round Body	70	12	VIO
PDO VIO 1 X 70-31MM J TYPE TC	PDX019334J	1		31 mm J Shape Taper Cut	70	12	VIO

Description	FG Code	Size		Needle Profile Shape, Type & Length	Suture Length (cm)	Pack Size	Suture Colour
FILAXYN (POLYDIOXANONE SUTURE)							
PDO VIO 2 X 90-45MM HC TC	PDX029251	2		45 mm 1/2 Circle Taper Cut	90	12	VIO
PDO VIO 2 X 90-36MM HC TC	PDX029252	2		36 mm 1/2 Circle Taper Cut	90	12	VIO
PDO VIO 2-0 X 70-36MM CU RC	PDX209632	2		36 mm 3/8 Circle Reverse Cutting	70	12	VIO
PDO VIO 2-0 X 90-40MM HC RB	PDX209151	2-0		40 mm 1/2 Circle Round Body	90	12	VIO
PDO VIO 2-0 X 90-25MM CU RC	PDX209794	2-0		25 mm 3/8 Circle Reverse Cutting	90	12	VIO
PDO VIO 2-0 X 70-26MM HC TC	PDX209184	2-0		26 mm 1/2 Circle Taper Cut	70	12	VIO
PDO VIO 2-0 X 70-26MM HC RB	PDX200333	2-0		26 mm 1/2 Circle Round Body	70	12	VIO
PDO VIO 2-0 X 90-30MM HC RB	PDX200334	2-0		30 mm 1/2 Circle Round Body	90	12	VIO
PDO VIO 2-0 X 70-35MM CU CT	PDX209631	2-0		35 mm 3/8 Circle Conventional Cutting	70	12	VIO
PDO VIO 2-0 X 70-27MM HC RC	PDX200968	2-0		27 mm 1/2 Circle Reverse Cutting	70	12	VIO
PDO VIO 2-0 X 90-26MM HC RC	PDX200664	2-0		26 mm 1/2 Circle Reverse Cutting	90	12	VIO
PDO VIO 2-0 X 70-21MM HC RB	PDX209220	2-0		21 mm 1/2 Circle Round Body	70	12	VIO
PDO VIO 2-0 X 70-31MM HC RB	PDX20 9133	2-0		31 mm 1/2 Circle Round Body	70	12	VIO
PDO VIO 2-0 X 90-36MM HC TC	PDX209380	2-0		36 mm 1/2 Circle Taper Cut	90	12	VIO
PDO VIO 2-0 X 70-24MM HC RB	PDX20 2333	2-0		24 mm 1/2 Circle Round Body	70	12	VIO
PDO VIO 2-0 X 70-26MM HC RB DN	PDX20 9333	2-0		26 mm 1/2 Circle Round Body Round Body	70	12	VIO
PDO VIO 2-0 X 70-36MM J TYPE RB	PDX207036J	2-0		36 mm J Shape Round Body	70	12	VIO
PDO VIO 2-0 X 90-45MM HC TC	PDX209381	2-0		45 mm 1/2 Circle Taper Cut	90	12	VIO

Description	FG Code	Size		Needle Profile Shape, Type & Length	Suture Length (cm)	Pack Size	Suture Colour
FILAXYN (POLYDIOXANONE SUTURE)							
PDO VIO 2-0 X 90-40MM HC TC	PDX209382	2-0		40 mm 1/2 Circle Taper Cut	90	12	VIO
PDO VIO 3-0 X 70-20MM HC RB	PDX30 9237	3-0		20 mm 1/2 Circle Round Body	70	12	VIO
PDO VIO 3-0 X 70-30MM HC RB DN	PDX30 9332DN	3-0		30 mm 1/2 Circle Round Body Round Body	70	12	VIO
PDO VIO 3-0 X 75-20MM HC RB DN	PDX30 9336DN	3-0		20 mm 1/2 Circle Round Body Round Body	75	12	VIO
PDO VIO 3-0 X 70-26MM HC RB	PDX300316	3-0		26 mm 1/2 Circle Round Body	70	12	VIO
PDO VIO 3-0 X 70-26MM HC RB H	PDX300332	3-0		26 mm 1/2 Circle Round Body	70	12	VIO
PDO VIO 3-0 X 75-27MM HC RB	PDX300333	3-0		27 mm 1/2 Circle Round Body	75	12	VIO
PDO VIO 3-0 X 70-27MM HC RC	PDX300968	3-0		27 mm 1/2 Circle Reverse Cutting	70	12	VIO
PDO VIO 3-0 X 45-22MM HC CT	PDX300460	3-0		22 mm 1/2 Circle Conventional Cutting	45	12	VIO
PDO VIO 3-0 X 70-25MM CU RC	PDX300442	3-0		25 mm 3/8 Circle Reverse Cutting	70	12	VIO
PDO VIO 3-0 X 90-31MM HC RB DN	PDX309134	3-0		31 mm 1/2 Circle Round Body Round Body	90	12	VIO
PDO VIO 3-0 X 70-36MM HC RC	PDX309152	3-0		36 mm 1/2 Circle Reverse Cutting	70	12	VIO
PDO VIO 3-0 X 70-31MM HC RB	PDX30 9132	3-0		31 mm 1/2 Circle Round Body	70	12	VIO
PDO VIO 3-0 X 70-20MM HC RB C1	PDX30 9116	3-0		20 mm 1/2 Circle Round Body	70	12	VIO
PDO VIO 3-0 X 45-19MM CU CT	PDX300497	3-0		19 mm 3/8 Circle Conventional Cutting	45	12	VIO
PDO VIO 3-0 X 75-30MM HC RB	PDX3013432	3-0		30 mm 1/2 Circle Round Body	75	12	VIO
PDO VIO 3-0 X 90-25MM HC TC	PDX309179	3-0		25 mm 1/2 Circle Taper Cut	90	12	VIO
PDO VIO 3-0 X 70-26MM CU CT	PDX309625	3-0		26 mm 3/8 Circle Conventional Cutting	70	12	VIO

Description	FG Code	Size		Needle Profile Shape, Type & Length	Suture Length (cm)	Pack Size	Suture Colour
FILAXYN (POLYDIOXANONE SUTURE)							
PDO VIO 3-0 X 70-22MM HC RB	PDX307322	3-0		22 mm 1/2 Circle Round Body	70	12	VIO
PDO VIO 4-0 X 90-16MM CU CT	PDX409867	4-0		16 mm 3/8 Circle Conventional Cutting	90	12	VIO
PDO VIO 4-0 X 90-25MM CU CT	PDX409870	4-0		25 mm 3/8 Circle Conventional Cutting	90	12	VIO
PDO VIO 4-0 X 70-20MM HC RB	PDX40 9304	4-0		20 mm 1/2 Circle Round Body	70	12	VIO
PDO VIO 4-0 X 90-17MM HC RB	PDX409109	4-0		17 mm 1/2 Circle Round Body	90	12	VIO
PDO VIO 4-0 X 90-31MM HC RB	PDX409131	4-0		31 mm 1/2 Circle Round Body	90	12	VIO
PDO VIO 4-0 X 70-16MM CU RB	PDX409077	4-0		16 mm 3/8 Circle Round Body	70	12	VIO
PDO VIO 4-0 X 70-16MM HC RC DN	PDX400507	4-0		16 mm 1/2 Circle Reverse Cutting Reverse Cutting	70	12	VIO
PDO VIO 4-0 X 90-19MM CU RC	PDX400422	4-0		19 mm 3/8 Circle Reverse Cutting	90	12	VIO
PDO VIO 4-0 X 70-17MM CU RB	PDX407417	4-0		17 mm 3/8 Circle Round Body	70	12	VIO
PDO VIO 4-0 X 70-26MM HC RB	PDX409426	4-0		26 mm 1/2 Circle Round Body	70	12	VIO
PDO VIO 4-0 X 45-13MM HC RB	PDX40 9102	4-0		13 mm 1/2 Circle Round Body	45	12	VIO
PDO VIO 4-0 X 90-17MM HC RB DN	PDX40 9109DN	4-0		17 mm 1/2 Circle Round Body Round Body	90	12	VIO
PDO VIO 4-0 X 70-20MM HC RB C1	PDX40 9115	4-0		20 mm 1/2 Circle Round Body	70	12	VIO
PDO VIO 4-0 X 45-19MM CU CT	PDX409950	4-0		19 mm 3/8 Circle Conventional Cutting	45	12	VIO
PDO VIO 4-0 X 75-26MM HC RC	PDX407526	4-0		26 mm 1/2 Circle Reverse Cutting	75	12	VIO
PDO VIO 4-0 X 70-60MM ST CT	PDX409716	4-0		60 mm Straight Conventional Cutting	70	12	VIO
PDO VIO 5-0 X 70-13MM HC RB DN	PDX50 9201DN	5-0		13 mm 1/2 Circle Round Body Round Body	70	12	VIO
PDO VIO 5-0 X 70-16MM HC RC DN	PDX500844	5-0		16 mm 1/2 Circle Reverse Cutting Reverse Cutting	70	12	VIO
PDO VIO 5-0 X 70-13MM HC RB	PDX509101	5-0		13 mm 1/2 Circle Round Body	70	12	VIO
PDO VIO 5-0 X 90-17MM HC RB	PDX500303	5-0		17 mm 1/2 Circle Round Body	90	12	VIO

Description	FG Code	Size		Needle Profile Shape, Type & Length	Suture Length (cm)	Pack Size	Suture Colour
FILAXYN (POLYDIOXANONE SUTURE)							
PDO VIO 5-0 X 70-13MM CU RB	PDX50 1013	5-0		13 mm 3/8 Circle Round Body	70	12	VIO
PDO VIO 5-0 X 45-13MM CU CT	PDX509863	5-0		13 mm 3/8 Circle Conventional Cutting	45	12	VIO
PDO VIO 5-0 X 70-17MM HC RB DN	PDX509108DN	5-0		17 mm 1/2 Circle Round Body Round Body	70	12	VIO
PDO VIO 5-0 X 90-17MM HC RB DN	PDX500303DN	5-0		17 mm 1/2 Circle Round Body Round Body	90	12	VIO
PDO VIO 5-0 X 70-16MM CU CT	PDX509733	5-0		16 mm 3/8 Circle Conventional Cutting	70	12	VIO
PDO VIO 5-0 X 45-19MM CU CT	PDX509869	5-0		19 mm 3/8 Circle Conventional Cutting	45	12	VIO
PDO VIO 5-0 X 70-16MM C RB	PDX5013489	5-0		16 mm 1/2 Circle Round Body	70	12	VIO
PDO VIO 6-0 X 90-13MM CU RB	PDX609096	6-0		13 mm 3/8 Circle Round Body	90	12	VIO
PDO VIO 6-0 X 90-8MM CU RB	PDX609091	6-0		8 mm 3/8 Circle Round Body	90	12	VIO
PDO VIO 6-0 X 45-13MM HC RB	PDX600432	6-0		13 mm 1/2 Circle Round Body	45	12	VIO
PDO VIO 6-0 X 70-16MM HC RC DN	PDX604323	6-0		16 mm 1/2 Circle Reverse Cutting Reverse Cutting	70	12	VIO
PDO VIO 6-0 X 90-9MM CU TC	PDX601712	6-0		9 mm 3/8 Circle Taper Cut	90	12	VIO
PDO VIO 6-0 X 45-8MM CU RB DN	PDX60 9091DN	6-0		8 mm 3/8 Circle Round Body Round Body	45	12	VIO
PDO VIO 6-0 X 45-11MM CU RB DN	PDX60 9093DN	6-0		11 mm 3/8 Circle Round Body Round Body	45	12	VIO
PDO VIO 6-0 X 75-13MM CU RB DN	PDX60 127DN	6-0		13 mm 3/8 Circle Round Body Round Body	75	12	VIO
PDO VIO 6-0 X 45-13MM CU CT	PDX609861	6-0		13 mm 3/8 Circle Conventional Cutting	45	12	VIO
PDO VIO 6-0 X 45-13MM HC SM	PDX607613	6-0		13 mm 1/2 Circle Spatula Micropoint	45	12	VIO
PDO VIO 6-0 X 45-11MM CU RC	PDX600487	6-0		11 mm 3/8 Circle Reverse Cutting	45	12	VIO
PDO VIO 7-0 X 60-9MM CU RB DN	PDX70 155DN	7-0		9 mm 3/8 Circle Round Body Round Body	60	12	VIO
PDO VIO 7-0 X 45-13MM CU RB DN	PDX709095DN	7-0		13 mm 3/8 Circle Round Body Round Body	45	12	VIO

Description	FG Code	Size		Needle Profile Shape, Type & Length	Suture Length (cm)	Pack Size	Suture Colour
FILAXYN (POLYDIOXANONE SUTURE)							
PDO VIO 7-0 X 70-9.3MM CU TC DN	PDX70 1711DN	7-0		9.3 mm 3/8 Circle Taper Cut Taper Cut	70	12	VIO
PDO VIO 3-0 X 70-26MM HC RB DN-3D	PDX301316Z	3-0		26 mm 1/2 Circle Round Body Round Body	70	36	VIO
PDO VIO 4-0 X 45-19MM CU RC	PDX400422SL	4-0		19 mm 3/8 Circle Reverse Cutting	45	12	VIO

MERISOFT PLAIN™

Plain Catgut Suture

Properties

The tensile strength is maintained for not more than 7 to 10 days post implantation, and its absorption is complete within 60 to 90 days.

Indications

MERISOFT PLAIN™ suture is indicated for use in general soft tissue approximation and/or ligation.

Advantages

- Highest level of collagen guarantees predictable absorption
- Accurate polishing ensures smooth tissue passage
- Packed in liquid to reduce memory & increase pliability
- Rigidly controlled spinning & drying to ensure higher tensile strength



Product Overview

Primary packaging	: Packed in an aluminum foil pouch
Secondary packaging	: Primary pack is loaded in a secondary pouch made of medical grade paper & laminated polyethylene film
Sizes	: USP 6/0 (metric 1) - USP 2 (metric 6)
Colour	: Yellow-tan
	
Sterilization	: Gamma radiation

Description	FG Code	Size		Needle Profile Shape, Type & Length	Suture Length (cm)	Pack Size	Suture Colour
MERISOFT PLAIN (PLAIN CATGUT SUTURE)							
PLAIN CATGUT 0 X 100-45MM HC RB	ME100429	0		45 mm 1/2 Circle Round Body	100	12	PLAIN
PLAIN CATGUT 0 X 75-30MM HC RB	ME104030	0		30 mm 1/2 Circle Round Body	75	12	PLAIN
PLAIN CATGUT 0 X 75-48MM HC RB	ME104048	0		48 mm 1/2 Circle Round Body	75	12	PLAIN
PLAIN CATGUT 1 X 75-30MM HC RB	ME0128126	1		30 mm 1/2 Circle Round Body	75	12	PLAIN
PLAIN CATGUT 1 X 75-40MM HC RB	ME010735	1		40 mm 1/2 Circle Round Body	75	12	PLAIN
PLAIN CATGUT 1 X 100-43MM HC RB	ME010427	1		43 mm 1/2 Circle Round Body	100	12	PLAIN
PLAIN CATGUT 1 X 100-48MM HC RB	ME010448	1		48 mm 1/2 Circle Round Body	100	12	PLAIN
PLAIN CATGUT 2-0 X 76-30MM CU RB	ME5033	2-0		30 mm 3/8 Circle Round Body	76	12	PLAIN
PLAIN CATGUT 2-0 X 100-37MM HC RB	ME204237	2-0		37 mm 1/2 Circle Round Body	100	12	PLAIN
PLAIN CATGUT 2-0 X 90-22MM HC RB	ME204222	2-0		22 mm 1/2 Circle Round Body	90	12	PLAIN
PLAIN CATGUT 2-0 X 90-26MM HC RB	ME204227	2-0		26 mm 1/2 Circle Round Body	90	12	PLAIN
PLAIN CATGUT 2-0 X 70-24MM CU RC	ME204224	2-0		24 mm 3/8 Circle Reverse Cutting	70	12	PLAIN
PLAIN CATGUT 2-0 X 90-50MM HC RB	ME207333	2-0		50 mm 1/2 Circle Round Body	50	12	PLAIN
PLAIN CATGUT 3-0 X 90-22MM HC RB	ME304322	3-0		22 mm 1/2 Circle Round Body	90	12	PLAIN
PLAIN CATGUT 3-0 X 90-17MM HC RB	ME304317	3-0		17 mm 1/2 Circle Round Body	90	12	PLAIN
PLAIN CATGUT 3-0 X 76-16MM CU CT	ME4285	3-0		16 mm 3/8 Circle Conventional Cutting	76	12	PLAIN
PLAIN CATGUT 3-0 X 75-26MM HC RB	ME304327SL	3-0		26 mm 1/2 Circle Round Body	75	12	PLAIN
PLAIN CATGUT 3-0 X 76-26MM HC RB	ME4236	3-0		26 mm 1/2 Circle Round Body	76	12	PLAIN
PLAIN CATGUT 4-0 X 76-16MM CU RC	ME4283	4-0		16 mm 3/8 Circle Reverse Cutting	76	12	PLAIN

Description	FG Code	Size		Needle Profile Shape, Type & Length	Suture Length (cm)	Pack Size	Suture Colour
MERISOFT PLAIN (PLAIN CATGUT SUTURE)							
PLAIN CATGUT 5-0 X 76-13MM CU RC	ME4282	5-0		13 mm 3/8 Circle Reverse Cutting	76	12	PLAIN
PLAIN CATGUT 0 X 152	SS2004	0		Non needled	152	12	PLAIN
PLAIN CATGUT 1 X 152	SS2005	1		Non needled	152	12	PLAIN
PLAIN CATGUT 2-0 X 152	SS2003	2-0		Non needled	152	12	PLAIN
PLAIN CATGUT 3-0 X 152	SS2002	3-0		Non needled	152	12	PLAIN

MERISOFT CHROMIC™

Chromic Catgut Suture

Properties

MERISOFT CHROMIC™ gut suture is treated with a chromium salt solution to resist the enzymatic action, prolonging the absorption time to 90 to 120 days.

Indications

MERISOFT CHROMIC™ suture is indicated for use in general soft tissue approximation and/or ligation.

Advantages


- Highest level of collagen guarantees predictable absorption
- Accurate polishing ensures smooth tissue passage
- Packed in liquid to reduce memory & increase pliability



Product Overview

Primary packaging	: Packed in an aluminum foil pouch
Secondary packaging	: Primary pack is loaded in a secondary pouch made of medical grade paper & laminated polyethylene film
Sizes	: USP 6/0 (metric 1) - USP 2 (metric 6)
Colour	: Brown dyed
	
Sterilization	: Gamma radiation

Description	FG Code	Size		Needle Profile Shape, Type & Length	Suture Length (cm)	Pack Size	Suture Colour
MERISOFT CHROMIC (CHROMIC CATGUT SUTURE)							
CHROMIC CATGUT 0 X 76-30MM CU RB	ME4217	0		30 mm 3/8 Circle Round Body	76	12	CHROMIC
CHROMIC CATGUT 0 X 76-40MM CU RB	ME4221	0		40 mm 3/8 Circle Round Body	76	12	CHROMIC
CHROMIC CATGUT 0 X 100-45MM HC RB	ME4229	0		45 mm 1/2 Circle Round Body	100	12	CHROMIC
CHROMIC CATGUT 0 X 76-30MM HC RB	ME4242	0		30 mm 1/2 Circle Round Body	76	12	CHROMIC
CHROMIC CATGUT 0 X 76-40MM HC RB TN	ME4246	0		40 mm 1/2 Circle Round Body	76	12	CHROMIC
CHROMIC CATGUT 0 X 76-40MM HC RB TH	ME4258	0		40 mm 1/2 Circle Round Body	76	12	CHROMIC
CHROMIC CATGUT 0 X 76-36MM HC RB	ME104036	0		36 mm 1/2 Circle Round Body	76	12	CHROMIC
CHROMIC CATGUT 0 X 100-40MM HC TC	ME4279	0		40 mm 1/2 Circle Taper Cut	100	12	CHROMIC
CHROMIC CATGUT 0 X 76-45MM ST RB	ME4205	0		45 mm Straight Round Body	76	12	CHROMIC
CHROMIC CATGUT 0 X 76-45MM CU RB	ME4225	0		45 mm 3/8 Circle Round Body	76	12	CHROMIC
CHROMIC CATGUT 0 X 76-45MM HC RB	ME4250	0		45 mm 1/2 Circle Round Body	76	12	CHROMIC
CHROMIC CATGUT 0 X 76-40MM HC TC	ME4257	0		40 mm 1/2 Circle Taper Cut	76	12	CHROMIC
CHROMIC CATGUT 0 X 76-40MM HC RC	ME104040	0		40 mm 1/2 Circle Reverse Cutting	76	12	CHROMIC
CHROMIC CATGUT 1 X 76-45MM CU RB	ME4226	1		45 mm 3/8 Circle Round Body	76	12	CHROMIC
CHROMIC CATGUT 1 X 100-45MM HC RB	ME4227	1		45 mm 1/2 Circle Round Body	100	12	CHROMIC
CHROMIC CATGUT 1 X 100-50MM HC RB	ME4248	1		50 mm 1/2 Circle Round Body	100	12	CHROMIC
CHROMIC CATGUT 1 X 76-40MM HC RB	ME4259	1		40 mm 1/2 Circle Round Body	76	12	CHROMIC

Description	FG Code	Size		Needle Profile Shape, Type & Length	Suture Length (cm)	Pack Size	Suture Colour
MERISOFT CHROMIC (CHROMIC CATGUT SUTURE)							
CHROMIC CATGUT 1 X 75-30MM HC RB	ME014995	1		30 mm 1/2 Circle Round Body	75	12	CHROMIC
CHROMIC CATGUT 1 X 100-40MM HC TRO	ME4239	1		40 mm 1/2 Circle Trocar	100	12	CHROMIC
CHROMIC CATGUT 1 X 76-63MM CU RB	Me4267	1		63 mm 3/8 Circle Round Body	76	12	CHROMIC
CHROMIC CATGUT 1 X 76-50MM HC RC	ME4273	1		50 mm 1/2 Circle Reverse Cutting	76	12	CHROMIC
CHROMIC CATGUT 1 X 100-40MM HC TC	ME4297	1		40 mm 1/2 Circle Taper Cut	100	12	CHROMIC
CHROMIC CATGUT 2 X 100-45MM HC RB	ME4228	2		45 mm 1/2 Circle Round Body	100	12	CHROMIC
CHROMIC CATGUT 2 X 100-50MM HC RB	ME024250	2		50 mm 1/2 Circle Round Body	100	12	CHROMIC
CHROMIC CATGUT 2 X 100-40MM HC TRO	ME4244	2		40 mm 1/2 Circle Trocar	100	12	CHROMIC
CHROMIC CATGUT 2 X 100-65MM HC RB	ME4278	2		65 mm 1/2 Circle Round Body	100	12	CHROMIC
CHROMIC CATGUT 2 X 75-40MM HC RC	ME024244	2		40 mm 1/2 Circle Reverse Cutting	75	12	CHROMIC
CHROMIC CATGUT 2-0 X 76-30MM CU RB	ME4216	2-0		30 mm 3/8 Circle Round Body	76	12	CHROMIC
CHROMIC CATGUT 2-0 X 76-30MM HC RB	ME4241	2-0		30 mm 1/2 Circle Round Body	76	12	CHROMIC
CHROMIC CATGUT 2-0 X 152-36MM HC RCRB DN	ME4277	2-0		36 mm 1/2 Circle Reverse Cutting Round Body	152	12	CHROMIC
CHROMIC CATGUT 2-0 X 75-40MM HC RB	ME200991	2-0		40 mm 1/2 Circle Round Body	75	12	CHROMIC
CHROMIC CATGUT 2-0 X 76-26MM HC RB	ME204225	2-0		26 mm 1/2 Circle Round Body	76	12	CHROMIC

Description	FG Code	Size		Needle Profile Shape, Type & Length	Suture Length (cm)	Pack Size	Suture Colour
MERISOFT CHROMIC (CHROMIC CATGUT SUTURE)							
CHROMIC CATGUT 2-0 X 76-40MM HC RB	ME4245	2-0		40 mm 1/2 Circle Round Body	76	12	CHROMIC
CHROMIC CATGUT 2-0 X 76-45MM HC RB	ME4249	2-0		45 mm 1/2 Circle Round Body	76	12	CHROMIC
CHROMIC CATGUT 2-0 X 76-30MM HC CT	ME4255	2-0		30 mm 1/2 Circle Conventional Cutting	76	12	CHROMIC
CHROMIC CATGUT 2-0 X 76-50MM ST RB	ME4208	2-0		50 mm Straight Round Body	76	12	CHROMIC
CHROMIC CATGUT 2-0 X 75-45MM HC CT	ME200863	2-0		45 mm 1/2 Circle Conventional Cutting	75	12	CHROMIC
CHROMIC CATGUT 3-0 X 76-22MM CU CT	ME4201	3-0		22 mm 3/8 Circle Conventional Cutting	76	12	CHROMIC
CHROMIC CATGUT 3-0 X 76-20MM HC RB	ME4237	3-0		20 mm 1/2 Circle Round Body	76	12	CHROMIC
CHROMIC CATGUT 3-0 X 75-26MM HC RB	ME300882	3-0		26 mm 1/2 Circle Round Body	75	12	CHROMIC
CHROMIC CATGUT 3-0 X 75-30MM HC RB	ME300172	3-0		30 mm 1/2 Circle Round Body	75	12	CHROMIC
CHROMIC CATGUT 3-0 X 75-26MM CU RC	ME304326	3-0		26 mm 3/8 Circle Reverse Cutting	75	12	CHROMIC
CHROMIC CATGUT 3-0 X 38-16MM FE CT	ME4265	3-0		16 mm 5/8 Circle Conventional Cutting	38	12	CHROMIC
CHROMIC CATGUT 3-0 X 76-16MM CU CT	ME4274	3-0		16 mm 3/8 Circle Conventional Cutting	76	12	CHROMIC
CHROMIC CATGUT 3-0 X 76-12MM HC RC	ME4286	3-0		12 mm 1/2 Circle Reverse Cutting	76	12	CHROMIC
CHROMIC CATGUT 3-0 X 76-26MM CU RC	ME4828	3-0		26 mm 3/8 Circle Reverse Cutting	76	12	CHROMIC
CHROMIC CATGUT 3-0 X 75-26MM HC RC	ME300886	3-0		26 mm 1/2 Circle Reverse Cutting	75	12	CHROMIC
CHROMIC CATGUT 4-0 X 76-16MM CU RB	ME4048	4-0		16 mm 3/8 Circle Round Body	76	12	CHROMIC
CHROMIC CATGUT 4-0 X 76-16MM CU RC	ME4287	4-0		16 mm 3/8 Circle Reverse Cutting	76	12	CHROMIC
CHROMIC CATGUT 4-0 X 76-17MM CU TC	ME5047	4-0		17 mm 3/8 Circle Taper Cut	76	12	CHROMIC
CHROMIC CATGUT 4-0 X 75-20MM HC RB	ME400181	4-0		20 mm 1/2 Circle Round Body	75	12	CHROMIC
CHROMIC CATGUT 4-0 X 75-17MM HC RB	ME400117	4-0		17 mm 1/2 Circle Round Body	75	12	CHROMIC

Description	FG Code	Size		Needle Profile Shape, Type & Length	Suture Length (cm)	Pack Size	Suture Colour
MERISOFT CHROMIC (CHROMIC CATGUT SUTURE)							
CHROMIC CATGUT 4-0 X 75-19MM CU RC	ME400419	4-0		19 mm 3/8 Circle Reverse Cutting	75	12	CHROMIC
CHROMIC CATGUT 4-0 X 38-12MM CU RC	ME4266	4-0		12 mm 3/8 Circle Reverse Cutting	38	12	CHROMIC
CHROMIC CATGUT 4-0 X 75-26MM HC RB	ME400121	4-0		26 mm 1/2 Circle Round Body	75	12	CHROMIC
CHROMIC CATGUT 5-0 X 90-12MM CU RC	ME500268	5-0		12 mm 3/8 Circle Reverse Cutting	90	12	CHROMIC
CHROMIC CATGUT 5-0 X 30-10MM FE TC	ME4263	5-0		10 mm 5/8 Circle Taper Cut	30	12	CHROMIC
CHROMIC CATGUT 5-0 X 76-13MM CU RC	ME4276	5-0		13 mm 3/8 Circle Reverse Cutting	76	12	CHROMIC
CHROMIC CATGUT 5-0 X 76-13MM CU RB	ME4289	5-0		13 mm 3/8 Circle Round Body	76	12	CHROMIC
CHROMIC CATGUT 1 X 75-60MM HC RB	ME016075	1		60 mm 1/2 Circle Round Body	75	12	CHROMIC
CHROMIC CATGUT 2 X 75-80MM HC RB	ME028075	2		80 mm 1/2 Circle Round Body	75	12	CHROMIC
CHROMIC CATGUT 0 X 152	SS2214	0		Non needed	152	12	CHROMIC
CHROMIC CATGUT 1 X 152	SS2215	1		Non needed	152	12	CHROMIC
CHROMIC CATGUT 2 X 152	SS2216	2		Non needed	152	12	CHROMIC
CHROMIC CATGUT 2-0 X 152	SS2213	2-0		Non needed	152	12	CHROMIC

FILASILK™

Silk Suture

Indications

- General surgery procedures
- Skin closure
- Oral surgery procedures
- Ophthalmic procedures
- Ligatures

Advantages

- Excellent handling properties
- Good knot security



Product Overview

Structure	: Braided
Colour	: Black
Origin	: Silkworm larval cocoons (Bombyx mori)
Composition	: Natural silk filaments
Sizes	: USP 8/0 (metric 0.4) - USP 2 (metric 5)
Sterilization	: Gamma radiation



Description	FG Code	Size		Needle Profile Shape, Type & Length	Suture Length (cm)	Pack Size	Suture Colour
FILASILK SILK SUTURE							
SILK BLK 1 X 75-30MM CU RC	SLK010689	1		30 mm 3/8 Circle Reverse Cutting	75	12	BLACK
SILK BLK 2 X 90-90MM CU RC	SLK020794	2		90 mm 3/8 Circle Reverse Cutting	90	12	BLACK
SILK BLK 2 X 90-40MM HC RB	SLK020936	2		40 mm 1/2 Circle Round Body	90	12	BLACK
SILK BLK 2 X 75-30MM HC RC	SLK025230	2		30 mm 1/2 Circle Reverse Cutting	75	12	BLACK
SILK BLK 2 X 75-40MM HC RC	SLK025340	2		40 mm 1/2 Circle Reverse Cutting	75	12	BLACK
SILK BLK 2 X 75-40MM HC RB	SLK020936SL	2		40 mm 1/2 Circle Round Body	75	12	BLACK
SILK BLK 2 X 100-90MM CU CT	SLK020792	2		90 mm 3/8 Circle Conventional Cutting	100	12	BLACK
SILK BLK 2 X 70-36MM CU RC	SLK021236	2		36 mm 3/8 Circle Reverse Cutting	70	12	BLACK
SILK BLK 2 X 70-39MM CU RC	SLK021239	2		39 mm 3/8 Circle Reverse Cutting	70	12	BLACK
SILK BLK 2 X 75-35MM HC RB	SLK028235	2		35 mm 1/2 Circle Round Body	75	12	BLACK
SILK BLK 1 X 50-90MM CU CT (2 suture per foil)	SLK012797	1		90 mm 3/8 Circle Conventional Cutting	50	12	BLACK
SILK BLK 2-0 X 90-20MM CU RC	SLK200367	2-0		20 mm 3/8 Circle Reverse Cutting	90	12	BLACK
SILK BLK 2-0 X 90-60MM ST CT	SLK200790	2-0		60 mm Straight Conventional Cutting	90	12	BLACK
SILK BLK 2-0 X 76-45MM CU RC	SLK20 5036	2-0		45 mm 3/8 Circle Reverse Cutting	76	12	BLACK
SILK BLK 2-0 X 76-50MM HC RC	SLK20 5052	2-0		50 mm 1/2 Circle Reverse Cutting	76	12	BLACK
SILK BLK 2-0 X 90-30MM HC RB	SLK20 5331	2-0		30 mm 1/2 Circle Round Body	90	12	BLACK
SILK BLK 2-0 X 76-25MM HC TC	SLK20 5670	2-0		25 mm 1/2 Circle Taper Cut	76	12	BLACK
SILK BLK 2-0 X 45-23MM HC RC	SLK200737	2-0		23 mm 1/2 Circle Reverse Cutting	45	12	BLACK
SILK BLK 2-0 X 75-22MM HC RB	SLK205222	2-0		22 mm 1/2 Circle Round Body	75	12	BLACK

Description	FG Code	Size		Needle Profile Shape, Type & Length	Suture Length (cm)	Pack Size	Suture Colour
FILASILK SILK SUTURE							
SILK BLK 2-0 X 75-26MM HC RB	SLK200327	2-0		26 mm 1/2 Circle Round Body	75	12	BLACK
SILK BLK 2-0 X 75-17MM HC RB	SLK205610	2-0		17 mm 1/2 Circle Round Body	75	12	BLACK
SILK BLK 2-0 X 75-25MM HC CT	SLK205225	2-0		25 mm 1/2 Circle Conventional Cutting	75	12	BLACK
SILK BLK 2-0 X 76-50MM ST RB	SLK20 5208	2-0		50 mm Straight Round Body	76	12	BLACK
SILK BLK 2-0 X 90-15MM CU RC	SLK200503	2-0		15 mm 3/8 Circle Reverse Cutting	90	12	BLACK
SILK BLK 2-0 X 75-22MM HC CT	SLK2052221	2-0		22 mm 1/2 Circle Conventional Cutting	75	12	BLACK
SILK BLK 2-0 X 76-60MM ST CT	SLK20 5290	2-0		60 mm Straight Conventional Cutting	76	12	BLACK
SILK BLK 2-0 X 100-90MM CU CT	SLK200791	2-0		90 mm 3/8 Circle Conventional Cutting	100	12	BLACK
SILK BLK 2-0 X 45-26MM CU RC	SLK200685	2-0		26 mm 3/8 Circle Reverse Cutting	45	12	BLACK
SILK BLK 2-0 X 75-30MM HC CT	SLK205231	2-0		30 mm 1/2 Circle Conventional Cutting	75	12	BLACK
SILK BLK 2-0 X 90-26MM HC RB	SLK200327XL	2-0		26 mm 1/2 Circle Round Body	90	12	BLACK
SILK BLK 2-0 X 75-26MM CU RC	SLK200685XL	2-0		26 mm 3/8 Circle Reverse Cutting	75	12	BLACK
SILK BLK 2-0 X 90-24MM CU RB	SLK205224	2-0		24 mm 3/8 Circle Round Body	90	12	BLACK
SILK BLK 2-0 X 90-26MM CU CT	SLK205226	2-0		26 mm 3/8 Circle Conventional Cutting	90	12	BLACK
SILK BLK 2-0 X 45-30MM HC RC	SLK200530	2-0		30 mm 1/2 Circle Reverse Cutting	45	12	BLACK
SILK BLK 2-0 X 75-30MM CU RC	SLK205230	2-0		30 mm 3/8 Circle Reverse Cutting	75	12	BLACK
SILK BLK 2-0 X 75-16MM CU CT	SLK200509	2-0		16 mm 3/8 Circle Conventional Cutting	75	12	BLACK
SILK BLK 2-0 X 76-25MM HC RB	SLK205671	2-0		25 mm 1/2 Circle Round Body	76	12	BLACK
SILK BLK 2-0 X 100-55MM ST CT	SLK200775	2-0		55 mm Straight Conventional Cutting	100	12	BLACK
SILK BLK 2-0 X 76-36MM CU RC	SLK205236	2-0		36 mm 3/8 Circle Reverse Cutting	76	12	BLACK
SILK BLK 2-0 X 76-26MM HC RB VB	SLK20 5065	2-0		26 mm 1/2 Circle Round Body	76	12	BLACK

Description	FG Code	Size		Needle Profile Shape, Type & Length	Suture Length (cm)	Pack Size	Suture Colour
FILASILK SILK SUTURE							
SILK BLK 2-0 X 90-17MM HC TC	SLK20 5611	2-0		17 mm 1/2 Circle Taper Cut	90	12	BLACK
SILK BLK 2-0 X 76-30MM HC RB VB	SLK20 533VB	2-0		30 mm 1/2 Circle Round Body	76	12	BLACK
SILK BLK 2-0 X 100-75MM ST CT	SLK2013102	2-0		75 mm Straight Conventional Cutting	100	12	BLACK
SILK BLK 2-0 X 90-35MM CU RC	SLK200667	2-0		35 mm 3/8 Circle Reverse Cutting	90	12	BLACK
SILK BLK 2-0 X 75-36MM CU RC	SLK205275	2-0		36 mm 3/8 Circle Reverse Cutting	75	12	BLACK
SILK BLK 2-0 X 75-35MM CU CT	SLK202429	2-0		35 mm 3/8 Circle Conventional Cutting	75	12	BLACK
SILK BLK 2-0 X 75-80MM CU CT	SLK200261	2-0		80 mm 3/8 Circle Conventional Cutting	75	12	BLACK
SILK BLK 2-0 X 75-30MM HC RB	SLK200333	2-0		30 mm 1/2 Circle Round Body	75	12	BLACK
SILK BLK 2-0 X 75-48MM HC CT	SLK200552	2-0		48 mm 1/2 Circle Conventional Cutting	75	12	BLACK
SILK BLK 2-0 X 75-60MM ST CT	SLK200790SL	2-0		60 mm Straight Conventional Cutting	75	12	BLACK
SILK BLK 2-0 X 45-24MM CU RC	SLK200695	2-0		24 mm 3/8 Circle Reverse Cutting	45	12	BLACK
SILK BLK 2-0 X 100-40MM CU CT	SLK203959	2-0		40 mm 3/8 Circle Conventional Cutting	100	12	BLACK
SILK BLK 2-0 X 45-30MM HC RC C1	SLK205227	2-0		30 mm 1/2 Circle Reverse Cutting	45	12	BLACK
SILK BLK 2-0 X 75-75MM ST CT	SLK200766	2-0		75 mm Straight Conventional Cutting	75	12	BLACK
SILK BLK 2-0 X 76-30MM HC RB	SLK20 5333	2-0		30 mm 1/2 Circle Round Body	76	12	BLACK
SILK BLK 2-0 X 90-35MM HC CT	SLK20 5335	2-0		35 mm 1/2 Circle Conventional Cutting	90	12	BLACK
SILK BLK 2-0 X 75-39MM CU RC	SLK200648	2-0		39 mm 3/8 Circle Reverse Cutting	75	12	BLACK
SILK BLK 2-0 X 100-60MM ST RC	SLK200723	2-0		60 mm Straight Reverse Cutting	100	12	BLACK
SILK BLK 2-0 X 75-40MM HC RB	SLK200584	2-0		40 mm 1/2 Circle Round Body	75	12	BLACK

SILK BLK 3-0 X 45-20MM HC RB	SLK30 5087SL	3-0		20 mm 1/2 Circle Round Body	45	12	BLACK
SILK BLK 3-0 X 45-26MM CU RC	SLK30 5028SL	3-0		26 mm 3/8 Circle Reverse Cutting	45	12	BLACK
SILK BLK 3-0 X 76-22MM HC RB VB	SLK30 5064	3-0		22 mm 1/2 Circle Round Body	76	12	BLACK
SILK BLK 3-0 X 90-35MM HC RB	SLK30 5100	3-0		35 mm 1/2 Circle Round Body	90	12	BLACK
SILK BLK 3-0 X 76-20MM HC RB VB	SLK30 587VB	3-0		20 mm 1/2 Circle Round Body	76	12	BLACK
SILK BLK 3-0 X 75-16MM CU CT	SLK3013086	3-0		16 mm 3/8 Circle Conventional Cutting	75	12	BLACK
SILK BLK 3-0 X 75-60MM ST CT	SLK3013089	3-0		60 mm Straight Conventional Cutting	75	12	BLACK
SILK BLK 3-0 X 45-24MM CU RC	SLK300684SL	3-0		24 mm 3/8 Circle Reverse Cutting	45	12	BLACK
SILK BLK 3-0 X 45-25MM HC RC	SLK305323	3-0		25 mm 1/2 Circle Reverse Cutting	45	12	BLACK
SILK BLK 3-0 X 45-22MM HC CT	SLK302508SL	3-0		22 mm 1/2 Circle Conventional Cutting	45	12	BLACK
SILK BLK 3-0 X 35-17MM HC RB (6 in 1)	SLK306540	3-0		17 mm 1/2 Circle Round Body	35	12	BLACK
SILK BLK 3-0 X 90-35MM CU RC	SLK300666	3-0		35 mm 3/8 Circle Reverse Cutting	90	12	BLACK
SILK BLK 3-0 X 76-16MM CU CT	SLK30 5002	3-0		16 mm 3/8 Circle Conventional Cutting	76	12	BLACK
SILK BLK 3-0 X 76-25MM HC RB	SLK30 5070	3-0		25 mm 1/2 Circle Round Body	76	12	BLACK
SILK BLK 3-0 X 75-17MM HC TC	SLK3054602	3-0		17 mm 1/2 Circle Taper Cut	75	12	BLACK
SILK BLK 3-0 X 45-15MM HC RC	SLK305345	3-0		15 mm 1/2 Circle Reverse Cutting	45	12	BLACK
SILK BLK 4-0 X 90-20MM HC RC	SLK400377	4-0		20 mm 1/2 Circle Reverse Cutting	90	12	BLACK
SILK BLK 4-0 X 90-16MM CU CT	SLK40 5000	4-0		16 mm 3/8 Circle Conventional Cutting	90	12	BLACK
SILK BLK 4-0 X 90-16MM HC RB	SLK40 5050	4-0		16 mm 1/2 Circle Round Body	90	12	BLACK
SILK BLK 4-0 X 76-16MM CU RB	SLK40 5082	4-0		16 mm 3/8 Circle Round Body	76	12	BLACK
SILK BLK 4-0 X 76-20MM HC RB	SLK40 5086	4-0		20 mm 1/2 Circle Round Body	76	12	BLACK
SILK BLK 4-0 X 75-24MM CU RC	SLK400629	4-0		24 mm 3/8 Circle Reverse Cutting	75	12	BLACK

Description	FG Code	Size	Needle Profile Shape, Type & Length	Suture Length (cm)	Pack Size	Suture Colour
FILASILK SILK SUTURE						
SILK BLK 1 X 45 (45X12)	SLS014512	1	Non needed	45	12	BLACK
SILK BLK 2 X 76 (76X2)	SLS02 216	2	Non needed	76	12	BLACK
SILK BLK 2 X 180	SLS02180	2	Non needed	180	12	BLACK
SILK BLK 2 X 270 (45X6)	SLS02364	2	Non needed	45	12	BLACK
SILK BLK 2 X 250	SLS02252	2	Non needed	250	12	BLACK
SILK BLK 2 X 45 (45X12)	SLS024512	2	Non needed	45	12	BLACK
SILK BLK 3 X 76 (76X2)	SLS03 217	3	Non needed	76	12	BLACK
SILK BLK 2-0 X 76 (76X2)	SLS20 213	2-0	Non needed	76	12	BLACK
SILK BLK 2-0 X 75 (75X12)	SLS20305	2-0	Non needed	75	12	BLACK
SILK BLK 2-0 X 180	SLS20180	2-0	Non needed	180	12	BLACK
SILK BLK 2-0 X 75 (75X5)	SLS20375	2-0	Non needed	75	12	BLACK
SILK BLK 2-0 X 35 (35X10)	SLS20 223	2-0	Non needed	35	12	BLACK
SILK BLK 2-0 X 45 (45X17)	SLS20203	2-0	Non needed	45	12	BLACK
SILK BLK 2-0 X 75 (75X10)	SLS201075	2-0	Non needed	75	12	BLACK
SILK BLK 2-0 X 250-2D	SLS20253Y	2-0	Non needed	250	24	BLACK
SILK BLK 2-0 X 150 (150X2)	SLS20605	2-0	Non needed	150	12	BLACK
SILK BLK 2-0 X 250	SLS20253	2-0	Non needed	250	12	BLACK
SILK BLK 2-0 X 45 (45X12)	SLS204512	2-0	Non needed	45	12	BLACK
SILK BLK 3-0 X 76 (76X2)	SLS30 212	3-0	Non needed	76	12	BLACK
SILK BLK 3-0 X 35 (35X10)	SLS30 222	3-0	Non needed	35	12	BLACK
SILK BLK 3-0 X 75 (75X10)	SLS30750	3-0	Non needed	75	12	BLACK
SILK BLK 3-0 X 180	SLS30180	3-0	Non needed	180	12	BLACK
SILK BLK 3-0 X 45 (45X17)	SLS30202	3-0	Non needed	45	12	BLACK
SILK BLK 3-0 X 45 (45X12)	SLS304512	3-0	Non needed	45	12	BLACK
SILK BLK 4-0 X 76 (76X2)	SLS40 211	4-0	Non needed	76	12	BLACK
SILK BLK 0 X 25MTR (6 EA PER BOX)	SLK10 824R	0	Non needed	25 Mtr	6	BLACK
SILK BLK 1 X 25MTR (6 EA PER BOX)	SLK01 825R	1	Non needed	25 Mtr	6	BLACK
SILK BLK 2 X 25MTR (6 EA PER BOX)	SLK02 826R	2	Non needed	25 Mtr	6	BLACK
SILK BLK 2-0 X 25MTR (6 EA PER BOX)	SLK20 823R	2-0	Non needed	25 Mtr	6	BLACK
SILK BLK 3-0 X 25MTR (6 EA PER BOX)	SLK30 822R	3-0	Non needed	25 Mtr	6	BLACK

MERISTEELTM

Steel Suture

Indications


- Orthopedic procedures
- Sternum closure

Advantages

- Exceptional tensile strength
- Excellent tissue compatibility



Product Overview

Structure	: Monofilament
Colour	: Silver/Steel
	
Chemical composition	: Stainless steel
Sizes	: USP 5 (metric 1) - USP 6 (metric 8)
Sterilization	: Gamma radiation

Description	FG Code	Size		Needle Profile Shape, Type & Length	Suture Length (cm)	Pack Size	Suture Colour
MERISTEEL (STEEL SUTURE)							
STEEL WIRE 0 X 45-17MM HC BP (45X2)	STC600	0		17 mm 1/2 Circle Blunt Point	45	12	-
STEEL WIRE 1 X 45-40MM HC TC (45X4)	STC660	1		40 mm 1/2 Circle Taper Cut	45	12	-
STEEL WIRE 1 X 45-40MM HC BP (45X4)	STC661	1		40 mm 1/2 Circle Blunt Point	45	12	-
STEEL WIRE 1 X 45-35MM HC RB BLUNT	STC9435	1		35 mm 1/2 Circle Blunt Point	45	12	-
STEEL WIRE 2 X 45-40MM HC TC (45X4)	STC650	2		40 mm 1/2 Circle Taper Cut	45	12	-
STEEL WIRE 2 X 45-40MM HC BP (45X4)	STC662	2		40 mm 1/2 Circle Blunt Point	45	12	-
STEEL WIRE 2 X 45-40MM HC RB	STC9540	2		40 mm 1/2 Circle Round Body	45	12	-
STEEL WIRE 4 X 45-40MM HC BP (45X4)	STC664	4		40 mm 1/2 Circle Blunt Point	45	12	-
STEEL WIRE 4 X 45-48MM HC CT (45X4)	STC652	4		48 mm 1/2 Circle Conventional Cutting	45	12	-
STEEL WIRE 4 X 75-48MM HC TC (75X2)	STC040651	4		48 mm 1/2 Circle Taper Cut	75	12	-
STEEL WIRE 5 X 75-55MM HC RC(75X2)	STC050945	5		55 mm 1/2 Circle Reverse Cutting	75	12	-
STEEL WIRE 5 X 45-48MM HC CT (45X4)	STC653	5		48 mm 1/2 Circle Conventional Cutting	45	12	-
STEEL WIRE 5 X 45-44MM HC RB BP (45X4)	STC9455	5		44 mm 1/2 Circle Blunt Point	45	6	-
STEEL WIRE 5 X 75-44MM HC RB BP (75X2)	STC9495	5		44 mm 1/2 Circle Blunt Point	75	6	-
STEEL WIRE 5 X 45-48MM HC RB BP (45X4)	STC9497	5		48 mm 1/2 Circle Blunt Point	45	12	-
STEEL WIRE 5 X 45-48MM HC CT (45X1)	STC9653	5		48 mm 1/2 Circle Conventional Cutting	45	12	-

Description	FG Code	Size		Needle Profile Shape, Type & Length	Suture Length (cm)	Pack Size	Suture Colour
MERISTEEL (STEEL SUTURE)							
STEEL WIRE 5 X 45-48MM HC TC (45X4)	STC05650	5		48 mm 1/2 Circle Taper Cut	45	12	-
STEEL WIRE 6 X 45-48MM HC CT (45X4)	STC654	6		48 mm 1/2 Circle Conventional Cutting	45	12	-
STEEL WIRE 6 X 45-44MM HC RB BP (45X4)	STC9456	6		44 mm 1/2 Circle Blunt Point	45	6	-
STEEL WIRE 6 X 75-44MM HC RB BP (75X2)	STC9494	6		44 mm 1/2 Circle Blunt Point	75	6	-
STEEL WIRE 6 X 45-48MM HC CT (45X1)	STC9654	6		48 mm 1/2 Circle Conventional Cutting	45	12	-
STEEL WIRE 6 X 45-48MM CU CT (45X6)	STC9654CU	6		48 mm 3/8 Circle Conventional Cutting	45	12	-
SS WIRE 6X45-48 HCCT (45X1)44 HCBP(45X3)	STC700	6		48/44 mm 1/2 Circle Conventional Cutting Blunt Point	45	12	-
STEEL WIRE 6 X 75-65MM HC RB BLUNT	STC9656	6		65 mm 1/2 Circle Blunt Point	75	12	-
STEEL WIRE 7 X 75-50MM HC RC (75X2)	STC070624	7		50 mm 1/2 Circle Reverse Cutting	75	12	-

Indications

- Skin closure (intracutaneous, subcutaneous)
- Plastic and reconstructive procedures
- Microsurgery procedures
- Ophthalmic procedures

Advantages

- Flexible, easy to handle and tie
- Smooth passage through tissue
- Excellent histocompatibility



Product Overview

Structure	: Monofilament
Colour	: Black
Chemical composition	: Polyamide 6 - 6.6
Sizes	: USP 10/0 (metric 0.2) - USP 2 (metric 5)
Sterilization	: Gamma radiation

Description	FG Code	Size		Needle Profile Shape, Type & Length	Suture Length (cm)	Pack Size	Suture Colour
FILAMIDE (POLYAMIDE SUTURE)							
POLYAMIDE BLK 1 X 100-30MM HC RB	NYL010500	1		30 mm 1/2 Circle Round Body	100	12	BLACK
POLYAMIDE BLK 1 X 100-48MM HC RB	NYL010749	1		48 mm 1/2 Circle Round Body	100	12	BLACK
POLYAMIDE BLK 1 X 100-48MM HC RC	NYL013148	1		48 mm 1/2 Circle Reverse Cutting	100	12	BLACK
POLYAMIDE BLK 1 X 200-45MM HC TC LOOP	NYL010758	1		45 mm 1/2 Circle Taper Cut	200	12	BLACK
POLYAMIDE BLK 1 X 150-50MM HC RB BP LOOP	NYL010729XL	1		50 mm 1/2 Circle Blunt Point	150	12	BLACK
POLYAMIDE BLK 1 X 200-48MM HC RB LOOP	NYL010768	1		48 mm 1/2 Circle Round Body	200	12	BLACK
POLYAMIDE BLK 2 X 90-48MM HC RC	NYL023259	2		48 mm 1/2 Circle Reverse Cutting	90	12	BLACK
POLYAMIDE BLK 2 X 75-45MM HC RB	NYL023245	2		45 mm 1/2 Circle Round Body	75	12	BLACK
POLYAMIDE BLK 2 X 90-90MM CU RC	NYL023290	2		90 mm 3/8 Circle Reverse Cutting	90	12	BLACK
POLYAMIDE BLK 2 X 75-40MM HC RC	NYL022407	2		40 mm 1/2 Circle Reverse Cutting	75	12	BLACK
POLYAMIDE BLK 2 X 70-90MM CU CT	NYL02 3398	2		90 mm 3/8 Circle Conventional Cutting	70	12	BLACK
POLYAMIDE BLK 2 X 75-35MM HC RB	NYL023235	2		35 mm 1/2 Circle Round Body	75	12	BLACK
POLYAMIDE BLK 2 X 45-78MM CU CT	NYL020290	2		78 mm 3/8 Circle Conventional Cutting	45	12	BLACK
POLYAMIDE BU 2 X 100-90MM CU CT	NYL020798BU	2		90 mm 3/8 Circle Conventional Cutting	100	12	BLUE
POLYAMIDE BLK 2 X 75-40MM HC RC EX	NYL02 3240	2		40 mm 1/2 Circle Reverse Cutting	75	12	BLACK

Description	FG Code	Size		Needle Profile Shape, Type & Length	Suture Length (cm)	Pack Size	Suture Colour
FILAMIDE (POLYAMIDE SUTURE)							
POLYAMIDE BLK 2-0 X 70-45MM CU RC	NYL20 3336	2-0		45 mm 3/8 Circle Reverse Cutting	70	12	BLACK
POLYAMIDE BLK 2-0 X 35-45MM CU RC	NYL20 3336SL	2-0		45 mm 3/8 Circle Reverse Cutting	35	12	BLACK
POLYAMIDE BLK 2-0 X 100-31MM HC RB	NYL200568	2-0		31 mm 1/2 Circle Round Body	100	12	BLACK
POLYAMIDE BLK 2-0 X 45-26MM CU RC	NYL200664	2-0		26 mm 3/8 Circle Reverse Cutting	45	12	BLACK
POLYAMIDE BLK 2-0 X 90-26MM CU RC	NYL201626	2-0		26 mm 3/8 Circle Reverse Cutting	90	12	BLACK
POLYAMIDE BLK 2-0 X 90-30MM CU RC	NYL203232	2-0		30 mm 3/8 Circle Reverse Cutting	90	12	BLACK
POLYAMIDE BLK 2-0 X 75-24MM CU CT	NYL203224	2-0		24 mm 3/8 Circle Conventional Cutting	75	12	BLACK
POLYAMIDE BLK 2-0 X 75-40MM HC RB	NYL205240	2-0		40 mm 1/2 Circle Round Body	75	12	BLACK
POLYAMIDE BLK 2-0 X 70-60MM ST CT	NYL20 3390	2-0		60 mm Straight Conventional Cutting	70	12	BLACK
POLYAMIDE BU 2-0 X 100-55MM ST CT	NYL2013506	2-0		55 mm Straight Conventional Cutting	100	12	BLUE
POLYAMIDE BU 2-0 X 90-39MM CU RC	NYL202751	2-0		39 mm 3/8 Circle Reverse Cutting	90	12	BLUE
POLYAMIDE BU 2-0 X 100-45MM CU RC	NYL205859	2-0		45 mm 3/8 Circle Reverse Cutting	100	12	BLUE
POLYAMIDE BU 2-0 X 75-26MM CU RC	NYL201626BU	2-0		26 mm 3/8 Circle Reverse Cutting	75	12	BLUE
POLYAMIDE BU 2-0 X 45-40MM CU CT	NYL200786	2-0		40 mm 3/8 Circle Conventional Cutting	45	12	BLUE
POLYAMIDE BLK 2-0 X 100-45MM CU RC	NYL20 3336XL	2-0		45 mm 3/8 Circle Reverse Cutting	100	12	BLACK
POLYAMIDE BU 2-0 X 100-90MM CU RC	NYL200795	2-0		90 mm 3/8 Circle Reverse Cutting	100	12	BLUE
POLYAMIDE BLK 2-0 X 100-45MM CU RC EX	NYL200736	2-0		45 mm 3/8 Circle Reverse Cutting	100	12	BLACK
POLYAMIDE BLK 2-0 X 70-40MM CU RC	NYL20 3333	2-0		40 mm 3/8 Circle Reverse Cutting	70	12	BLACK
POLYAMIDE BLK 2-0 X 100-60MM CU RC	NYL200728	2-0		60 mm 3/8 Circle Reverse Cutting	100	12	BLACK
POLYAMIDE BU 2-0 X 100-90MM CU CT	NYL200796	2-0		90 mm 3/8 Circle Conventional Cutting	100	12	BLUE

Description	FG Code	Size		Needle Profile Shape, Type & Length	Suture Length (cm)	Pack Size	Suture Colour
FILAMIDE (POLYAMIDE SUTURE)							
POLYAMIDE BLK 2-0 X 90-30MM CU RC-3D	NYL203232Z	2-0		30 mm 3/8 Circle Reverse Cutting	90	36	BLACK
POLYAMIDE BLK 2-0 X 90-26MM HC RC	NYL209026	2-0		26 mm 1/2 Circle Reverse Cutting	90	12	BLACK
POLYAMIDE BLK 2-0 X 70-30MM CU RB	NYL207030	2-0		30 mm 3/8 Circle Round Body	70	12	BLACK
POLYAMIDE BLK 2-0 X 90-55MM ST CT	NYL200776	2-0		55 mm Straight Conventional Cutting	90	12	BLACK
POLYAMIDE BLK 2-0 X 90-30MM HC RC	NYL203682	2-0		30 mm 1/2 Circle Reverse Cutting	90	12	BLACK
POLYAMIDE BLK 2-0 X 90-60MM ST CT	NYL200628	2-0		60 mm Straight Conventional Cutting	90	12	Black
POLYAMIDE BLK 2-0 X 75-19MM CU RC	NYL2033210	2-0		19 mm 3/8 Circle Reverse Cutting	75	12	BLACK
POLYAMIDE BLK 2-0 X 45-24MM CU RC	NYL2093368	2-0		24 mm 3/8 Circle Reverse Cutting	45	12	BLACK
POLYAMIDE BLK 3-0 X 90-25MM CU RB	NYL306631	3-0		25 mm 3/8 Circle Round Body	90	12	BLACK
POLYAMIDE BLK 3-0 X 90-35MM CU RC	NYL301627	3-0		35 mm 3/8 Circle Reverse Cutting	90	12	BLACK
POLYAMIDE BLK 3-0 X 70-26MM CU RC	NYL30 3328	3-0		26 mm 3/8 Circle Reverse Cutting	70	12	BLACK
POLYAMIDE BLK 3-0 X 35-26MM CU RC	NYL30 3328SL	3-0		26 mm 3/8 Circle Reverse Cutting	35	12	BLACK
POLYAMIDE BLK 3-0 X 70-16MM CU CT	NYL30 3321	3-0		16 mm 3/8 Circle Conventional Cutting	70	12	BLACK
POLYAMIDE BLK 3-0 X 45-19MM CU RC	NYL301621	3-0		19 mm 3/8 Circle Reverse Cutting	45	12	BLACK
POLYAMIDE BLK 3-0 X 90-26MM CU RC	NYL301625	3-0		26 mm 3/8 Circle Reverse Cutting	90	12	BLACK
POLYAMIDE BLK 3-0 X 90-19MM CU RC	NYL301622	3-0		19 mm 3/8 Circle Reverse Cutting	90	12	BLACK
POLYAMIDE BLK 3-0 X 75-35MM HC RB	NYL303335	3-0		35 mm 1/2 Circle Round Body	75	12	BLACK
POLYAMIDE BLK 3-0 X 75-25MM CU CT	NYL303326	3-0		25 mm 3/8 Circle Conventional Cutting	75	12	BLACK
POLYAMIDE BLK 3-0 X 75-40MM HC RB	NYL305340	3-0		40 mm 1/2 Circle Round Body	75	12	BLACK
POLYAMIDE BLK 3-0 X 70-60MM ST CT	NYL30 3388	3-0		60 mm Straight Conventional Cutting	70	12	BLACK

Description	FG Code	Size		Needle Profile Shape, Type & Length	Suture Length (cm)	Pack Size	Suture Colour
FILAMIDE (POLYAMIDE SUTURE)							
POLYAMIDE BU 3-0 X 45-19MM CU RC	NYL3013498	3-0		19 mm 3/8 Circle Reverse Cutting	45	12	BLUE
POLYAMIDE BU 3-0 X 75-26MM CU RC	NYL301685	3-0		26 mm 3/8 Circle Reverse Cutting	75	12	BLUE
POLYAMIDE BU 3-0 X 75-60MM ST CT	NYL306441	3-0		60 mm Straight Conventional Cutting	75	12	BLUE
POLYAMIDE BU 3-0 X 45-26MM CU RC	NYL301625BU	3-0		26 mm 3/8 Circle Reverse Cutting	45	12	BLUE
POLYAMIDE BU 3-0 X 45-24MM CU RC	NYL301857	3-0		24 mm 3/8 Circle Reverse Cutting	45	12	BLUE
POLYAMIDE BLK 3-0 X 100-30MM CU RC-3D	NYL303330Z	3-0		30 mm 3/8 Circle Reverse Cutting	100	36	BLACK
POLYAMIDE BLK 3-0 X 90-22MM HC RB	NYL309022	3-0		22 mm 1/2 Circle Round Body	90	12	BLACK
POLYAMIDE BLK 3-0 X 45-26MM CU RC	NYL301625SL	3-0		26 mm 3/8 Circle Reverse Cutting	45	12	BLACK
POLYAMIDE BLK 3-0 X 75-24MM CU RC	NYL305669	3-0		24 mm 3/8 Circle Reverse Cutting	75	12	BLACK
POLYAMIDE BLK 3-0 X 75-60MM ST CT	NYL300627	3-0		60 mm Straight Conventional Cutting	75	12	Black
POLYAMIDE BU 3-0 X 100-26MM CU CT	NYL3013500	3-0		26 mm 3/8 Circle Conventional Cutting	100	12	Blue
POLYAMIDE BLK 3-0 X 100-30MM CU RC	NYL303330	3-0		30 mm 3/8 Circle Reverse Cutting	100	12	Black
POLYAMIDE BLK 3-0 X 75-20MM HC RC	NYL303370	3-0		20 mm 1/2 Circle Reverse Cutting	75	12	Black
POLYAMIDE BLK 4-0 X 90-16MM CU CT	NYL40 3319	4-0		16 mm 3/8 Circle Conventional Cutting	90	12	BLACK
POLYAMIDE BLK 4-0 X 45-11MM CU RC	NYL401612	4-0		11 mm 3/8 Circle Reverse Cutting	45	12	BLACK
POLYAMIDE BLK 4-0 X 75-19MM CU RC	NYL401619	4-0		19 mm 3/8 Circle Reverse Cutting	75	12	BLACK
POLYAMIDE BLK 4-0 X 45-19MM CU RC	NYL401620	4-0		19 mm 3/8 Circle Reverse Cutting	45	12	BLACK
POLYAMIDE BLK 4-0 X 45-16MM CU RC	NYL401861	4-0		16 mm 3/8 Circle Reverse Cutting	45	12	BLACK
POLYAMIDE BLK 4-0 X 90-19MM CU RC	NYL401618	4-0		19 mm 3/8 Circle Reverse Cutting	90	12	BLACK
POLYAMIDE BLK 4-0 X 75-19MM CU CT	NYL403419	4-0		19 mm 3/8 Circle Conventional Cutting	75	12	BLACK
POLYAMIDE BLK 4-0 X 75-28MM HC RB	NYL405428	4-0		28 mm 1/2 Circle Round Body	75	12	BLACK

Description	FG Code	Size		Needle Profile Shape, Type & Length	Suture Length (cm)	Pack Size	Suture Colour
FILAMIDE (POLYAMIDE SUTURE)							
POLYAMIDE BU 5-0 X 45-16MM CU CT	NYL5013486	5-0		16 mm 3/8 Circle Conventional Cutting	45	12	BLUE
POLYAMIDE BLK 5-0 X 45-16MM CU CT	NYL501854SL	5-0		16 mm 3/8 Circle Conventional Cutting	45	12	BLACK
POLYAMIDE BLK 5-0 X 70-12MM CU RC	NYL50 3317	5-0		12 mm 3/8 Circle Reverse Cutting	70	12	BLACK
POLYAMIDE BLK 5-0 X 45-15MM CU CT-3D	NYL5018951Z	5-0		15 mm 3/8 Circle Conventional Cutting	45	36	BLACK
POLYAMIDE BLK 5-0 X 75-16MM CU RC	NYL501616	5-0		16 mm 3/8 Circle Reverse Cutting	75	12	BLACK
POLYAMIDE BLK 5-0 X 90-13MM CU RB DN	NYL509013DN	5-0		13 mm 3/8 Circle Round Body Round Body	90	12	BLACK
POLYAMIDE BU 5-0 X 45-26MM CU CT	NYL5013487	5-0		26 mm 3/8 Circle Conventional Cutting	45	12	BLUE
POLYAMIDE BLK 5-0 X 45-16MM CU RC	NYL505668	5-0		16 mm 3/8 Circle Reverse Cutting	45	12	BLACK
POLYAMIDE BLK 5-0 X 45-8MM HC RC	NYL501765	5-0		8 mm 1/2 Circle Reverse Cutting	45	12	BLACK
POLYAMIDE BLK 6-0 X 70-10MM CU RC	NYL60 3320	6-0		10 mm 3/8 Circle Reverse Cutting	70	12	BLACK
POLYAMIDE BLK 6-0 X 45-16MM CU RC	NYL601614	6-0		16 mm 3/8 Circle Reverse Cutting	45	12	BLACK
POLYAMIDE BLK 6-0 X 45-13MM CU RC	NYL601851	6-0		13 mm 3/8 Circle Reverse Cutting	45	12	BLACK
POLYAMIDE BLK 6-0 X 90-16MM CU CT	NYL601615	6-0		16 mm 3/8 Circle Conventional Cutting	90	12	BLACK
POLYAMIDE BLK 6-0 X 70-13MM CU CT	NYL60 3351	6-0		13 mm 3/8 Circle Conventional Cutting	70	12	BLACK
POLYAMIDE BLK 6-0 X 70-26MM CU CT	NYL60 3331	6-0		26 mm 3/8 Circle Conventional Cutting	70	12	BLACK
POLYAMIDE BLK 6-0 X 75-9MM CU RB	NYL607509	6-0		9 mm 3/8 Circle Round Body	75	12	BLACK
POLYAMIDE BLK 6-0 X 45-15MM CU CT	NYL603615	6-0		15 mm 3/8 Circle Conventional Cutting	45	12	BLACK
POLYAMIDE BLK 6-0 X 45-11MM CU RC	NYL601610	6-0		11 mm 3/8 Circle Reverse Cutting	45	12	Black
POLYAMIDE BLK 7-0 X 60-9MM CU RB	NYL706009	7-0		9 mm 3/8 Circle Round Body	60	12	BLACK
POLYAMIDE BLK 8-0 X 90-09MM CU RB	NYL808170	8-0		9 mm 3/8 Circle Round Body	90	12	BLACK
POLYAMIDE BLK 8-0 X 38-6MM CU RB	NYL80 3708	8-0		6 mm 3/8 Circle Round Body	38	12	BLACK
POLYAMIDE BLK 8-0 X 38-6MM CU SM	NYL80 3322	8-0		6 mm 3/8 Circle Spatula Micropoint	38	12	BLACK

Description	FG Code	Size		Needle Profile Shape, Type & Length	Suture Length (cm)	Pack Size	Suture Colour
FILAMIDE (POLYAMIDE SUTURE)							
POLYAMIDE BLK 8-0 X 38-6MM CU SM DN	NYL80 3324DN	8-0		6 mm 3/8 Circle Spatula Micropoint Spatula Micropoint	38	12	BLACK
POLYAMIDE BLK 8-0 X 38-6MM CU RC	NYL80 3343	8-0		6 mm 3/8 Circle Reverse Cutting	38	12	BLACK
POLYAMIDE BLK 8-0 X 13-5MM CU RB	NYL802812	8-0		5 mm 3/8 Circle Round Body	13	12	BLACK
POLYAMIDE BLK 8-0 X 45-9MM CU RB DN	NYL808170DN	8-0		9 mm 3/8 Circle Round Body Round Body	45	12	BLACK
POLYAMIDE BLK 8-0 X 15-6.5MM CU RB	NYL80 2808	8-0		6.5 mm 3/8 Circle Round Body	15	12	BLACK
POLYAMIDE BLK 9-0 X 38-06MM CU RB SN	NYL90 3709	9-0		6 mm 3/8 Circle Round Body	38	12	BLACK
POLYAMIDE BLK 9-0 X 13-05MM CU RB SN	NYL903905	9-0		5 mm 3/8 Circle Round Body	13	12	BLACK
POLYAMIDE BLK 9-0 X 13-05MM CU RB	NYL90 2829	9-0		5 mm 3/8 Circle Round Body	13	12	BLACK
POLYAMIDE BLK 9-0 X 38-6MM CU SM DN	NYL90 3714DN	9-0		6 mm 3/8 Circle Spatula Micropoint Spatula Micropoint	38	12	BLACK
POLYAMIDE BLK 9-0 X 38-6MM CU SM	NYL90 3715	9-0		6 mm 3/8 Circle Spatula Micropoint	38	12	BLACK
POLYAMIDE BLK 9-0 X 25-6.5MM CU RB DN	NYL90 2800DN	9-0		6.5 mm 3/8 Circle Round Body Round Body	25	12	BLACK
POLYAMIDE BLK 9-0 X 30-6.5MM CU SM DN	NYL901717	9-0		6.5 mm 3/8 Circle Spatula Micropoint Spatula Micropoint	30	12	BLACK
POLYAMIDE BLK 1 X 150-48MM HC RB LOOP - 3D	NYL013348Z	1		48 mm 1/2 Circle Round Body	150	36	BLACK
POLYAMIDE BLK 3-0 X 90-19MM CU RC - 3D	NYL301622Z	3-0		19 mm 3/8 Circle Reverse Cutting	90	36	BLACK
POLYAMIDE BLK 4-0 X 90-13MM CU RC - 3D	NYL401613Z	4-0		13 mm 3/8 Circle Reverse Cutting	90	36	BLACK
POLYAMIDE BLK 5-0 X 75-16MM CU RC - 3D	NYL501616Z	5-0		16 mm 3/8 Circle Reverse Cutting	75	36	BLACK

Description	FG Code	Size	Needle Profile Shape, Type & Length	Suture Length (cm)	Pack Size	Suture Colour
FILAMIDE (POLYAMIDE SUTURE)						
POLYAMIDE BLK 0 X 76 (76X2)	NYS10 904	0	Non needled	76	12	BLACK
POLYAMIDE BLK 1 X 76 (76X2)	NYS01 905	1	Non needled	76	12	BLACK
POLYAMIDE BLK 1 X 750 (75X10)	NYS011075	1	Non needled	75	12	BLACK
POLYAMIDE BLK 1 X 35 (35X5)	NYS01 705	1	Non needled	35	12	BLACK
POLYAMIDE BLK 1 X 100 (100X2)	NYS01185	1	Non needled	100	12	BLACK
POLYAMIDE BLK 2 X 76 (76X2)	NYS02 906	2	Non needled	76	12	BLACK
POLYAMIDE BLK 2-0 X 76 (76X2)	NYS20 903	2-0	Non needled	76	12	BLACK
POLYAMIDE BLK 2-0 X 35 (35X5)	NYS20 703	2-0	Non needled	35	12	BLACK

FILAPROPTM

Polypropylene Suture

Indications

- Skin closure (intracutaneous, subcutaneous)
- Cardiovascular procedures
- Neurological procedures
- Ophthalmic procedures

Advantages

- Smooth passage through tissue
- Minimal acute inflammatory reaction in tissues
- Unsurpassed strength
- Optimal elasticity and elongation properties
- Available in Less Memory Tray (LMT)




Product Overview

Structure	: Monofilament
Colour	: Blue
Chemical composition	: Polypropylene
Sizes	: USP 10/0 (metric 0.2) - USP 2 (metric 5)
Sterilization	: Ethylene oxide

Description	FG Code	Size		Needle Profile Shape, Type & Length	Suture Length (cm)	Pack Size	Suture Colour
FILAPROP (POLYPROPYLENE SUTURE)							
POLYPROPYLENE BU 0 X 90-30MM CU RC	PPL102815	0		30 mm 3/8 Circle Reverse Cutting	90	12	BLUE
POLYPROPYLENE BU 0 X 70-30MM HC RB	PPL10 830	0		30 mm 1/2 Circle Round Body	70	12	BLUE
POLYPROPYLENE BU 0 X 90-30MM HC RB	PPL10 846	0		30 mm 1/2 Circle Round Body	90	12	BLUE
POLYPROPYLENE BU 0 X 90-30MM HC TC	PPL10 894	0		30 mm 1/2 Circle Taper Cut	90	12	BLUE
POLYPROPYLENE BU 0 X 75-40MM HC RB	PPL108434	0		40 mm 1/2 Circle Round Body	75	12	BLUE
POLYPROPYLENE BU 0 X 75-37MM HC RB	PPL108424	0		37 mm 1/2 Circle Round Body	75	12	BLUE
POLYPROPYLENE BU 0 X 75-30MM HC RB	PPL100830	0		30 mm 1/2 Circle Round Body	75	12	BLUE
POLYPROPYLENE BU 0 X 75-24MM CU CT	PPL108024	0		24 mm 3/8 Circle Conventional Cutting	75	12	BLUE
POLYPROPYLENE BU 0 X 75-26MM HC RB	PPL108026	0		26 mm 1/2 Circle Round Body	75	12	BLUE
POLYPROPYLENE BU 0 X 75-45MM HC RC	PPL108045	0		45 mm 1/2 Circle Reverse Cutting	75	12	BLUE
POLYPROPYLENE BU 0 X 75-45MM HC RB	PPL108075	0		45 mm 1/2 Circle Round Body	75	12	BLUE
POLYPROPYLENE BU 0 X 90-26MM HC RB DN	PPL10 824MIPH	0		26 mm 1/2 Circle Round Body Round Body	90	1	BLUE
POLYPROPYLENE BU 0 X 100-40MM HC RB	PPL108434XL	0		40 mm 1/2 Circle Round Body	100	12	BLUE
POLYPROPYLENE BU 0 X 100-31MM J TYPE TC	PPL10982J	0		31 mm J Shape Taper Cut	100	12	BLUE
POLYPROPYLENE BU 0 X 100-31MM HC RB	PPL108430	0		31 mm 1/2 Circle Round Body	100	12	BLUE
POLYPROPYLENE BU 0 X 90-40MM HC RC	PPL108470	0		40 mm 1/2 Circle Reverse Cutting	90	12	BLUE

Description	FG Code	Size		Needle Profile Shape, Type & Length	Suture Length (cm)	Pack Size	Suture Colour
FILAPROP (POLYPROPYLENE SUTURE)							
POLYPROPYLENE BU 0 X 70-30MM HC RB HB	PPL10 842	0		30 mm 1/2 Circle Round Body	70	12	BLUE
POLYPROPYLENE BU 0 X 100-26MM HC RB DN	PPL108026XDN	0		26 mm 1/2 Circle Round Body Round Body	100	12	BLUE
POLYPROPYLENE BU 0 X 75-37MM HC TC	PPL108037	0		37 mm 1/2 Circle Taper Cut	75	12	BLUE
POLYPROPYLENE BU 0 X 75-38MM CU CT	PPL108038	0		38 mm 3/8 Circle Conventional Cutting	75	12	BLUE
POLYPROPYLENE BU 1 X 90-45MM HC TC	PPL018450	1		45 mm 1/2 Circle Taper Cut	90	12	BLUE
POLYPROPYLENE BU 1 X 90-50MM HC RC	PPL018475	1		50 mm 1/2 Circle Reverse Cutting	90	12	BLUE
POLYPROPYLENE BU 1 X 150-48MM HC RB LOOP	PPL01 834	1		48 mm 1/2 Circle Round Body	150	12	BLUE
POLYPROPYLENE BU 1 X 70-30MM HC RB	PPL01 840	1		30 mm 1/2 Circle Round Body	70	12	BLUE
POLYPROPYLENE BU 1 X 70-40MM HC RB	PPL01 843	1		40 mm 1/2 Circle Round Body	70	12	BLUE
POLYPROPYLENE BU 1 X 100-45MM HC RC	PPL01 883	1		45 mm 1/2 Circle Reverse Cutting	100	12	BLUE
POLYPROPYLENE BU 1 X 90-90MM CU CT	PPL012989	1		90 mm 3/8 Circle Conventional Cutting	90	12	BLUE
POLYPROPYLENE BU 1 X 100-40MM HC RC	PPL01884	1		40 mm 1/2 Circle Reverse Cutting	100	12	BLUE
POLYPROPYLENE BU 1 X 100-36MM HC RB	PPL01842	1		36 mm 1/2 Circle Round Body	100	12	BLUE
POLYPROPYLENE BU 1 X 90-90MM CU RC	PPL018190	1		90 mm 3/8 Circle Reverse Cutting	90	12	BLUE
POLYPROPYLENE BU 1 X 75-30MM HCRB TRODN	PPL018130	1		30 mm 1/2 Circle Round Body TRO	75	12	BLUE
POLYPROPYLENE BU 1 X 100-40MM HC RB	PPL01 843DS	1		40 mm 1/2 Circle Round Body	100	12	BLUE

Description	FG Code	Size		Needle Profile Shape, Type & Length	Suture Length (cm)	Pack Size	Suture Colour
FILAPROP (POLYPROPYLENE SUTURE)							
POLYPROPYLENE BU 1 X 100-40MM HC RB BP	PPL018999	1		40 mm 1/2 Circle Blunt Point	100	12	BLUE
POLYPROPYLENE BU 1 X 75-40MM CU RC	PPL018140	1		40 mm 3/8 Circle Reverse Cutting	75	12	BLUE
POLYPROPYLENE BU 1 X 100-40MM HC RB EX	PPL018435	1		40 mm 1/2 Circle Round Body	100	12	BLUE
POLYPROPYLENE BU 2 X 75-45MM HC RC	PPL028245	2		45 mm 1/2 Circle Reverse Cutting	75	12	BLUE
POLYPROPYLENE BU 2 X 75-45MM HC RB	PPL028275	2		45 mm 1/2 Circle Round Body	75	12	BLUE
POLYPROPYLENE BU 10-0 X 13-06MM CU RB DN	PPL1002795	10-0		6 mm 3/8 Circle Round Body Round Body	13	12	BLUE
POLYPROPYLENE BU 10-0 X 23-16MM ST CT DN	PPL1001713DN	10-0		16 mm Straight Conventional Cutting Conventional Cutting	23	12	BLUE
POLYPROPYLENE BU 10-0 X 13-5MM CU RB DN	PPL1002794	10-0		5 mm 3/8 Circle Round Body Round Body	13	12	BLUE
POLYPROPYLENE BU 10-0 X 90-6.3MM CU SM	PPL1001777	10-0		6.3 mm 3/8 Circle Spatula Micropoint	90	12	BLUE
POLYPROPYLENE BU 10-0 X 45-06MM CU SM	PPL1001713E	10-0		6 mm 3/8 Circle Spatula Micropoint	45	12	BLUE
POLYPROPYLENE BU 2-0 X 90-40MM CU RC	PPL208689	2-0		40 mm 3/8 Circle Reverse Cutting	90	12	BLUE
POLYPROPYLENE BU 2-0 X 90-90MM CU RC	PPL200987	2-0		90 mm 3/8 Circle Reverse Cutting	90	12	BLUE
POLYPROPYLENE BU 2-0 X 90-40MM HC TC	PPL208852	2-0		40 mm 1/2 Circle Taper Cut	90	12	BLUE
POLYPROPYLENE BU 2-0 X 90-30MM HC RB	PPL20 844	2-0		30 mm 1/2 Circle Round Body	90	12	BLUE
POLYPROPYLENE BU 2-0 X 90-26MM HC TC DN	PPL20 8977	2-0		26 mm 1/2 Circle Taper Cut Taper Cut	90	12	BLUE
POLYPROPYLENE BU 2-0 X 90-17.5MM HCTC DN	PPL20 8937	2-0		17.5 mm 1/2 Circle Taper Cut Taper Cut	90	12	BLUE
POLYPROPYLENE BU 2-0 X 90-26MM HC RB	PPL208523	2-0		26 mm 1/2 Circle Round Body	90	12	BLUE

Description	FG Code	Size		Needle Profile Shape, Type & Length	Suture Length (cm)	Pack Size	Suture Colour
FILAPROP (POLYPROPYLENE SUTURE)							
POLYPROPYLENE BU 2-0 X 75-26MM HC RB	PPL208833	2-0		26 mm 1/2 Circle Round Body	75	12	BLUE
POLYPROPYLENE BU 2-0 X 75-24MM CU RC	PPL208685	2-0		24 mm 3/8 Circle Reverse Cutting	75	12	BLUE
POLYPROPYLENE BU 2-0 X 75-26MM HC CT	PPL208978	2-0		26 mm 1/2 Circle Conventional Cutting	75	12	BLUE
POLYPROPYLENE BU 2-0 X 90-60MM ST CT	PPL200623	2-0		60 mm Straight Conventional Cutting	90	12	BLUE
POLYPROPYLENE BU 2-0 X 100-24MM CU RC	PPL208685XL	2-0		24 mm 3/8 Circle Reverse Cutting	100	12	BLUE
POLYPROPYLENE BU 2-0 X 75-26MM HC RB DN	PPL208834	2-0		26 mm 1/2 Circle Round Body Round Body	75	12	BLUE
POLYPROPYLENE BU 2-0 X 90-36MM HC RB	PPL208236	2-0		36 mm 1/2 Circle Round Body	90	12	BLUE
POLYPROPYLENE BU 2-0 X 90-70MM ST CT	PPL200753	2-0		70 mm Straight Conventional Cutting	90	12	BLUE
POLYPROPYLENE BU 2-0 X 90-31MM HC RB VB	PPL20295	2-0		31 mm 1/2 Circle Round Body	90	12	BLUE
POLYPROPYLENE BU 2-0 X 75-45MM HC RC	PPL208245	2-0		45 mm 1/2 Circle Reverse Cutting	75	12	BLUE
POLYPROPYLENE BU 2-0 X 75-45MM HC RB	PPL208275	2-0		45 mm 1/2 Circle Round Body	75	12	BLUE
POLYPROPYLENE BU 2-0 X 70-30MM CU RB	PPL208230	2-0		30 mm 3/8 Circle Round Body	70	12	BLUE
POLYPROPYLENE BU 2-0 X 70-25MM HC TC	PPL20 807	2-0		25 mm 1/2 Circle Taper Cut	70	12	BLUE
POLYPROPYLENE BU 2-0 X 75-20MM HC RB	PPL2085200	2-0		20 mm 1/2 Circle Round Body	75	12	BLUE
POLYPROPYLENE BU 2-0 X 75-30MM CU RC	PPL207697	2-0		30 mm 3/8 Circle Reverse Cutting	75	12	BLUE
POLYPROPYLENE BU 2-0 X 75-60MM ST CT	PPL200631	2-0		60 mm Straight Conventional Cutting	75	12	BLUE
POLYPROPYLENE BU 2-0 X 75-31MM HC RB	PPL208231	2-0		31 mm 1/2 Circle Round Body	75	12	BLUE
POLYPROPYLENE BU 2-0 X 75-30MM HC RB	PPL20844SL	2-0		30 mm 1/2 Circle Round Body	75	12	BLUE

Description	FG Code	Size		Needle Profile Shape, Type & Length	Suture Length (cm)	Pack Size	Suture Colour
FILAPROP (POLYPROPYLENE SUTURE)							
POLYPROPYLENE BU 2-0 X 90-31MM HC RB DN	PPL208526	2-0		31 mm 1/2 Circle Round Body Round Body	90	12	BLUE
POLYPROPYLENE BU 2-0 X 90-31MM HC TC DN	PPL208850	2-0		31 mm 1/2 Circle Taper Cut Taper Cut	90	12	BLUE
POLYPROPYLENE BU 2-0 X 45-26MM CU RC	PPL208026	2-0		26 mm 3/8 Circle Reverse Cutting	45	12	BLUE
POLYPROPYLENE BU 2-0 X 90-40MM HC TC DN	PPL208852DN	2-0		40 mm 1/2 Circle Taper Cut Taper Cut	90	12	BLUE
POLYPROPYLENE BU 2-0 X 100-75MM ST CT	PPL200752	2-0		75 mm Straight Conventional Cutting	100	12	BLUE
POLYPROPYLENE BU 2-0 X 90-55MM ST CT	PPL208571	2-0		55 mm Straight Conventional Cutting	90	12	BLUE
POLYPROPYLENE BU 2-0 X 90-40MM HC RB DN	PPL20 2040DN	2-0		40 mm 1/2 Circle Round Body Round Body	90	12	BLUE
POLYPROPYLENE BU 2-0 X 70-30MM HC RB	PPL20 841	2-0		30 mm 1/2 Circle Round Body	70	12	BLUE
POLYPROPYLENE BU 2-0 X 90-90MM CU CT	PPL200988	2-0		90 mm 3/8 Circle Conventional Cutting	90	12	BLUE
POLYPROPYLENE BU 2-0 X 90-30MM HC RB-3D	PPL20844Z	2-0		30 mm 1/2 Circle Round Body	90	36	BLUE
POLYPROPYLENE BU 2-0 X 90-26MM HCTCDN-3D	PPL208977Z	2-0		26 mm 1/2 Circle Taper Cut Taper Cut	90	36	BLUE
POLYPROPYLENE BU 2-0 X90-31MM HCRBVB-3D	PPL20295Z	2-0		31 mm 1/2 Circle Round Body	90	36	BLUE
POLYPROPYLENE BU 2-0 X 90-26MM HC TC	PPL208926	2-0		26 mm 1/2 Circle Taper Cut	90	12	BLUE
POLYPROPYLENE BU 2-0 X 90-30MM HC RC	PPL208931	2-0		30 mm 1/2 Circle Reverse Cutting	90	12	BLUE

Description	FG Code	Size		Needle Profile Shape, Type & Length	Suture Length (cm)	Pack Size	Suture Colour
FILAPROP (POLYPROPYLENE SUTURE)							
POLYPROPYLENE BU 2-0 X 75-70MM ST RB DN	PPL208400	2-0		70 mm Straight Round Body Round Body	75	12	BLUE
POLYPROPYLENE BU 2-0 X 90-30MM HC RB DN	PPL20844DN	2-0		30 mm 1/2 Circle Round Body Round Body	90	12	BLUE
POLYPROPYLENE BU 2-0 X 90-17MM HC TC	PPL208938	2-0		17 mm 1/2 Circle Taper Cut	90	12	BLUE
POLYPROPYLENE BU 2-0 X 100-40MM CU CT	PPL208690	2-0		40 mm 3/8 Circle Conventional Cutting	100	12	BLUE
POLYPROPYLENE BU 3-0 X 70-22MM CU CT	PPL30 018	3-0		22 mm 3/8 Circle Conventional Cutting	70	12	BLUE
POLYPROPYLENE BU 3-0 X 90-25MM CU CT	PPL30 800	3-0		25 mm 3/8 Circle Conventional Cutting	90	12	BLUE
POLYPROPYLENE BU 3-0 X 90-25MM HC RB	PPL30 825	3-0		25 mm 1/2 Circle Round Body	90	12	BLUE
POLYPROPYLENE BU 3-0 X 70-25MM HC RB	PPL30 838	3-0		25 mm 1/2 Circle Round Body	70	12	BLUE
POLYPROPYLENE BU 3-0 X 70-20MM HC RB	PPL30 887	3-0		20 mm 1/2 Circle Round Body	70	12	BLUE
POLYPROPYLENE BU 3-0 X 70-60MM ST CT	PPL30 888	3-0		60 mm Straight Conventional Cutting	70	12	BLUE
POLYPROPYLENE BU 3-0 X 90-26MM HC RB DN	PPL30 8522	3-0		26 mm 1/2 Circle Round Body Round Body	90	12	BLUE
POLYPROPYLENE BU 3-0 X 90-17.5MM HCTC DN	PPL30 8936	3-0		17.5 mm 1/2 Circle Taper Cut Taper Cut	90	12	BLUE
POLYPROPYLENE BU 3-0 X 90-17MM HC RB	PPL308558	3-0		17 mm 1/2 Circle Round Body	90	12	BLUE
POLYPROPYLENE BU 3-0 X 90-16MM HC TC DN	PPL308667	3-0		16 mm 1/2 Circle Taper Cut Taper Cut	90	12	BLUE
POLYPROPYLENE BU 3-0 X 75-22MM HC RB	PPL308770	3-0		22 mm 1/2 Circle Round Body	75	12	BLUE
POLYPROPYLENE BU 3-0 X 90-26MM HC TC DN	PPL308976	3-0		26 mm 1/2 Circle Taper Cut Taper Cut	90	12	BLUE
POLYPROPYLENE BU 3-0 X 45-19MM CU RC	PPL308665	3-0		19 mm 3/8 Circle Reverse Cutting	45	12	BLUE

Description	FG Code	Size		Needle Profile Shape, Type & Length	Suture Length (cm)	Pack Size	Suture Colour
FILAPROP (POLYPROPYLENE SUTURE)							
POLYPROPYLENE BU 3-0 X 45-24MM CU RC	PPL308684	3-0		24 mm 3/8 Circle Reverse Cutting	45	12	BLUE
POLYPROPYLENE BU 3-0 X 75-26MM HC RB	PPL308832	3-0		26 mm 1/2 Circle Round Body	75	12	BLUE
POLYPROPYLENE BU 3-0 X 75-19MM CU RC	PPL308687	3-0		19 mm 3/8 Circle Reverse Cutting	75	12	BLUE
POLYPROPYLENE BU 3-0 X 90-16MM HC RB	PPL308559	3-0		16 mm 1/2 Circle Round Body	90	12	BLUE
POLYPROPYLENE BU 3-0 X 90-25MM HC TC	PPL308977	3-0		25 mm 1/2 Circle Taper Cut	90	12	BLUE
POLYPROPYLENE BU 3-0 X 90-30MM HC RB	PPL308525	3-0		30 mm 1/2 Circle Round Body	90	12	BLUE
POLYPROPYLENE BU 3-0 X 90-16MM CU RC	PPL308549	3-0		16 mm 3/8 Circle Reverse Cutting	90	12	BLUE
POLYPROPYLENE BU 3-0 X 90-26MM CU RC	PPL308685	3-0		26 mm 3/8 Circle Reverse Cutting	90	12	BLUE
POLYPROPYLENE BU 3-0 X 90-17MM HC TC	PPL308317	3-0		17 mm 1/2 Circle Taper Cut	90	12	BLUE
POLYPROPYLENE BU 3-0 X 90-36MM HC RB	PPL308336	3-0		36 mm 1/2 Circle Round Body	90	12	BLUE
POLYPROPYLENE BU 3-0 X 75-22MM HC RB DN	PPL308771	3-0		22 mm 1/2 Circle Round Body Round Body	75	12	BLUE
POLYPROPYLENE BU 3-0 X 75-60MM ST CT	PPL3013413	3-0		60 mm Straight Conventional Cutting	75	12	BLUE
POLYPROPYLENE BU 3-0 X 75-16MM CU CT	PPL3013453	3-0		16 mm 3/8 Circle Conventional Cutting	75	12	BLUE
POLYPROPYLENE BU 3-0 X 120-26MM HC RB	PPL308534	3-0		26 mm 1/2 Circle Round Body	120	12	BLUE
POLYPROPYLENE BU 3-0 X 75-24MM CU RC	PPL307694	3-0		24 mm 3/8 Circle Reverse Cutting	75	12	BLUE
POLYPROPYLENE BU 3-0 X 90-36MM HC RB DN	PPL307499	3-0		36 mm 1/2 Circle Round Body Round Body	90	12	BLUE
POLYPROPYLENE BU 3-0 X 120-26MM HC TC DN	PPL307778	3-0		26 mm 1/2 Circle Taper Cut Taper Cut	120	12	BLUE
POLYPROPYLENE BU 3-0 X 45-26MM CU RC	PPL308025	3-0		26 mm 3/8 Circle Reverse Cutting	45	12	BLUE


Description	FG Code	Size		Needle Profile Shape, Type & Length	Suture Length (cm)	Pack Size	Suture Colour
FILAPROP (POLYPROPYLENE SUTURE)							
POLYPROPYLENE BU 3-0 X 90-31MM HC RB DN	PPL308525DN	3-0		31 mm 1/2 Circle Round Body Round Body	90	12	BLUE
POLYPROPYLENE BU 3-0 X 45-22MM HC CT	PPL308731	3-0		22 mm 1/2 Circle Conventional Cutting	45	12	BLUE
POLYPROPYLENE BU 3-0 X 45-19MM CU CT	PPL308878	3-0		19 mm 3/8 Circle Conventional Cutting	45	12	BLUE
POLYPROPYLENE BU 3-0 X 90-26MM HC RB DNV	PPL308355	3-0		26 mm 1/2 Circle Round Body Round Body	90	12	BLUE
POLYPROPYLENE BU 3-0 X 90-31MM HC TC DN	PPL308849DN	3-0		31 mm 1/2 Circle Taper Cut Taper Cut	90	12	BLUE
POLYPROPYLENE BU 3-0 X 90-40MM HC TC DN	PPL308851	3-0		40 mm 1/2 Circle Taper Cut Taper Cut	90	12	BLUE
POLYPROPYLENE BU 3-0 X90-26MM HCRB TRODN	PPL308833	3-0		26 mm 1/2 Circle Round Body TRO	90	12	BLUE
POLYPROPYLENE BU 3-0 X 90-20MM HC RBTC DN	PPL309557	3-0		20 mm 1/2 Circle Round Body Taper Cut	90	12	BLUE
POLYPROPYLENE BU 3-0 X 90-17MM HC RB DN	PPL308558DN	3-0		17 mm 1/2 Circle Round Body Round Body	90	12	BLUE
POLYPROPYLENE BU 3-0 X 90-16MM CU RC-3D	PPL308549Z	3-0		16 mm 3/8 Circle Reverse Cutting	90	36	BLUE
POLYPROPYLENE BU 3-0 X 90-60MM ST CT	PPL300621	3-0		60 mm Straight Conventional Cutting	90	12	BLUE
POLYPROPYLENE BU 3-0 X 45-26MM CU CT EX	PPL308021	3-0		26 mm 3/8 Circle Conventional Cutting	45	12	BLUE
POLYPROPYLENE BU 3-0 X 75-36MM CU CT	PPL300626	3-0		36 mm 3/8 Circle Conventional Cutting	75	12	BLUE
POLYPROPYLENE BU 4-0 X 90-19MM CU RC	PPL408015	4-0		19 mm 3/8 Circle Reverse Cutting	90	12	BLUE
POLYPROPYLENE BU 4-0 X 90-16MM CU CT	PPL408875	4-0		16 mm 3/8 Circle Conventional Cutting	90	12	BLUE

Description	FG Code	Size		Needle Profile Shape, Type & Length	Suture Length (cm)	Pack Size	Suture Colour
FILAPROP (POLYPROPYLENE SUTURE)							
POLYPROPYLENE BU 4-0 X 90-25MM CU CT	PPL408879	4-0		25 mm 3/8 Circle Conventional Cutting	90	12	BLUE
POLYPROPYLENE BU 4-0 X 70-16MM HC RB	PPL40 849	4-0		16 mm 1/2 Circle Round Body	70	12	BLUE
POLYPROPYLENE BU 4-0 X 70-16MM CU CT	PPL40 870	4-0		16 mm 3/8 Circle Conventional Cutting	70	12	BLUE
POLYPROPYLENE BU 4-0 X90-25MM HC TC DNC1	PPL40 867	4-0		25 mm 1/2 Circle Taper Cut Taper Cut	90	12	BLUE
POLYPROPYLENE BU 4-0 X 90-17.5MM HCTC DN	PPL40 8935	4-0		17.5 mm 1/2 Circle Taper Cut Taper Cut	90	12	BLUE
POLYPROPYLENE BU 4-0 X 90-26MM HC RB DN	PPL40 8521	4-0		26 mm 1/2 Circle Round Body Round Body	90	12	BLUE
POLYPROPYLENE BU 4-0 X 90-16MM HC RB DN	PPL40 8557	4-0		16 mm 1/2 Circle Round Body Round Body	90	12	BLUE
POLYPROPYLENE BU 4-0 X 90-26MM HC RB	PPL408355	4-0		26 mm 1/2 Circle Round Body	90	12	BLUE
POLYPROPYLENE BU 4-0 X 90-17MM HC RB	PPL408558	4-0		17 mm 1/2 Circle Round Body	90	12	BLUE
POLYPROPYLENE BU 4-0 X 75-13MM HC CT DN	PPL408666	4-0		13 mm 1/2 Circle Conventional Cutting Conventional Cutting	75	12	BLUE
POLYPROPYLENE BU 4-0 X 90-20MM HC RB	PPL408761	4-0		20 mm 1/2 Circle Round Body	90	12	BLUE
POLYPROPYLENE BU 4-0 X 90-20MM HC CT DN	PPL408840	4-0		20 mm 1/2 Circle Conventional Cutting Conventional Cutting	90	12	BLUE
POLYPROPYLENE BU 4-0 X 90-26MM HC CT DN	PPL408846	4-0		26 mm 1/2 Circle Conventional Cutting Conventional Cutting	90	12	BLUE
POLYPROPYLENE BU 4-0 X 45-19MM CU RC	PPL408682	4-0		19 mm 3/8 Circle Reverse Cutting	45	12	BLUE
POLYPROPYLENE BU 4-0 X 45-20MM CU RC	PPL408683	4-0		20 mm 3/8 Circle Reverse Cutting	45	12	BLUE
POLYPROPYLENE BU 4-0 X 90-13MM HC RB DN	PPL408204	4-0		13 mm 1/2 Circle Round Body Round Body	90	12	BLUE
POLYPROPYLENE BU 4-0 X 90-18MM CU RB DN	PPL408861	4-0		18 mm 3/8 Circle Round Body Round Body	90	12	BLUE

Description	FG Code	Size		Needle Profile Shape, Type & Length	Suture Length (cm)	Pack Size	Suture Colour
FILAPROP (POLYPROPYLENE SUTURE)							
POLYPROPYLENE BU 4-0 X 90-24MM CU RC	PPL408684	4-0		24 mm 3/8 Circle Reverse Cutting	90	12	BLUE
POLYPROPYLENE BU 4-0 X 70-19MM CU RC	PPL408682XL	4-0		19 mm 3/8 Circle Reverse Cutting	70	12	BLUE
POLYPROPYLENE BU 4-0 X90-25MM HC TC	PPL408425	4-0		25 mm 1/2 Circle Taper Cut	90	12	BLUE
POLYPROPYLENE BU 4-0 X 90-16MM CU RC	PPL408416	4-0		16 mm 3/8 Circle Reverse Cutting	90	12	BLUE
POLYPROPYLENE BU 4-0 X 75-22MM HC RB DN	PPL40 8526DN	4-0		22 mm 1/2 Circle Round Body Round Body	75	12	BLUE
POLYPROPYLENE BU 4-0 X 70-17MM HC TC	PPL40 835	4-0		17 mm 1/2 Circle Taper Cut	70	12	BLUE
POLYPROPYLENE BU 4-0 X 45-16MM CU RC	PPL408884	4-0		16 mm 3/8 Circle Reverse Cutting	45	12	BLUE
POLYPROPYLENE BU 4-0 X 45-16MM CU CT	PPL408007	4-0		16 mm 3/8 Circle Conventional Cutting	45	12	BLUE
POLYPROPYLENE BU 4-0 X 45-26MM CU CT	PPL408020SL	4-0		26 mm 3/8 Circle Conventional Cutting	45	12	BLUE
POLYPROPYLENE BU 4-0 X90-26MM HC TC DNC2	PPL40 8845	4-0		26 mm 1/2 Circle Taper Cut Taper Cut	90	12	BLUE
POLYPROPYLENE BU 4-0 X 45-19MM CU CT	PPL408011	4-0		19 mm 3/8 Circle Conventional Cutting	45	12	BLUE
POLYPROPYLENE BU 4-0 X 90-17MM HC TC	PPL40 857	4-0		17 mm 1/2 Circle Taper Cut	90	12	BLUE
POLYPROPYLENE BU 4-0 X 45-12MM CU RC	PPL40 879	4-0		12 mm 3/8 Circle Reverse Cutting	45	12	BLUE
POLYPROPYLENE BU 4-0 X 90-16MM HC TC DN	PPL408417	4-0		16 mm 1/2 Circle Taper Cut Taper Cut	90	12	BLUE
POLYPROPYLENE BU 4-0 X 90-17MM HC TC EX	PPL400835	4-0		17 mm 1/2 Circle Taper Cut	90	12	BLUE
POLYPROPYLENE BU 4-0 X 70-17MM HC TC C1	PPL40 899	4-0		17 mm 1/2 Circle Taper Cut	70	12	BLUE
POLYPROPYLENE BU 4-0 X 90-26MM HC RB EX	PPL408831	4-0		26 mm 1/2 Circle Round Body	90	12	BLUE
POLYPROPYLENE BU 4-0 X 90-20MM HC TC DN	PPL408841DN	4-0		20 mm 1/2 Circle Taper Cut Taper Cut	90	12	BLUE
POLYPROPYLENE BU 4-0 X 70-16MM CU CT-3D	PPL40870Z	4-0		16 mm 3/8 Circle Conventional Cutting	70	36	BLUE

Description	FG Code	Size		Needle Profile Shape, Type & Length	Suture Length (cm)	Pack Size	Suture Colour
FILAPROP (POLYPROPYLENE SUTURE)							
POLYPROPYLENE BU 4-0 X 90-16MM HCRBDN-3D	PPL408557Z	4-0		16 mm 1/2 Circle Round Body Round Body	90	36	BLUE
POLYPROPYLENE BU 4-0 X 45-20MM CU RC-3D	PPL408683Z	4-0		20 mm 3/8 Circle Reverse Cutting	45	36	BLUE
POLYPROPYLENE BU 4-0 X90-26MM HCRB EX-3D	PPL408831Z	4-0		26 mm 1/2 Circle Round Body	90	36	BLUE
POLYPROPYLENE BU 4-0 X 90-20MM HC RB DN	PPL408761DN	4-0		20 mm 1/2 Circle Round Body Round Body	90	12	BLUE
POLYPROPYLENE BU 4-0 X 90-22MM HC RB	PPL408762	4-0		22 mm 1/2 Circle Round Body	90	12	BLUE
POLYPROPYLENE BU 4-0 X 60-13MM HC RB DN	PPL40 8204SL	4-0		13 mm 1/2 Circle Round Body Round Body	60	12	BLUE
POLYPROPYLENE BU 4-0 X 90-17MM HC RB DN	PPL408558DN	4-0		17 mm 1/2 Circle Round Body Round Body	90	12	BLUE
POLYPROPYLENE BU 4-0 X 75-13MM HC TC DN	PPL408665	4-0		13 mm 1/2 Circle Taper Cut Taper Cut	75	12	BLUE
POLYPROPYLENE BU 4-0 X 90-17MM HC RBDNVB	PPL408329	4-0		17 mm 1/2 Circle Round Body Round Body	90	12	BLUE
POLYPROPYLENE BU 4-0 X 45-13MM CU RC	PPL405742	4-0		13 mm 3/8 Circle Reverse Cutting	45	12	BLUE
POLYPROPYLENE BU 5-0 X 90-19MM CU CT	PPL508010	5-0		19 mm 3/8 Circle Conventional Cutting	90	12	BLUE
POLYPROPYLENE BU 5-0 X 90-16MM CU RC	PPL508882	5-0		16 mm 3/8 Circle Reverse Cutting	90	12	BLUE
POLYPROPYLENE BU 5-0 X 90-19MM CU RC	PPL508014	5-0		19 mm 3/8 Circle Reverse Cutting	90	12	BLUE
POLYPROPYLENE BU 5-0 X 70-16MM CU RB	PPL50 881	5-0		16 mm 3/8 Circle Round Body	70	12	BLUE
POLYPROPYLENE BU 5-0 X 70-16MM CU RB DN	PPL50 882DN	5-0		16 mm 3/8 Circle Round Body Round Body	70	12	BLUE
POLYPROPYLENE BU 5-0 X 90-12MM CU RC	PPL50 889	5-0		12 mm 3/8 Circle Reverse Cutting	90	12	BLUE
POLYPROPYLENE BU 5-0 X 90-16MM HCRB DNC2	PPL50 856	5-0		16 mm 1/2 Circle Round Body Round Body	90	12	BLUE
POLYPROPYLENE BU 5-0 X 75-13MM CU RB DN	PPL50 8721	5-0		13 mm 3/8 Circle Round Body Round Body	75	12	BLUE

Description	FG Code	Size		Needle Profile Shape, Type & Length	Suture Length (cm)	Pack Size	Suture Colour
FILAPROP (POLYPROPYLENE SUTURE)							
POLYPROPYLENE BU 5-0 X 60-10MM HC RB DN	PPL50 8803	5-0		10 mm 1/2 Circle Round Body Round Body	60	12	BLUE
POLYPROPYLENE BU 5-0 X 60-13MM HC RB DN	PPL50 8205	5-0		13 mm 1/2 Circle Round Body Round Body	60	12	BLUE
POLYPROPYLENE BU 5-0 X 75-13MM HC RB DN	PPL50 8710	5-0		13 mm 1/2 Circle Round Body Round Body	75	12	BLUE
POLYPROPYLENE BU 5-0 X 75-11MM HC RB DVB	PPL50 8802BX	5-0		11 mm 1/2 Circle Round Body Round Body	75	12	BLUE
POLYPROPYLENE BU 5-0 X 60-09MM CU TC DN	PPL508663	5-0		9 mm 3/8 Circle Taper Cut Taper Cut	60	12	BLUE
POLYPROPYLENE BU 5-0 X 60-13MM HC TC DN	PPL508664	5-0		13 mm 1/2 Circle Taper Cut Taper Cut	60	12	BLUE
POLYPROPYLENE BU 5-0 X 75-13MM CU TC DN	PPL508720	5-0		13 mm 3/8 Circle Taper Cut Taper Cut	75	12	BLUE
POLYPROPYLENE BU 5-0 X 60-11MM HC TC DN	PPL508804	5-0		11 mm 1/2 Circle Taper Cut Taper Cut	60	12	BLUE
POLYPROPYLENE BU 5-0 X 75-13MM HC TC DN	PPL508816	5-0		13 mm 1/2 Circle Taper Cut Taper Cut	75	12	BLUE
POLYPROPYLENE BU 5-0 X 75-13MM CU RB	PPL508890	5-0		13 mm 3/8 Circle Round Body	75	12	BLUE
POLYPROPYLENE BU 5-0 X 90-17MM HC RB	PPL508557	5-0		17 mm 1/2 Circle Round Body	90	12	BLUE
POLYPROPYLENE BU 5-0 X 45-13MM CU RC	PPL508698	5-0		13 mm 3/8 Circle Reverse Cutting	45	12	BLUE
POLYPROPYLENE BU 5-0 X 90-17MM HC RB DN	PPL508756	5-0		17 mm 1/2 Circle Round Body Round Body	90	12	BLUE
POLYPROPYLENE BU 5-0 X 90-13MM CU TC	PPL508722	5-0		13 mm 3/8 Circle Taper Cut	90	12	BLUE
POLYPROPYLENE BU 5-0 X 90-16MM CU CT	PPL508006	5-0		16 mm 3/8 Circle Conventional Cutting	90	12	BLUE
POLYPROPYLENE BU 5-0 X 75-26MM HC RB DN	PPL508526	5-0		26 mm 1/2 Circle Round Body Round Body	75	12	BLUE
POLYPROPYLENE BU 5-0 X 75-17MM HC TC DN	PPL508517	5-0		17 mm 1/2 Circle Taper Cut Taper Cut	75	12	BLUE

Description	FG Code	Size		Needle Profile Shape, Type & Length	Suture Length (cm)	Pack Size	Suture Colour
FILAPROP (POLYPROPYLENE SUTURE)							
POLYPROPYLENE BU 5-0 X 70-15MM CU CT	PPL50 826	5-0		15 mm 3/8 Circle Conventional Cutting	70	12	BLUE
POLYPROPYLENE BU 5-0 X 90-12MM CU TC DN	PPL505969	5-0		12 mm 3/8 Circle Taper Cut Taper Cut	90	12	BLUE
POLYPROPYLENE BU 5-0 X 45-19MM CU RC	PPL508686	5-0		19 mm 3/8 Circle Reverse Cutting	45	12	BLUE
POLYPROPYLENE BU 5-0 X 45-13MM CU RB DN	PPL508717	5-0		13 mm 3/8 Circle Round Body Round Body	45	12	BLUE
POLYPROPYLENE BU 5-0 X 90-16MM HC TC DN	PPL508830	5-0		16 mm 1/2 Circle Taper Cut Taper Cut	90	12	BLUE
POLYPROPYLENE BU 5-0 X 45-16MM CU CT	PPL508006SL	5-0		16 mm 3/8 Circle Conventional Cutting	45	12	BLUE
POLYPROPYLENE BU 5-0 X 45-16MM CU RC	PPL508882SL	5-0		16 mm 3/8 Circle Reverse Cutting	45	12	BLUE
POLYPROPYLENE BU 5-0 X 90-16MM HCRB DNC1	PPL50 8556	5-0		16 mm 1/2 Circle Round Body Round Body	90	12	BLUE
POLYPROPYLENE BU 5-0 X 70-12MM CU RB DN	PPL50 827DN	5-0		12 mm 3/8 Circle Round Body Round Body	70	12	BLUE
POLYPROPYLENE BU 5-0 X 70-13MM CU RC	PPL50 862	5-0		13 mm 3/8 Circle Reverse Cutting	70	12	BLUE
POLYPROPYLENE BU 5-0 X75-13MM CU RB DNVB	PPL508321	5-0		13 mm 3/8 Circle Round Body Round Body	75	12	BLUE
POLYPROPYLENE BU 5-0 X90-16MM HC RB TC DN	PPL508757	5-0		16 mm 1/2 Circle Round Body Taper Cut	90	12	BLUE
POLYPROPYLENE BU 5-0 X 75-11MM HC RB DN	PPL50 8802	5-0		11 mm 1/2 Circle Round Body Round Body	75	12	BLUE
POLYPROPYLENE BU 5-0 X 75-13MM HC RB DNV	PPL50 8710BL	5-0		13 mm 1/2 Circle Round Body Round Body	75	12	BLUE
POLYPROPYLENE BU 5-0 X 70-16MM CU RB-3D	PPL50881Z	5-0		16 mm 3/8 Circle Round Body	70	36	BLUE
POLYPROPYLENE BU 5-0X90-16MM HCRBDNC1-3D	PPL508556Z	5-0		16 mm 1/2 Circle Round Body Round Body	90	36	BLUE
POLYPROPYLENE BU 5-0 X 75-13MM HC RB	PPL508711	5-0		13 mm 1/2 Circle Round Body	75	12	BLUE

Description	FG Code	Size		Needle Profile Shape, Type & Length	Suture Length (cm)	Pack Size	Suture Colour
FILAPROP (POLYPROPYLENE SUTURE)							
POLYPROPYLENE BU 5-0 X 90-13MM CUTC DNVB	PPL508720BL	5-0		13 mm 3/8 Circle Taper Cut Taper Cut	90	12	BLUE
POLYPROPYLENE BU 6-0 X 70-15MM CU CT	PPL60 823	6-0		15 mm 3/8 Circle Conventional Cutting	70	12	BLUE
POLYPROPYLENE BU 6-0 X 45-12MM CU RC	PPL60 878	6-0		12 mm 3/8 Circle Reverse Cutting	45	12	BLUE
POLYPROPYLENE BU 6-0 X 70-10MM CU RB DN	PPL60 895DN	6-0		10 mm 3/8 Circle Round Body Round Body	70	12	BLUE
POLYPROPYLENE BU 6-0 X 60-13MM CU RB DN	PPL60 818	6-0		13 mm 3/8 Circle Round Body Round Body	60	12	BLUE
POLYPROPYLENE BU 6-0 X 70-13MM CU RB DN	PPL60 829	6-0		13 mm 3/8 Circle Round Body Round Body	70	12	BLUE
POLYPROPYLENE BU 6-0 X 60-9MM CU RB DNC1	PPL60 872	6-0		9 mm 3/8 Circle Round Body Round Body	60	12	BLUE
POLYPROPYLENE BU 6-0 X 70-8MM CU RB DN	PPL60 893	6-0		8 mm 3/8 Circle Round Body Round Body	70	12	BLUE
POLYPROPYLENE BU 6-0 X 60-13MM CURB DNTN	PPL60 8726	6-0		13 mm 3/8 Circle Round Body Round Body	60	12	BLUE
POLYPROPYLENE BU 6-0 X 60-9MM CU RB DNC2	PPL60 8805	6-0		9 mm 3/8 Circle Round Body Round Body	60	12	BLUE
POLYPROPYLENE BU 6-0 X 60-13MM CURB DNVB	PPL60 8307BL	6-0		13 mm 3/8 Circle Round Body Round Body	60	12	BLUE
POLYPROPYLENE BU 6-0 X 60-11MM CU RB DN	PPL60 8597	6-0		11 mm 3/8 Circle Round Body Round Body	60	12	BLUE
POLYPROPYLENE BU 6-0 X 60-13MM CURB DNTH	PPL60 8718	6-0		13 mm 3/8 Circle Round Body Round Body	60	12	BLUE
POLYPROPYLENE BU 6-0 X 60-13MM CU RB DVB	PPL60 8726BL	6-0		13 mm 3/8 Circle Round Body Round Body	60	12	BLUE
POLYPROPYLENE BU 6-0 X 75-13MM CU RB DVB	PPL60 8726BX	6-0		13 mm 3/8 Circle Round Body Round Body	75	12	BLUE
POLYPROPYLENE BU 6-0 X 75-13MM CU RB DN	PPL60 8726XL	6-0		13 mm 3/8 Circle Round Body Round Body	75	12	BLUE
POLYPROPYLENE BU 6-0 X 60-10MM CU RB DN	PPL60 8804	6-0		10 mm 3/8 Circle Round Body Round Body	60	12	BLUE

Description	FG Code	Size		Needle Profile Shape, Type & Length	Suture Length (cm)	Pack Size	Suture Colour
FILAPROP (POLYPROPYLENE SUTURE)							
POLYPROPYLENE BU 6-0 X 75-10MM CU RB DVB	PPL60 8804BX	6-0		10 mm 3/8 Circle Round Body Round Body	75	12	BLUE
POLYPROPYLENE BU 6-0 X 60-11MM CU TC DN	PPL60 8802	6-0		11 mm 3/8 Circle Taper Cut Taper Cut	60	12	BLUE
POLYPROPYLENE BU 6-0 X 75-13MM CU RB DNV	PPL608706	6-0		13 mm 3/8 Circle Round Body Round Body	75	12	BLUE
POLYPROPYLENE BU 6-0 X 60-13MM CU TC DN	PPL608707	6-0		13 mm 3/8 Circle Taper Cut Taper Cut	60	12	BLUE
POLYPROPYLENE BU 6-0 X 60-09MM CU TC DN	PPL608807	6-0		9 mm 3/8 Circle Taper Cut Taper Cut	60	12	BLUE
POLYPROPYLENE BU 6-0 X 75-13MM CU TC DN	PPL608814	6-0		13 mm 3/8 Circle Taper Cut Taper Cut	75	12	BLUE
POLYPROPYLENE BU 6-0 X 60-13MM HC TC DN	PPL608815	6-0		13 mm 1/2 Circle Taper Cut Taper Cut	60	12	BLUE
POLYPROPYLENE BU 6-0 X 75-13MM HC RB DN	PPL608714	6-0		13 mm 1/2 Circle Round Body Round Body	75	12	BLUE
POLYPROPYLENE BU 6-0 X 90-09MM CU TC	PPL608806	6-0		9 mm 3/8 Circle Taper Cut	90	12	BLUE
POLYPROPYLENE BU 6-0 X 75-10MM HC RB DN	PPL608610	6-0		10 mm 1/2 Circle Round Body Round Body	75	12	BLUE
POLYPROPYLENE BU 6-0 X 75-16MM CU RC	PPL608616	6-0		16 mm 3/8 Circle Reverse Cutting	75	12	BLUE
POLYPROPYLENE BU 6-0 X 60-13MM HC RB DVB	PPL608860	6-0		13 mm 1/2 Circle Round Body Round Body	60	12	BLUE
POLYPROPYLENE BU 6-0 X 75-13MM CU RB	PPL608889	6-0		13 mm 3/8 Circle Round Body	75	12	BLUE
POLYPROPYLENE BU 6-0 X 75-13MM CU TC	PPL607975	6-0		13 mm 3/8 Circle Taper Cut	75	12	BLUE
POLYPROPYLENE BU 6-0 X 75-9.3MM CU RB DN	PPL607406	6-0		9.3 mm 3/8 Circle Round Body Round Body	75	12	BLUE
POLYPROPYLENE BU 6-0 X 45-16MM CU TC	PPL608005	6-0		16 mm 3/8 Circle Taper Cut	45	12	BLUE
POLYPROPYLENE BU 6-0 X 45-11MM CU RC	PPL608868	6-0		11 mm 3/8 Circle Reverse Cutting	45	12	BLUE
POLYPROPYLENE BU 6-0 X 60-9.3MM CU RB DN	PPL608712	6-0		9.3 mm 3/8 Circle Round Body Round Body	60	12	BLUE

Description	FG Code	Size		Needle Profile Shape, Type & Length	Suture Length (cm)	Pack Size	Suture Colour
FILAPROP (POLYPROPYLENE SUTURE)							
POLYPROPYLENE BU 6-0 X 70-13MM CU CT	PPL60 817	6-0		13 mm 3/8 Circle Conventional Cutting	70	12	BLUE
POLYPROPYLENE BU 6-0 X60-9MM CURB DNVB	PPL608305	6-0		9 mm 3/8 Circle Round Body Round Body	60	12	BLUE
POLYPROPYLENE BU 6-0 X 45-13MM CU RC	PPL608871	6-0		13 mm 3/8 Circle Reverse Cutting	45	12	BLUE
POLYPROPYLENE BU 6-0 X 75-11MM CU TC DN	PPL608802XL	6-0		11 mm 3/8 Circle Taper Cut Taper Cut	75	12	BLUE
POLYPROPYLENE BU 6-0 X 75-12MM CU RB DN	PPL608776	6-0		12 mm 3/8 Circle Round Body Round Body	75	12	BLUE
POLYPROPYLENE BU 6-0 X60-13MM CURB DNVB	PPL608307	6-0		13 mm 3/8 Circle Round Body Round Body	60	12	BLUE
POLYPROPYLENE BU 6-0 X 60-10MM HC RB DN	PPL60 8207	6-0		10 mm 1/2 Circle Round Body Round Body	60	12	BLUE
POLYPROPYLENE BU 6-0 X 60-08MM CU RB DNV	PPL60 8305BL	6-0		8 mm 3/8 Circle Round Body Round Body	60	12	BLUE
POLYPROPYLENE BU 6-0 X 75-10MM HCRBDN C1	PPL60 8906	6-0		10 mm 1/2 Circle Round Body Round Body	75	12	BLUE
POLYPROPYLENE BU 6-0 X 45-13MM CURB DN	PPL60 8726SL	6-0		13 mm 3/8 Circle Round Body Round Body	45	12	BLUE
POLYPROPYLENE BU 6-0 X 75-9MM CU RB DN	PPL608907	6-0		9 mm 3/8 Circle Round Body Round Body	75	12	BLUE
POLYPROPYLENE BU 6-0 X 90-16MM CU CT	PPL608873	6-0		16 mm 3/8 Circle Conventional Cutting	90	12	BLUE
POLYPROPYLENE BU 6-0 X 45-10MM CU RC	PPL6013433	6-0		10 mm 3/8 Circle Reverse Cutting	45	12	BLUE
POLYPROPYLENE BU 6-0 X 45-11MM CU SM	PPL608697	6-0		11 mm 3/8 Circle Spatula Micropoint	45	12	BLUE
POLYPROPYLENE BU 6-0 X 75-9MM CU RB DNVB	PPL608305XL	6-0		9 mm 3/8 Circle Round Body Round Body	75	12	BLUE
POLYPROPYLENE BU 6-0 X 75-9MM CU TC DN	PPL608807XL	6-0		9 mm 3/8 Circle Taper Cut Taper Cut	75	12	BLUE
POLYPROPYLENE BU 6-0 X 75-11MM CU RB DVB	PPL608597BXL	6-0		11 mm 3/8 Circle Round Body Round Body	75	12	BLUE
POLYPROPYLENE BU 6-0 X 75-13MM CU TC DVB	PPL608814BL	6-0		13 mm 3/8 Circle Taper Cut Taper Cut	75	12	BLUE

Description	FG Code	Size	Needle Profile Shape, Type & Length	Suture Length (cm)	Pack Size	Suture Colour
-------------	---------	------	---	--------------------------	--------------	------------------

FILAPROP (POLYPROPYLENE SUTURE)

POLYPROPYLENE BU 6-0 X 75-11MM CUTC DNVB	PPL608802BXL	6-0		11 mm 3/8 Circle Taper Cut Taper Cut	75	12	BLUE
POLYPROPYLENE BU 7-0 X 60-9MM CU RB DNC1	PPL70 8702	7-0		9 mm 3/8 Circle Round Body Round Body	60	12	BLUE
POLYPROPYLENE BU 7-0 X 60-9MM CU TC DN	PPL70 8704	7-0		9 mm 3/8 Circle Taper Cut Taper Cut	60	12	BLUE
POLYPROPYLENE BU 2 X 7-0X60-8MM CU RB DN	PPL70 866ML	7-0		8 mm 3/8 Circle Round Body Round Body	60	12	BLUE
POLYPROPYLENE BU 7-0 X 60-8MM CU TC DN	PPL70 804	7-0		8 mm 3/8 Circle Taper Cut Taper Cut	60	12	BLUE
POLYPROPYLENE BU 7-0 X 60-6MM CU RB DN	PPL70 811	7-0		6 mm 3/8 Circle Round Body Round Body	60	12	BLUE
POLYPROPYLENE BU 7-0 X 60-8MM CU RB DNN1	PPL70 866	7-0		8 mm 3/8 Circle Round Body Round Body	60	12	BLUE
POLYPROPYLENE BU 7-0 X 60-9MM CU RB DNC2	PPL70 873	7-0		9 mm 3/8 Circle Round Body Round Body	60	12	BLUE
POLYPROPYLENE BU 7-0 X 45-9MM CU RB DN	PPL70 8701	7-0		9 mm 3/8 Circle Round Body Round Body	45	12	BLUE
POLYPROPYLENE BU 7-0 X 60-8MM CU RB DN	PPL70 8735	7-0		8 mm 3/8 Circle Round Body Round Body	60	12	BLUE
POLYPROPYLENE BU 7-0 X 60-8MM CU RB DNN2	PPL70 8078	7-0		8 mm 3/8 Circle Round Body Round Body	60	12	BLUE
POLYPROPYLENE BU 7-0 X 60-10MM CU RB DN	PPL70 8070	7-0		10 mm 3/8 Circle Round Body Round Body	60	12	BLUE
POLYPROPYLENE BU 7-0 X 60-8MM CU RB DNVB	PPL70 8078BL	7-0		8 mm 3/8 Circle Round Body Round Body	60	12	BLUE
POLYPROPYLENE BU 7-0 X 60-09MM CU RB DN	PPL708304	7-0		9 mm 3/8 Circle Round Body Round Body	60	12	BLUE
POLYPROPYLENE BU 7-0 X 60-13MM CU RB DN	PPL708725	7-0		13 mm 3/8 Circle Round Body Round Body	60	12	BLUE
POLYPROPYLENE BU 7-0 X 60-13MM CU TC DN	PPL708813	7-0		13 mm 3/8 Circle Taper Cut Taper Cut	60	12	BLUE
POLYPROPYLENE BU 7-0 X 60-9MM CU RB DN VB	PPL70 8304BL	7-0		9 mm 3/8 Circle Round Body Round Body	60	12	BLUE

Description	FG Code	Size		Needle Profile Shape, Type & Length	Suture Length (cm)	Pack Size	Suture Colour
FILAPROP (POLYPROPYLENE SUTURE)							
POLYPROPYLENE BU 7-0 X 75-8MM CU RB DNVB	PPL708078BXL	7-0		8 mm 3/8 Circle Round Body Round Body	75	12	BLUE
POLYPROPYLENE BU 7-0 X 60-11MM CU TC DN	PPL70 8700	7-0		11 mm 3/8 Circle Taper Cut Taper Cut	60	12	BLUE
POLYPROPYLENE BU 7-0 X 60-9MM CU RB DNVB	PPL708021	7-0		9 mm 3/8 Circle Round Body Round Body	60	12	BLUE
POLYPROPYLENE BU 7-0 X 75-9MM CU TC DN	PPL708704XL	7-0		9 mm 3/8 Circle Taper Cut Taper Cut	75	12	BLUE
POLYPROPYLENE BU 7-0 X 75-9MM CU RB	PPL707402	7-0		9 mm 3/8 Circle Round Body	75	12	BLUE
POLYPROPYLENE BU 7-0 X 75-9MM CU RB DN	PPL708701XL	7-0		9 mm 3/8 Circle Round Body Round Body	75	12	BLUE
POLYPROPYLENE BU 7-0 X 75-9MM CU RB DN VB	PPL707222	7-0		9 mm 3/8 Circle Round Body Round Body	75	12	BLUE
POLYPROPYLENE BU 7-0 X 45-9MM CU RB DNVB	PPL707225	7-0		9 mm 3/8 Circle Round Body Round Body	45	12	BLUE
POLYPROPYLENE BU 7-0 X 90-11MM CU TC DN	PPL708801	7-0		11 mm 3/8 Circle Taper Cut Taper Cut	90	12	BLUE
POLYPROPYLENE BU 7-0 X 60-10MM CU RB DN C1	PPL70 2070	7-0		10 mm 3/8 Circle Round Body Round Body	60	12	BLUE
POLYPROPYLENE BU 7-0 X 60-6MM CU RB DN C1	PPL70 2076	7-0		6 mm 3/8 Circle Round Body Round Body	60	12	BLUE
POLYPROPYLENE BU 7-0 X 75-8MM CURBDNVBC1	PPL70 2078BX	7-0		8 mm 3/8 Circle Round Body Round Body	75	12	BLUE
POLYPROPYLENE BU 7-0 X 60-9MM CU RB DN	PPL70 2079	7-0		9 mm 3/8 Circle Round Body Round Body	60	12	BLUE
POLYPROPYLENE BU 7-0 X 75-9MM CU RBDN C1	PPL70 2079XL	7-0		9 mm 3/8 Circle Round Body Round Body	75	12	BLUE
POLYPROPYLENE BU 7-0 X 60-9MM CU TC DNVB	PPL70 8121BL	7-0		9 mm 3/8 Circle Taper Cut Taper Cut	60	12	BLUE
POLYPROPYLENE BU 7-0 X 60-6.5MM CU RB DN	PPL708207	7-0		6.5 mm 3/8 Circle Round Body Round Body	60	12	BLUE
POLYPROPYLENE BU 7-0 X 60-09MM CU TC DN	PPL708705	7-0		9 mm 3/8 Circle Taper Cut Taper Cut	60	12	BLUE

Description	FG Code	Size		Needle Profile Shape, Type & Length	Suture Length (cm)	Pack Size	Suture Colour
FILAPROP (POLYPROPYLENE SUTURE)							
POLYPROPYLENE BU 7-0 X 75-9.3MM CU TC DN	PPL707237	7-0		9.3 mm 3/8 Circle Taper Cut Taper Cut	75	12	BLUE
POLYPROPYLENE BU 7-0 X 60-9.5MM CU RB DN	PPL70 8704RB	7-0		9.5 mm 3/8 Circle Round Body Round Body	60	12	BLUE
POLYPROPYLENE BU 7-0 X 75-9MM CU TC DNVB	PPL707237BL	7-0		9 mm 3/8 Circle Taper Cut Taper Cut	75	12	BLUE
POLYPROPYLENE BU 7-0 X 75-6MM CU RB DN	PPL702076XL	7-0		6 mm 3/8 Circle Round Body Round Body	75	12	BLUE
POLYPROPYLENE BU 2 X 8-0X60-8MM CU RB DN	PPL80 876ML	8-0		8 mm 3/8 Circle Round Body Round Body	60	12	BLUE
POLYPROPYLENE BU 8-0 X 60-6MM CU RB DNTN	PPL80 812	8-0		6 mm 3/8 Circle Round Body Round Body	60	12	BLUE
POLYPROPYLENE BU 8-0 X 60-9MM CU RB DN	PPL80 813	8-0		9 mm 3/8 Circle Round Body Round Body	60	12	BLUE
POLYPROPYLENE BU 8-0 X 60-8MM CU RB DNTN	PPL80 876	8-0		8 mm 3/8 Circle Round Body Round Body	60	12	BLUE
POLYPROPYLENE BU 8-0 X 60-6.5MM CU RB DN	PPL80 8732	8-0		6.5 mm 3/8 Circle Round Body Round Body	60	12	BLUE
POLYPROPYLENE BU 8-0 X 60-8MM CU RB DNTH	PPL80 8741	8-0		8 mm 3/8 Circle Round Body Round Body	60	12	BLUE
POLYPROPYLENE BU 8-0 X 60-8MM CU TC VBDN	PPL80 8101BL	8-0		8 mm 3/8 Circle Taper Cut Taper Cut	60	12	BLUE
POLYPROPYLENE BU 8-0 X 60-8MM CU RB DN	PPL80 8088	8-0		8 mm 3/8 Circle Round Body Round Body	60	12	BLUE
POLYPROPYLENE BU 8-0 X 60-6MM CU RB DNTH	PPL80 8086XF	8-0		6 mm 3/8 Circle Round Body Round Body	60	12	BLUE
POLYPROPYLENE BU 8-0 X 13-06MM CU RB	PPL802775	8-0		6 mm 3/8 Circle Round Body	13	12	BLUE
POLYPROPYLENE BU 8-0 X 75-8MM CU RB DNTH	PPL808741XL	8-0		8 mm 3/8 Circle Round Body Round Body	75	12	BLUE
POLYPROPYLENE BU 8-0 X 45-6MM CU RB DN	PPL802777	8-0		6 mm 3/8 Circle Round Body Round Body	45	12	BLUE
POLYPROPYLENE BU 8-0 X 60-9MM CU TC DNV	PPL808703	8-0		9 mm 3/8 Circle Taper Cut Taper Cut	60	12	BLUE

Description	FG Code	Size		Needle Profile Shape, Type & Length	Suture Length (cm)	Pack Size	Suture Colour
FILAPROP (POLYPROPYLENE SUTURE)							
POLYPROPYLENE BU 8-0 X 75-6.5MM CU RB DN	PPL808732XL	8-0		6.5 mm 3/8 Circle Round Body Round Body	75	12	BLUE
POLYPROPYLENE BU 8-0 X 60-6.5MM CURBDNEX	PPL805662	8-0		6.5 mm 3/8 Circle Round Body Round Body	60	12	BLUE
POLYPROPYLENE BU 8-0 X 75-6.5MM CURB DNV	PPL808732BXL	8-0		6.5 mm 3/8 Circle Round Body Round Body	75	12	BLUE
POLYPROPYLENE BU 9-0 X 13-5.1MM CU RB DN	PPL902783	9-0		5 mm 3/8 Circle Round Body Round Body	13	12	BLUE
POLYPROPYLENE BU 9-0 X 90-5MM CU RB DN	PPL902780	9-0		5 mm 3/8 Circle Round Body Round Body	90	12	BLUE
POLYPROPYLENE BU 0 X 100-40MM HC RC	PPL108470XL	0		40 mm 1/2 Circle Reverse Cutting	100	12	BLUE
POLYPROPYLENE BU 1 X 100-48MM HC RC	PPL018576	1		48 mm 1/2 Circle Reverse Cutting	100	12	BLUE
POLYPROPYLENE BU 4-0 X90-17MM HC RB DNTP	PPG40 846	4-0		17 mm 1/2 Circle Round Body Round Body	90	12	BLUE
POLYPROPYLENE BU 4-0X90-17.5MM HCTC DTP3	PPG40 8935	4-0		17.5 mm 1/2 Circle Taper Cut Taper Cut	90	12	BLUE
POLYPROPYLENE BU 5-0 X 75-13MM HCRB DNTP	PPG50 810	5-0		13 mm 1/2 Circle Round Body Round Body	75	12	BLUE
POLYPROPYLENE BU 5-0 X60-10MM HCRB DNTP3	PPG50 8803	5-0		10 mm 1/2 Circle Round Body Round Body	60	12	BLUE
POLYPROPYLENE BU 5-0 X 75-13MM HC RB DN TP3	PPG50 8710	5-0		13 mm 1/2 Circle Round Body Round Body	75	12	BLUE
POLYPROPYLENE BU 6-0 X60-10MM HCRB DNTP3	PPG60 8207	6-0		10 mm 1/2 Circle Round Body Round Body	60	12	BLUE
POLYPROPYLENE BU 6-0X60-10MM HCRB DNTP3	PPG60 8260SR	6-0		10 mm 1/2 Circle Round Body Round Body	60	12	BLUE
POLYPROPYLENE BU 7-0X60-09MM HCRB DNTP3	PPG70 8270SR	7-0		9 mm 1/2 Circle Round Body Round Body	60	12	BLUE

MERICRON XL™

Polyester Suture

Indications

- Cardiovascular procedures
- Ophthalmic procedures
- Neurological procedures


Advantages

- Smoother suture surface
- Less resistance to knot sliding
- Better strength
- Pliable structure



Product Overview



Structure	: Braided
Colour	: Green, White
Chemical composition	: Polyethylene terephthalate
Coating	: Polybutylate
Sizes	: USP 6/0 (metric 0.7) - USP 5 (metric 7)
Sterilization	: Gamma Radiation





Description	FG Code	Size		Needle Profile Shape, Type & Length	Suture Length (cm)	Pack Size	Suture Colour
MERICRON XL (POLYESTER SUTURE)							
POLYESTER GR 0 X 90-50MM HC TC	Me100979	0		50 mm 1/2 Circle Taper Cut	90	12	GREEN
POLYESTER GR 0 X 90-26MM HC TC DN	ME6978	0		26 mm 1/2 Circle Taper Cut Taper Cut	90	12	GREEN
POLYESTER GR 0 X 75-22MM HC RC	ME517G	0		22 mm 1/2 Circle Reverse Cutting	75	12	GREEN
POLYESTER GR 0 X 90-36MM HC RB	ME844G	0		36 mm 1/2 Circle Round Body	90	12	GREEN
POLYESTER GR 0 X 75-30MM HC RB	ME106030	0		30 mm 1/2 Circle Round Body	75	12	GREEN
POLYESTER GR 0 X 75-26MM HC RB	ME10361	0		26 mm 1/2 Circle Round Body	75	12	GREEN
POLYESTER GR 0 X 75-31MM HC RB	ME100975	0		31 mm 1/2 Circle Round Body	75	12	GREEN
POLYESTER GR 0 X 90-36MM HC RC	ME1032049	0		36 mm 1/2 Circle Reverse Cutting	90	12	GREEN
POLYESTER GR 0 X 90-30MM CU RC	ME106680	0		30 mm 3/8 Circle Reverse Cutting	90	12	GREEN
POLYESTER GR 0 X 75-36MM CU RC	ME107562	0		36 mm 3/8 Circle Reverse Cutting	75	12	GREEN
POLYESTER GREEN 1 X 90-30MM HC RC	ME0132047	1		30 mm 1/2 Circle Reverse Cutting	90	12	GREEN
POLYESTER GREEN 1 X 75-22MM HC RC	ME010518	1		22 mm 1/2 Circle Reverse Cutting	75	12	GREEN
POLYESTER GREEN X 90-37MM HC RB	ME016137	1		37 mm 1/2 Circle Round Body	90	12	GREEN
POLYESTER GREEN 1 X 90-30MM HC RB	ME016130	1		30 mm 1/2 Circle Round Body	90	12	GREEN
POLYESTER GREEN 1 X 75-31MM J TYPE RB H	ME010969J	1		31 mm J Shape Round Body	75	12	GREEN
POLYESTER GREEN 1 X 90-30MM HC RB-3D	ME016130Z	1		30 mm 1/2 Circle Round Body	90	36	GREEN

Description	FG Code	Size		Needle Profile Shape, Type & Length	Suture Length (cm)	Pack Size	Suture Colour
MERICRON XL (POLYESTER SUTURE)							
POLYESTER GREEN 1 X 90-36MM HC TC	ME014407	1		36 mm 1/2 Circle Taper Cut	90	12	GREEN
POLYESTER GR 2 X 100-45MM HC TC	ME643	2		45 mm 1/2 Circle Taper Cut	100	12	GREEN
POLYESTER GR 2 X 75-45MM HC TC (4 in 1)	ME4843	2		45 mm 1/2 Circle Taper Cut	75	6	GREEN
POLYESTER GR 2 X 75-45MM HC TC	ME4842	2		45 mm 1/2 Circle Taper Cut	75	12	GREEN
POLYESTER GR 2 X 90-37MM HC RB	ME026237	2		37 mm 1/2 Circle Round Body	90	12	GREEN
POLYESTER GR 5 X 75-55MM HC TC	ME646	5		55 mm 1/2 Circle Taper Cut	75	12	GREEN
POLYESTER GR 5 X 75-55MM HC TC (4 in 1)	ME4846	5		55 mm 1/2 Circle Taper Cut	75	6	GREEN
POLYESTER GR 2-0 X 90-17MM HC TC DN C2	ME617	2-0		17 mm 1/2 Circle Taper Cut Taper Cut	90	12	GREEN
POLYESTER WH 2-0 X 76-17MM HC TC	ME673W	2-0		17 mm 1/2 Circle Taper Cut	76	12	WHITE
POLYESTER WH 2-0 X 90-25MM HC TC DN	ME687W	2-0		25 mm 1/2 Circle Taper Cut Taper Cut	90	12	WHITE
POLYESTER GR 2-0 X 90-17MM HC TC DN C1	ME6937	2-0		17 mm 1/2 Circle Taper Cut Taper Cut	90	12	GREEN
POLYESTER GR 2-0 X 90-26MM HC TC DN	ME6977	2-0		26 mm 1/2 Circle Taper Cut Taper Cut	90	12	GREEN
POLYESTER WH 2-0 X 90-26MM HC TC DN	ME6987	2-0		26 mm 1/2 Circle Taper Cut Taper Cut	90	12	WHITE
POLYESTER WH 2-0 X 90-17MM HC TC DN	ME617W	2-0		17 mm 1/2 Circle Taper Cut Taper Cut	90	12	WHITE

Description	FG Code	Size		Needle Profile Shape, Type & Length	Suture Length (cm)	Pack Size	Suture Colour
MERICRON XL (POLYESTER SUTURE)							
POLYESTER GREEN 2-0 X 90-17MM HC RB DN	ME206760	2-0		17 mm 1/2 Circle Round Body Round Body	90	12	GREEN
POLYESTER GREEN 2-0 X 90-26MM HC RB DN	ME206767	2-0		26 mm 1/2 Circle Round Body Round Body	90	12	GREEN
POLYESTER GREEN 2-0 X 90-26MM CU TC DN	ME206997	2-0		26 mm 3/8 Circle Taper Cut Taper Cut	90	12	GREEN
POLYESTER GREEN 2-0 X 75-31MM HC RB	ME200932	2-0		31 mm 1/2 Circle Round Body	75	12	GREEN
POLYESTER GREEN 2-0 X 75-22MM HC RB	ME206833	2-0		22 mm 1/2 Circle Round Body	75	12	GREEN
POLYESTER GREEN 2-0 X 75-26MM HC RB	ME200411	2-0		26 mm 1/2 Circle Round Body	75	12	GREEN
POLYESTER XL WH 2-0 X 75-26MM HC RB	ME200412	2-0		26 mm 1/2 Circle Round Body	75	12	WHITE
POLYESTER WHGR 2-0 X 75-26MM HC TC DN5	ME201075	2-0		26 mm 1/2 Circle Taper Cut Taper Cut	75	6	GREEN/ WHITE
POLYESTER WHGR 2-0 X 75-17MM HC TC DN5	ME201017	2-0		17 mm 1/2 Circle Taper Cut Taper Cut	75	6	GREEN/ WHITE
POLYESTER WHITE 2-0 X 90-16MM HC RB DN	ME208051	2-0		16 mm 1/2 Circle Round Body Round Body	90	12	WHITE
POLYESTER WHGR 2-0 X 75-26MM HC TC DN8	ME201356	2-0		26 mm 1/2 Circle Taper Cut Taper Cut	75	6	GREEN/ WHITE
POLYESTER GR 2-0 X 75-36MM HC RB	ME206236	2-0		36 mm 1/2 Circle Round Body	75	12	GREEN
POLYESTER GR 2-0 X 90-25MM CU TC	ME200095	2-0		25 mm 3/8 Circle Taper Cut	90	12	GREEN
POLYESTER WH 2-0 X 75-17MM HC RB	ME2043122	2-0		17 mm 1/2 Circle Round Body	75	12	WHITE
POLYESTER GR 2-0 X 90-24MM CU CT	ME206224	2-0		24 mm 3/8 Circle Conventional Cutting	90	12	GREEN
POLYESTER GR 2-0 X 75-30MM HC RC	ME206230	2-0		30 mm 1/2 Circle Reverse Cutting	75	12	GREEN

Description	FG Code	Size		Needle Profile Shape, Type & Length	Suture Length (cm)	Pack Size	Suture Colour
MERICRON XL (POLYESTER SUTURE)							
POLYESTER GR 4-0 X 76-8MM QC SM	ME682	4-0		8 mm 1/4 Circle Spatula Micropoint	76	12	GREEN
POLYESTER GREEN 4-0 X 75-26MM HC RB	ME406831	4-0		26 mm 1/2 Circle Round Body	75	12	GREEN
POLYESTER GR 4-0 X 90-17MM HC TC DN	ME4013395	4-0		17 mm 1/2 Circle Taper Cut Taper Cut	90	12	GREEN
POLYESTER GREEN 4-0 X 75-20MM HC RB	ME406761	4-0		20 mm 1/2 Circle Round Body	75	12	GREEN
POLYESTER GREEN 4-0 X 75-13MM CU RB DN	ME6891G	4-0		13 mm 3/8 Circle Round Body Round Body	75	12	GREEN
POLYESTER GREEN 4-0 X 90-17.5MM HC RC DN	ME6935G	4-0		17.5 mm 1/2 Circle Reverse Cutting Reverse Cutting	90	12	GREEN
POLYESTER GR 4-0 X 75-20MM HC RB DN	ME406761DN	4-0		20 mm 1/2 Circle Round Body Round Body	75	12	GREEN
POLYESTER GR 5-0 X 45-8MM QC SM DN	ME683	5-0		8 mm 1/4 Circle Spatula Micropoint Spatula Micropoint	45	12	GREEN
POLYESTER GREEN 5-0 X 60-17MM HC RB DN	ME506757	5-0		17 mm 1/2 Circle Round Body Round Body	60	12	GREEN
POLYESTER GREEN 5-0 X 45-8MM HC SM DN	ME500883DN	5-0		8 mm 1/2 Circle Spatula Micropoint Spatula Micropoint	45	12	GREEN
POLYESTER GREEN 5-0 X 75-13MM HC RB DN	ME6889	5-0		13 mm 1/2 Circle Round Body Round Body	75	12	GREEN
POLYESTER GREEN 5-0 X 75-13MM CU RB DN	ME6890G	5-0		13 mm 3/8 Circle Round Body Round Body	75	12	GREEN
POLYESTER WH 5-0 X 45-8MM QC SM DN	ME50893DN	5-0		8 mm 1/4 Circle Spatula Micropoint Spatula Micropoint	45	12	WHITE
POLYESTER WH 5-0 X 45-8MM HC SM DN	ME50883DN	5-0		8 mm 1/2 Circle Spatula Micropoint Spatula Micropoint	45	12	WHITE
POLYESTER GR 6-0 X 75-11MM CU RB DN	ME6597	6-0		11 mm 3/8 Circle Round Body Round Body	75	12	GREEN
POLYESTER GREEN 1 X 90-36MM HC RC - 3D	ME012520Z	1		36 mm 1/2 Circle Reverse Cutting	90	36	GREEN

Description	FG Code	Size		Needle Profile Shape, Type & Length	Suture Length (cm)	Pack Size	Suture Colour
MERICRON XL (POLYESTER SUTURE)							
POLYESTER GREEN 3-0 X 75-17MM HC TC-3D	ME306935Z	3-0		17 mm 3/8 Circle Reverse Cutting	75	36	GREEN
POLYESTER GREEN 4-0 X 45-13MM CU RC-3D	ME407573Z	4-0		13 mm 3/8 Circle Reverse Cutting	45	36	GREEN
POLYESTER GREEN 4-0 X 45-19MM CU RC-3D	ME407142Z	4-0		19 mm 3/8 Circle Reverse Cutting	45	36	GREEN
POLYESTER WHGR 2-0X75-25MM HC TC DN TP6	MME6577X	2-0		25 mm 1/2 Circle Taper Cut Taper Cut	75	6	GREEN/ WHITE
POLYESTER WHGR 2-0X75-17MM HC TC DN TP6	MME6555X	2-0		17 mm 1/2 Circle Taper Cut Taper Cut	75	6	GREEN/ WHITE
POLYESTER WHGR 2-0X75-25MM HC TC DN TP5	MME6577	2-0		25 mm 1/2 Circle Taper Cut Taper Cut	75	6	GREEN/ WHITE
POLYESTER WHGR 2-0X75-17MM HC TC DN TP5	MME6555	2-0		17 mm 1/2 Circle Taper Cut Taper Cut	75	6	GREEN/ WHITE
POLYESTER GR 2-0 X 90-17MM HC TC DN TP	PME617	2-0		17 mm 1/2 Circle Taper Cut Taper Cut	90	12	GREEN
POLYESTER WH 2-0 X 90-17MM HC TC DN TP	PME617W	2-0		17 mm 1/2 Circle Taper Cut Taper Cut	90	12	WHITE
POLYESTER GR 2-0 X 90-25MM HC TC DN TP	PME687	2-0		25 mm 1/2 Circle Taper Cut Taper Cut	90	12	GREEN
POLYESTER WH 2-0 X 90-25MM HC TC DN TP	PME687W	2-0		25 mm 1/2 Circle Taper Cut Taper Cut	90	12	WHITE
POLYSTER XL GR 2-0 X 75-26MM HC TC DN TP	ME202613	2-0		26 mm 1/2 Circle Taper Cut Taper Cut	75	12	GREEN
POLYSTER XL GR 2-0 X 75-17MM HC TC DN TP	ME200553	2-0		17 mm 1/2 Circle Taper Cut Taper Cut	75	12	GREEN
POLYESTER XL WHGR 2-0X75-17MM HCTCDN TP3	MME6555MP	2-0		17 mm 1/2 Circle Taper Cut Taper Cut	75	6	GREEN/ WHITE

Description	FG Code	Size		Needle Profile Shape, Type & Length	Suture Length (cm)	Pack Size	Suture Colour
MERICRON XL (POLYESTER SUTURE)							
POLYESTER XL WHGR 2-0X75-25MM HCTCDN TP3	MME6577MP	2-0		25 mm 1/2 Circle Taper Cut Taper Cut	75	6	GREEN/ WHITE
POLYESTER XL GR 2-0 X 90-17MM HC TCDN P3	PME617MP	2-0		17 mm 1/2 Circle Taper Cut Taper Cut	90	12	GREEN
POLYESTER XL WH 2-0 X 90-17MM HC TCDN P3	PME617WMP	2-0		17 mm 1/2 Circle Taper Cut Taper Cut	90	12	WHITE
POLYESTER GR 2-0 X 90-17MM HC TC DN TP6	PME6937	2-0		17 mm 1/2 Circle Taper Cut Taper Cut	90	12	GREEN
POLYESTER GR 2-0 X 90-26MM HC TC DN TP6	PME6977G	2-0		26 mm 1/2 Circle Taper Cut Taper Cut	90	12	GREEN
POLYESTER WH 2-0 X 90-26MM HC TC DN TP6	PME6987W	2-0		26 mm 1/2 Circle Taper Cut Taper Cut	90	12	WHITE
POLYESTER XL WH/GR 2-0X90-25MM HC TC DNTP5	SM10B77	2-0		25 mm 1/2 Circle Taper Cut Taper Cut	90	6	GREEN/ WHITE
POLYESTER GR 5-0 X 75-13MM HC RB DN TP3	PME6889G	5-0		13 mm 1/2 Circle Round Body Round Body	75	12	GREEN

Meril Endo Surgery Pvt. Ltd

Third Floor E1-E3, Meril Park, Survey Number 135/2/B & 172/2
Muktanand Marg, Chala, Vapi 396 191. Gujarat. India.

T : +91 260 305 2100

F : +91 260 305 2125

E : askinfo@merillife.com

Meril life sciences Pvt. Ltd.

301, A - Wing, Business Square, Chakala,
Andheri Kurla Road, Andheri East, Mumbai 400 093

T : +91 22 39350700

F : +91 22 39350777

E : askinfo@merillife.com

Meril, inc.

2436 Emrick Boulevard, Bethlehem, PA - 18020

P : (610)500-2080

F : (610)317-1672

Meril South America

Doc Med LTDA, Al. dos Tupiniquins, 1079 -
Cep: 04077-003 - Moema., Sao Paulo. Brazil.

T : +55 11 3624 5935

F : +55 11 3624 5936

Meril GmbH.

Bornheimer Strasse 135-137, D-53119 Bonn., Germany.

T : +49 228 7100 4000

F : +49 228 7100 4001

Meril Tibbi Cihazlar

Meril Tibbi Cihazlar İmalat ve Ticaret A.Ş. İçerenköy Mah.
Çetinkaya Sok. Prestij Plaza No:28, Kat:4 Ataşehir, 34752,
İstanbul / Turkey

T : +90 216 641 44 24

F : +90 216 641 44 25

Meril China Co. Ltd.

2301b,23f, Lixin Plaza, no 90, South Hubin Road,Xiamen,
China.

T : 0086-592-5368505

F : 0086-592-5368519

Meril SA Pty Ltd

102, 104, S101 And S102, Boulevard West Office Park,
142 Western Service Road, Erf 813 Woodmead Extension 17,
Sandton, Johannesburg – 2191 South Africa

T : +27 11 465-2049

F : +27 86 471 7941

Meril Medical LLC.

Nauchnyi Proezd 19, Moscow , Russia – 117 246.

Office - +7 495 772 7643

EU representative.

Obelis S.A., Bd, General Wahis 53, 1030, Brussels, Belgium.

T : +32 2 732 5954

F : +32 2 732 6003

E : mail@obelis.net

Meril Bangladesh Pvt. Ltd.

Union Heights - 01, 55 - 2, Bir Uttam Qazi Nuruzzaman Sarak,
West Panthapath Dhaka - 1205 Bangladesh

T : +88 02 9122714 / 24

