

Manufacturer Declaration Letter

We hereby clarify some functions of YA-D5-5 electric hospital bed as follow.

1. Patient control / Functions - sedentary positioning.

YES, It means both back board and knee board adjust to some degrees, it will be available for patient sedentary positioning.

Meanwhile, we confirm IV pole MK-IS08C for hospital bed is height adjustable type.

2. Nurse control / Patient control lock - mandatory.

YES, there is manual key in remote hand control.

3. Mattress / thickness - model MK-M01, minimum 14 cm;

YES, mattress foam thickness can be customized.



Signed & Stamp:

A handwritten signature in black ink is written over a red circular stamp. The stamp contains the company name in English and Chinese, along with the registration number 44009978502E.

Name: Eric CC

Title: Sales Manager

Dated: 30-SEP-2024