

duraporter™

The **duraporter™** translucent transport box is **ideal for holding and storing** a variety of laboratory equipment and samples, including:

- Hazardous/sensitive research samples
- Clinical and Phlebotomy samples
- Blood bags and/or tissue/organ samples
- Public Health samples
- Laboratory instruments

Across the globe we are experiencing a heightened focus on disease control within hospitals, field hospitals, sea ports, airports, clinical and research laboratories, etc. Help these establishments safely handle, store and transport potentially infectious specimens within the facility with the **duraporter™** by Heathrow Scientific.

The **duraporter™** is a durable, leak-proof container that is easy to load, securely locks and is ready for transport. Available in clear, translucent red and yellow, the DuraPorter™ provides needed visibility for content identification and protection during sample transport.

The DuraPorter™ is in compliance with U.S. standards, OSHA 29 CFR 1910.1030.



Load



Lock



Transport

A silicone water-tight seal and three clasps provide secure closure

Heathrow Scientific®

HS120052 DuraPorter shown in use with a 72 place OneRack® test tube rack (order separately)

**Provides
Needed
Visibility and
Protection**
During Sample
Transport

The **duraporter**[™] Advantage

- Durable polycarbonate construction highly resistant to impact
- Available in three colors for use/owner identification
- Silicone water-tight seal and three lid clasps provide secure closure
- Carry handle folds neatly into lid, allowing the units to be stacked
- For laboratory testing, will accommodate 13 and 16 mm tubes in a 36 or 72-place OneRack[®] or similar size 13/16 mm tube rack
- Use as one large storage space, or insert separators to create smaller compartments
- The unit can be autoclaved between uses (leave open and unlocked with the silicone seal in place)



The **duraporter**[™] lid opens a full 180° to expedite loading and emptying

The entire unit can also be placed in a freezer or refrigerator if the enclosed samples need to be kept chilled.



Included separators can also be used to create inner compartments. Or, use without racks or separators as one large storage space for larger samples or supplies



Part No.	Color	Dimensions (mm)
HS120052	Clear with blue handles	380 x 196 x 160
HS120078	Yellow with yellow handles	380 x 196 x 160
HS120077	Red with red handles	380 x 196 x 160

Note: Always engage all three clamps before transporting samples. Do not use if clamps or housing are broken or cracked. Not for radioactive materials.

Ask for these and other innovative Heathrow Scientific products:

one rack[™]
test tube rack

36 or 72-Place Test Tube Racks

Choose autoclavable fiberglass-reinforced polyoxymethylene (Delrin[®]) or polypropylene.

Both versions have molded coordinates for tube orientation and a flat labeling area for easier identification. **8 racks per case.**



White ●	Blue ●	Green ●	Magenta ●	Description	Tubes	Wells	L x W x H, mm
HS27511A HS27512A	HS27511B HS27512B	HS27511C HS27512C	HS27511D HS27512D	Full-Size Delrin [®] Racks	13 mm 16 mm	72 72	104 x 202 x 59 127 x 250 x 70
HS27501A HS27502A	HS27501B HS27502B	HS27501C HS27502C	HS27501D HS27502D	Half-Size Delrin [®] Racks	13 mm 16 mm	36 36	104 x 104 x 59 127 x 127 x 70
HS27551A HS27552A	HS27551B HS27552B	HS27551C HS27552C	HS27551D HS27552D	Full-Size Polypropylene Racks	13 mm 16 mm	72 72	104 x 202 x 59 127 x 250 x 70