



## NON VACUUM BLOOD COLLECTION TUBE

## **Clot Activator**

REF	Additive	Cap Colour	Tube Size (mm)	Draw Volume (mL)	Tube Materia	Pack/Carton	Label
RJ1904103-5P	Clot Activator	🔴 Red	13x75	3.5mL	PET	100PCS/1200PCS	Paper/Transparent
RJ1904103-5G	Clot Activator	🔴 Red	13x75	3.5mL	GLASS	100PCS/1200PCS	Paper/Transparent
RJ1904104-0P	Clot Activator	Red	13x75	4mL	PET	100PCS/1200PCS	Paper/Transparent
RJ1904104-0G	Clot Activator	🔴 Red	13x75	4mL	GLASS	100PCS/1200PCS	Paper/Transparent
RJ1904205-0P	Clot Activator	🔴 Red	13×100	5mL	PET	100PCS/1200PCS	Paper/Transparent
RJ1904205-0G	Clot Activator	Red	13×100	5mL	GLASS	100PCS/1200PCS	Paper/Transparent
RJ1904206-0P	Clot Activator	🔴 Red	13x100	6mL	PET	100PCS/1200PCS	Paper/Transparent
RJ1904206-0G	Clot Activator	Red	13×100	6mL	GLASS	100PCS/1200PCS	Paper/Transparent
RJ1904207-0P	Clot Activator	🔴 Red	13×100	7mL	PET	100PCS/1200PCS	Paper/Transparent
RJ1904207-0G	Clot Activator	🔴 Red	13×100	7mL	GLASS	100PCS/1200PCS	Paper/Transparent
RJ1904308-0P	Clot Activator	🔴 Red	16×100	8mL	PET	100PCS/1200PCS	Paper/Transparent
RJ1904308-0G	Clot Activator	🔴 Red	16×100	8mL	GLASS	100PCS/1200PCS	Paper/Transparent
RJ1904309-0P	Clot Activator	🔴 Red	16×100	9mL	PET	100PCS/1200PCS	Paper/Transparent
RJ1904309-0G	Clot Activator	Red	16×100	9mL	GLASS	100PCS/1200PCS	Paper/Transparent
RJ1904310-0P	Clot Activator	🔴 Red	16×100	10mL	PET	100PCS/1200PCS	Paper/Transparent
RJ1904310-0G	Clot Activator	🔴 Red	16×100	10mL	GLASS	100PCS/1200PCS	Paper/Transparent

## Gel & Clot Activator

REF	Additive	Cap Colour	Tube Size (mm)	Draw Volume (mL)	Tube Materia	Pack/Carton	Label
RJ1905103-5P	Gel & Clot Activator	😑 Yellow	13x75	3.5mL	PET	100PCS/1200PCS	Paper/Transparent
RJ1905103-5G	Gel & Clot Activator	😑 Yellow	13x75	3.5mL	GLASS	100PCS/1200PCS	Paper/Transparent
RJ1905104-0P	Gel & Clot Activator	😑 Yellow	13x75	4mL	PET	100PCS/1200PCS	Paper/Transparent
RJ1905104-0G	Gel & Clot Activator	😑 Yellow	13x75	4mL	GLASS	100PCS/1200PCS	Paper/Transparent
RJ1905205-0P	Gel & Clot Activator	😑 Yellow	13×100	5mL	PET	100PCS/1200PCS	Paper/Transparent
RJ1905205-0G	Gel & Clot Activator	😑 Yellow	13×100	5mL	GLASS	100PCS/1200PCS	Paper/Transparent
RJ1905206-0P	Gel & Clot Activator	😑 Yellow	13×100	6mL	PET	100PCS/1200PCS	Paper/Transparent
RJ1905206-0G	Gel & Clot Activator	Yellow	13×100	6mL	GLASS	100PCS/1200PCS	Paper/Transparent
RJ1905207-0P	Gel & Clot Activator	😑 Yellow	13×100	7mL	PET	100PCS/1200PCS	Paper/Transparent
RJ1905207-0G	Gel & Clot Activator	😑 Yellow	13×100	7mL	GLASS	100PCS/1200PCS	Paper/Transparent
RJ1905308-0P	Gel & Clot Activator	😑 Yellow	16×100	8mL	PET	100PCS/1200PCS	Paper/Transparent
RJ1905308-0G	Gel & Clot Activator	😑 Yellow	16×100	8mL	GLASS	100PCS/1200PCS	Paper/Transparent
RJ1905309-0P	Gel & Clot Activator	😑 Yellow	16×100	9mL	PET	100PCS/1200PCS	Paper/Transparent
RJ1905309-0G	Gel & Clot Activator	😑 Yellow	16×100	9mL	GLASS	100PCS/1200PCS	Paper/Transparent
RJ1905310-0P	Gel & Clot Activator	😑 Yellow	16×100	10mL	PET	100PCS/1200PCS	Paper/Transparent
RJ1905310-0G	Gel & Clot Activator	😑 Yellow	16x100	10mL	GLASS	100PCS/1200PCS	Paper/Transparent