



The purpose of this guide is to show our NIV mask range with examples of breathing systems and accessories that can be used with each type of mask. Always refer to the user guide supplied with your NIV device for the specific set-up required.

Vented NIV masks with an anti-asphyxiation valve

Vented masks have integrated holes or leak ports in the mask to remove exhaled CO₂. Vented masks can also be referred to as ported masks or masks with CO₂ leak port. These masks have an anti-asphyxiation valve that opens in the case of an unexpected ventilator failure, or breathing system disconnection, allowing the patient to breathe spontaneously from the atmosphere.

The masks shown on this page are vented masks, with an anti-asphyxiation valve. This type of mask can be used with single limb breathing systems which do not have a built-in CO₂ leak port.

Examples of Intersurgical breathing systems that are appropriate for use with this type of mask are shown on page 3, opposite together with filter accessories that may be required.

VariFit NIV vented mask with anti-asphyxiation valve

For use with single limb breathing systems that do not have a CO₂ leak port.


Make an enquiry

Video available at www.intersurgical.com



Clear elbow ☐



Code	Description	Colour coded elbow	Box Qty.
2260000	VariFit NIV vented mask with anti-asphyxiation valve, small adult	<input type="radio"/>	8
 2261000	VariFit NIV vented mask with anti-asphyxiation valve, medium adult	<input type="radio"/>	8
2262000	VariFit NIV vented mask with anti-asphyxiation valve, large adult	<input type="radio"/>	8
2277000	VariFit NIV vented mask with anti-asphyxiation valve, extra large adult	<input type="radio"/>	8


FaceFit NIV Ported Mask - with anti-asphyxiation valve and CO₂ leak port

For use with single limb breathing systems that do not have a CO₂ leak port.

Make an enquiry

Video available at www.intersurgical.com



Code	Description	Box Qty.
2252000	FaceFit, NIV Ported Mask, with anti-asphyxiation valve and CO ₂ leak port, large adult	12
 2251000	FaceFit, NIV Ported Mask, with anti-asphyxiation valve and CO ₂ leak port, medium adult	12
2250000	FaceFit, NIV Ported Mask, with anti-asphyxiation valve and CO ₂ leak port, small adult	12