

## GIMA ULTRASOUND GEL



Multilingual label: GB, FR, IT, ES, PT

| GIMA code | ULTRASOUND GEL                                  | Minimum order |
|-----------|-------------------------------------------------|---------------|
| 33269     | Sterile ultrasound gel - sachet - 20ml - transp | box of 48     |
| 33271     | Ultrasound gel - sachet - 20ml - transparent    | box of 200    |
| 33270     | Ultrasound gel - tube 60 ml - transparent       | box of 12     |
| 33272     | Ultrasound gel - tube 250 ml - blue             | box of 40     |
| 33273     | Ultrasound gel - tube 250 ml - transparent      | box of 40     |
| 33288     | Ultrasound gel - bottle 1 l - blue              | box of 18     |
| 33289     | Ultrasound gel - bottle 1 l - transparent       | box of 18     |
| 33276     | Ultrasound gel - tank 5 l - blue                | box of 2*     |
| 33277     | Ultrasound gel - tank 5 l - transparent         | box of 2*     |
| 33286     | Ultrasound gel - bag 5 l - blue                 | box of 2**    |
| 33287     | Ultrasound gel - bag 5 l - transparent          | box of 2**    |
| 33283     | Pump for tank 5 l and bag 5 l                   |               |

\*Includes 2 empty bottles 250 ml \*\*Includes 1 empty bottle 250 ml and 1 nozzle

### GIMA ULTRASOUND GEL

Hypoallergenic, water soluble gel suitable for ultrasonography, doppler and defibrillator applications. It ensures the ultrasound waves come to the equipment display in a clear and uninterrupted way. Salt, formaldehyde, fragrance free.



Sterile

## GEL CLEANER

### • 33281 GEL PROBE CLEANER 250 ml

It removes the gel left over the probe, makes the image clearer and extends the probe life. Thanks to the active agents contained, it prevents the bacteria to pass from patient to patient, and it ensures hygiene. Alcohol-free.



33281



## COSMETIC GEL

33252

33253



33256

CE  
FDA

## TECARTHERAPY CREAM

### • 33251 TECARTHERAPY CREAM 1 l

It allows to reduce the contact impedance between the probe and the skin. The high conductivity ensures the best results both in the capacitive phase and in the resistive phase of the treatment. Preservatives and dyes used are non-toxic. Latex free.

33251



## EEG PASTE

### • 33262 EEG PASTE - 400 g

Used with EEG electrodes in EEG application for brain mapping. Thanks to the active substances in its content, it allows the electrode to cling to the strand of hair easily and provides conductivity. It is hypoallergenic, does not irritate skin, soluble in water and can be cleaned easily. Does not contain oil and oily substances. It has no toxic effects and odour.



33262

### KONIX COSMETIC LASER & CAVITATION GEL

Laser & Cavitation gel reduces the pain that patients experience during depilation and makes the device last longer by cooling the device head and window. Thanks to the vacuum technology, its active substances and its high viscosity preserve probe from damages. It is not affected by the body salt during the use, ensures an easy application.

**Technical properties:** it is hypoallergenic, does not irritate skin, is water soluble and easily cleanable.

Oily substances, formaldehyde and salt free. No toxic effects and odourless.

| GIMA code | COSMETIC LASER & CAVITATION GEL | Minimum order |
|-----------|---------------------------------|---------------|
| 33253     | Cosmetic gel - tube 250 ml      | box of 40     |
| 33252     | Cosmetic gel - bottle 1 l       | box of 18     |
| 33256     | Cosmetic gel - bag 5 l          | box of 2**    |
| 33283     | Pump for bag 5 l                |               |

\*\*Includes 1 empty bottle 250 ml and 1 nozzle for easy refill