


Fastsorb

Surgical Suture, Polyglycolicacid (PGA)
Fast Absorbable, Synthetic, Braided, Undyed

Surgical Sutures
Absorbable



Boz®

Definition

Synthetic absorbable, coated, short term wound support and short term mass absorption.

Material

Polyglycolic Acid (PGA)

Coating

Calcium Stearate and polycaprolactone

Structure

Braided

Colour

Beige (undyed)

Size

EP : 0.4 to 8

USP : 8/0 to 6

Tensile Strength Retention

7 Days 75 %

Mass Absorption

42 days

Sterilization Method

Gamma Irradiation

Main Indication

Gynaecology (episiotomies)

Urology

Ophtalmic surgery (conjunctive suturing)

Odontology (oral mucosa)






















Pediatric surgery

Skin closer - intracutaneous, subcutaneous.

Ligatures

















Fastorb

Surgical Suture, Polyglycolicacid (PGA) Fast Absorbable, Synthetic, Braided, Undyed

Needle Length and description			Product No	USP EP	Length (cm)	Unit Pcs	Color
12 mm							
3/8 Reverse Cutting 330 Micron			FS1545	6/0 0,7	75	12	Undyed
13 mm							
1/2 Round Body, Taper Point 330 Micron			FS1563	6/0 0,7	75	12	Undyed
1/2 Round Body, Taper Point 380 Micron			FS1886	5/0 1,0	75	12	Undyed
3/8 Reverse Cutting 330 Micron			FS5912	6/0 0,7	75	12	Undyed
3/8 Reverse Cutting 380 Micron			FS1929	5/0 1,0	75	12	Undyed
3/8 Reverse Cutting Premium 380 Micron			FS1941	5/0 1,0	45	12	Undyed
16 mm							
1/2 Round Body, Taper Point			FS1989	5/0 1,0	75	12	Undyed
			FS2503	4/0 1,5	75	12	Undyed
3/8 Reverse Cutting			FS1673	6/0 0,7	75	12	Undyed
			FS2019	5/0 1,0	75	12	Undyed
			FS2539	4/0 1,5	75	12	Undyed
3/8 Inside Cutting			FS2557	4/0 1,5	75	12	Undyed
			FS3038	3/0 2,0	75	12	Undyed
19 mm							
3/8 Reverse Cutting			FS2639	4/0 1,5	75	12	Undyed
			FS3108	3/0 2,0	75	12	Undyed
20 mm							
1/2 Round Body, Taper Point			FS2228	5/0 1,0	75	12	Undyed






Fastorb

Surgical Suture, Polyglycolicacid (PGA) Fast Absorbable, Synthetic, Braided, Undyed

Needle Length and description			Product No	USP EP	Length (cm)	Unit Pcs	Color
20 mm							
1/2 Round Body, Taper Point			FS2683	4/0 1,5	75	12	Undyed
			FS3141	3/0 2,0	75	12	Undyed
3/8 Reverse Cutting			FS2250	5/0 1,0	75	12	Undyed
22 mm							
1/2 Round Body, Taper Point			FS2745	4/0 1,5	75	12	Undyed
25 mm							
1/2 Round Body, Taper Point			FS2787	4/0 1,5	75	12	Undyed
			FS3251	3/0 2,0	75	12	Undyed
3/8 Reverse Cutting			FS3273	3/0 2,0	75	12	Undyed
			FS3756	2/0 3,0	75	12	Undyed
26 mm							
1/2 Round Body, Taper Point			FS3785	2/0 3,0	75	12	Undyed
			FS4214	0 3,5	75	12	Undyed
3/8 Reverse Cutting			FS3821	2/0 3,0	75	12	Undyed
			FS4236	0 3,5	75	12	Undyed
			FS3345	3/0 2,0	75	12	Undyed
30 mm							
1/2 Round Body, Taper Point			FS3368	3/0 2,0	75	12	Undyed
			FS3857	2/0 3,0	75	12	Undyed
			FS4243	0 3,5	75	12	Undyed
			FS4607	1 4,0	75	12	Undyed







Fastorb

Surgical Suture, Polyglycolicacid (PGA) Fast Absorbable, Synthetic, Braided, Undyed

Needle Length and description		Product No	USP EP	Length (cm)	Unit Pcs	Color
30 mm						
3/8 Reverse Cutting		FS2885	4/0 1,5	75	12	Undyed
		FS3396	3/0 2,0	75	12	Undyed
		FS4628	1 4,0	75	12	Undyed
35 mm						
1/2 Round Body, Taper Point		FS3915	2/0 3,0	75	12	Undyed
		FS4299	0 3,5	75	12	Undyed
		FS4305	0 3,5	90	12	Undyed
		FS4643	1 4,0	75	12	Undyed
		FS4648	1 4,0	90	12	Undyed
3/8 Reverse Cutting		FS3422	3/0 2,0	75	12	Undyed
		FS3932	2/0 3,0	75	12	Undyed
37 mm						
1/2 Round Body, Taper Point		FS3954	2/0 3,0	90	12	Undyed
40 mm						
1/2 Round Body, Taper Point		FS3981	2/0 3,0	75	12	Undyed
		FS3988	2/0 3,0	90	12	Undyed
		FS4367	0 3,5	75	12	Undyed
		FS4374	0 3,5	90	12	Undyed

Fastsorb

Surgical Suture, Polyglycolicacid (PGA) Fast Absorbable, Synthetic, Braided, Undyed

Needle Length and description	Product No	USP EP	Length (cm)	Unit Pcs	Color
40 mm					
1/2 Round Body, Taper Point 	FS4712	1 4,0	75	12	Undyed
	FS4718	1 4,0	90	12	Undyed
	FS5001	2 5,0	75	12	Undyed
3/8 Reverse Cutting 	FS4005	2/0 3,0	75	12	Undyed
45 mm					
1/2 Round Body, Taper Point 	FS4421	0 3,5	75	12	Undyed
	FS4771	1 4,0	75	12	Undyed
	FS5034	2 5,0	75	12	Undyed
Straight Reverse Cutting 	FS2911	4/0 1,5	75	12	Undyed
50 mm					
1/2 Round Body, Taper Point 	FS4037	2/0 3,0	75	12	Undyed
	FS4453	0 3,5	75	12	Undyed
	FS4811	1 4,0	75	12	Undyed
	FS5065	2 5,0	75	12	Undyed
65 mm					
1/2 Round Body, Taper Point 	FS4864	1 4,0	100	12	Undyed

Fastsorb

Surgical Suture, Polyglycolicacid (PGA) Fast Absorbable, Synthetic, Braided, Undyed

Non-Needle Suture	Product No	USP EP	Length (cm)	Unit Pcs	Color
150 cm					
FASTSORB Non-Needle Suture	FS1373	6/0 0,7	150	12	Undyed
	FS2325	4/0 1,5	150	12	Undyed
	FS2934	3/0 2,0	150	12	Undyed
	FS3510	2/0 3,0	150	12	Undyed
	FS4113	0 3,5	150	12	Undyed
	FS4531	1 4,0	150	12	Undyed
	FS4907	2 5,0	150	12	Undyed