

Anti Decubitus Bedsore Alternating Pressure Bubble Mattress (QDC-303)



• Product Description

Anti Decubitus Bedsore Alternating Pressure Bubble Mattress QDC-303 is designed for the treatment and prevention of pressure sore stage I in short-term and domiciliary therapies. The mattress in PVC is composed of 130 7cm-high bubble cells, particularly comfortable, it is fixed to the bed by extra flaps on both top and bottom side. Very easy to install, maintain and use thanks to the possibility of adjust comfort range. The pump can be hung to the end of the bed by means of two hooks.

Technical Parameters

Technical Parameters	
Bubble Mattress	
Dimensions without flaps (L×W×H)	200×90×7cm
Extension Flaps	Head & Foot Side
Mattress material	Odourless Medical-Grade PVC (Rohs / Flame Resistance / EN-71)
Material thickness	0.32 mm
Cold-resistamt	-30°C
Number of cells	130
Weight support	135 kg
Pump	
Electrical	AC 110V/220V 50Hz/60Hz
Pressure range	40-100 mmHg
Air output	7-8 L/min

Guangdong Yuehua Medical Instrument Factory Co., Ltd.

Add: Rongsheng Science and Technology Zone, Daxue Road, Shantou, Guang;long, China.

Tel: +86 0754-82525836