


Product Data Sheet



CT Gantry Cover

Ref.:	Product name:	Dimensions:	Weight:	Packaging dimensions:
1231-1	CT Gantry Cover C33	1190 x 1140 x 10 mm	700 g	750 x 450 x 60mm
Intended use		Reusable hygienic CT gantry protection during invasive and non-invasive procedures to keep the CT gantry protected from liquids. Covers the lower half of the CT gantry. May only be used by medical personnel. Non-sterile and not to be used in the sterile field..		
Directives & classification		CE according to MDR – Regulation 2017/745 Class 1 product		
Applied quality assurance		SN EN ISO 13485:2016		
Instructions for use		Enclosed to the product		
Restrictions of use		Do not use the product for purposes other than those described.		
General safety notes		Non sterile. Do not sterilize. Dispose of any product that is cracked or appears damaged Caution: Federal law restricts this device to sale by or on the order of a licensed health-care professional (Rx only).		
Materials		Material TPU PA LDPE Cardboard	Component Foil Velcro packaging packaging	
Marking		CE		
Storage		Store in dry place / room temperature		

Figures are subject to usual tolerances. All information to the best of our current knowledge. They do not guarantee a warranted quality and are no foundation for a contractual relationship. Errors excepted.