

Item Code 910004

INTERNATIONAL  
PATENT

Volume Accuracy  
 $\pm 2,5\%$  or 50 ml



# FlowMIR™

Disposable Turbine Flowmeter  
MIR exclusive product  
International patent

*Each turbine, which includes a cardboard mouthpiece, has been individually factory tested with a computerized system*

*For each patient, after the Spirometry test, both turbine and mouthpiece are thrown away*



FlowMIR Dispensers (10 and 60 pcs)



Comfortable  
Packaging



Singularly Tested  
and Packed



Eliminates staff  
clean-up time



No Sterilisation



No Calibration



No Cross  
Contamination



No Antibacterial  
Filter



Always 100%  
accurate and hygienic



ATS/ERS Compliant  
for accuracy



Not affected by  
Vapour Condensation



Not affected by  
Ambient conditions



The Best Sensor  
for Spirometry

## Consumables / Accessories



### Paper roll for printer

Code 910350

MIR Model Spirolab

Standard thermal paper roll for printer.



### Plastic noseclip

Code 910320

MIR Model All Spiro

Plastic noseclip.



### Calibration syringe 3L

Code 919000

MIR Model All Spiro

3L calibration syringe.



### Oximetry finger sensor

Code 919024

MIR Model All devices with oximetry option

Reusable, ADULT.



### Oximetry finger sensor

Code 939005

MIR Model All devices with oximetry option

Reusable, PAEDIATRIC 5-45 kg.