
## Flextube™ breathing systems for active humidification with heated wire

Available with single or dual heated wires and an auto fill humidification chamber if required. All Intersurgical active humidification breathing systems are suitable for use with the majority of humidification base units.

Silver Knight™ breathing systems include anti-microbial protection designed to help slow down the proliferation of micro-organisms and reduce cross infection within the Intensive Care environment. On heated wire Silver Knight breathing systems, both the green inspiratory and the lilac expiratory limbs include the anti-microbial additive for added protection against healthcare acquired infections.

Please see page 5 for a guide to our tube types.



## **Dual heated wire**

Code	Description	Box Oty.	
2026310 (S*) 🚳	22mm Flextube dual heated wire breathing system with auto-fill chamber and extra limb. ≥ 1.6m	7 (6*)	
2026131	22mm Flextube - Silver Knight dual heated wire breathing system with auto-fill chamber and extra limb. ≥ 1.6m	7	•
4525850	15mm Flextube dual heated wire breathing system with extra limb. ≥ 1.6m	15	
6392810	10mm Flextube dual heated wire breathing system with auto-fill chamber and extra limb. ≥ 1.6m	7	

## Single heated wire

Code	Description	Box Oty.	
2025310 (S*)	22mm Flextube single heated wire breathing system with water trap, auto-fill chamber and extra limb. ≥ 1.6m	7 (5*)	
2025131	22mm Flextube - Silver Knight single heated wire breathing system with water trap, auto-fill chamber and extra limb. ≥ 1.6m	7	(A)
4504810	15mm Flextube single heated wire breathing system with water trap, auto-fill chamber and extra limb. ≥ 1.2m	7	
6091810	10mm Flextube single heated wire breathing system with neo water trap, auto-fill chamber and extra limb. ≥ 1.6m	20	

## Single limb heated wire

Code	Description	Box Oty.
2060310	22mm Flextube single heated wire limb with auto-fill chamber and extra limb. ≥ 1.6m	7
4560810 <i>NEW</i>	15mm Flextube single heated wire limb with auto-fill chamber and extra limb, ≥ 1.2m	7

(S\*) Add an S to the seven digit code number for the sterile version of this product eg. 2026310S (sterile box quantity is shown in brackets).

