Republica Moldova mun. Chişinău, MD-2005

str. Albișoara 64/2 tel/fax.: (373-22) 54-73-73 tel./fax: (373-22) 02-44-00



Rechizitele bancare:
Cod fiscal/1003600117582
Cod TVA: /0205086
BC "Moldova-Agroindbank"SA
filiala M.Eminescu
cod: AGRNMD2X864
Cod IBAN:MD14AG000000225184801542

URYXXON® Relax

URYXXON® Relax – automated urine analysis at the point of care



The URYXXON® Relax provides dependable urine status results to detect early stages of many diseases such as diabetes, kidney disease and urinary tract infections. Instrument-read results have long proven to be advantageous for both, busy health care professionals and patients. URYXXON®Relax readings eliminate the subjectivity of visual colour interpretation. The comprehensive interface options and the optimized printouts minimize risks associated with manual transcriptions. Reliable results can be obtained virtually immediate at the point of care.

The $\mathsf{URYXXON}^{\ensuremath{\mbox{\scriptsize B}}}$ Relax makes urine analysis simpler and more reliable.

Technical data and ordering information

reclinical data and ordering information	
Test strips	URYXXON®Stick 10
Capacity	50 strips/h
Instrument memory	200 patient test results including name or patient ID
Interface	User: Touch screen display, alphanumeric input, password protection, Computer: USB interface for connection to PC alternatively RS232 interface for connection to PC and PS/2 interface for connection of keyboard and/or barcode reader
Power requirements	110-240 V AC, automatic, Battery powered operation (optional) with 6 AA batteries
Dimensions	Depth: 20 cm (8 inches) Width: 16 cm (6 inches) Height: 7.5 cm (3 inches)
Weight	Weight: 710 g (1.90 lb) (without batteries and power supply)
Operation	Temperature range: 5-40 °C (41 °F-104 °F), Humidity range: 20-80% relative humidity, noncondensing, Calibration: automatic, self-calibrating
Reference	93088

