





ST-SU07

Safety Lancet

Ergonomic design offers a more secure hold when to activate the device

Tri-bevel needle tip ensures the sharpness and

The spring buffer the puncture making the procedure a more comfortable experience

The complete range of needle sizes from 18G to 32G ensure the accurate

Blood samples for all popular meters

The needle retracts immediately after use Sterilization: Gamma

Size: 21G, 23G, 26G





ST-SU12

Twist Lancets

Conventional twist-off cap design Sterilized by gamma-radiation Smooth tri-bevel needle tip for comfortable sampling experience Versitile with most lancing devices

Size: 21G, 23G, 26G, 28G, 30G Packing: 100pcs/Box, 20000pcs/ctn





ST-SU08

Pressure Activated Safety Lancets

Needle is fully shielded before and after use Pre-loaded, pressure activated design Ensures a good blood sample every time
Ultra sharp tri-bevel needle or bi-bevel blade,
high speed penetration Variety of needle specifications to suit most

capillary blood requirements Sterilization: Gamma Size: 21G, 23G, 26G



ST-SU13

Flat Lancets (Exclusive to Roche Diagnostics)

5 twist lancets (28G) are used with Roche Accu-Chek Softclix Lancing Device and the needle tip is siliconized for easy and smooth motion during lancing

Size: 28G, 30G

Packing: 100pcs/Box, 20000pcs/ctn





Button Activated Safety Lancets

Needle tip is sealed with cap through injection

Pre-loaded, button operated design Single-use only, sterilised by gamma irradiation,

cannot be reused High quality, ultra sharp tri-bevel needle, high speed penetration

Variety of needle specifications to suit most capillary

blood requirements Sterilization: Gamma Size: 21G, 23G, 26G



ST-SU14

Disposable Electrosurgical Pencil

170mm-ABS pencil appearance 7 cm stainless steel tip 2.36mm-diameter of the tip 3m-PVC medical cable

Maximal output current: 1A Frequency confine: 0-1.0MHZ

1pc/blister bag



ST-SU10

Lancing Device (Releasable)



ST-SU11

Lancing Device

The user can install the lancet with the device easily. And the device will work with the lancet

Smoothly without any blockage.

The lancet won't fall on its own when to fire the device. Retracting automatically after firing, the launcher is safe and reliable.







ST-SU15

Surgical Hand Brush

With iodine or without iodine Packing: 1/mold bag, 20 bags/box, 10boxes/ctn

ST-SU16

Surgical Hand Brush

Material: PP+Nylon



www.sintrue.com