

HMS Plus hemostasis management system

HANS PLUS

1111



Ordering Information

HMS Plus Test Cartridges and Liquid Controls

Disposable cartridges undergo strict manufacturing standards and quality control to provide accurate, dependable performance. Cartridge room temperature and/or refrigerated shelf life are stamped on each box.

Heparin Assay Cartridges and Controls

| Cartridges: Each b | Cartridges: Each box contains 9 cartridges with syringes and blunt tip needles | | | Controls: Each box contains 10 vials of control and deionized water | |
|--------------------|--|---------------|---------------|---|---|
| Catalog # | Description | # of Channels | Heparin Level | Control Catalog # | Description - (10 vials of control and 10 vials of deionized water for reconstitution) |
| 304-01POR | Red | 4 | 0.0-0.9 mg/kg | 206 01000 | Red/Yellow |
| 304-02POR | Yellow | 4 | 0.0-1.5 mg/kg | 306-01POR | |
| 304-03POR | Tan | 4 | 1.5-3.0 mg/kg | 306-02POR | Tan/Silver |
| 304-04POR | Silver | 4 | 2.0-3.5 mg/kg | | |
| 304-05POR | Blue | 4 | 2.5-4.0 mg/kg | 306-03POR | Blue/Gold |
| 304-06POR | Green | 4 | 3.5-5.0 mg/kg | 306-04POR | Green/White |
| 304-07POR | Orange | 6 | 0.0-2.5 mg/kg | 306-05POR | Orange |
| 304-08POR | Gold | 6 | 1.5-4.0 mg/kg | 306-03POR | Blue/Gold |
| 304-09POR | White | 6 | 2.5-5.0 mg/kg | 306-04POR | Green/White |
| 304-10POR | Purple | 4 | 4.5-6.0 mg/kg | 206 00000 | Duran la /Dia ala |
| 304-11POR | Black | 6 | 3.5-6.0 mg/kg | 306-09POR | Purple/Black |

Heparin Dose Response Cartridge

| Cartridges: Each box contains 9 cartridges, syringes, and blunt tip needles | | | |
|--|-------------|------------------|--|
| Catalog # | Description | # of Channels | |
| 304-20POR | HDR | 6 | |

High Range ACT

| Cartridges: Each box contains 18 cartridges and 9 syringes, and blunt tip needles | | Control: Each box contains 15 vials of control and deionized water | | |
|---|-------------|--|----------------------|--|
| Catalog # | Description | # of Channels | Control Catalog # | Description |
| 304-30 | HRACT | 2 | 550-07 | CLOTtrac HR Coagulation Control |
| | | | 550-08 | CLOTtrac HR Abnormal Coagulation Control |
| | | | 550-13 | CLOTtrac HR Control Pak Coagulation Control (1 box each of 550-07 and 550-08) |

HMS Plus Instruments

Each instrument includes an internal printer, one-year warranty, operating manual, and case analysis pad.

| Catalog # | | Description | |
|-----------|----------|--------------|---------|
| 30514 | HMS PLUS | 100-120 volt | |
| 30522 | HMS PLUS | 200-240 volt | English |
| 30515 | HMS PLUS | 200-240 volt | German |
| 30517 | HMS PLUS | 200-240 volt | French |
| 30518 | HMS PLUS | 200-240 volt | Italian |
| 30524 | HMS PLUS | 200-240 volt | Spanish |
| 30527 | HMS PLUS | 200-240 volt | Dutch |

HEPline Kit

Kit contains 5 tubes of increasing amounts of heparin concentrations designed to give final concentrations of 1, 2, 3, 4, and 5 U/ml when fresh whole blood is added to a final volume of 5 ml.

| Catalog # | Description | |
|-----------|-------------|--|
| 313-50 | HEPline Kit | |

Accessories

| Catalog # | Description | Quantity |
|-------------|------------------------------------|-------------------|
| HMSPLUSSC | Bar Code Scanner | 1 |
| HMSPUSSCRS | Bar Code Scanner - European Union | 1 |
| HMSPLUSSCYY | Bar Code Scanner - Japan | 1 |
| HMSPLUSEDM | External Data Manager | 1 |
| 31351 | Electronic Quality Control | 1 |
| 300-01 | 3cc Monoject Syringes | 100 per Box |
| 300-02 | Blunt Needles, 1-7/16", 19 GA | 100 per Box |
| 300-04 | Thermal Printer Paper | 5 Rolls per Box |
| 300-05 | HMS Case Analysis Pad | 50 Sheets per Pad |
| 300-10 | Temperature Verification Cartridge | 1 |
| 313-18 | QA Records Packet | 1 |
| 31506 | Salvage Reservoir Cups | 100 per Box |
| 30032 | HMS PLUS Custom Cart | 1 |

Manuals

Manuals can be ordered through the CardioVascular Service Department.

| Part Number | Description | |
|-------------|-------------------|--|
| 86506001 | Operator's Manual | |