

Breathing System - Silicone Breathing Circuit (Ventilation)

* Heated Wire Adapter port



* Heated Wire Adapter port connected to heater wire cable



Heated wire breathing system accessories

Reusable heated wire adapter

	002-018-1120	1.2m, 22M for tube, 22F for chamber, 7.6F for temperature probe, LEMO 2 pin
	002-018-1160	1.6m, 22M for tube, 22F for chamber, 7.6F for temperature probe, LEMO 2 pin

Reusable humidifier chamber

	002-023-001	PC, Reusable, 22M-22F, Max. capacity 700ml; Compress volumn - Full water 500ml; Compress volumn - Empty 200ml
--	-------------	---

Reusable humidifier chamber

	002-023-002	PC, Reusable, 22M-15F, Max. capacity 360ml, Compress volumn - Full water 215ml, Compress volumn - Empty 145ml
--	-------------	---

Temperature probe

	002-025-002	Temperature probe, 1 to 2
--	-------------	---------------------------

Adaptor cable heater wire

	002-024-003	Single heater wire cable, 1 to 2, LEMO 2 pin and LEMO 3 pin (for single-use heated wire)
--	-------------	--

Adaptor cable heater wire

	002-024-004	Reusable single heater wire cable, 1 to 1, LEMO 2 pin (for reusable heated wire)
--	-------------	--

Draw wire

	002-037-001	2 meter
--	-------------	---------

Adaptor cable heater wire

	002-024-005	Reusable single heater wire cable, 1 to 1, LEMO 3 pin (for single-use heated wire)
--	-------------	--



For every child,
the right to **SURVIVE**

