


Anaesthetic Face Masks



The air cushion with soft and good adhesion is closely matched with different facial contours to reduce gas leakage during assisted breathing. Through the transparent mask, it is easy to observe vomiting and aspiration as soon as possible, and is suitable for surgical anaesthesia.

The air cushion valve can change air cushion pressure, which is convenient for quick adjustment. The specifications are 0#~5#, and the sizes are various, which can be matched with the faces of patients of different ages. 2#~5# interfaces have an inside diameter of 22mm, with fixing rings, and the 0# and 1# interfaces have an outside diameter of 15mm, without fixing rings.

1#



2#



3#



4#



5#



Anaesthetic face masks

Code	Description	Package
EM01-101	Anaesthetic face mask, infant 1#	50pcs/Box
EM01-102	Anaesthetic face mask, pediatric 2#	50pcs/Box
EM01-103	Anaesthetic face mask, large pediatric 3#	50pcs/Box
EM01-104	Anaesthetic face mask, small adult 4#	50pcs/Box
EM01-105	Anaesthetic face mask, adult 5#	50pcs/Box