- Alternating Therapy kit mattress and pump
- F Kit matelas anti-escarres et compresseur à cycle alterné
- Kit colchón antiescaras y compresor de aire alternante
- Kit colchão anti-escaras e compressor alterno

Interchangeable cells mattress and pump with pressure regulation







Hooks to hang the compressor to the bed

Transport function

with cap

Suitable for prevention and treatment of decubitus sores STAGE III. The first three elements are static for pillow function; the remainder is divided into 2 sections with 2-way alternating operation; complete with padded cover in biocompatible, waterproof, perspiring Nylon/PU, tested according to European fire reaction standards and equipped with zip on all four sides for easy removal and cleaning. Mattress fixed to the bed with supplied straps and equipped with grommet eyelets. Compressor is equipped with manual pressure regulation, visual alarm and sound for low pressure and interchangeable air filter. The STATIC function blocks the alternating system and effectively makes it a normal mattress. Mattress/compressor connection through fast connectors with cap, in case of patient transport.

Technical features

PIUMA UP 3 SERIES		
Product codes	LAD470	LAD471
SPECIFICATIONS		
Compressor		
Power	220V / 50 Hz	
Pressure	55~120mmHg	
Power cable	5 m	
Air release	7.5/8.5LPM	
Cycle time	12 min	
Fuse	T1A 250V	
Noise	30-35 dB	
Mattress		
Material	· ·	or LAD470 mattress) or LAD471 mattress)
Elements no.	17 (of which 8	with micro-holes)
CPR valve	Deflation time	around 15 sec.

	W:	D:	1 H:		
Compressor	30 cm	13,5 cm	11 cm	-	
Mattress while functioning	196 cm	86 cm	13 cm	160 kg	
Spare parts			Accessories		
Cod I AD722	Cod. LAR700 Nylon + PU co	over	Cod. LAR702 Bielastic, antimic waterproof, persp verified according European fire rea	piring and g to the	