

MW2A



## Male Catheterization/Enema Simulator

MW52



## Female Catheterization/Enema Simulator II



MW2A



MW52

### High quality simulators for urethral catheterization, perineal care and enema practice

Achieve better training results with our soft and realistic catheterization simulator that can also be used on SP or full body manikins.

#### FEATURES

- 1 | Soft and flexible genitalia allow proper cleaning and procedural preparation
- 2 | Successful insertion confirmation by fluid outflow
- 3 | Enema insertion in lateral position
- 4 | Detachable external genitalia

#### SKILLS

- |                                 |                                  |
|---------------------------------|----------------------------------|
| Intermittent catheterization    | Cleansing and aseptic technique  |
| Indwelling catheterization      | Catheter insertion and placement |
| Perineal care                   | Catheter Removal                 |
| Enema                           | Manual bladder compression       |
| Localization of urethral meatus |                                  |

Compact sets consist of a genitalia unit and accessories, for hybrid training and table top demonstration.

**MW2C** 11381-199-B

Male Genital Unit only



**MW2D** 11433-299-B

Female Genital Unit only



Enema Training



Pass of catheter can be observed through transparent bladder reservoir



### DESCRIPTIONS

#### SPECIFICATIONS

Size: 34 x 33 x 24cm / 13.5 x 13 x 9.5 inch  
Weight: 3.5kg / 7.7lbs(M) 3.4kg / 7.5lbs(F)

#### SET INCLUDES

- 1 lower torso manikin
- 1 genital organ unit
- 3 valves (one is incorporated in the assembled unit)
- 1 irrigation bottle

#### MATERIALS

Soft resin  
Latex free

- 1 rectum tube
- 3 irrigation drainage tubes
- 1 lubricant
- 1 storage case
- 1 instruction manual

#### REPLACEMENT PARTS

- |              |                                         |
|--------------|-----------------------------------------|
| 11381-100-01 | clip (small) for urethral tube for MW2A |
| 11381-100-02 | clip (large) for MW2A/MW2B/MW2C         |
| 11433-000-01 | clip (large) for MW52                   |
| 11381-100-03 | urinary bladder reservoir               |
| 11381-100-04 | rectum unit                             |

- |           |                             |
|-----------|-----------------------------|
| 11381-300 | 6 replacement valves        |
| 11433-010 | 1 replacement female skin   |
| 11433-020 | 2 replacement female meatus |
| 11381-800 | lubricant                   |
| 11381-410 | 1 replacement male skin     |
| 11381-510 | 2 replacement male meatus   |

#### RECOMMENDED DEVICES

\* Catheters are not included.