

UNIVERSAL SUTURES

"HEALING BEYOND COMFORT"




UNISUR
Lifecare Pvt. Ltd.


Sutures
Mesh
Bone Wax

AN ISO & CE CERTIFIED SURGICAL SUTURES AND MESH MANUFACTURING COMPANY
We Manufacture complete range of Surgical Sutures, Hernia Mesh & Bonewax

Needle Shape



Needle Type:




International Coding Pattern

International Coding Pattern											
Code 1	Product	Code 2	USP Size	Code 3	Code 4	Needle Curvature	Code 5	Needle Type	Code 6		
1	Braided & Coated Polyglycolic Acid Undyed	A	3	SUTURE LENGTH in cm	C	3/8 Circle	Z	Round Body / Taper Point	NEEDLE LENGTH in mm		
2	Braided & Coated Polyglycolic Acid	B	2		D	1/2 Circle	Y	Cutting			
3	Monofilament Polyamide (Nylon)	C	1		A	Straight	X	Reverse Cutting			
4	Catgut Chromic	D	0		E	5/8 Circle	K	Taper Cutting			
5	Black Braided Silk	E	2-0					T		Trocar Point	
6	Braided Polyester	F	3-0					R		Blunt Point	
7	Catgut Plain	G	4-0					S		Spatula	
8	Monofilament Polypropylene	H	5-0								
9	Monofilament Polydioxanone	I	6-0								
10	Braided & Coated Polyglactin 910	J	7-0								
11	Monofilament Polyglecaprone 25	K	8-0								
12	Braided & Coated Polyglactin 910 Undyed										

Example:

Braided & Coated Polyglycolic Acid	Size 1	75cm	1/2 Circle	Round Body	40mm
------------------------------------	--------	------	------------	------------	------



UNILON

Monofilament Polyamide

Non Absorbable Surgical Suture U.S.P.



UNILON – Monofilament Polyamide (Nylon) is a non absorbable sterile surgical suture, composed of the long-chain aliphatic homo polymers Nylon 6 and Nylon 6,6.

UNILON Monofilament Polyamide Sutures are dyed black or blue to enhance visibility in tissue.

In-house needle manufacturing provides you the better needles with optimum sharpness which finally delivers smooth riding on tissues during suturing.

Indications:

Monofilament Polyamide Suture is indicated for use in general soft tissue approximation and/or ligation, including use in cardiovascular, ophthalmic and neurological procedures.

- The best material for skin approximation
- Being a monofilament with smooth surface never support bacterial growth followed by capillary action
- Uniform elasticity supports small movements results cosmetic effect
- Ultra smooth surface results negligible pain

Best Used By:

- General Surgeons,
- Gynaecologists,
- Orthopaedics

Ultimate Suture For:

- Skin Closure
- Hard Tissues approximation
- Casualty stitching

Additional Information:

Suture Characteristics	:	Non Absorbable Surgical Suture
Type	:	Monofilament
Material	:	Polyamide (Nylon)
Color	:	Black
Absorption	:	Non Absorbable
USP Range	:	6/0 – 2
Sterilization	:	EO (Ethylene Oxide)
Shelf Life	:	5 Years
Product Characteristics	:	Excellent tensile strength, Excellent knot tying, Minimum tissue reaction.
Packaging	:	In food grade folder packed in a medical grade pouch, with or without needle
Needle Type	:	Round bodied, cutting edge, taper point, straight, blunt point
MFG. LIC No	:	MFG/MD/2020/000109
OEM. LIC No	:	N. Code: KA/DEVICE/MFG/MD/2020/000109

UNILON

Monofilament Polyamide

Non Absorbable Surgical Suture U.S.P.



Code No.	SUTURE		NEEDLE DESCRIPTION
	USP SIZE	LENGTH	
N 3303	2-0	75 cm	1/2 Circle Reverse Cutting 38mm
N 3304	0	75 cm	1/2 Circle Reverse Cutting 38mm
N 3304L	0	100 cm	1/2 Circle Reverse Cutting 38mm
N 3305	1	75 cm	1/2 Circle Reverse Cutting 38mm
N 3306	2	75 cm	1/2 Circle Reverse Cutting 38mm
N 3315	5-0	35 cm	3/8 Circle Cutting 16mm
N 3315RB	5-0	75 cm	3/8 Circle Round Bodied 16mm
N 3316	5-0	75 cm	3/8 Circle Round Bodied 12mm
N 3317	5-0	70 cm	3/8 Circle Reverse Cutting 12mm
N 3318	4-0	70 cm	3/8 Circle Cutting 16mm
N 3319RB	4-0	75 cm	1/2 Circle Round Bodied 19mm
N 3320N	6-0	75 cm	3/8 Circle Reverse Cutting 13mm
N 3320RB	6-0	75 cm	3/8 Circle Round Bodied 13mm
N 3321	3-0	70 cm	3/8 Circle Cutting 16mm
N 3326N	4-0	38 cm	3/8 Circle Cutting 10mm
N 3328	3-0	70 cm	3/8 Circle Reverse Cutting 26mm
N 3328T	3-0	70 cm	3/8 Circle Reverse Cutting 26mm
N 3328A	3-0	35 cm	3/8 Circle Reverse Cutting 26mm
N 3328AT	3-0	35 cm	3/8 Circle Reverse Cutting 26mm
N 3331	6-0	70 cm	3/8 Circle Cutting 26mm
N 3336	2-0	70 cm	3/8 Circle Reverse Cutting 45mm
N 3336T	2-0	70 cm	3/8 Circle Reverse Cutting 45mm
N 3336A	2-0	35 cm	3/8 Circle Reverse Cutting 45mm
N 3336AT	2-0	35 cm	3/8 Circle Reverse Cutting 45mm
N 3336DK	2-0	110 cm	3/8 Circle Reverse Cutting 45mm
N 3338	1	100 cm	1/2 Circle Reverse Cutting 45mm Heavy
N 3338A	1	45 cm	1/2 Circle Reverse Cutting 45mm Heavy
N 3340	0 (Loop)	150 cm	1/2 Circle Round Bodied 40mm
N 3341	2-0	75 cm	1/2 Circle Round Bodied 30mm
N 3346	0	100 cm	1/2 Circle Round Bodied 40mm Heavy

UNILON

Monofilament Polyamide
Non Absorbable Surgical Suture U.S.P.



Code No.	SUTURE		NEEDLE DESCRIPTION	
	USP SIZE	LENGTH		
3E75DZ26	Lot 171	2-0	75 cm	1/2 Circle Round Body 26cm
3E45CX26	Lot 172	2-0	45 cm	3/8 Circle Reverse Cutting 26cm
3F45CZ15	Lot 173	3-0	45 cm	3/8 Circle Round Body 15cm
3F75DZ16	Lot 174	3-0	75 cm	1/2 Circle Round Body 16cm
3G45CZ15	Lot 175	4-0	45 cm	3/8 Circle Round Body 15cm
3G75DZ22	Lot 176	4-0	75 cm	1/2 Circle Round Body 22cm
3H35DZ16	Lot 177	5-0	35 cm	1/2 Circle Round Body 16cm
3H45CX15	Lot 178	5-0	45 cm	3/8 Circle Reverse Cutting 15cm