



X-RAY

Weight & Size overview

Front aprons

Article code	Sizes	Weight 0.35 mm Leq	Weight 0.5 mm Leq	Chest size	Length
Front-S	Small	2.7 kg	3.6 kg	Up to 90	100
Front-M	Medium	3.1 kg	4.1 kg	Up to 100	100
Front-L	Large	3.4 kg	4.5 kg	Up to 110	100
Front-XL	Extra-large	3.7 kg	4.9 kg	Up to 120	105

UFront aprons

Article code	Sizes	Weight 0.35 mm Leq	Weight 0.5 mm Leq	Chest size	Length
UFront-S	Small	2.9 kg	3.9 kg	Up to 90	110*
UFront-M	Medium	3.3 kg	4.4 kg	Up to 100	115*
UFront-L	Large	3.6 kg	4.8 kg	Up to 110	120*
UFront-XL	Extra-large	3.9 kg	5.2 kg	Up to 120	120*

* Longer than our standard front apron in able to work with it while sitting down.

Coat aprons

Article code	Sizes male	Weight 0.5/0.25 mm Leq	Weight 0.7/0.35 mm Leq	Chest size	Hip	Length
Coat-A	Small	4.0 kg	5.6 kg	Up to 90	95-100	100
Coat-B	Medium	4.3 kg	6.0 kg	Up to 110	110-115	100
Coat-C	Large	4.4 kg	6.2 g	Up to 115	115-120	100
Coat-CXL	Extra-large	5.0 kg	7.0 kg	Up to 120	120-125	105
Article code	Sizes female	Weight 0.5/0.25 mm Leq	Weight 0.7/0.35 mm Leq	Chest size	Hip	Length
Coat-D	Small	3.8 kg	5.3 kg	Up to 85	90-95	100
Coat-E	Medium	4.1 kg	5.7 kg	Up to 95	95-100	100
Coat-F	Large	4.2 kg	5.9 kg	Up to 105	100-110	100
Coat-FXL	Extra-large	4.8 kg	6.7 kg	Up to 115	115-120	105

Vest & Kilt aprons

Article code	Sizes male	Weight 0.5/0.25 mm Leq	Weight 0.7/0.35 mm Leq	Chest size	Hip	Vest length	Kilt length
VK-A	Small	4.7 kg	6.5 kg	Up to 90	95-100	57	59
VK-B	Medium	6.3 kg	8.8 kg	Up to 110	110-115	62	64
VK-C	Large	6.8 kg	9.5 kg	Up to 115	115-120	65	67
VK-CXL	Extra-large	7.4 kg	10.3 kg	Up to 120	120-125	68	70
Article code	Sizes female	Weight 0.5/0.25 mm Leq	Weight 0.7/0.35 mm Leq	Chest size	Hip	Vest length	Kilt length
VK-D	Small	4.5 kg	6.3 kg	Up to 85	90-95	51	55
VK-E	Medium	4.8 kg	6.7 kg	Up to 95	95-100	55	58
VK-F	Large	5.4 kg	7.5 kg	Up to 105	100-110	56	61
VK-FXL	Extra-large	6.2 kg	8.6 kg	Up to 115	115-120	57	63

0.5mm Leq = 6kg/m2

Weight according to EN 61331-3:2014, this list can only be used as a reference, all listed weights have a tolerance of 7%. W/S = 0.25/0.25mm Leq = 3kg/m2 / W/S = 0.35/0.35mm Leq = 4,2kg/m2



Thyroid shields

3 different models, adjustable size using magnet closure

0.5 mm Leq

Article code	Weight*
VIPE	250 grams
TM-3	340 grams
TS-4	320 grams

Gonad Protection

0.5 mm Leq

Article code	Weight*
GP-S	600 grams
GP-M	1.0 kg
GP-L	1.4 kg
GP-Belt	750 grams

Head Protection

0.5 mm Leq

Article code	Weight*
HP-S	250 grams
HP-M	300 grams
HP-L/LX	350 grams
HP-XXL	400 grams

0.5mm Leq = 6kg/m2

Weight according to EN 61331-3:2014, this list can only be used as a reference, all listed weights have a tolerance of 7%. W/S = 0.25/0.25mm Leq = 3kg/m2 / W/S = 0.35/0.35mm Leq = 4,2kg/m2