

AB-0583-T

21003762-
İNG

02-21

TEST RESULTS

TEST METHOD: EN 13795-1:2019

SURGICAL CLOTHING AND DRAPES –REQUIREMENTS AND TEST METHODS

ANNEX 1: SURGICAL CLOTHING AND DRAPES;

TENSILE STRENGTH; EN 29073-3:1996

Instron 5969 (Load: 50 kN), Strip Method.
Speed: 100 mm/min $\pm$ 10, Gauge length 200 mm.
Pre-load was not applied. Without wetting samples.
The average results are given for weft and warp direction of five samples
Performed in the conditioned room (20 $\pm$ 2°C-65% $\pm$ 4).

Critical area

Dry ;

RESULT

Weft

67,4 N

Warp

142,4 N

REQUIREMENT

$\geq$ 20N (Dry)

$\geq$ 20N (Dry)

TENSILE STRENGTH; EN 29073-3:1996

Instron 5969 (Load: 50 kN), Strip Method.
Speed: 100 mm/min $\pm$ 10, Gauge length 200 mm.
Pre-load was not applied. Without wetting samples.
The average results are given for weft and warp direction of five samples
Performed in the conditioned room (20 $\pm$ 2°C-65% $\pm$ 4).

Less critical area

Dry ;

RESULT

Weft

50,2 N

Warp

88,1 N

REQUIREMENT

$\geq$ 20N (Dry)

$\geq$ 20N (Dry)

**ZENMED TIBBİ CİHAZLAR
VE SAĞLIK ÜRÜNLERİ SAN. VE TİC. LTD.ŞTİ.**
Mahmutbey Mah. 2567.Sk. No:1 İş Yeri No:11,21
Bağcılar / İSTANBUL / TÜRKİYE
Tel: 0212 445 45 25 Tic.Sic.No: 794704
Güneşli Y.D.: 997 067 1942