

DATE	28.03.2025
DOC.NO	FRM-KLT-BYM-069
PAGE NO	1
REV.NO	
REV.DATE	

TECHNICAL DATA SHEET

PRODUCT: Large surgical set for interventions with increased risk of fluids, disposable

Description of Product:

Large surgical set for interventions with increased risk of fluids, disposable			
1	1. Mayo dressing for instrument table 145x80cm (+/-10cm) – 1 pc., Waterproof, absorbent, absorbent area minimum 65x85cm (+/-5cm). Pictogram for guidance.	80x145 cm	3
2	2. Sheet with adhesive border and absorbent field 90x75cm (+/-5cm) – 2 pcs., Non-woven material minimum 4 layers, not less than 42 g/m2, Poly layered, waterproof, absorbent field area minimum 90x40cm (+/-5cm), with increased absorption level. Adhesive area applied directly to the edge of the dressing, impregnated in the fabric, without wrinkles, 100% adhesion capacity. Double-sided glue is not allowed. latex-free, rosin-free. Pictogram for guidance.	90x75 cm	2
3	3. Sterile absorbent cellulose wipes 20x25cm – 4 pcs.	20x25 cm	4
4	4. Surgical tape minimum 10x50cm (±5cm) – 1 pc. waterproof, absorbent, with glue on the entire surface, free of rosin. For firm fixation of tubes and instruments during surgery.	10x50 cm	1
5	5. Sheet with adhesive border and absorbent field 175x175cm (+/-10cm) – 1 pc., Non-woven material minimum 4 layers, not less than 42 g/m2, Poly laminate, waterproof, absorbent field area minimum 60x40cm, with increased absorption level. Adhesive area applied directly to the edge of the field, impregnated in the fabric, without creases, 100% adhesion capacity. Double-sided glue will not be allowed. latex-free, rosin-free. Pictogram for guidance.	175x175 cm	1
6	6. Sterile sheet with adhesive border and absorbent field 240x150cm (+/-10cm) – 1 pc., Non-woven material minimum 4 layers, not less than 42 g/m2, Poly layered, waterproof, absorbent, absorbent field area minimum 60x40cm, with increased absorption level. With "Tube Holders" (holes for fixing tubes and wires). Adhesive area applied directly to the edge of the field, impregnated in the fabric, without creases, 100% adhesion capacity. Double-sided glue will not be allowed. latex-free, rosin-free. Pictogram for guidance.	240x150 cm	1
7	7. Absorbent patient sheet 150x190cm (+/-10cm) – 1 pc., absorption area minimum 75x190cm (+/-5cm). 2 layers, waterproof, absorbent.	150x190 cm	1

Raw Materials: Triflex / PE + NW / Cellulose

Product Colour: Medical Blue

Reference Code: SP-02003-89

Weight in Grams: 72 gsm Triflex / 28 gsm(SS)+65 mic.(PE) / 60 gsm (Cellulose)


Package: Flat Pouch

Product: Medical Device Directive MDR 2017/745 ECC ANNEX XI A ALL VERSIONS. / Products and materials used do not contain metal.
Product Shelf Life - 5 Years

Product Materials

Unit / Size

PROPERTIES



Tolerances: +/- 2% cm

Measurement: cm

Package Information

The products in the sterilization bag are double-packaged to reduce all risks during transportation. Double packaged products are put into Baymed's standard sized carton; dimensions are as follows:
Height = 44 cm; Length = 40 cm ve Width = 60 cm.

Preparation Date

QUALITY CONTROL APPROVAL