

To: IP Oficiul de Gestionare a Programelor de Asistentă Externă

Tender ID: ocde-b3wdp1-MD-1672306960054

Zoetermeer, February 6th 2023

Statement

We, **CANON MEDICAL SYSTEMS EUROPE BV**, Zilverstraat 1, 2718 RP Zoetermeer, The Netherlands, as Authorized European Representative of **Canon Medical Systems Corporation**, 1385, Shimoishigami, Otawara-Shi, Tochigi Pref. 324, Japan,, who are official manufacturers of computer tomography system *Canon Aquilion One Genesis Edition (TSX-305A)*, having factories Canon Corporation, Medical System Division (Nasu Works), 1385, Shimoishigami, Otawara-Shi, Tochigi Pref. 324, Japan, with regard to bid submission via electronic tendering portal by our authorized distributor in Republic of Moldova - Medexcom-Teh S.R.L. do hereby, as supplement to the product data and brochures submitted on the portal would like to confirm the following technical parameters of our equipment offered – computed tomography system Canon Aquilion One Genesis Edition (TSX-305A):

- Total width of the detector, z-axis, mm: 160
- Reconstructed image slice thickness mm: 0.5, 1, 2, 3, 4, 5, 7, 8, 10 (selectable, maximum 10mm). Please see below the extract from the Canon tender questionnaire for Aquilion One Genesis Edition, Model TSX-305A, document nr. MSNCT0055EAB - which is a company confidential document and should be treated as such.
- Necessary quality assurance phantoms including the phantom holders to check the image quality and calibration of the CT scanner are included as standard with all Canon CT scanners.
- Expected x-ray tube life (warranty): 24 months without any limitations.
- All software applications on the post-processing and acquisition workstations are provided with a life-time licence.

Extract from the Canon tender questionnaire for Aquilion One Genesis Edition, Model TSX-305A,
document nr. MSNCT0055EAB



No. MSNCT0055EAB

15-14 Slice Sensitivity profile


15-14-1 State measurement method for Slice Sensitivity profile : Ramp in Catphan CTP-401

15-14-2 Slice Sensitivity profile at FWHM

Acq. slice thickness	Reconstructed slice thickness								
	0.5 mm	1.0 mm	2.0 mm	3.0 mm	4.0 mm	5.0 mm	7.0mm	8.0 mm	10.0 mm

CANON MEDICAL SYSTEMS EUROPE BV

Zilverstraat 1, 2718 RP Zoetermeer, The Netherlands



CANON MEDICAL SYSTEMS EUROPE B.V.
Zilverstraat 1, 2718 RP Zoetermeer
THE NETHERLANDS
TEL +31 (0)79 368 9222

P. Bishop

Power of Attorney