



Consistent. Automated. Brilliant.

Automated photo documentation of face and body

The FotoFinder ATBM master technology is a benchmark for Total Body Dermoscopy for the examination and follow-up of patients with multiple lesions. With the DX version for Dermaesthetics, the ATBM master can also be used in aesthetics, whether consultation, planning or photo documentation, whether facial, body or vein treatment.

ATBM
MASTER



DX



The automated imaging process with FotoFinder **Guided Photography** delivers consistent before-and-after photos in no time at all, impressively visualizing the treatment progress and success.

Brilliant images, personalized treatment plans and convincing before-and-after reports can increase patient satisfaction and compliance.

FotoFinder **ATBM master DX** is available as a stand-alone solution for aesthetics or as an extension of **ATBM master** for Total Body Dermoscopy.



Dr. Maurice Adatto

Dermatologic Surgeon,
Medical Director Skinpulse
Derm & Laser Centre, Geneva,
Switzerland.

Taking excellent quality pictures is mandatory when performing aesthetic treatments, especially in the field of body shaping. Patients need to see before-and-after pictures, as they often forget how they looked like before.

The FotoFinder system provides this excellence through quality pictures with very high definition, standardization of lighting and exposure, making the before-and-after comparison very accurate.

The system has a very user-friendly interface, making photography an easy task for everybody, also without having any specific knowledge.

ATBM master DX

**Automated photo documentation
of face and body**

Made in Germany

Simple and delegable

The **Guided Photography** Templates assist to quickly capture full face or body images. Choose from a large variety of photography templates that display step-by-step onscreen instructions. FotoFinder's virtual patient shows your staff treatment-specific patient poses and body area to be photographed. Every photo is automatically labeled and stored in your database for effortless retrieval.

Studio-quality photos

The computer-controlled **PolFlash XE** xenon flash enables both cross-polarized, reflection-free images and non-polarized photos with studio lighting. The integrated zoom motor automatically adjusts the focal length depending on which **Guided Photography** template is selected.

Reproducible images in the blink of an eye

The **Ghost feature** displays a transparent baseline photo on top of the live image, making it easy to align your patient at follow-ups. Onscreen instructions, consistent camera settings, **Laser Liner** and **FotoFinder Ghost** make follow ups from the right angle and same distance a cinch.



Truly flexible with the FotoFinder Universe iPad app

The FotoFinder **Universe iPad App** offers access to your images from anywhere in your office, for viewing, editing and comparison. Quick snapshots can be taken and are automatically archived in your FotoFinder database. Create treatment plans on your iPad during consultation and review the results with your clients to reinforce treatment success.

Personalized treatment plans

The **Proposer** is an impressive client consultation tool for effective, time-saving customization. Simply highlight treatment areas and drag and drop your recommended procedures on your client's photo. Add from your personalized product database your recommended skin care products and treatment protocols and present your proposal on a professional report for an exciting patient experience.

Compelling Marketing and Reports

FotoFinder's time saving editing tools allow you to quickly create impressive side-by-side comparisons for your online marketing activities. Choose from a variety of compelling print layouts, which feature your company logo and address for a polished and professional look. Improve your social media presence and get more referrals of satisfied patients.